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© 2011 Illumina, Inc. All rights reserved.
Illumina, illuminaDx, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, GoldenGate Indexing, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CPro, GenomeStudio, Genetic Energy, HiSeq, HiScan, TruSeq, Eco, MiSeq and Nextera are registered trademarks or trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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
The screenshot shows the DesignStudio interface for TruSeq Custom Amplicon. The main panel displays a table of targets with columns for Target, Chr:Start-Stop, Length, Amplicons, Avoid SNP, SNP Count, Coverage, Score, and Added. The table lists several targets, including CLCC1_Exon_202118240, CLCC1_Exon_202167440, CLCC1_Exon_202184940, CLCC1_Exon_202411140, TP53_Exon_188701140, and RFX4_Exon_190713840. The interface also shows a summary of 55 targets, 117 amplicons, and 0 gaps, with a cumulative target size of 18,004 bp and 88% coverage. A 'WARNINGS' section indicates 4 low scoring targets.

TARGET [?]	CHR:START-STOP	LENGTH	AMPLICONS [?]	AVOID SNP [?]	SNP COUNT [?]	COVERAGE [?]	SCORE [?]	ADDED
CLCC1_Exon_202118240	1 :109,484,027 - 109,484,120	94	1	Yes	0	100 %	95 %	09/30/2011
CLCC1_Exon_202167440	1 :109,490,233 - 109,490,340	108	1	Yes	0	100 %	95 %	09/30/2011
CLCC1_Exon_202184940	1 :109,505,983 - 109,506,111	129	1	Yes	0	100 %	95 %	09/30/2011
CLCC1_Exon_202411140	1 :109,482,258 - 109,482,404	147	1	Yes	0	100 %	95 %	09/30/2011
TP53_Exon_188701140	17 : 7,573,927 - 7,574,033	107	1	Yes	4	100 %	95 %	09/30/2011
TP53_Exon_188709140	17 : 7,590,695 - 7,590,863	169	1	Yes	0	100 %	95 %	09/30/2011
RFX4_Exon_190713840	3 :186,515,324 - 186,515,403	80	1	Yes	0	100 %	95 %	09/30/2011

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Myllumina (iCom)への登録手順

Create an Account

Personal Information * Required fields appear in bold

Salutation
Select a Value ▾

* First Name Middle Name * Last Name

* Job Title

Why register?

- Place orders and track order status
- View and redeem quotes
- Design and order Illumina custom products
- Maintain favorite lists
- Manage subscriptions
- Register for upcoming webinars and view technical bulletins

Your username will be the same as the email address you enter.
Please use your work email to expedite account verification.

Log In Information

* Email Address * Confirm Email

* Password * Confirm Password

Customer Information

Customer Number

Your customer number is the number Illumina uses to identify your billing information. You may have received this in a welcome packet after you placed your first order or during your training. If you do not know your customer number it can be found on your invoice, order confirmation, or packing slip. If you do not have a recent invoice or packing slip, contact your institution's purchasing department or Illumina's Customer Service Team.

Security Validation

* Type the two words:

fine 2つの単語を入力します: no CAPTCHA™ stop spam. read books.

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

illumina® TruSeq Custom Amplicon

Create Project Manage Targets Review Design

Links

- TruSeq Custom Amplicon User Guide
- UCSC Browser
- HUGO

Webinar

- TruSeq Custom Enrichment: Learn How to Customize Content with DesignStudio

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TARGET REGION [?]	CHR:START-STOP	TARGET [?]	AMPLICONS [?]	AVOID SNP [?]	ADDED (UTC) ↓
<input type="checkbox"/> test1	11 : 1,000,000 - 1,010,000	7 EX	0	Yes	02/22/2012

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No

Links

Edit Project Manage Targets Review Design

Manage Targets

+ NEW TARGET REGION

START DESIGN

Coordinate Gene

GENE NAME

KRAS 50 50

KRAS 12:25,358,180 - 25,403,854 + Bases Downstream

LABELS (SEPARATE MULTIPLE BY COMMA OR TABS, MAX 25 (0-9A-Z))

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FILTERS

All Page x DELETE SNP ON SNP OFF EXPORT

View: Regions | Targets

TARGET REGION [?]	CHR:START-STOP	TARGET [?]	AMPLICONS [?]	AVOID SNP [?]	ADDED (UTC)	INFO
<input type="checkbox"/> KRAS ±50	12 : 25,358,130 - 25,403,904	6 EX	0	Yes	02/22/2012	

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EST AMPLICONS: 62

START DESIGN


TARGET REGION [?]	CHR:START-STOP	TARGET [?]	AMPLICONS [?]	AVOID SNP [?]	ADDED (UTC) ↓
chr1_A	1 : 1,000,000 - 1,000,500	FR	0	Yes	02/16/2012
BRAF ±50	7 : 140,433,763 - 140,624,614	18 EX	0	Yes	02/16/2012
KRAS ±50	12 : 25,358,130 - 25,403,904	6 EX	0	Yes	02/16/2012

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Target Proximity Notification

Please note that targets within 250 bp of one another will be merged into one contiguous target to improve amplicon design. Do you want to Start Design now?

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


<input type="checkbox"/> All <input type="checkbox"/> Page <input checked="" type="checkbox"/> DELETE <input type="checkbox"/> SNP ON <input type="checkbox"/> SNP OFF <input type="button" value="EXPORT"/>		View: Regions Targets				
TARGET REGION [?]	CHR:START-STOP	TARGET [?]	AMPLICONS [?]	AVOID SNP [?]	ADDED (UTC) ↓	INFO
<input type="checkbox"/> KRAS ±50	12 : 25,358,130 - 25,403,904	6 EX	0	Yes	02/22/2012	D
<input type="checkbox"/> BRAF ±50	7 : 140,433,763 - 140,624,614	18 EX	0	Yes	02/22/2012	
<input checked="" type="checkbox"/> KRAS ±50	12 : 25,358,130 - 25,403,904	6 EX	0	Yes	02/22/2012	D

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 **Your Amplicon Design for "test" is Ready for Review**

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Congratulations!

Your project test has completed design and is ready for your review. Please follow this link to review the designs for your TruSeq Custom Amplicon Kit in DesignStudio:

<http://designstudio.illumina.com/truseqca/project/4531>

Project Summary
test

Number Of Targets: 5
Number Of Attempted Amplicons: 16
Number Of Gaps: 0
Total Gap Distance: 0
Cumulative Target (BP): 2264
Redundancy: 0
Coverage: 100 %

Sincerely,
The Illumina DesignStudio Team

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The screenshot shows the DesignStudio interface for a TruSeq Custom Amplicon project. The left sidebar contains a 'DESIGN SUMMARY' section with 'Orderable Amplicons: 117' circled in red. The main area displays a 'DESIGN WARNINGS' section with a '?' icon circled in red. Below this is a table of amplicons with columns for 'AMPLICON', 'STRAND', 'AVOID SNP', 'SNP COUNT', 'SCORE', and 'INFO'. A table callout provides detailed information for one amplicon, including coordinates, strand, and SNP data. Red callouts point to various elements: '設計できたアンプリコン' (Designable amplicons), '設計できなかったアンプリコン' (Non-designable amplicons), 'デザインのサマリー' (Design summary), 'LOW SCORING TARGETS: 4' (circled in red), 'ストランド情報' (Strand information), 'SNP情報' (SNP information), 'プローブスコア' (Probe score), 'ターゲット領域の情報' (Target region information), 'UCSCヘルク' (UCSC link), and 'DESIGN WARNINGS [?] の? をクリックすると' (Clicking the ? in DESIGN WARNINGS [?] will display target regions according to the type of notification).

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This displays target regions according to the type of notification.
[learn more ...](#)

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Region, Target, and Amplicon Details

DesignStudio enables you to look at high- and low-level target details. To view low-level details for a target, ensure that you are on the Manage Targets page, and then click the link for the target in the Target Region column of the table in the bottom half of the window (which contains the target regions you have added). When you are looking at a specific target region, you can drill down to a specific probe by clicking a specific target. Finally, when you are looking at individual probes, you can toggle the View to see either probes or gaps.

Table: Regions and Targets Table

Field	Description
ADDED (UTC)	The Universal Time Code indicating when the target was added
AMPLICONS	The estimated number of amplicons that will be designed for the region or target
AVOID SNP	Indicates whether SNPs are being avoided for the target region and/or amplicon
CHR: START-STOP	The chromosome, followed by the start position and then the stop position in bp
COVERAGE	The total number of non-overlapping bases covered by the attempted amplicons divided by the total number of bases in design
INFO	Lists error codes for low scoring and duplicate probe targets. Hover over an error code to see what it means. See the Error Codes section for details.
LENGTH	The total number of bases for each target where the software will design amplicons. The length will help estimate the number of amplicons required to cover the region. Length is calculated in the following manner: (Target Coordinate Stop - Target Coordinate Start) + 1 - Target Length.
SCORE	This score reflects the predicted success rate of amplicons given Illumina's evaluation of the amplicon probes.

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Reset Designed Items

Once your amplicon design is complete, you cannot return to Manage Targets without removing all existing amplicon designs. Do you want to proceed?

It is recommended that you export your design data prior to returning.

Save a copy of my current project state:

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The screenshot shows the Illumina website interface. At the top left is the Illumina logo. To the right, it says "Hello, Kensuke. (Not Kensuke?)". Below this are contact numbers "800.809.4566", "MyIllumina", and a "Tools" dropdown menu. A navigation bar contains links for "APPLICATIONS", "SYSTEMS", "SERVICES", "SCIENCE", "SUPPORT", and "COMPANY", along with a search bar. On the left side, there is a vertical menu with links for "Home", "Account", "Support Webinars", "Support Bulletins", and "Downloads". The main content area features a "Welcome Back, Kensuke" message and a link to "To download the DMAP decode file client click here."

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illumina®

Log in to get personalized account information. Quick Order View Cart (0)

800.809.4566 MyIllumina Tools

APPLICATIONS SYSTEMS SERVICES SCIENCE SUPPORT COMPANY

Support / Sequencing / Kits and Reagents / TruSeq Custom Amplicon

TruSeq Custom Amplicon Support

Illumina product support pages provide fast, convenient access to product information, requirements and compatibility, training, downloadable resources, and tips and techniques identified by other researchers and Illumina scientists. We are constantly building out these pages to add more resources, so if you don't see what you need right now, come back soon.

Catalog IDs: FC-130-1001, FC-130-1003, FC-130-1006, FC-130-1007

Download Reagent Prep Guide

What's New

- Experiment Manager User Guide 02/01/2012
- Illumina Experiment Manager 01/30/2012
- TruSeq Custom Amplicon Library Prep Reagents MSDS 12/20/2011

Support

- Contents & Storage
- Input Requirements
- Documentation & Literature
- Training
- Questions & Answers
- Best Practices

Features:

- Amplify hundreds of genomic targets in a single reaction and sequence up to 96 samples in a single MiSeq run.
- Generate up to 384 amplicons across 96 samples within a single plate, with less than 3 hours hands on time.
- Easily create and manage projects using Illumina's online DesignStudio software.
- Perform variant calling and analysis across all samples using simple on-instrument, automated analysis software.
- Get the convenience of a fully integrated DNA-to-data solution including online probe design and ordering, assay, sequencing, automated data analysis and offline software for reviewing results.

TruSeq Custom Amplicon Calculator

This convenient tool will help you calculate your supported sample multiplexing level based on the number of targeted regions and desired coverage. A customer login is required.

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Description	Document Type	File Info
Experiment Manager User Guide	User Guide	PDF (< 1 MB)
TruSeq Custom Amplicon Index Adapters MSDS	MSDS	PDF (170 KB)
TruSeq Custom Amplicon Library Prep Tracking Form	Lab Tracking Form	PDF (< 1 MB)
TruSeq Custom Amplicon Library Prep Reagents MSDS	MSDS	PDF (237 KB)
TruSeq Custom Amplicon Library Preparation Experienced User Card	Experienced User Card	PDF (< 1 MB)
TruSeq Custom Amplicon Library Preparation Guide	Sample Preparation Guide	PDF (2 MB)

Description	Document Type	File Info
TruSeq Custom Amplicon Data Sheet	Data Sheet	PDF (977 KB)

Documentation & Literature の
「TruSeq Custom Amplicon Library Preparation
Guide」

Training の
「Sample Preparation: TruSeq Custom Amplicon
(TSCA)」


Download Not Available

The file you requested either doesn't exist or you don't have permission to view it.

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

Equipment

Pre-PCR

Equipment	Supplier
37° Incubator	Forced Air Oven, VWR International or comparable
Heat Block, 96-well	SciGene, Hybex Microsample Incubator for PCR plate
Tabletop Centrifuge	General lab supplier (Plate centrifuge that attains designated speeds of protocol)

NOTE
Use a dedicated set of pipettes, pipette tips, vortexer, and centrifuge for pre-amplification steps.

Post-PCR

Equipment	Supplier
Magnetic Stand-96	DynaMag™-96 Side Skirted
Post-PCR Plate Shaker	Q Instruments, BioShake IQ, Part#1808-1806
Tabletop Centrifuge	General lab supplier (Plate centrifuge that attains designated speeds of protocol)
Gel Electrophoresis Supplies and Apparatus	General lab supplier
[Optional] Bioanalyzer System	Agilent Technologies
Heat Block for 1.5mL centrifuge tubes	General lab supplier

NOTE
Use a dedicated set of pipettes, pipette tips, vortexer, heat block, and centrifuge for post-amplification steps.

Additional Components

Consumable	Catalog #
TruSeq Index Plate Fixture Kit	FC-130-1005
TruSeq Custom Amplicon Filter Plate	FC-130-1006
TruSeq Index Plate Fixture and Collar Kit	FC-130-1007

User-Supplied Consumables

Quantity	Consumable	Supplier
As needed	10 N NaOH (prepare from tablets or use a standard solution)	General lab supplier
3	96-well skirted PCR plates, 0.2 ml, polypropylene	Bio-Rad, Part # MSP-9601
3	96-well storage plates, 0.8 ml (MIDI plates)	Fisher Scientific, Part # AB-0859
5 ml	Agencourt AMPure XP, 60 ml kit	Beckman Coulter, Part # A63881/A63880
3	Adhesive aluminum foil seal	Beckman Coulter, Part # 538619
As needed	Conical tubes, 15 ml	General lab supplier
2	Eppendorf microcentrifuge tubes (screw top recommended)	General lab supplier
40 ml	Ethanol, 200 proof for molecular biology	General lab supplier
1	Microseal 'A' adhesive seals	Bio-Rad, Part # MSA-5001
As needed	Microseal 'B' adhesive seals	Bio-Rad, Part # MSB-1001
2	PCR Eight-Tube Strips	General lab supplier
As needed	Solution basin, PVC, non-sterile (trough)	Labcor, Part# 730-001
As needed	Agarose gel electrophoresis apparatus	General Lab Supplier
As needed	Agarose gel (2% or 4%)	General Lab Supplier
As needed	DNA molecular weight markers	General Lab Supplier
As needed	Ice bucket	General Lab Supplier

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
SampleSheet

[Header]		
Investigator Name	Demo	
Project Name	Demo	
Experiment Name	Demo	
Date	5/31/2011	
Workflow	Resequencing	
[Reads]		
151		
151		
[Manifests]		
A	ManifestNameHere	
[Settings]		
FilterPCRDuplicates	0	
[Data]		
Sample_ID	Sample_Name	GenomeFolder
PhiX	PhiX	PhiX

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