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## 16S Metagenomic Sequencing Library Preparation

*Preparing 16S Ribosomal RNA Gene Amplicons for the Illumina MiSeq System*

Introduction	2
16S Library Preparation Workflow	5
Amplicon PCR	6
PCR Clean-Up	8
Index PCR	10
PCR Clean-Up 2	13
[Optional] Validate Library	15
Library Quantification, Normalization, and Pooling	16
Library Denaturing and MiSeq Sample Loading	17
MiSeq Reporter Metagenomics Workflow	20
Supporting Information	21

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These represent only some of the acceptable combinations. Alternatively, please check the real sequences of each index in the table above to make sure each base position will have signal in both color channels for the index read:

Good				Bad			
Index 1		Index 2		Index 1		Index 2	
705	GGACTCCT	503	TATCCTCT	705	GGACTCCT	502	CTCTCTAT
706	TAGGCATG	503	TATCCTCT	706	TAGGCATG	502	CTCTCTAT
701	TAAGCGA	504	AGAGTAGA	701	TAAGCGA	503	TATCCTCT
702	CGTACTAG	504	AGAGTAGA	702	CGTACTAG	503	TATCCTCT
	√√√√√√√√		√√√√√√√√		√√√√√√√√		√√√√xxxx

√=signal in both color  
x=signal missing in one color channel

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**Table 5** Libraries Pooled: 6 or fewer; Sequencing Workflow: Single Index

Plex	Index 1 (i7) Selection	Index 2 (i5) Selection
1-plex (no pooling)	Any Index 1 adapter	Any Index 2 adapter
2-plex	<ul style="list-style-type: none"><li>• [option 1] N702 and N701</li><li>• [option 2] N702 and N704</li></ul>	
3-plex	<ul style="list-style-type: none"><li>• [option 1] N701, N702, and N704</li><li>• [option 2] N703, N705, and N706</li></ul>	
4- or 5-plex	<ul style="list-style-type: none"><li>• [option 1] N701, N702, N704, and any other Index 1 adapter</li><li>• [option 2] N703, N705, N706, and any other Index 1 adapter</li></ul>	
6-plex	N701, N702, N703, N704, N705, and N706	

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The screenshot displays the Illumina Sequencing Analysis Viewer 1.7.32 interface. The 'Analysis' tab is selected and highlighted with a red box. The 'Run Folder' field is empty, with 'Browse' and 'Refresh' buttons to its right. Below the 'Analysis' tab, the 'Status' section shows 'Extracted: 0', 'Called: 0', and 'Scored: 0'. The main interface is divided into four panels:

- Flowcell Chart:** A vertical bar chart showing intensity across different lanes. A color scale on the right ranges from 0.20 (blue) to 0.80 (yellow). Controls for 'Intensity', 'Top Surface', 'Cycle 1', 'Base A', and 'Fix Scale' are on the left.
- Data By Cycle:** A line graph showing intensity over cycles. A dropdown menu is open, with '% Base' selected and highlighted by a red box. A yellow box with the text '% Base' is overlaid on the graph. The x-axis ranges from 0 to 100, and the y-axis from 0 to 80.
- QScore Distribution:** A line graph showing the distribution of Q-scores. The x-axis ranges from 0 to 40, and the y-axis from 0 to 80. Controls for 'All Lanes', 'Both Surfaces', 'All Reads', and 'All Cycles' are on the left.
- Data By Lane:** A line graph showing data across lanes. The x-axis is labeled 1 through 8, and the y-axis from 0 to 80. Controls for 'Clusters' and 'Both Surfaces' are on the left.
- Qscore Heatmap:** A heatmap showing Q-scores across lanes. Controls for 'All Lanes' and 'Both Surfaces' are on the left.

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The image illustrates how to identify low diversity samples in sequencing data. It compares high diversity (stable base percentages) with low diversity (highly fluctuating base percentages) across different lanes and cycles. The central panel highlights the '% Base' metric in the 'Data By Cycle' view, which is used to detect these fluctuations.

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