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© 2013 Illumina, Inc. All rights reserved.
Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners. Illumina, Inc. and the Genetic Energy streaming bases design are trademarks or Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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The screenshot displays the BaseSpace web interface. At the top, there is a navigation bar with the BaseSpace logo, several icons (bar chart, person, folder, folder, microscope, building, question mark), a search icon, a user profile icon, and the Illumina logo. Below the navigation bar, the main content area is divided into several sections:

- NOTIFICATIONS:** A list of notifications, including "BaseSpace Updates: NextSeq 500 and More" (6 hours ago), "Tiffany Malek invited you to the Project Lasso_test" (Jan 10, 2014, Accepted), and "Wishing you Happy Holidays with a Festive Competition!" (Dec 21, 2013).
- BaseSpace iCredits:** A section showing 0 iCredits and an "Add More" button.
- TIP OF THE DAY:** A tip stating "Did you know that you can build native apps on this environment" with a "Learn more" link.
- RUNS:** A list of completed runs: "PhixRun2", "PhiX_Run1", and "RNA test run", each with a green progress bar labeled "Complete".
- PROJECTS:** A list of projects: "Tattooine", "BaseSpace_Demo_Project", "Our Run", and "Lasso_test".
- ANALYSES:** A list of analyses: "BWA Whole Geno...", "BWA Enrichment...", "Illumina's Uploader", and "Illumina's Uploader".

At the bottom of the page, there is a footer with links for "Help", "FAQ", "Developers", "Terms", and "Blog". A vertical "contact us" button is located on the right side of the interface.

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The image shows the BaseSpace navigation bar with six icons: Dashboard, Prep, Runs, Projects, Apps, and Public Data. Each icon is accompanied by a callout box with a description in Japanese. The callout boxes are color-coded: Dashboard (blue), Prep (teal), Runs (green), Projects (light green), Apps (medium green), and Public Data (light green). The background shows a blurred view of the BaseSpace interface with a search bar and the Illumina logo.

BaseSpace Dashboard Prep Runs Projects Apps Public Data ? illumina

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BaseSpace® Dashboard Prep Runs Projects Apps **Public Data** Help ? Eri Kibukawa illumina®

> [HiSeq 2500: TruSeq Stranded mRNA LT \(SEQC: UHR & Brain\)](#)
RNA-Seq

> [HiSeq 2000: TruSeq Stranded Total RNA \(MAQC\)](#)
RNA-Seq Differential Expression

▼ [NextSeq 500: RNA-Seq \(8plex\)](#)
NextSeq 500 data generated for reference human brain RNA and universal human reference RNA (UHRR). Libraries were prepared using TruSeq stranded mRNA or TruSeq Total RNA with Ribo-Zero Gold reagent kits.

Run	NextSeq 500: RNA-Seq (8plex) (49.41 GB)	Import
Project	NextSeq 500: RNA-Seq (8plex) 0	Import

Research Areas

- Cancer Research
- Genetic Disease
- Complex Disease
- Microbial Research

Categories

- Exome
- Resequencing
- Small RNA
- Targeted Sequencing
- RNA-Seq**
- De Novo Assembly
- Chr-Seq
- Gene Fusion Detection
- Methyl-Seq
- Metagenomics
- Tumor Normal
- Variant Analysis
- Differential Expression
- Quality

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The screenshot shows the 'SHARING SETTINGS' dialog for a project named '2X151 Rhodobacter Resequencing'. Under the 'INVITE A COLLABORATOR' section, there is a text input field for 'Email address' and a larger text area for 'Add a personal message (optional)'. A green 'Invite' button is located below these fields. In the 'COLLABORATORS (1)' section, the email 'esmith@illumina.com' is listed with a 'Read Only' permission level. A yellow box highlights a dropdown arrow and a close 'x' icon next to the permission level. At the bottom, there are 'Cancel' and 'Save Settings' buttons.

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The screenshot shows the 'Sharing by link' dialog for the same project. It prompts the user to 'Copy the below URL and send it to anyone you want to share with.' A blue highlighted text box contains the URL 'https://basespace.illumina.com/s/lhVA72PqVhPL'. A red 'Deactivate' button is located at the bottom right of the dialog.

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Genomatix GePS
RNA-Seq Pathway Analysis

Omixon HLA Typing
HLA Typing

Open-Source PicardSpace
Alignment QC

Broad IGV
Genome Visualization

GeneTalk Variant
Analyzer
GENETALK GMBH

iPathwayGuide
ADVAITA BIO

Elsevier Genomics Data
Publishing

PathSEQ Virome
PATHGEN DX PTE
LTD

DNASTAR SeqMan Ngen
Assembly

DeepChek®-HIV
ABL S.A.

Strand Elastic Genome Browser
Advanced Genome Visualization

Biomatters Molecular Profiler
Polyomic Analysis

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

TOPHAT ALIGNMENT BETA

Illumina

App Session Name: TopHat Alignment Beta 01/15/2014 11:41:09

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
App Session Name: Cufflinks Assembly & DE Beta 01

Save Results To: ①いち、

TopHat Alignments Selection Criteria

Reference Genome: Homo sapiens/hg19 (RefSeq)

Stranded:

Options

Novel Transcript Assembly:

Control Group

Group Label: control

TopHat Alignment App Result (s):

Adjust transcript assembly for samples without polyA selection:

Comparison Group

Group Label: comp

TopHat Alignment App Result (s):


Adjust transcript assembly for samples without polyA selection:

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



- Targeted Regions:
- Nextera Rapid Capture Exome v1.1
 - Nextera Rapid Capture Exome v1.1**
 - Nextera Rapid Capture Exome v1.2
 - Nextera Rapid Capture Expanded Exome v1.0
 - TruSight Autism v1.0
 - TruSight Cancer v1.0
 - TruSight Cardiomyopathy v1.0
 - TruSight Exome v1.0
 - TruSight Inherited Disease v1.0
 - TruSight One v1.0

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 PDF Summary Report

SAMPLE INFORMATION ⁱ

Total PF Reads	Percent Q30	Percent Duplicate Paired Reads	Fragment Length Median	Fragment Length Standard Deviation
65,262,238	82.3%	6.5%	227 bp	80 bp

ENRICHMENT SUMMARY ⁱ

Total Length of Targeted Reference	Padding Size
44,296,481 bp	150 bp

Read Level Enrichment ⁱ

Total Aligned Reads	Percent Aligned Reads	Target Aligned Reads	Read Enrichment	Padded Target Aligned Reads	Padded Read Enrichment
62,912,089	96.4%	39,472,997	62.7%	42,236,377	67.1%

Base Level Enrichment ⁱ

Total Aligned Bases	Target Aligned Bases	Bases Enrichment	Padded Target Aligned Bases	Padded Base Enrichment
8,778,180,848	4,378,065,923	49.9%	5,837,362,539	66.5%

SMALL VARIANTS SUMMARY ⁱ

	SNVs	Insertions	Deletions
Total Passing	31,080	989	1,185
Percent Found in dbSNP	99.0%	88.0%	78.3%

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Consequence

Show only variants that are:

- Missense
- Frameshift
- Stop gained
- Stop lost
- Initiator codon
- In-frame insertion
- In-frame deletion
- Splice

Polyphen "damaging"

SIFT "deleterious"

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Gene Filter

Disease

Include List

Exclude List

Min Variant Alleles

Show only variants:

- Inside genes
- Inside coding regions
- In conserved regions

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Variant Frequency

- Global Frequency <
- American Pop Frequency <
- Asian Pop Frequency <
- African Pop Frequency <
- European Pop Frequency <
- EVS Frequency <

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Next is Now

Thank You.

© 2014 Illumina, Inc. All rights reserved.
Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

illumina®