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**OPEN** Long-range PCR in next-generation sequencing: comparison of six enzymes and evaluation on the MiSeq sequencer

SUBJECT AREAS:  
NEXT-GENERATION SEQUENCING  
BIOINFORMATICS

Haiying Jia<sup>1,2</sup>, Yunfei Guo<sup>2,3</sup>, Weiwei Zhao<sup>4</sup> & Kai Wang<sup>2,3,5</sup>

Received  
14 May 2014  
Accepted  
24 June 2014  
Published  
18 July 2014


<sup>1</sup>The First Affiliated Hospital, Jinan University, Guangzhou, Guangdong 510632, China, <sup>2</sup>Zilkha Neurogenetic Institute, Keck School of Medicine, University of Southern California, Los Angeles, CA 90089, USA, <sup>3</sup>Department of Preventive Medicine, Keck School of Medicine, University of Southern California, Los Angeles, CA 90089, USA, <sup>4</sup>Kingmed Diagnostics, Guangzhou, Guangdong 510330, China, <sup>5</sup>Department of Psychiatry, Keck School of Medicine, University of Southern California, Los Angeles, CA 90089, USA.

Correspondence and requests for materials should be addressed to K.W. (kaiwang@usc.edu)

Long-range PCR remains a flexible, fast, efficient and cost-effective choice for sequencing candidate genomic regions in a small number of samples, especially when combined with next-generation sequencing (NGS) platforms. Several long-range DNA polymerases are advertised as being able to amplify up to 15 kb or longer genomic DNA. However, their real-world performance characteristics and their suitability for NGS remain unclear. We evaluated six long-range DNA polymerases (Invitrogen SequelPrep, Invitrogen AccuPrime, TaKaRa PrimeSTAR GXL, TaKaRa LA Taq Hot Start, KAPA Long Range HotStart and QIAGEN LongRange PCR Polymerase) to amplify three amplicons, with sizes of 12.9 kb, 9.7 kb, and 5.8 kb, respectively. Subsequently, we used the PrimeSTAR enzyme to amplify entire BRCA1 (83.2 kb) and BRCA2 (84.2 kb) genes from nine subjects and sequenced them on an Illumina MiSeq sequencer. We found that the TaKaRa PrimeSTAR GXL DNA polymerase can amplify almost all amplicons with different sizes and Tm values under identical PCR conditions. Other enzymes require alteration of PCR conditions to obtain optimal performance. From the MiSeq run, we identified multiple intronic and exonic single-nucleotide variations (SNVs), including one mutation (c.5946delT in BRCA2) in a positive control. Our study provided useful results for sequencing research focused on large genomic regions.

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# Nexteraの原理

参考文献 Rapid, low-input, low-bias construction of shotgun fragment libraries by high-density in vitro transposition. Genome Biology 2010 Dec 8;11(12):R119.

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#	Sample ID	Sample Name	Chr	Position	Score	Variant Type	Call	Frequency	Depth	Filter
1	1471-8	NA18507	chr7	117235055	3070	SNP	T->T			

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Application Note: Sequencing illumina®

## Analysis of Skeletal Remains Using Deep mtDNA Sequencing

MISeq® system aids in identification of skeletal remains by sequencing mtDNA hypervariable regions.

**Introduction**

Forensic researchers face significant challenges in identifying unknown remains found in mass graves, such as those located in the war-torn countries of Vietnam, Kosovo, and Somalia. Dental records and circumstantial evidence are often inconclusive. DNA collected from skeletal samples is typically highly degraded, contaminated, and present in small quantities. The poor quality of these samples limits the success of conventional short tandem repeat (STR) typing techniques. In contrast, mitochondrial DNA (mtDNA) sequencing is quite useful in analyzing these samples, offering the following advantages:

- It is abundant and present in hundreds of copies per cell compared to nuclear DNA.
- It has a low mutation rate.
- It is maternally inherited and does not undergo sexual recombination, enabling maternal relationships to be established across generations.

For years, forensic laboratories have used capillary electrophoresis (CE)-based Sanger sequencing to perform mtDNA analysis. While the full mtDNA genome is approximately 16 kb, the hypervariable I and II (HVI and HVII) regions are usually the targets of mtDNA sequencing for human identification. These regions are located within a 1.1 kb region called the "D-loop". Capable of providing deeper coverage of specific regions of interest, next-generation sequencing (NGS) systems enable mtDNA sequencing to be performed at a high throughput, using simpler workflows.

The Illumina MISeq system is a fully integrated NGS instrument, providing a rapid workflow and streamlined sample preparation that is ideal for sequencing small genomes and amplicon samples, such as the D-loop region of mtDNA. This application note describes the use of the MISeq system for deep sequencing of the HVI and HVII regions in mtDNA and how this technique is being used to aid the identification of human skeletal remains found in mass graves in Vietnam (Figure 1).

**mtDNA Sequencing of Bone Fragments**


Forensic scientists used the MISeq system to sequence bone fragments uncovered at a mass grave site in Vietnam after DNA analysis of the samples using CE-based sequencing failed. The process began with the extraction of total DNA from blood and bone samples using established methods. Four individual amplicons spanning the HVI and HVII regions were amplified from mtDNA using PCR primers specific to the D-loop/Control region<sup>12</sup> (Figure 2). Library preparation was performed using the Nextera® XT DNA Sample Preparation Kit. The samples were then sequenced on the MISeq system using

**Figure 1: Recovery of Human Remains**

MISeq technology is helping to identify human skeletal remains such as these, recovered in Vietnam. There are an estimated 300,000 unidentified human skeletal remains lying in mass graves across Vietnam.

**Figure 2: NGS mtDNA Sequencing Workflow**



Specific amplicons spanning the HVI and HVII regions were sequenced on the MISeq system using a 2 x 150 bp run protocol.

[http://applications.illumina.com/content/dam/illumina-marketing/documents/products/appnotes/appnote\\_vietnam\\_mtDNA.pdf](http://applications.illumina.com/content/dam/illumina-marketing/documents/products/appnotes/appnote_vietnam_mtDNA.pdf)

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- 1) mtDNA and its role in ancestry. *Genebase Tutorials*. Retrieved October 1, 2014, from <http://www.genebase.com/learning/article/17>
- 2) *Journal of Human Genetics* **56**, 689-694

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## Human mtDNA Genome For the Illumina Sequencing Platform

Introduction	2
Workflow	3
Sample Sheet	4
First PCR Amplification	6
Quantity	9
Tagment Amplified DNA	11
Limited-Cycle PCR Amplification	13
Limited-Cycle PCR Clean-Up	17
Bead-based Normalization	19
Pool and Load from Bead-based Normalization	23
Bioanalyzer-based Normalization	25
Pool and Load from Bioanalyzer-based Normalization	27
Analysis	29
Supporting Information	36


**IMPORTANT NOTICE** This document provides information for an application for Illumina technology that has been demonstrated internally and may be of interest to customers. This information is provided as-is and is not an Illumina product and is not accompanied by any rights or warranties. Customers using or adapting this information should obtain any licenses required and materials from authorized vendors. Illumina products mentioned herein are for research use only unless marked otherwise. While customer feedback is welcomed, this application is not supported by Illumina Technical Support and Field Application Scientists.

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<sup>1</sup>Stawski, H., B. J. Bintz, E. S. Burnside, and M. Wilson. 2013. Preparing Whole Genome Human Mitochondrial DNA Libraries for Next Generation Sequencing (NGS) Using Illumina Nextera XT. Poster presentation at the 65th Annual American Academy of Forensic Sciences Conference. In: Proceedings of the American Academy of Forensic Sciences. Washington, D.C. [www.aafs.org/sites/default/files/pdf/ProceedingsWashingtonDC2013.pdf](http://www.aafs.org/sites/default/files/pdf/ProceedingsWashingtonDC2013.pdf)

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## Human mtDNA D-loop Hypervariable Region

For the Illumina Sequencing Platform

Introduction	2
Workflow	3
Sample Sheet	4
First PCR Amplification	6
First PCR Clean Up	10
Quantify	12
Tagment Amplified DNA	14
Limited-Cycle PCR Amplification	16
Limited-Cycle PCR Clean-Up	20
Bead-based Normalization	23
Pool and Load from Bead-based Normalization	27
Bioanalyzer-based Normalization	29
Pool and Load from Bioanalyzer-based Normalization	31
Analysis	33
Supporting Information	40

**IMPORTANT NOTICE** This document provides information for an application for Illumina technology that has been demonstrated internally and may be of interest to customers. This information is provided as-is and is not an Illumina product and is not accompanied by any rights or warranties. Customers using or adapting this information should obtain any licenses required and materials from authorized vendors. Illumina products mentioned herein are for research use only unless marked otherwise. While customer feedback is welcomed, this application is not supported by Illumina Technical Support and Field Application Scientists.

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