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© 2013 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners. The Genetic Energy streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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
http://support.illumina.com/sequencing/sequencing_software/experiment_manager.html

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| | | | | | | | | |
|-------------------|-----------------------------------|-------------|--------------|-------------|-------------|----------|----------------|-------------|
| [Header] | | | | | | | | |
| IEMFileVersion | | | | | | | 4 | |
| Investigator Name | Isabelle | | | | | | | |
| Experiment Name | HiSeq X ten run | | | | | | | |
| Date | | | | 3/14/2014 | | | | |
| Workflow | GenerateFASTQ | | | | | | | |
| Application | HiSeq FASTQ Only | | | | | | | |
| Assay | TruSeq HT | | | | | | | |
| Description | none | | | | | | | |
| Chemistry | Default | | | | | | | |
| [Reads] | | | | | | | | |
| | 151 | | | | | | | |
| | 151 | | | | | | | |
| [Settings] | | | | | | | | |
| ReverseComplement | | | | | | | 0 | |
| Adapter | AGATCGGAAGAGCACACGTCTGAACTCCAGTCA | | | | | | | |
| AdapterRead2 | AGATCGGAAGAGCGTCGTGTAGGGAAAGAGTGT | | | | | | | |
| [Data] | | | | | | | | |
| Lane | Sample_ID | Sample_Name | Sample_Plate | Sample_Well | I7_Index_ID | index | Sample_Project | Description |
| | 1 sample_ID1 | test1 | | | D701 | ATTACTCG | Hxten | |
| | 2 sample_ID2 | test2 | | | D712 | AGCGATAG | Hxten | |
| | 3 sample_ID3 | test3 | | | D710 | TCCGCGAA | Hxten | |
| | 4 sample_ID4 | test4 | | | D708 | TAATGCGC | Hxten | |

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- ▶ http://support.illumina.com/sequencing/sequencing_instruments/nextseq-500.ilmn

illumina[®] Log in to get personalized account information. Quick Order View Cart

Contact Us MyIllumina Tools

APPLICATIONS SYSTEMS INFORMATICS CLINICAL SERVICES SCIENCE SUPPORT Search

COMPANY

Support » Sequencing » Sequencing Instruments » NextSeq 500 | Follow us:

NextSeq 500 Support

- Overview
- Site Prep/Lab Environment
- Requirements & Compatibility
- Supported Kits
- Downloads
- Documentation & Literature
- Training
- Questions & Answers
- Services & Warranties
- Bulletins
- Webinars
- Product Ordering Information

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NextSeq 500

Latest Updates

- [Third-Party Analysis Software and Utilities Tech Note](#) 06/16/2014
- [Illumina Two-Channel SBS Sequencing Technology](#) 03/08/2014
- [NextSeq 500 Sequencing System Brochure](#) 01/12/2014

User Guides

- [NextSeq 500 System User Guide \(15046563 D\)](#)

Current Consumables

NextSeq 500 High Output Kit—Available in three sizes (300 cycles, 150 cycles, and 75 cycles), the NextSeq 500 High Output Kit includes one high-output flow cell, the high-output reagent cartridge, and the buffer cartridge.

NextSeq 500 Mid Output Kit—Available in two sizes (300 cycles and 150 cycles), the NextSeq 500 Mid Output Kit includes one mid-output flow cell, the mid-output reagent cartridge, and the buffer cartridge.

Current Software

NextSeq Control Software (NCS) v1.2—The control software interface guides you through the steps to load the flow cell and reagents before beginning a sequencing run. During the run, NCS

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- ▶ http://support.illumina.com/sequencing/sequencing_instruments/nextseq-500/downloads.ilmn

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- ▶ http://supportres.illumina.com/documents/documentation/software_documentation/bcl2fastq/bcl2fastq2-user-guide-15051736-b.pdf

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Thank You!

© 2014 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

illumina[®]