# TruSight One

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© 2014 Illumina, Inc. All rights reserved. Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, GAIx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDX, MiSeqFGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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
The screenshot shows the BaseSpace Public Data interface. The top navigation bar includes the BaseSpace logo and several icons, with a red arrow pointing to the public data icon. Below the navigation bar, there is a search bar and a list of data categories. The 'TruSight One Sequencing Panel MiSeq Trio Data' section is highlighted with a red box, and a red arrow points to the 'Import' button for the 'Project' row in the table below it.

Public Data

Search

> [Nextera Rapid Capture Exome Data Sets](#)

Targeted Sequencing Exome

> [MiSeq Cancer Demo](#)

Targeted Sequencing Cancer Research

▼ [TruSight One Sequencing Panel MiSeq Trio Data](#)

This example dataset is a family trio sequenced on a MiSeq using a MiSeq Reagent Kit, v3. This trio was obtained from the Coriell Institute's NIGMS Human Genetic Cell Repository. These samples are from the CEPH/Utah Pedigree 1463. [View accompanying blog post](#)

Targeted Sequencing

Run	TruSight One Sequencing Panel Dataset (11.52 GB)	Import
Project	TruSight One Sequencing Panel Dataset ()	Import

> [HiSeq 2500: Nextera Rapid Capture Exome \(CEPH Trio\)](#)

Resequencing Exome Targeted Sequencing

> [MiSeq: TruSeq Amplicon \(96plex\)](#)

Targeted Sequencing

Research Areas

- Cancer Research
- Genetic Disease
- Complex Disease
- Microbial Research

Categories


- Exome
- Resequencing
- Small RNA
- Targeted Sequencing

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


All Posts Segments

Search

**Blog @ Illumina**  
Real scientists. Real commentary.

### A Bird With No Wings?

Sep 16, 2014

The kiwi bird, native to New Zealand, doesn't have wings and it doesn't have a tail. With nostrils at the end of its beak, it has a great sense of smell, which makes up for its poor sense of sight. ...

[Read more...](#)


[Subscribe to feed](#)

[Subscribe to Illumina's Blog by Email](#)


[Tweets by @illumina](#)

### Looking Behind the Curtain at What Controls How You Think and Move

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
VCF Import Options

Please select which variants you would like to import from your VCF file(s)

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Variants in exons

Padding  bases

Variants in genes specified by gene list

Browse...

Variants in regions specified by BED file

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
The screenshot shows the Illumina VariantStudio interface. The 'Current Sample' dropdown menu is set to 'One\_12880', indicated by a red arrow. The 'Gene View' table displays the following data:

Gene	Variant	Count	Position
AGRN	C>C/T	1	976580
AGRN	T>C/C	1	977330
AGRN	A>G/G	1	981931
AGRN	T>C/C	1	982994
AGRN	T>C/C	1	984302
AGRN	G>G/A	1	986732
AGRN	C>T/T	1	987200
AGRN	C>T/T	1	990280

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Population Frequency  
Cross Sample Subtraction  
Family Based  
Custom  
Classification

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

### Autosomal Recessive Transmission Workflow

There are two possibilities for recessive transmission.

- 1 A single gene contains a variant that is:
  - Heterozygous (0/1) in the mother
  - Heterozygous (0/1) in the father
  - Homozygous (1/1) in the affected children

Figure 27 Autosomal Recessive Transmission Logic #1

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



**ClinVar**  
ClinVar aggregates information about sequence variation and its relationship to human health.

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Variant

Pass Filter

- Quality > 0
- Read Depth > 0
- Alt Variant Freq > 0

Show only variants:

- Inside genes
- In conserved regions
- Only variants without dbSNP ID
- Only variants with Cosmic annotation
  - where matches mutant allele
  - where not matches mutant allele
- Only variants with ClinVar annotation
  - where matches mutant allele
  - where not matches mutant allele

Apply Filters ==>

Gene	Variant	Chr	Coordinate	Classification	Type	Genot...	Exonic	Filters	Quality	GQX	Inherited From	Alt Var Freq
N/A	A>T/T	1	324822	...	snv	hom	no	PASS	2757.34	99	father	ε
N/A	A>T/T	1	187610852	...	snv	hom	no	PASS	3138.56	99	both	
N/A	A>G/G	1	187611003	...	snv	hom	no	PASS	2335.34	99	both	
N/A	T>G/G	1	202956340	...	snv	hom	no	PASS	2735.14	99	both	
N/A	T>C/C	1	219706428	...	snv	hom	no	PASS	1436.15	99	both	
N/A	A>G/G	2	3488502	...	snv	hom	no	PASS	1097.04	81	both	
N/A	G>C/C	2	10959475	...	snv	hom	no	PASS	1067.83	81	both	
N/A	A>G/G	2	131672950	...	snv	hom	no	PASS	2419.38	99	both	
N/A	C>T/T	2	131673316	...	snv	hom	no	PASS	431.14	36	both	
N/A	T>C/C	2	132229720	...	snv	hom	no	PASS	984.72	84	both	
N/A	C>G/G	3	136739406	...	snv	hom	no	PASS	588.87	45	both	
N/A	C>T/T	3	195467853	...	snv	hom	no	PASS	856.13	66	both	
N/A	A>G/G	3	5091384	...	snv	hom	no	PASS	1481.95	99	both	
N/A	T>C/C	4	5991384	...	snv	hom	no	PASS	5862.59	99	both	ε
N/A	T>C/C	4	5991384	...	snv	hom	no	PASS	1753.22	99	both	
N/A	A>G/G	9	10959475	...	snv	hom	no	PASS	1097.04	81	both	ε
N/A	A>G/G	9	10959475	...	snv	hom	no	PASS	1097.04	81	both	
N/A	A>G/G	9	10959475	...	snv	hom	no	PASS	1097.04	81	both	
N/A	G>A/A	6	32223632	...	snv	hom	no	PASS	1407.03	99	both	

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ASN.1 terms	ClinVar and VCF
0 – unknown	Uncertain significance
1 – untested	not provided (includes the cases where data are not available or unknown)
2 - non-pathogenic	Benign
3 - probable-non-pathogenic	Likely benign
4 - probable-pathogenic	Likely pathogenic
5 – pathogenic	Pathogenic
6 - drug-response	drug response
7 – histocompatibility	histocompatibility
255 - other	other
	confers sensitivity
	risk factor
	association
	protective

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Gene View

Gene	Variant	Chr	Coordinate	ClinVar Allele Type	ClinVar Significance	ClinVar Disease Name
					pathogenic	
A2M	T>C/C	12	9232268	germline	other non-patho...	ALPHA-2-MACROGLOBULIN_POLYMORPHI.
ASPM	G>G/T	1	197070442		non-pathogenic	Primary_autosomal_recessive_microcephal.
ASPM	T>T/C	1	197070697		non-pathogenic	Primary_autosomal_recessive_microcephal.
ASPM	G>G/A	1	197070707		non-pathogenic	Primary_autosomal_recessive_microcephal.
ASPM	C>C/T	1	197070706		non-pathogenic	Primary_autosomal_recessive_microcephal.
ASPM	T>T/C	1	197072420		non-pathogenic	Primary_autosomal_recessive_microcephal.
ASPM	T>T/C	1	197072420		non-pathogenic	Primary_autosomal_recessive_microcephal.
ASPM	A>A/G	1	197097700		non-pathogenic	Primary_autosomal_recessive_microcephal.
ATP7B	C>C/T	13	52508979	germline	non-pathogenic	Wilson's_disease
ATP7B	G>G/A	13	52511606	germline	non-pathogenic	Wilson's_disease
ATP7B	A>A/G	13	52515354	germline	probable-non-pat...	Wilson's_disease
ATP7B	C>C/T	13	52520507	germline	non-pathogenic	Wilson's_disease

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The image displays two overlapping windows of the 'フィルター エディタ' (Filter Editor) application. The top window shows a filter rule: 'And' followed by '[ClinVar Significance] Contains pathogenic'. A red arrow points to the word 'Contains'. The bottom window shows the same filter rule, but the operator has been changed to 'Equals', and the word 'Equals' is highlighted with a dashed border. A red box highlights the '適用(A)' (Apply) button at the bottom right of the bottom window. A red callout box with white text and a red arrow points to the 'Contains' operator in the top window, containing the text: '“Contains”をクリックして “Equals”に変更'.

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And

[ClinVar Significance] Contains pathogenic

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And

[ClinVar Significance] Equals pathogenic

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Gene ▲	Variant	Chr	Coordinate	ClinVar Significance <sup>📶</sup>	ClinVar Disease Name	Inherited From
DHODH	G>G/A	16	72050942	pathogenic	Miller_syndrome	father
DHODH	C>C/T	16	72057435	pathogenic	Miller_syndrome	mother

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The screenshot displays a software interface for filtering genomic data. On the left, there is a panel titled "Unaffected Siblings" with an empty text box and an "Edit Unaffected Siblings" button. Below this is a checkbox labeled "Use only passing variants in relatives" which is currently unchecked. Further down are two buttons: "Apply Filters =>" and "Clear Filters". At the bottom left of this panel, a button labeled "Filter History" is highlighted with a red rectangular box. To the right of the filter panel is a table with a list of gene symbols: CEP4, BARC, CDH1, EXO1, HBG1, SH3B, NPHS, MYO, FBP1, FUT6, ITGB, SPAT, AQP2, and PTC1. Below the table, there are navigation icons and the text "レコード 265 /". At the bottom of the interface, a status bar shows "Sample: One\_12880\_homo Genes, Variants: (3034, 11538) -> (445, 1976)".

Unaffected Siblings

Edit Unaffected Siblings

Use only passing variants in relatives

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
Filter History  
Sample: One\_12880\_homo Genes, Variants: (3034, 11538) -> (445, 1976)

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**OMIM**

OMIM is a comprehensive, authoritative compendium of human genes and genetic phenotypes that is freely available and updated daily. OMIM is authored and edited at the McKusick-Nathans Institute of Genetic Medicine, Johns Hopkins University School of Medicine, under the direction of Dr. Ada Hamosh. Its official home is [omim.org](http://omim.org).

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Number of Entries in OMIM (Updated October 10th, 2014) :

Prefix	Autosomal	X Linked	Y Linked	Mitochondrial	Totals
* Gene description	13,963	683	48	35	14,729
+ Gene and phenotype, combined	88	2	0	2	92
# Phenotype description, molecular basis known	3,930	287	4	28	4,249
% Phenotype description or locus, molecular basis unknown	1,542	134	5	0	1,681
Other, mainly phenotypes with suspected mendelian basis	1,730	113	2	0	1,845
Totals	21,253	1,219	59	65	22,596

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The screenshot displays a table with columns for PubMed, GeneReviews, and Disease. The 'Disease' column is highlighted with a red box. Below the table, the 'Genes' tab is selected and highlighted with a red box and a red arrow. A red arrow also points from the text 'Genes' タブを選択' to the 'Genes' tab.

PubMed	GeneReviews	Disease
<a href="#">A2M</a>		Alzheimer's disease; Alzheimer disease, type 3; Alzheimer disease, type 4; Alz
<a href="#">ABCA12</a>	NBK1420	Autosomal recessive congenital ichthyosis 4B; Autosomal recessive congenita
<a href="#">ABCA3</a>		Disorder of lung; Pulmonary alveolar proteinosis; Respiratory distress; Congen
<a href="#">ABCA7</a>		
<a href="#">ABCB4</a>		Cholecystitis; Cholestasis of pregnancy; Jaundice, familial obstructive, of infan
<a href="#">ABCC9</a>	NBK1309	Atrial fibrillation familial; Primary dilated cardiomyopathy; Arrhythmogenic right
<a href="#">ABCG8</a>	NBK131810	Sitosterolemia; Gallbladder disease 4
<a href="#">ABO</a>		
<a href="#">ACACB</a>		
<a href="#">ACBD5</a>		
<a href="#">ACO2</a>		Infantile cerebellar-retinal degeneration
<a href="#">ACSF3</a>		Combined malonic and methylmalonic aciduria
<a href="#">ACVR1C</a>		
<a href="#">ADAMTS16</a>		
<a href="#">ADAMTS17</a>		Weill-Marchesani-like syndrome
<a href="#">ADAMTS2</a>		Ehlers-Danlos syndrome; Dermatosparaxis
<a href="#">ADARB1</a>		
<a href="#">ADCY9</a>		

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The screenshot shows the Illumina VariantStudio interface. The 'Filters' panel on the left has a red box around the 'Gene' dropdown menu and a red arrow pointing to the 'Disease' input field, which contains the text 'Miller'. Below the input field are several filter options: 'Include List', 'Exclude List', 'Min Variant Alleles' (set to 2), and 'where custom gene annotation contains:'. The 'Gene View' table on the right displays a list of genes associated with the 'Miller' disease.


External Variants	Ambig... Variants	Custom Gene Annotation	PubMed	GeneReview
1	0		<a href="#">A2M</a>	
1	0		<a href="#">ABCA12</a>	NBK1420
1	0		<a href="#">ABCA3</a>	
1	5		<a href="#">ABCA7</a>	
2	0		<a href="#">ABCB4</a>	
1	0		<a href="#">ABCC9</a>	NBK1309
2	1		<a href="#">ABCG8</a>	NBK131810
1	0		<a href="#">ABO</a>	
1	0		<a href="#">ACACB</a>	
1	0		<a href="#">ACBD5</a>	
1	0		<a href="#">ACO2</a>	

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
Gene ▲	Variant	Chr	Coordinate	ClinVar Significance <sup>?</sup>	ClinVar Disease Name	Inherited From
DHODH	G>G/A	16	72050942	pathogenic	Miller_syndrome	father
DHODH	C>C/T	16	72057435	pathogenic	Miller_syndrome	mother

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The screenshot shows the Illumina VariantStudio interface. The 'Filters' panel on the left has 'Population Frequency' selected in the 'Consequence' dropdown. The 'Global Frequency' is set to 100, 'Asian Pop Frequency' is set to 1, and 'EVS Frequency' is set to 100. The 'Apply Filters =>' button is highlighted with a red box. A red arrow points from the 'Asian Pop Frequency' dropdown to the text 'アジア集団1%以下'. The main table displays a list of variants with columns for Gene, Variant, Chr, Coordinate, Classification, Type, Genot..., Exonic, Filters, Quality, and GQX. The table shows 1976 variants before filtering and 679 variants after filtering. A dark grey box at the bottom of the screenshot contains the text '1,976 => 679 バリエーション'.

Gene	Variant	Chr	Coordinate	Classification	Type	Genot...	Exonic	Filters	Quality	GQX
N/A	A>T/T	1	324822	...	snv	hom	no	PASS	2757.34	9%
N/A	A>G/G	1	262956340	...	snv	hom	no	PASS	3138.56	9%
N/A	T>G/G	1	262956340	...	snv	hom	no	PASS	2735.14	9%
N/A	T>C/C	1	219706428	...	snv	hom	no	PASS	1436.15	9%
N/A	A>G/G	2	3488502	...	snv	hom	no	PASS	1097.04	8:
N/A	G>C/C	2	10959475	...	snv	hom	no	PASS	1067.83	8:
N/A	A>G/G	2	131672950	...	snv	hom	no	PASS	2419.38	9%
N/A	C>T/T	2	131673316	...	snv	hom	no	PASS	431.14	3%
N/A	T>C/C	2	132229720	...	snv	hom	no	PASS	984.72	8%
N/A	C>G/G	3	136739406	...	snv	hom	no	PASS	588.87	4%
N/A	C>T/T	3	195467853	...	snv	hom	no	PASS	856.13	6%
N/A	G>A/A	4	5990339	...	snv	hom	no	PASS	1481.95	9%
N/A	A>G/G	4	5990791	...	snv	hom	no	PASS	5862.59	9%
N/A	T>C/C	4	5991384	...	snv	hom	no	PASS	1753.22	9%

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Gene	Variant	Chr	Coordinate	Classification	Type	Genot...	Exonic	Filters	Quality	GQX
N/A									2757.34	9%
N/A								PASS	3138.56	9%
N/A								PASS	2335.34	9%
N/A								PASS	984.72	8%
N/A								PASS	856.13	6%
N/A								PASS	1705.04	5%
N/A								PASS	1171.69	6%
N/A								PASS	2212.65	9%
N/A	A>G/G	14	106610469	...	snv	hom	no	PASS	2889.52	9%
N/A	C>T/T	14	106610509	...	snv	hom	no	PASS	2252.86	9%
N/A	T>A/A	14	106610535	...	snv	hom	no	PASS	1835.49	8%
N/A	A>G/G	14	106791119	...	snv	hom	no	PASS	2861.8	9%
N/A	C>G/G	14	106815944	...	snv	hom	no	PASS	1078.98	8%
N/A	A>G/G	15	22482858	...	snv	hom	no	PASS	1779.39	9%
N/A	A>G/G	15	22483277	...	snv	hom	no	PASS	790.92	6%

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
- ▶ VEP v2.8 (Ensembl)
  - <http://asia.ensembl.org/info/docs/tools/vep/index.html>
  - McLaren W, Pritchard B, Rios D, Chen Y, Flicek P, Cunningham F. Deriving the consequences of genomic variants with the Ensembl API and SNP Effect Predictor. **Bioinformatics** 26(16):2069-70(2010)
- ▶ PolyPhen
  - <http://genetics.bwh.harvard.edu/pph2/>
  - Adzhubei IA, Schmidt S, Peshkin L, Ramensky VE, Gerasimova A, Bork P, Kondrashov AS, Sunyaev SR. A method and server for predicting damaging missense mutations. **Nat Methods** 7(4):248-249 (2010).
- ▶ SIFT
  - <http://sift.jcvi.org/>
  - Kumar P, Henikoff S, Ng PC. Predicting the effects of coding non-synonymous variants on protein function using the SIFT algorithm. **Nat Protoc.** 2009;4(7):1073-81.
- ▶ ClinVar Version September 5, 2013 (NCBI)
  - [www.ncbi.nlm.nih.gov/clinvar/](http://www.ncbi.nlm.nih.gov/clinvar/)
  - Landrum MJ *et al.* ClinVar: public archive of relationships among sequence variation and human phenotype. **Nucleic Acids Res.** 2014 Jan 1;42(1):D980-5.

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- ▶ TruSight One User Guideなど  
[http://support.illumina.com/sequencing/sequencing\\_kits/trusight\\_one\\_kit/documentation.html](http://support.illumina.com/sequencing/sequencing_kits/trusight_one_kit/documentation.html)
- ▶ TruSight One Data Sheetなど  
[http://support.illumina.com/content/dam/illumina-marketing/documents/products/datasheets/datasheet\\_trusight\\_one\\_panel.pdf](http://support.illumina.com/content/dam/illumina-marketing/documents/products/datasheets/datasheet_trusight_one_panel.pdf)
- ▶ TruSight One 遺伝子一覧  
[http://products.illumina.com/content/dam/illumina-marketing/documents/products/gene\\_lists/gene\\_list\\_trusight\\_one.zip](http://products.illumina.com/content/dam/illumina-marketing/documents/products/gene_lists/gene_list_trusight_one.zip)
- ▶ TruSight One BaseSpace Public data  
<https://basespace.illumina.com/datacentral>

The screenshot shows the BaseSpace Public Data interface. The navigation bar includes Dashboard, Prep, Runs, Projects, Apps, Public Data, and Help. The main content area is titled 'Public Data' and features a search bar. Below the search bar, a list of data sets is displayed, each with a title and associated tags. Two entries are highlighted with red boxes:

- > **TruSight One Sequencing Panel MiSeq Trio Data**  
Targeted Sequencing
- > **NextSeq 500: TruSight One (CEPH Trio replicates)**  
Targeted Sequencing

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