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© 2012 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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Fox GE et al. (1977) PNAS

Classification of methanogenic bacteria by 16S ribosomal RNA characterization
(comparative oligonucleotide cataloging/phylogeny/molecular evolution)

GEORGE E. FOX*[†], LINDA J. MAGRUM*, WILLIAM E. BALCH[‡], RALPH S. WOLFE[‡], AND CARL R. WOESE*[‡]

Departments of *Genetics and Development and [†]Microbiology, University of Illinois, Urbana, Illinois 61801

Communicated by H. A. Barker, August 10, 1977



ABSTRACT The 16S ribosomal RNAs from 10 species of methanogenic bacteria have been characterized in terms of the oligonucleotides produced by T₁ RNase digestion. Comparative analysis of these data reveals the methanogens to constitute a distinct phylogenetic group containing two major divisions. These organisms appear to be only distantly related to typical bacteria.

residues represented by hexamers and larger in catalog A and in catalog B and their overlap of common sequences, respectively. The association coefficient, S_{AB} , so defined provides what is generally an underestimate of the true degree of homology between two catalogs because related but nonidentical oligomers are not considered. The matrix of S_{AB} values for each binary comparison among the members of a given set of

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1)

Open

The ISME Journal (2012), 1–5
© 2012 International Society for Microbial Ecology. All rights reserved. 1712-0202/12
www.ismej.com/content/12

SHORT COMMUNICATION

Selection of primers for optimal taxonomic classification of environmental 16S rRNA gene sequences

David AW Soergel¹, Neelendu Dey^{2,3,4}, Rob Knight³ and Steven E Brenner¹

¹Department of Plant and Microbial Biology, University of California, Berkeley, CA, USA; ²Division of Gastroenterology, Department of Medicine, University of California, San Francisco, CA, USA and ³Howard Hughes Medical Institute and Department of Chemistry and ⁴Biochemistry, University of Colorado, Boulder, CO, USA

Microbial community profiling using 16S rRNA gene sequences requires accurate taxonomy as algorithms. Universally primers target conserved sequences and amplify sequences from many taxa, but they provide variable coverage of different environments, and regions of the rRNA gene differ in taxonomic informativeness—especially when high-throughput short-read sequencing technologies (for example, 454 and Illumina) are used. We introduce a new evaluation procedure that provides an improved measure of expected taxonomic precision when classifying environmental sequence reads from a given primer. Applying this measure to thousands of combinations of primers and read lengths, simulating single-ended and paired-ended sequencing, reveals that these choices greatly affect taxonomic informativeness. The most informative sequence region may differ by environment, partly due to variable coverage of different environments in reference databases. Using our Rax method of classifying paired-end reads, we found that paired-end sequencing provides substantial benefit in some environments including human gut, but not in others. Optimal primer choices for short reads totaling 96 nt provides 82–100% of the confident genus classifications available from longer reads. The ISME Journal advance online publication: 12 January 2012; doi:10.1089/ismej.2011.208

Subject Category: microbial population and community ecology

Keywords: 16S ribosomal RNA; taxonomy; phylogeny; classification; bacteria; sequencing

Variation in 16S ribosomal gene sequences has been used since the mid-1980s to characterize microbial diversity (Stahl *et al.*, 1984). Interest in sequence-based surveys of environmental microbes has exploded in recent years with the availability of sequencing technologies that produce ever-larger data sets at ever-decreasing cost; in particular, the Illumina platform is attractive because of throughput, despite its short reads (Sogin *et al.*, 2006; Lazarovic *et al.*, 2009; Claesson *et al.*, 2012; Caporaso *et al.*, 2011; Dugan and Ochman, 2012). Here, we examine the reliability of assignment of novel sequences to known taxa under thousands of simulated conditions, varying primer choice, read length and environment.

Previous work on taxonomic classification of environmental 16S rRNA gene sequences has focused on whether reference sequences matching a given query share taxonomic annotations (Jonsson *et al.*, 2002; Desnais *et al.*, 2006; Sogin *et al.*, 2006; Wang *et al.*, 2007). Validations of taxonomic classifiers have typically compared a limited range of primers, read lengths and environments (Stundquist *et al.*, 2007; Huse *et al.*, 2008; Liu *et al.*, 2008; Wu *et al.*, 2008; Hays *et al.*, 2009). Reference databases contain many sequences derived from some environments and few associated with others, however (Supplementary Figure S1), leading to substantial variation in classification quality. In addition, the use of ‘best-ones-out’ cross-validation at the sequence level (Stundquist *et al.*, 2007; Wang *et al.*, 2007; Liu *et al.*, 2008; Wu *et al.*, 2008)—where a single sequence with a known annotation is held out from a reference database and classified using the remainder—is problematic: reference sequences matching best-of-sets sequences are likely to originate from the same sample, because natural environments contain ‘microdiverse’ clusters of closely related strains (Achilles *et al.*, 2004). We addressed these issues by simulating taxonomic reads from slight large environmental data

Correspondence: DAW Soergel, Current address: Department of Computer Science, University of Massachusetts, 140 Governors Drive, Amherst, MA 01003-0004, USA. E-mail: soergel@cs.umass.edu

Current address: Center for Genome Sciences and Systems Biology, Washington University School of Medicine, Saint Louis, MO 63108, USA.

Current address: Division of Gastroenterology, Department of Internal Medicine, Washington University School of Medicine, Saint Louis, MO 63108, USA.

Received 20 March 2011; revised 21 November 2011; accepted 3 December 2011

2)

Published online 28 August 2012

Nucleic Acids Research, 2013, Vol. 41, No. 1, e1
doi:10.1093/nar/gks808

Evaluation of general 16S ribosomal RNA gene PCR primers for classical and next-generation sequencing-based diversity studies

Anna Klindworth^{1,2}, Elmar Pruesse^{1,2}, Timmy Schweer¹, Jörg Peplies³, Christian Quast¹, Matthias Horn¹ and Frank Oliver Glöckner^{1,2,4}*

¹Max Planck Institute for Marine Microbiology, Microbial Genomics and Bioinformatics Research Group, Celsiusstr. 1, 28359 Bremen, ²Jacobs University Bremen, School of Engineering and Sciences, Campus 1, 28759 Bremen, ³Ribicon GmbH, D-28359 Bremen, Germany and ⁴Department of Microbial Ecology, University of Vienna, Althanstr. 14, 1090 Vienna, Austria

Received December 12, 2011; Accepted July 31, 2012

ABSTRACT

16S ribosomal RNA gene (rDNA) amplicon analysis remains the standard approach for the cultivation-independent investigation of microbial diversity. The accuracy of these analyses depends strongly on the choice of primers. The overall coverage and phylum spectrum of 175 primers and 512 primer pairs were evaluated *in silico* with respect to the SILVA 16S/18S rDNA non-redundant reference dataset (SSURF 108 NR). Based on this evaluation a selection of ‘best available’ primer pairs for Bacteria and Archaea for three amplicon size classes (100–400, 400–1000, >1000 bp) is provided. The most promising bacterial primer pair (S-D-Bact-0341-b-S-17/S-D-Bact-0785-a-A-21), with an amplicon size of 464 bp, was experimentally evaluated by comparing the taxonomic distribution of the 16S rDNA amplicons with 16S rDNA fragments from directly sequenced metagenomes. The results of this study may be used as a guideline for selecting primer pairs with the best overall coverage and phylum spectrum for specific applications, therefore reducing the bias in PCR-based microbial diversity studies.

INTRODUCTION

Understanding microbial diversity has been the ambition of scientists for decades. Because diversity analysis by cultivation is problematic for a significant fraction of Bacteria and Archaea, culture-independent surveys have been developed. In the past, the most commonly used approach was cloning and sequencing of the 16S ribosomal RNA gene (rDNA) using conserved broad-range

PCR primers (1). With the advent of massive parallel sequencing technologies, direct sequencing of PCR amplicons became feasible (2–4). In 2006, Roche’s 454 GS 20 pyrosequencing (5) became the first high-throughput sequencing technology to be successfully applied for large scale biodiversity analysis and was key to uncovering the ‘rare biosphere’ (6). The continuous development of the technology, offering read lengths of up to 1000 bp nowadays, further improved throughput and resolution of 16S rDNA sequencing (7). Since then, additional high-throughput sequencing technologies have become commercially available. The attractiveness of Illumina (8) lies in the reduced per base costs and comparatively high sequencing depth (9), despite having short read lengths. While the major advantage of Ion Torrent (10) are its relatively low cost and rapid sequencing speed. Furthermore, Pacific Bioscience (PacBio) now employs the ‘single-molecule real-time’ (SMRT) sequencing technology, designed to achieve average read lengths of more than 3000 bp (11). For a detailed review of sequencing technologies please refer to Loman *et al.* (12). There is no doubt that the rapid development of sequencing technologies has opened a new dimension in biodiversity analysis, but the diversity of technologies also adds complexity to the experimental design of a study. The most critical step for accurate rDNA amplicon analysis is still the choice of primers (4,13). Using suboptimal primers, or more precisely, primer pairs, can lead to under-representation (14) or selection against single species or even whole groups (15–17). Using inappropriate primers consequently leads to questionable biological conclusions (17–19).

In this study, 175 broad range 16S rDNA primers and 512 primer pairs were investigated *in silico* with respect to overall coverage and phylum spectrum for Bacteria and Archaea as well as amplicon length. Primer sequences were compared with all 376 437 16S/18S rDNA sequences

*To whom correspondence should be addressed. Tel: +49 421 2028970; Fax: +49 421 2028580; Email: fog@mpi-bremen.de

© The Author(s) 2012. Published by Oxford University Press. This is an Open Access article distributed under the terms of the Creative Commons Attribution License (<http://creativecommons.org/licenses/by-nc/2.0/>), which permits non-commercial reuse, distribution, and reproduction in any medium, provided the original work is properly cited. For commercial re-use, please contact journals.permissions@oup.com.

1) David AW Soergel, et al., Selection of primers for optimal taxonomic classification of environmental 16S rRNA gene sequences. *The ISME Journal* (2012), 1–5

2) Klindworth A, et al., Evaluation of general 16S ribosomal RNA gene PCR primers for classical and next-generation sequencing-based diversity studies. *Nucleic Acids Res* (2013), 41(1)

16S Metagenomics Sequencing Library Preparation

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
Thanks to Bo-Young Hong and Patricia Diaz, Univ. Conn. Health Center.
830 sequences from the Human Oral Microbiome Database classified
using mothur/RDP.

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
Analysis | Imaging | Summary | **Tile Status** | TruSeq Controls | Indexing

Run Summary

Level	Yield Total (G)	Projected Total Yield (G)	Yield Perfect (G)	Yield <=3 errors (G)	Aligned (%)	% Perfect [Num Cycles]	% <=3 errors [Num Cycles]	Error Rate (%)	Intensity Cycle 1	% Intensity Cycle 20	% >= Q30
Read 1	6.2	6.2	0.6	0.8	13.59	67.4 [300]	95.3 [300]	0.92	25	118.1	84.1
Read 2 (I)	0.1	0.1	0.0	0.0	0.00	0.0 [5]	0.0 [5]	0.00	518	0.0	81.6
Read 3	6.2	6.2	0.5	0.8	13.52	58.9 [300]	91.6 [300]	1.06	38	91.2	71.2
Total	12.5	12.5	1.1	1.6	13.55	63.2	93.5	0.99	194	104.6	77.7

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- ▶ Weinstock Lab
 - George Weinstock
 - Erica Sodergren
 - Brandi Herter
 - Kathie Mihindukulasuriya
- ▶ TGI Illumina Production
 - Matt Cordes

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Rank	Species	Total # Sequences in sample	Total % Sequences in Sample	Notes
1	<i>Calothrix parietina</i>	106,004	38.11%	Blue-green algae
2	<i>Rickettsia</i> sp.	15,631	5.62%	Tick, flea, lice borne
3	<i>Acinetobacter</i> sp.	5,764	2.07%	Soil bacteria
4	<i>Symploca atlantica</i>	4,917	1.77%	Algae
5	<i>Proteobacteria</i>	4,109	1.48%	Variable (some pathogen known)
6	<i>Thiomonas thermosulfata</i>	3,847	1.38%	Extremophile
8	<i>Blautia</i> sp.	3,705	1.33%	Gut bacteria
9	<i>Blautia coccooides</i>	3,330	1.20%	Gut Bacteria
10	<i>Pedobacter</i>	2,881	1.04%	Soil Bacteria
Other Species of note	<i>Mycoplasma agassizii</i>	138	0.05%	

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Rank	Species	Total # Sequences in sample	Total % sequences in sample	Notes
1	<i>Mycoplasma agassizii</i>	31,924	18.2%	Etiologic agent*
2	<i>Myroides odoratus</i>	11,654	6.6%	
3	<i>Flavobacteriaceae sp.</i>	11,481	6.5%	
4	<i>Chelonobacter sp.</i>	11,384	6.5%	Associated diseased tortoise **
5	<i>Pedobacter sp.</i>	7,894	4.5%	Soil bacteria
6	<i>Flavobacterium swingsii</i>	5,492	3.1%	
7	<i>Proteus penneri</i>	5,263	3.0%	Invasive pathogen
8	<i>Pseudomonas brenneri</i>	4,538	2.6%	
9	<i>Pseudomonas sp.</i>	4,193	2.4%	Water borne biofilms
10	<i>Pseudomonas marginalis</i>	4,160	2.4%	

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Rank	Species	Total # Sequences in Sample	Total % sequences in sample	Notes
1	<i>Chelonobacter sp.</i>	141,734	57.6%	Associated with disease tortoise **
2	<i>Chelonobacter oris</i>	37,883	15.4%	
3	<i>Mycoplasma agassizii</i>	7,639	3.1%	Etiologic agent *
4	<i>Granulicatella adiacens</i>	5,825	2.4%	Normal human mucosal membranes
5	<i>Flavobacteriaceae</i>	5,443	2.2%	Water borne
6	<i>Myroides odoratus</i>	4,855	2.0%	Hm nosocomial infection
7	<i>Pedobacter sp.</i>	4,108	1.7%	Soil bacteria
8	<i>Deinococcus sp.</i>	2,051	0.8%	Soil, water
9	<i>Flavobacterium swingsii</i>	1,800	0.7%	Water borne
10	<i>Proteus penneri</i>	1,625	0.7%	Invasive pathogen

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