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© 2012 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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

Status Icon	Status Name	Description
	Status OK	No change. System is normal.
	Attention	Important information. Action is recommended.
	Warning	Warnings do not stop a run. However, some warnings require action before proceeding.
	Error	Errors usually stop a run and generally require action before proceeding with the run.

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Resolve Run Setup Errors

If any checks in the pre-run check fail, a red icon **X** appears next to the item. A message appears on the screen that describes the error and how to correct it.

Error	Action
X Flow Rate Measured	<p>The flow rate check screen opens. Use the drop-down list or on-screen keyboard to enter the following:</p> <ul style="list-style-type: none">• Solution: PR2• Volume: 250• Aspirate Rate: 2500• Dispense Rate: 2500 <p>Select Pump. If the error persists, set the volume to pump 500 µl PR2 and repeat the process. When fluids have been pumped, select Restart Check.</p> <p>When the pre-run check is successful, the Start Run button becomes active.</p> <p>If the flow check fails again, reseal the flow cell to make sure that flow is not interrupted due to misalignment. Inspect the flow cell gasket for lint or irregularities.</p>
X Free Disk Space	<p>If disk space is low, a message appears indicating how much disk space is required. Use the Manage Files feature to clear the required space from the instrument computer.</p>
X Network Connection	<p>Make sure that the network cable is plugged into the instrument.</p> <p>If the network connection is not restored, select Reboot on the Manage</p>

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The screenshot shows the Illumina website homepage. At the top left is the Illumina logo. To the right are links for "お問い合わせ" (Contact Us), "MyIllumina", and "Tools". Below this is a horizontal navigation menu with categories: "アプリケーション" (Applications), "システム" (Systems), "インフォマティクス" (Informatics), "臨床研究" (Clinical Research), "受託サービス" (Contract Services), "サイエンス" (Science), "サポート" (Support), and "カンパニー" (Company). A search bar is located on the right side of the navigation menu. Below the navigation menu are links for "Subscribe" and "Follow us:". The main content area features a large orange headline: "研究者による最新発表からシステムの使い方まで" (From the latest publications by researchers to the system's usage). Below the headline is the text: "ウェビナーを活用しオンラインで情報を収集" (Utilize webinars to collect information online). The background of the main content area shows a woman in a white lab coat holding a tablet, with various scientific charts and graphs displayed on the screen. On the right side of the page, there is a "製品検索" (Product Search) section with a dropdown menu for "システムに適したキットを探す" (Find kits suitable for the system). Below this are links for "NGSサンプル調製キット検索" (NGS sample preparation kit search), "新製品 & 人気の製品" (New products & popular products), and "すべての製品をみる" (View all products). At the bottom right, there is a "ログイン" (Login) section with input fields for "Username/Email" and "Password".

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illumina® Log in to get personalized account information. Quick Order View Cart

Contact Us MyIllumina Tools

APPLICATIONS SYSTEMS INFORMATICS CLINICAL SERVICES SCIENCE SUPPORT COMPANY Search

Support » Sequencing » Sequencing Instruments » MiSeq

MiSeq Support

- Overview
- Site Prep/Lab Environment
- Requirements & Compatibility
- Supported Kits
- Downloads
- Documentation & Literature
- Training
- Questions & Answers
- Troubleshooting
- Services & Warranties
- Bulletins
- Webinars
- Product Ordering Information

MiSeq

Latest Updates

- De Novo Assembly of Bacterial Genomes Application Note 05/15/2014
- Interview with Dr. James Hadfield 01/20/2014
- Low-Diversity Sequencing on the Illumina MiSeq Platform 06/20/2013

User Guides

- MiSeq System User Guide (15027617 M)

Troubleshooting

MiSeq Online Troubleshooting

Please log in with your MyIllumina customer account information to access the troubleshooting guide.

* Username/Email:

* Password: [Forgot your Password?](#)

LOGIN

By clicking on the "Login" button you indicate that you agree with our [Terms and Conditions](#).

Demo

See how to troubleshoot common MiSeq issues with this interactive tool. View the demo.

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illumina [Contact Us](#) [MyIllumina](#) [Tools](#)


APPLICATIONS SYSTEMS INFORMATICS CLINICAL SERVICES SCIENCE SUPPORT COMPANY

Support » Sequencing » Sequencing Instruments » MiSeq » **Troubleshooting**

Troubleshooting

使用方法をビデオで解説!
MiSeq Online Troubleshooting Demo

MiSeq Online Troubleshooting Demo

See how to troubleshoot common MiSeq issues with this interactive tool. View the demo.

[View Demo >](#)

Online Troubleshooting

[MiSeq Home](#) [Feedback](#) [Contact Technical Support](#)

Instructions Mouse over and click on the system image to troubleshoot.	MiSeq Would you like to troubleshoot a MCS Error or another MiSeq problem? <input type="radio"/> MiSeq Control Software Error <input checked="" type="radio"/> Other MiSeq Issues
--	---

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The screenshot shows the MiSeq Online Troubleshooting page. At the top, there is a navigation menu with categories: APPLICATIONS, SYSTEMS, INFORMATICS, CLINICAL, SERVICES, SCIENCE, SUPPORT, and COMPANY. A search bar is located to the right of the menu. Below the menu, a breadcrumb trail reads: Support » Sequencing » Sequencing Instruments » MiSeq » Troubleshooting. The main heading is "Troubleshooting". A central section titled "MiSeq Online Troubleshooting Demo" features a video player thumbnail. Two red-bordered boxes are overlaid on the page: the first box contains the text "MiSeq Control Software(MCS)でのエラーについてはトラブルシュートに進むため、MiSeq Control Software Errorを選択" and the second box contains "MiSeq Reporter(MSR)やComputer, ランのクオリティについてはOther MiSeq Issuesを選択". Below the demo section is a "View Demo »" link. At the bottom, there is a footer with links for "MiSeq Home", "Feedback", and "Contact Technical Support". The "Instructions" section on the left states: "Mouse over and click on the system image to troubleshoot." The "MiSeq MCS Error" section contains a question: "Would you like to troubleshoot a MCS Error or another MiSeq problem?" and two radio button options: "MiSeq Control Software Error" (which is selected and highlighted with a yellow box) and "Other MiSeq Issues" (also highlighted with a yellow box).

APPLICATIONS SYSTEMS INFORMATICS CLINICAL SERVICES SCIENCE SUPPORT COMPANY Search

Support » Sequencing » Sequencing Instruments » MiSeq » Troubleshooting

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View Demo »

MiSeq Home Feedback Contact Technical Support

Instructions
Mouse over and click on the system image to troubleshoot.

MiSeq MCS Error
Would you like to troubleshoot a MCS Error or another MiSeq problem?

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The screenshot displays the MiSeq Online Troubleshooting interface. At the top, there are navigation links: [MiSeq Home](#), [Feedback](#), and [Contact Technical Support](#). The main content area is divided into three columns. The left column is titled "MCS Error" and contains a sub-section "MCS Error". The middle column is also titled "MCS Error" and contains a sub-section "MCS Error" with a dropdown arrow. Below this, it asks "What MCS error message do you see?" and lists 15 radio button options. The first option, "Unable to measure flow rate/flow check failed", is selected and highlighted with a red circle and a yellow rounded rectangle. The right column is titled "MiSeq MCS Error" and contains the text: "The current question has no associated content. Please contact Illumina Technical Support using the button at the top of the screen if additional assistance is required to answer the question." A red-bordered box on the right contains the Japanese text: "MCS Errorのリストから、該当するエラーを選択する".

MiSeq Home Feedback Contact Technical Support

MCS Error

MCS Error

MCS Error

MCS Error

What MCS error message do you see?

- Unable to measure flow rate/flow check failed
- Sample sheet will not load
- PR2 bottle not detected
- RFID not recognized
- Object reference not set to an instance of an object
- No usable signal found in the images; clustering may have failed.
- Z-motor outside soft limits.
- Best focus too near edge of range
- Failed to get temperature
- Focus Mirror error
- FPGA error/FPGA Timeout
- Lane pump error
- Reagent valve error
- RFID cannot write error
- SIPDOWN error.
- Y motor error

MiSeq MCS Error

The current question has no associated content.

Please contact Illumina Technical Support using the button at the top of the screen if additional assistance is required to answer the question.

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MiSeq Online Troubleshooting

MiSeq Control Software Error

The screenshot displays the MiSeq Online Troubleshooting interface. At the top, there are navigation links: "MiSeq Home", "Feedback", and "Contact Technical Support". The main content area is divided into three columns:

- MCS Error**: A sidebar menu with the text "MCS Error".
- MCS Error**: The main content area showing an error message: "The MiSeq cannot complete a flow check, or measure the flow rate during the flow check. This is often caused by the flowcell not being aligned with the flowcell manifold. It can also be caused by blockages in the flowcell or reagent lines, or by communications issues in the instrument." Below the message are two radio button options: "Proceed to Troubleshooting" (highlighted with a yellow box) and "Return to Symptoms". A "back" button is also visible.
- MiSeq MCS Error**: A sidebar menu with the text "The current question has no associated content." Below this, there are three red-bordered boxes containing Japanese text:
 - エラーの内容についてと、一般的にこのエラーが生じる原因が表示される。
 - このままトラブルシューティングに進む場合 Proceed to Trouble shooting
 - エラーの選択リストに戻る場合 Return to Symptoms

MiSeq Online Troubleshooting

MiSeq Control Software Error

MiSeq Home Feedback Contact Technical Support

MCS Error

MCS Error

MCS Error

A volume test checks for obstructions in the fluid lines of the MiSeq. This test requires a used flowcell and a wash tray and wash bottle filled with laboratory-grade water. Full instructions for performing the test are shown at right.

Did the volume test fail?

Yes

No

[back](#)

Perform a Volume Test

An obstruction in the fluidics lines can cause poor reagent delivery and affect sequencing results. If an obstruction in the fluidics lines is suspected, perform a volume test.

A volume test checks the health of the fluidics system by estimating the volume between two bubbles as they pass by the sensors. To perform a volume test, the wash tray and wash bottle must be loaded with laboratory-grade water and a used flow cell must be in place. Follow the onscreen prompts to perform the test.

- 1 Make sure that a used flow cell is loaded on the instrument.
- 2 From the Manage Instrument screen, select **System Check**.
- 3 Select **Conduct Volume Test**, and then select **Next**.

System Check Screen

Please select the system check options to perform.

Select All

Motion System	Optics	Other
<input type="checkbox"/> Y Stage Test	<input type="checkbox"/> LED Power Test	<input type="checkbox"/> Prime Reagent Lines
<input type="checkbox"/> M3 Mirror Test	<input type="checkbox"/> 1p/18 Test	<input checked="" type="checkbox"/> Conduct Volume Test

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MiSeq Home Feedback Contact Technical Support

MCS Error

MCS Error

MCS Error

If the volume test failed, this may mean that there is an obstruction in the fluid lines. Please run a full maintenance wash as shown in the instructions to the right. This will require a used flowcell.

Did this resolve the problem?

Yes

No

« back

Perform a Maintenance Wash

Perform a maintenance wash every 30 days to ensure optimal performance. A maintenance wash includes a series of three wash steps using a wash solution of laboratory-grade water mixed with Tween 20. Allow approximately 90 minutes to complete the wash.

User-Supplied Consumables

- Tween 20 (Sigma-Aldrich, catalog # P7949)
- Laboratory-grade water

Part # 15027617 Rev. L

MiSeq Online Troubleshooting

MiSeq Control Software Error

MiSeq Home Feedback Contact Technical Support

MCS Error

MCS Error

MCS Error

Please Contact Illumina Technical Support

The next step in the troubleshooting process requires input from Illumina Technical Support. Please choose "Help me contact Illumina Technical Support" below to access contact information for Tech Support and a code to allow representatives to access troubleshooting steps completed so far.

[Help me contact Technical Support.](#)

[back](#)

Perform a Maintenance Wash

Perform a maintenance wash every 30 days to ensure optimal performance. maintenance wash includes a series of three wash steps using a wash solution of laboratory-grade water mixed with Tween 20. Allow approximately 90 minutes to complete the wash.

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MiSeq Online Troubleshooting

Contact Technical Support

Contact Illumina Technical Support

Do you want to leave the Illumina Online Troubleshooting Tool and contact Technical Support?

Yes

No, Return To Troubleshooting

Contacting Illumina Technical Support

1. Write down your troubleshooting session ID: **18020**. Technical Support uses this ID to continue your troubleshooting session.
2. Contact Technical Support:
 - ◆ phone click [here](#).
 - ◆ email click [here](#).
4. Your troubleshooting session history is available as a PDF at the link below.

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
illumina®		
This online troubleshooting investigation (#18320) was initiated on July 08, 2014 06:51 PM and completed on July 08, 2014 06:57 PM EST.		
Date and Time	Question	Answer
07/08/2014 06:51	System: MiSeq	
07/08/2014 06:51	Would you like to troubleshoot a MCS Error or another MiSeq problem?	mcs_error
07/08/2014 06:51	What MCS error message do you see?	flow_check_failed
07/08/2014 06:51	The MiSeq cannot complete a flow check, or measure the flow rate during the flow check. This is often caused by the flowcell not being aligned with the flowcell manifold. It can also be caused by blockages in the flowcell or reagent lines, or by communications issues in the instrument.	proceed
07/08/2014 06:51	A volume test checks for obstructions in the fluid lines of the MiSeq. This test requires a used flowcell and a wash tray and wash bottle filled with laboratory-grade water. Full instructions for performing the test are shown at right. Did the volume test fail?	yes
07/08/2014 06:55	If the volume test failed, this may mean that there is an obstruction in the fluid lines. Please run a full maintenance wash as shown in the instructions to the right. This will require a used flowcell. Did this resolve the problem?	no
07/08/2014 06:57	The next step in the troubleshooting process requires input from Illumina Technical Support. Please choose "Help me contact Illumina Technical Support" below to access contact information	Help me contact Technical Support.

MiSeq Online Troubleshooting

Other MiSeq Issues

Instructions

Mouse over and click on the system image to troubleshoot.

MiSeq MCS Error

Would you like to troubleshoot a MCS Error or another MiSeq problem?

- MiSeq Control Software Error
- Other MiSeq Issues

Instructions

Mouse over and click on the system image to troubleshoot.

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
MiSeq Online Troubleshooting

Feedback and Contact Technical Support

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The screenshot displays the MiSeq online troubleshooting interface. At the top, there are three buttons: "MiSeq", "Feedback", and "Contact Technical Support". The "Feedback" and "Contact Technical Support" buttons are highlighted with yellow boxes. A large yellow arrow points from the "Feedback" button to a "Feedback" dialog box that is open in the foreground. The dialog box contains the following text:

Feedback [close](#)

Thank you for your feedback!

[Cancel](#) [Submit Feedback](#)

We use your feedback to improve the Illumina online troubleshooting tool.


For all other issues, contact Illumina Technical Support.

Category:

Summary:

The dialog box also features a rich text editor toolbar with various icons for text formatting, alignment, and insertion.

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Illumina MiSeq

BaseSpace Options | Load Flow Cell | Load Reagents | Review | **Pre-Run Check** | Sequence | Post-Run Wash

Pre-Run Check is successful. Select Start Run to begin sequencing.

- ✓ Network Connection Active
- ✓ Flow Cell Cartridge Loaded
- ✓ Waste Bottle Loaded
- ✓ Reagent Cartridge Loaded
- ✓ Chiller Door Closed
- ✓ Primary Analysis Ready
- ✓ Sample Sheet Present
- ✓ Free Disk Space (20%)
- ✓ PR2 Bottle Loaded
- ✓ Flow Cell Door Closed
- ✓ **Flow Rate Measured**

Reagent ID: MS2000033-150V3
PR2 Bottle ID: MS2000045-00PR2
Flow Cell Barcode: 000000000-A4UBM

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<p>MCS Error</p> <p>MCS Error</p>	<p>MCS Error</p> <p>The MiSeq cannot complete a flow check, or measure the flow rate during the flow check. This is often caused by the flowcell not being aligned with the flowcell manifold. It can also be caused by blockages in the flowcell or reagent lines, or by communications issues in the instrument.</p> <p><input checked="" type="radio"/> Proceed to Troubleshooting</p> <p><input type="radio"/> Return to Symptoms</p> <p>⏪ back</p>	<p>MiSeq MCS Error</p> <p>The current question has no associated content.</p> <p>Please contact Illumina Technical Support using the button at the top of the screen if additional assistance is required to answer the question.</p>
--	--	--

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
⏪ back

Perform a Maintenance Wash

Perform a maintenance wash every 30 days to ensure optimal performance. The maintenance wash includes a series of three wash steps using a solution of laboratory-grade water mixed with Tween 20. Allow approximately 90 minutes to complete the wash.

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The image shows a screenshot of an Illumina instrument's user interface. On the left, there is a panel titled "MCS Error" with a sub-header "MCS Error". To its right is a larger panel, also titled "MCS Error", which contains a yellow highlighted box with the following text: "Please Contact Illumina Technical Support". Below this, it explains that the next step in the troubleshooting process requires input from Illumina Technical Support and provides a button labeled "Help me contact Technical Support." To the right of this is another panel titled "Perform a Maintenance Wash", which provides instructions on how to perform a maintenance wash every 30 days and lists "User-Supplied Consumables" such as Tween 20 and Laboratory-grade water. At the bottom right of this panel, it says "Part # 15027817 Rev. L". A red-bordered box is overlaid on the bottom left of the screenshot, containing the Japanese text "テクニカルサポートまでお問い合わせください。".

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[Help me contact Technical Support.](#)

Perform a Maintenance Wash

Perform a maintenance wash every 30 days to ensure optimal perform
The maintenance wash includes a series of three wash steps using a solution of laboratory-grade water mixed with Tween 20. Allow approxir
90 minutes to complete the wash.

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
A volume test checks the health of the fluidics system by estimating the volume between two bubbles as they pass by the sensors. To perform volume test, the wash tray and wash bottle must be loaded with laboratory-grade water and a used flow cell must be in place. Follow the onscreen prompts to perform the test.

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System Check Screen

The screenshot shows the 'System Check Screen' on the MiSeq instrument. It displays a list of system check options with checkboxes. Under the 'Other' category, the 'Conduct Volume Test' option is checked. A green checkmark icon is visible in the top right corner of the screen.

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MCS Error

Reset the flowcell:

1. Remove the flowcell from the stage.
2. Clean the flowcell and stage, making sure there is no salt or dried reagent near the gaskets or on the stage.
3. Clean the underside of the latch where it comes in contact with the gaskets.
4. Inspect the gaskets on the underside of the flowcell to make sure they are in place.
5. Set the flowcell back into place and close the latch.

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
A volume test checks the health of the fluidics system by estimating the volume between two bubbles as they pass by the sensors. To perform volume test, the wash tray and wash bottle must be loaded with laboratory-grade water and a used flow cell must be in place. Follow the onscreen prompts to perform the test.

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System Check Screen

The screenshot shows the 'System Check Screen' with a title bar 'Illumina MiSeq'. Below the title bar, it says 'Please select the system check options to perform.' There are three columns of checkboxes: 'Motion System' (Y Stage Test, M3 Mirror Test), 'Optics' (LED Power Test, Tip/Tip Test), and 'Other' (Prime Reagent Lines, Conduct Volume Test). The 'Conduct Volume Test' checkbox is checked. A green checkmark icon is visible in the top right corner of the screen.

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MCS Error

MCS Error

Turning the instrument off and on again often clears communication issues .

Note that if a run is in progress, it will not be able to be resumed. If this is the case you may wish to wait until the run is complete or consult Technical Support.

Contacting Technical Support is also advisable if the issue appears frequently.

MiSeq Powercycle Instructions:

1. Go to the Start menu and select Shut Down.
2. Once the screen has gone black, turn off the power to the MiSeq using the switch on the back.
3. Wait five minutes.
4. Turn on the MiSeq using the switch on the back.

Did powercycling the MiSeq resolve the problem?

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
MCS Error

Reseat the flowcell:

1. Remove the flowcell from the stage.
2. Clean the flowcell and stage, making sure there is no salt or dried reagent near the gaskets or on the stage.
3. Clean the underside of the latch where it comes in contact with the gaskets.
4. Inspect the gaskets on the underside of the flowcell to make sure they are in place.
5. Set the flowcell back into place and close the latch.

Try a different flowcell. Clean the stage and underside of the latch before loading the new flowcell.

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
A volume test checks the health of the fluidics system by estimating the volume between two bubbles as they pass by the sensors. To perform volume test, the wash tray and wash bottle must be loaded with laboratory-grade water and a used flow cell must be in place. Follow the onscreen prompts to perform the test.

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
A volume test checks the health of the fluidics system by estimating the volume between two bubbles as they pass by the sensors. To perform volume test, the wash tray and wash bottle must be loaded with laboratory-grade water and a used flow cell must be in place. Follow the onscreen prompts to perform the test.

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System Check Screen

The screenshot shows the 'System Check Screen' with a green checkmark in the top right corner. Below the title, it says 'Please select the system check options to perform.' There are three columns of options: 'Motion System' with 'Y Stage Test' and 'M3 Motor Test', 'Optics' with 'LED Power Test' and 'Light Test', and 'Other' with 'Reagent Lines' and 'Conduct Volume Test'. The 'Conduct Volume Test' option is checked with a blue checkmark.

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
Thank You!

Dear Customer,

We hope you have found the Online Troubleshooter useful!

For your convenience, a PDF of your troubleshooting session is available. Please click an option below to access the PDF, or if more assistance is required please click "Contact Tech Support" to receive a code that Technical Support can use to access your troubleshooting session.

Problem Addressed. Click for
contact Technical Support.

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A volume test checks the health of the fluidics system by estimating the volume between two bubbles as they pass by the sensors. To perform volume test, the wash tray and wash bottle must be loaded with laboratory-grade water and a used flow cell must be in place. Follow the onscreen prompts to perform the test.

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
Motion System	Optics	Other
<input type="checkbox"/> Select All		
<input type="checkbox"/> Y Stage Test	<input type="checkbox"/> LED Power Test	<input type="checkbox"/> Reagent Lines
<input type="checkbox"/> M3 Motor Test	<input type="checkbox"/> Light Test	<input checked="" type="checkbox"/> Conduct Volume Test

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Failed To Read PR2 RFID

Failed to read the RFID on the PR2. If this is the first failure, please attempt to retry. If you repeatedly receive this message with this PR2, please select 'Get Code' below to report the issue and receive a temporary RFID bypass code. The bypass code will allow you to enter the PR2 barcode using the on-screen keyboard.

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Reagent Compartment Symptoms

- Chiller Temperature
- Leak
- Noises
- PR2**
- Reagent / wash cartridge
- Visibly broken part

MiSeq

Run Quality
Software
Computer

The image shows a MiSeq sequencing machine. A yellow rounded rectangle highlights the reagent compartment area, which includes a reagent tray, a reagent bottle, and a wash cartridge. To the right of the machine, a bracket groups the terms 'Run Quality', 'Software', and 'Computer'.

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Reagent Compartment Symptoms

- Chiller Temperature
- Leak
- Noises
- PR2**
- Reagent / wash cartridge
- Visibly broken part

PR2

PR2 Problem

What problem do you see with PR2?

- RFID errors
- Visibly broken part
- Green liquid in PR2
- PR2 bottle not recognized

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The image shows a troubleshooting interface for a PR2 reagent compartment. On the left, a sidebar lists symptoms: Chiller Temperature, Leak, Noises, PR2 (highlighted), Reagent / wash cartridge, and Visibly broken part. The main panel is titled 'PR2' and contains a 'Bypass Code' section with instructions to obtain a code from the user guide. Below this, a question asks 'Does the bypass code work?' with two radio button options: 'Yes, bypass code works' and 'No, bypass code does not work'. The 'No' option is selected and circled in red. A 'back' button is visible below the options.

Resolve RFID Read Failure

If the system cannot read the RFID of a consumable, you can obtain a temporary bypass code from the Illumina website. A temporary bypass code expires in seven days.

- 1 Always select **Retry** before proceeding. If the RFID failed a second time, select **Get Code**.
- 2 From a computer with internet access, go to my.illumina.com and log in to your MyIllumina account.
- 3 From the MyIllumina page, click **Account**. In the Resources column, click **MiSeq Self-Service**.
- 4 On the MiSeq Self-Service page, enter the **MiSeq serial number**.
- 5 From the Type of Override Code drop-down list, select **RFID Override**.

MiSeq Self-Service Page

MyIllumina / Account / MiSeq Self Service

Self Service for MiSeq

Home | MiSeq Serial Number | Note: The MiSeq serial number can be found under the title "Instrument Name" in the About menu.

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PR2

New Reagents

Try loading new reagents affected by the RFID. Does this resolve the issue?

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Cause: Reagent RFID

Reagent RFID

- Problem Successfully Addressed!
- Problem still exists. Help me contact Technical Support.

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PR2

Power Cycle

Please power cycle your instrument following the instructions in the right hand panel. Does it fix the problem?

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PR2

Cause: Communication Problem

This is an intermittent communication problem

-  Problem Successfully Addressed!
-  Problem still exists. Help me contact Technical Support.

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<p>Reagent Compartment Symptoms</p> <ul style="list-style-type: none">Chiller TemperatureLeakNoisesPR2Reagent / wash cartridgeVisibly broken part	<p>PR2</p> <p>Please Contact Illumina Technical Support</p> <p>The next step in the troubleshooting process requires input from Illumina Technical Support. Please choose "Help me contact Illumina Technical Support" below to access contact information for Tech Support and a code to allow representatives to access troubleshooting steps completed so far.</p> <p> Help me contact Technical Support.</p>	<p>MiSeq MCS Error</p> <p>The current question has no associated content.</p> <p>Please contact Illumina Technical Support using the button at the top of the screen if additional assistance is required to answer the question.</p>
--	---	--

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Directory: D:\Illumina\MiSeqOutput Browse...

Bundle Logs

Run QC Files Thumbnail Images Run Logs Maintenance Logs


Date	Reagent Barcode	Workflow	R1	I1	R2	I2	Username
<input type="radio"/> 03/20/14	MS2006246-600V3	GenerateFASTQ	251	0	251	0	sbsuser
<input type="radio"/> 04/18/14	AA1234567-12345	Resequencing	5	0	0	0	sbsuser
<input type="radio"/> 05/01/14	MS2624960-300V2	Amplicon	151	8	151	8	sbsuser
<input type="radio"/> 05/15/14	MS2644485-300V2	Metagenomics	151	12	151	0	sbsuser
<input type="radio"/> 05/23/14	MS2006116-150V3	GenerateFASTQ	38	6	0	0	sbsuser

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The screenshot displays the Illumina MiSeq software interface. At the top, the window title is "Illumina MiSeq". Below the title bar, there is a "Manage Files" section with a green checkmark icon and a blue question mark icon. The interface is divided into several tabs: "Runs", "Sample Sheets", "Manifests", "Genomes", "Recipes", and "Bundle Logs". The "Bundle Logs" tab is active, showing a directory path of "D:\Illumina\MiSeqOutput" and a "Browse..." button. Below this, there are several radio buttons: "Bundle Logs" (highlighted with a yellow box), "Run QC Files" (highlighted with a yellow box), "Thumbnail Images", "Run Logs", and "Maintenance Logs". A table below lists several runs with columns for Date, Reagent Barcode, Workflow, R1, I1, R2, I2, and Username. The first row is selected with a checkmark in the "Date" column.

Date	Reagent Barcode	Workflow	R1	I1	R2	I2	Username
<input checked="" type="checkbox"/> 03/20/14	MS2006246-600V3	GenerateFASTQ	251	0	251	0	sbsuser
<input type="checkbox"/> 04/18/14	AA1234567-12345	Resequencing	5	0	0	0	sbsuser
<input type="checkbox"/> 05/01/14	MS2624960-300V2	Amplicon	151	8	151	8	sbsuser
<input type="checkbox"/> 05/15/14	MS2644485-300V2	Metagenomics	151	12	151	0	sbsuser
<input type="checkbox"/> 05/23/14	MS2006116-150V3	GenerateFASTQ	38	6	0	0	sbsuser

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