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


OK

Gene

The screenshot shows the BaseSpace ONSITE analysis report for the gene NEUROD1. It includes sections for 'Most Correlated Tissues', 'Most Correlated Diseases', 'Most Correlated Compounds', and 'Most Correlated Gene Perturbations'.

- Most Correlated Tissues:**
  - Cerebellar hemisphere
  - Cerebellar vermis
  - Cerebellum
  - Cerebellum peduncles
  - Pons
- Most Correlated Diseases:**
  - Diabetes mellitus type 1
  - Diabetes mellitus
  - Disorder of endocrine pancreas
  - Allergic disorder
  - Measles
- Most Correlated Compounds:**
  - Chronic Acid
  - Thioctic Acid
  - Tetanus Toxin
  - Sumatriptan
  - Cholera Toxin
- Most Correlated Gene Perturbations:**
  - NEUROD1
  - SETDB1
  - FBX1
  - SLC11A1
  - HR

Test ID	Gene	Locus	Status	log <sub>2</sub> (control FPKM)	log <sub>2</sub> (comparison FPKM)	log <sub>2</sub> (Ratio)	q Value	Significant
XLOC_000987	-	chr1:152020810-152021644	OK	-10	2.99	-12.99	0.164383	x
XLOC_001014	S100A9	chr1:153330329-153333503	OK	4.61	14.4	-9.79	0.582102	x
XLOC_001015	-	chr1:153359119-153359585	OK	-10	2.23	-12.23	0.22233	x
XLOC_001017	S100A1	chr1:153591275-153618799	OK	7.09	-10	17.09	0.264328	x
XLOC_001018	CHTOP	chr1:153591275-153618799	OK	4.99	4.67	0.32	0.954448	x
XLOC_001019	SNAPIN	chr1:153631120-153643504	OK	4.61	4.11	0.5	0.957225	x

<http://res.illumina.com/documents/products/technotes/technote-basespace-rna-seq.pdf>

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<http://seqanswers.com/wiki/Software/list>

<http://seqanswers.com/wiki/RNA-Seq>

The bioinformatics applications assigned the Biological domain **RNA-Seq** (topic\_3170 @) are tabulated below.

**Definition:**

- A topic concerning high-throughput sequencing of cDNA to measure the RNA content (transcriptome) of a sample, for example, to investigate how different alleles of a gene are expressed, detect post-transcriptional mutations or identify gene fusions.

**Synonyms:**

- WTSS
- Small RNA-Seq
- Whole transcriptome shotgun sequencing
- RNA-seq
- Small RNA-seq

Query returned 46 results.

	Biological domain	Bioinformatics method	Input format	Output format
ArrayExpressHTS	<b>RNA-Seq</b> RNA-Seq Quantitation		FASTQ	
Avadis NGS	ChIP-Seq DNA-Seq <b>RNA-Seq</b> Small RNA Pathway analysis	Alignment Quality Control Sequence analysis Visualization Biological Contextualization	SAM BAM BED ELAND FASTA FASTQ	
Chipster	ChIP-Seq <b>RNA-Seq</b> MiRNA-Seq MeDIP-Seq	QC Filtering Trimming Mapping Peak calling Motif detection Differential expression Pathway analysis Methylation analysis Genomic region matching Genome browser	FASTQ SAM BAM BED GTF	FASTQ SAM BAM BED GTF
	Genomics Whole Genome Resequencing	Mapping Assembly Alignment Colorspace	FASTA	FASTA FASTQ GFF

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[http://res.illumina.com/documents/products/datasheets/datasheet\\_rnaseq\\_analysis.pdf](http://res.illumina.com/documents/products/datasheets/datasheet_rnaseq_analysis.pdf)

Table 1: RNA-Seq Analysis Tools

Tool or Suite	Description	Availability	Link
Galaxy	Free form-based access to Bowtie, TopHat, and Cufflinks. Web browser client.	Academic/Open Source	galaxy.psu.edu
GenePattern	Free form-based access to Bowtie, TopHat, and Cufflinks. Local client.	Academic/Closed Source	www.broadinstitute.org/cancer/software/genepattern
Partek	Advanced statistics and interactive visualization for microarray and sequencing data.	Commercial	www.partek.com
CLC Bio	Software analyzing and visualizing sequencing data	Commercial	www.clcbio.com
GeoSpiza	Cloud-based analysis for microarray and sequencing data	Commercial	www.geospiza.com
GenomeQuest	Software for sequence data management	Commercial	www.genomequest.com
Avadis NGS	Software for sequence data analysis and management	Commercial	www.avadis-ngs.com
Ingenuity IPA	Software for biological pathway analysis	Commercial	www.ingenuity.com/products/pathways_analysis.html

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

BaseSpace®  
Genomics Cloud Computing

illuminatm

Sign up Log in

Now Available on the BaseSpace AppStore

RNA-Seq Workflow

Filters

|log2(ratio)|  
0.0 45.0

Significant  
Choose a value...

Status  
OK

Gene

Coverage (Normalized)

Normalized Position Along Transcript (5' to 3')

log2(comparison FPKM)

log2(control FPKM)

TopHat Alignment

Cufflinks Assembly & Differential Expression

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The screenshot shows the BaseSpace Cloud interface. The top navigation bar includes 'BaseSpace', 'Dashboard', 'Prep', 'Runs', 'Projects', 'Apps', 'Public Data' (circled in red), and 'Help'. Below the navigation bar, there are three data entries:

- > HiSeq 2500: TruSeq Stranded mRNA LT (SEQC: UHR & Brain)  
RNA-Seq
- > HiSeq 2000: TruSeq Stranded Total RNA (MAQC)  
RNA-Seq | Differential Expression
- ▼ NextSeq 500: RNA-Seq (8plex)  
NextSeq 500 data generated for reference human brain RNA and universal human reference RNA (UHRR). Libraries were prepared using TruSeq stranded mRNA or TruSeq Total RNA with Ribo-Zero Gold reagent kits.  
RNA-Seq

Below the 'NextSeq 500: RNA-Seq (8plex)' entry, there is a table with two rows and two columns:

Run	NextSeq 500: RNA-Seq (8plex) (49.41 GB)	Import
Project	NextSeq 500: RNA-Seq (8plex) 0	Import

The 'Import' button for the Project row is circled in red. A white arrow points from this button to the text below.

On the right side of the interface, there are two sections:

- Research Areas**
  - Cancer Research
  - Genetic Disease
  - Complex Disease
  - Microbial Research
- Categories**
  - Exome
  - Resequencing
  - Small RNA
  - Targeted Sequencing
  - De Novo Assembly
  - RNA-Seq (circled in red)
  - Gene Fusion Detection
  - ChIP-Seq
  - Methyl-Seq
  - Metagenomics
  - Tumor Normal
  - Variant Analysis
  - Differential Expression
  - Quality

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# Cufflink Assembly & DE 設定画面

App Session Name: Cufflinks Assembly & DE 05/16/2

Save Results To: Select Project(s):  
Transcriptome Demo

## TopHat Alignments Selection Criteria

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Adjust transcript assembly for samples without polyA selection:

## Comparison Group

Group Label: comparison

TopHat Alignment App Result(s):

Select App Result(s):

mRNA-Brain-C4	x
mRNA-Brain-C6	x

Adjust transcript assembly for samples without polyA selection:

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# Cufflink Assembly & DE 実行結果

Assembly <sup>i</sup>

	Control	Comparison	Merged
Gene Count	45,635	52,884	62,231
Transcript Count	87,392	96,262	119,726
Link to gene models	<a href="#">GTF result</a>	<a href="#">GTF result</a>	<a href="#">GTF result</a>
Relation to reference transcripts			
Equal (=)	44,934	45,437	46,093
Potentially novel (l)	18,903	20,155	31,453
Unknown, Intergenic (u)	21,848	28,622	39,115
Overlap with opposite-strand exon (x)	1,361	1,674	2,459
Other	346	374	606

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Differential Expression <sup>i</sup>

Gene Count	62,225
$\Delta$ Gene Count	27,052
Transcript Count	119,696
$\Delta$ Transcript Count	27,379
CuffDiff results	<a href="#">differential gene expression, differential transcript expression</a>

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
## Differential Expression Gene Browser

### Filters

**||log<sub>2</sub>(ratio)|**  
 2.0  23.0


### Significant

true ▾

### Status

OK ▾

### Gene

Save Plot as SVG

### Differential Expression Gene Browser

The Differential Expression Gene Browser can be used to filter and plot the differential expression results dynamically (Table 5). The data can be sorted by clicking any column heading and can be saved as either an SVG graphic or CSV table.

Table 5: Differential Expression Gene Browser Filter Options

Filter	Description
Log Ratio Cutoff	Allows the differentially expressed gene table to be filtered based on the observed log ratio of differential expression.
Significance	Filters for results achieving statistical significance at a q-value < 0.05.
Status	Filters based on the reliability of the result. If the result passes all stringency filters, the status is returned as "OK". If the result is flagged due to one of several factors (such as insufficient read depth), a failure mode of NOTEST, LOWDATA, or FAIL is returned.

Additional information is available at <http://cufflinks.cbcb.umd.edu/faq.html#notest>.

Test ID	Gene	Locus	Status	log <sub>2</sub> (Brain FPKM)	log <sub>2</sub> (UHR FPKM)	log <sub>2</sub> (Ratio)	q Value	Significant
A1CF	A1CF	chr10:52559168-52645435	OK	-3.860	1.410	-5.270	0.000	✓
A2ML1	A2ML1	chr12:89:					0.000	✓
AACSP1	AACSP1	chr5:178:					0.000	✓
AAK1	AAK1	chr2:698:					0.000	✓

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[www.illumina.com/content/dam/illumina-marketing/documents/products/technotes/technote-basespace-rna-seq.pdf](http://www.illumina.com/content/dam/illumina-marketing/documents/products/technotes/technote-basespace-rna-seq.pdf)

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
[www.illumina.com/informatics/research/sequencing-data-analysis-management/rna-seq-data-analysis.html](http://www.illumina.com/informatics/research/sequencing-data-analysis-management/rna-seq-data-analysis.html)

**RNA-Seq BaseSpace Apps: A Guided Tour**

See step-by-step instructions on how to navigate through the data analysis.

[View Video](#)

## RNA-Seq BaseSpace Apps: A Guided Tour

Test ID	Gene	Locus	Status	log <sub>2</sub> (uhr FPKM)	log <sub>2</sub> (brain FPKM)	log <sub>2</sub> (Ratio)	q Value	Significant
XLOC_000010	-	chr1:797247-799101	OK	-10.000	-0.740	-9.260	0.000	✓
XLOC_000011	-	chr1:800366-801218	OK	-10.000	-0.630	-9.370	0.000	✓
XLOC_000014	-	chr1:844861-845337	OK	-0.540	-10.000	9.460	0.001	✓
XLOC_000028	-	chr1:1100832-1101478	OK	-0.920	-10.000	9.080	0.000	✓
XLOC_000043	-	chr1:1314123-1314431	OK	0.040	-10.000	10.040	0.003	✓
XLOC_000062	GABRD	chr1:1950767-1962192	OK	-1.350	6.950	-8.300	0.000	✓
XLOC_000068	PLCH2	chr1:2398901-2439211	OK	-0.690	4.460	-5.150	0.000	✓
XLOC_000069	-	chr1:2462986-2463443	OK	-10.000	-0.700	-9.300	0.002	✓
XLOC_000070	-	chr1:2469531-2469979	OK	-10.000	-0.870	-9.130	0.002	✓
XLOC_000071	-	chr1:2470381-2470659	OK	-10.000	0.630	-10.630	0.001	✓
XLOC_000072	-	chr1:2472405-2473064	OK	-10.000	-0.980	-9.020	0.000	✓
XLOC_000073	-	chr1:2475563-2477232	OK	-10.000	-0.540	-9.460	0.000	✓
XLOC_000078	-	chr1:2495662-2500317	OK	-0.370	-10.000	9.630	0.000	✓
XLOC_000085	-	chr1:4036939-	OK	3.320	-4.210	-7.530	0.030	✓

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The screenshot shows the BaseSpace interface for a 'TopHat Alignment' job. The 'Analysis Info' tab is selected and highlighted with a red circle. The job details are as follows:

Name	TopHat Alignment 10/07/2014 1:57:42
Application	TopHat Alignment   Version: 1.0.0
Date started	Tuesday, October 7th 2014, 8:58:04 pm
Date completed	Wednesday, October 8th 2014, 7:15:18 pm
Duration	22 hours 17 minutes 14 seconds
Session Type	Multi-Node
Size	81.95 GB
Status	Complete (8 Nodes Complete)

At the bottom of the page, there is a 'Log Files' link and a note: 'Please view the Multi-Node details page to see logs for this analysis.'

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
mRNA ***	Differential Expression	10-20M Reads *
	Allele Specific Expression	50 – 100M Reads *
	Splice variation **	50 - 100M Reads *
	Complete Annotation	100M - 1B Reads *
	Transcript Based Assembly	50 – 200M Reads *

76bp~ x PE

- Based on human sized transcriptomes
- \*\* Also applies to RNA fusion transcripts in cancer
- \*\*\* Applies to poly A-selected libraries
  - Ribo-Zero, high quality RNA libraries: requires ~ 2X more reads
  - Ribo-Zero, FFPE libraries: requires ~ 4X more reads

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**License Info**  
[Privacy](#)  
[EULA](#)

**Description**

The TopHat Alignment workflow performs the following functions

- Read mapping using the TopHat 2 aligner
- FPKM estimation of reference genes and transcripts using Cufflinks 2
- Variant calling (SNVs and small indels) with the Isaac Variant caller
- Optional fusion calling with TopHat-Fusion

Available reference genomes include

- Homo sapiens UCSC hg19 (RefSeq & Gencode gene annotations)
- Mus musculus UCSC mm10 (RefSeq gene annotations)
- Rattus norvegicus UCSC m5 (RefSeq gene annotations)

Current limitations:

- Minimum and maximum of 100,000 and 400,000,000 reads per-sample, respectively
- Maximum of 2,000,000,000 read across all samples in a single analysis
- Minimum Read length 35bp and Maximum Read Length 125bp
- Paired-end reads are required for fusion detection

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




High Reproducibility  
with Low Inputs

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

# Upcoming Webinar

- ▶ **"Next generation tools for gene expression profiling and gene fusion detection in FFPE tumor samples"**
  - Gary P. Schroth Ph.D., Distinguished Scientist, Illumina
- ▶ **10月 13日 (火) 13:00~**
- ▶ **Register NOW!**
  - <https://illumina.webex.com/illumina/onstage/g.php?2f9bc840e5ddbfc22fcb31d25>

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- Accurate whole human genome sequencing using reversible terminator chemistry.  
Nature 456: 53-59 [PMID: 18987734]

## ▶ RNA-Seq, RPKM

- Mapping and quantifying mammalian transcriptomes by RNA-Seq  
Nature Methods, Volume 5, 621 – 628 [PMID: 18516045]

## ▶ BaseSpace RNA-Seq アプリ

- TopHat: discovering splice junctions with RNA-Seq.  
Bioinformatics. 2009 May 1;25(9):1105-11 [PMID: 19289445]
- Cufflinks/Cuffdiff  
Transcript assembly and quantification by RNA-Seq reveals unannotated transcripts and isoform switching during cell differentiation. Nat Biotechnol. 2010 May;28(5):511-5 [PMID: 20436464]
- Cuffdiff2 Differential analysis of gene regulation at transcript resolution with RNA-seq  
Nature Biotechnology 31, 46–53 (2013) [PMID: 23222703]
- TopHat-Fusion: an algorithm for discovery of novel fusion transcripts.  
Genome Biol. 2011 Aug 11;12(8):R72. [PMID: 21835007]
- STAR: ultrafast universal RNA-seq aligner.  
Bioinformatics\_ 2013 Jan 1;29(1):15-21. [PMID: 23104886]
- DESeq  
Differential expression analysis for sequence count data.  
Genome Biol. 2010;11(10):R106. [PMID: 20979621]  
DESeq2: [www.bioconductor.org/packages/2.13/bioc/html/DESeq2.html](http://www.bioconductor.org/packages/2.13/bioc/html/DESeq2.html) <http://www.ncbi.nlm.nih.gov/pubmed>

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<http://www.illumina.com/landing/basespace-core-apps-for-rna-sequencing/>  
<http://res.illumina.com/documents/products/technotes/technote-basespace-rna-seq.pdf>  
[http://support.illumina.com/help/BaseSpace\\_App\\_RNAseq\\_help/RNAseq\\_Apps\\_Help.htm](http://support.illumina.com/help/BaseSpace_App_RNAseq_help/RNAseq_Apps_Help.htm)  
<http://www.illumina.com/applications/sequencing/rna.ilmn>

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