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© 2015 Illumina, Inc. All rights reserved.
Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDX, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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
BaseSpaceとは

Seamless upload, monitoring, storage, analysis, and collaboration

The Future: Cloud Solutions will be everywhere and specifically valuable to researchers

BaseSpaceとは

The Future: Cloud Solutions will be everywhere and specifically valuable to researchers

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Download

BaseSpace utilizes a desktop client that allows for fast, reliable downloads transferred securely over SSL.

Install the BaseSpace downloader

The downloader requires a one-time installation. Click to install the downloader on this machine. If the downloader is already installed, you can skip this step.

[Install the Downloader](#) [Don't show me this again](#)

Select files to be downloaded

- All FASTQ files for this run ← アクティブにならない
- SAV file [Learn more about SAV](#)
- No analysis files are available for this run

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BaseSpace® Dashboard Prep Runs Projects Apps Public Data Help

Projects : PhiX 2x250 : Samples List : PhiX

PhiX

Launch App

Sample name	PhiX
Date created	08/20/2013 4:36:41
Genome	Phix - Illumina
Paired end	Paired
Number of Clusters	18,306,599
Number of Reads	36,613,198
Read 1 length	251
Read 2 length	251
Size	7.20 GB

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
- Confirming that the file name conforms to Illumina standards:
SampleName_SampleNumber_Lane_Read_FlowCellIndex.fastq.gz
SampleName_S#_L00#_R#_00#.fastq.gz

https://support.illumina.com/content/dam/illumina-support/documents/documentation/software_documentation/basespace/basespace-user-guide-15044182-e.pdf

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BaseSpace appsとは

<イルミナコアアプリ>

				
16S Metagenomics	TopHat Alignment	Cufflinks Assembly & DE	RNA Express	Variant Studio
				
BWA Enrichment	Isaac Enrichment	Broad IGV	TruSeq Amplicon	Amplicon DS
				
BWA WGS	Isaac WGS	Tumor Normal	Long Read Assembly	Long Read Phasing

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TopHat Alignment v1.0.0

Illumina, Inc.

Please review the validation errors below

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[Privacy](#)
[EULA](#)

Description

The TopHat Alignment workflow performs the following functions

- Read mapping using the TopHat 2 aligner
- FPKM estimation of reference genes and transcripts using Cufflinks 2
- Variant calling (SNVs and small indels) with the Isaac Variant caller
- Optional fusion calling with TopHat-Fusion

Available reference genomes include

- Homo sapiens UCSC hg19 (RefSeq & Gencode gene annotations)
- Mus musculus UCSC mm10 (RefSeq gene annotations)
- Rattus norvegicus UCSC rn5 (RefSeq gene annotations)

Current limitations:

- Minimum and maximum of 100,000 and 400,000,000 reads per-sample, respectively
- Maximum of 2,000,000,000 read across all samples in a single analysis
- Minimum Read length 35bp and Maximum Read Length 125bp
- Paired-end reads are required for fusion detection

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
BaseSpace® HelpCenter [Sign in](#)

Welcome to our HelpCenter

Learn how to unlock the potential of Illumina's next-generation informatics ecosystem.

- New to BaseSpace?**
Get started and Explore BaseSpace >
- Articles**
Everything in one place; All our articles >
- BaseSpace Developers**
Start building the next generation of sequencing apps >

What's New

- Nov 2015 | [New and Updated Illumina Core Apps](#)
- Oct 2015 | [New BaseSpace Features and Tiers](#)
- Oct 2015 | [Interpret Protein-Coding and Splicing Variants with MutationForecaster](#)


Popular Articles

- [How to share data using Email option](#)
- [How to transfer ownership](#)
- [How to delete/restore data](#)
- [How to combine samples](#)

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The screenshot shows the BaseSpace Blog homepage. The header includes the logo 'BaseSpace Blog by illumina', navigation links for 'About' and 'Press', and a 'Subscribe to RSS' button. The main content area features a 'Featured' section with a post titled 'New BaseSpace Features and Tiers' by PrithwishPal. The post text discusses the company's commitment to simplifying the informatics experience and the addition of new features and subscription tiers. A search bar and a 'RECENT POSTS' section are visible on the right side of the page.

BaseSpace Blog
by illumina

About Press  Subscribe to RSS

Featured

by PrithwishPal

in Announcement, Features, Releases

Leave a comment

New BaseSpace Features and Tiers

We launched BaseSpace in 2011 with a commitment towards 'simplifying the informatics experience', and we have delivered via frequent feature releases, feature improvements, and a growing ecosystem of Informatics Applications (Apps). Currently we have more than 60 BaseSpace Apps covering a wide range of bioinformatics analysis needs from alignment and variant calling to metagenomics, and even proteomics. Users can also run their own analysis pipelines on BaseSpace as private Apps. With the recent NeoPrep integration, users can set up a library prep, monitor sequencing while uploading data directly from the instrument, and analyze the data for a truly end-to-end experience, all within the secure BaseSpace cloud environment.

In the last several years, Illumina has continued to innovate, expanding its sequencing systems portfolio to offer a wide range of solutions from the focused power of the MiSeq to the population power of the HiSeq X. The rapidly growing adoption of these systems and their varied applications have resulted in increasingly diverse customer needs in terms of bioinformatics support. To serve these diverse needs, we are excited to announce the addition of a number of new features to the BaseSpace platform. These features have been packaged into two new subscription levels based on use cases: Professional and Enterprise.

The basic BaseSpace account will be available for free, and will include 1TB storage space. With the basic subscription, users can monitor sequencing runs, stream sequencing data directly to the cloud, store and share data, and access all Apps. The early access Linux-based command-line interface, [BaseMount](#), will continue to be

SEARCH

RECENT POSTS

- New and Updated Illumina Core Apps
- New BaseSpace Features and Tiers
- Interpret Protein-Coding and Splicing Variants with MutationForecaster
- VCAT 2.3 with Haplotype Compare and FASTQ Toolkit 2.0
- BaseMount: Directly Linking NGS Data to R Packages for RNA-Seq Differential Expression Analyses

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


bcl conversionとは

Analysis Type	Software	Outputs
<p style="color: blue; font-size: 1.2em;">Sequencing</p>	 <p style="text-align: center;">ICS/RTA</p>	 <p style="text-align: center;">Images/TIFF files</p>
<p style="color: green; font-size: 1.2em;">Primary Analysis</p>	 <p style="text-align: center;">ICS/RTA</p>	<div style="display: flex; justify-content: space-around;"> <div data-bbox="1304 635 1516 806">  </div> <div data-bbox="1545 635 1767 806">  </div> </div> <p style="text-align: center;"> Intensities Base Calling bclファイル </p>
<p style="color: purple; font-size: 1.2em;">Secondary Analysis</p>	<div style="display: flex; justify-content: space-around; align-items: center;"> <div data-bbox="691 928 975 1106">  </div> <div data-bbox="1014 978 1197 1149">  <p>HiSeq Analysis Software</p> </div> </div> <div style="display: flex; justify-content: space-around; align-items: center; margin-top: 20px;"> <div data-bbox="753 1120 927 1256">  <p>CASAVA</p> </div> </div>	 <p style="text-align: center;"> Alignments and Variant Detection </p>

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bcl2fastq2 Conversion Software

Guide

<code>--minimum-trimmed-read-length</code>	<p>Minimum read length after adapter trimming. bcl2fastq trims the adapter from the read down to the value of this parameter. If there is more adapter match below this value, then those bases are masked, not trimmed (replaced by N rather than removed).</p> <p><u>Default: 35</u></p>
<code>--mask-short-adapter-reads</code>	<p>This option applies when a read is trimmed to below the length specified by the <code>--minimum-trimmed-read-length</code> option (default of 35). These parameters specify the following behavior:</p> <p>If the number of bases left after adapter trimming is less than <code>--minimum-trimmed-read-length</code>, force the read length to be equal to <code>--minimum-trimmed-read-length</code> by masking adapter bases (replace with Ns) that fall below this length.</p> <p>If the number of ACGT bases left after this process falls below <code>--mask-short-adapter-reads</code>, mask all bases, resulting in a read with <code>--minimum-trimmed-read-length</code> number of Ns.</p> <p><u>Default: 22</u></p>

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bcl2fastq2 Conversion Software Guide

`--minimum-trimmed-read-length`

Minimum read length after adapter trimming. bcl2fastq trims the adapter from the read down to the value of this parameter. If there is more adapter match below this value, then those bases are masked, not trimmed (replaced by N rather than removed).

Default: 35

`--mask-short-adapter-reads`

This option applies when a read is trimmed to below the length specified by the `--minimum-trimmed-read-length` option (default of 35). These parameters specify the following behavior:

If the number of bases left after adapter trimming is less than `--minimum-trimmed-`

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If the number of ACGT bases left after this process falls below `--mask-short-adapter-reads`, mask all bases, resulting in a read with `--minimum-trimmed-read-length` number of Ns.

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Examples)

- Failed to read number of clusters from .../s_1_1101.bcl.gz
- Failed to read BCL file .../s_4_1103.bcl.gz
- Error: 2015-Jul-09 08:46:18: Invalid argument: /usr/.../BclReader.cpp(115)
: Throw in function unsigned int casava::alignment::BclReader::readClusterCount()'

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https://support.illumina.com/downloads/basespace_core_apps_user_guides.html

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