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© 2014 Illumina, Inc. All rights reserved. Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeqFGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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② Experiment Managerを選択

The screenshot shows the Illumina website's navigation menu. The 'サポート' (Support) tab is selected. Below the navigation bar, the breadcrumb trail reads 'サポート >> シーケンス >> ソフトウェア'. The main content area is a grid of software tools categorized by 'STAGE' and 'SOFTWARE'. The 'EXPERIMENTAL DESIGN' stage is highlighted, and 'Experiment Manager' is circled in red with a mouse cursor pointing to it.

STAGE	EXPERIMENTAL DESIGN	CLUSTER GENERATION	INSTRUMENT CONTROL IMAGE ANALYSIS	BASE CALLING Q-SCORING	FILE CONVERSION	DEMULTIPLEXING	ALIGNMENT	VARIANT CALLING VARIANT COUNTING	VISUALIZATION
SOFTWARE	Design Studio								
	Experiment Manager								
	blueFuse Workflow Manager								
		cBot Control							
			MiSeq Control/RTA						
			NextSeq Control/RTA2						
				HiSeq Control/RTA					
				GA Sequencing Control RTA					
				Sequencing Analysis Viewer					
					Off-Line Basecaller				
						Analysis Visual Controller			
							MiSeq Reporter		
								BaseSpace	
								BaseSpace Onsite	
								CASAVA	
							HiSeq Analysis Software		
								iGenomes	

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illumina®

Log in to get personalized account information. Quick Order MyIllumina Tools

Contact Us MyIllumina Tools

APPLICATIONS SYSTEMS INFORMATICS CLINICAL SERVICES SCIENCE SUPPORT COMPANY

Support » Sequencing » Sequencing Software » Experiment Manager

Experiment Manager Support

- Overview
- Computing Requirements
- Downloads**
- Documentation & Literature
- Training

Experiment Manager

Latest Updates

- [Illumina Experiment Manager Guide \(15031335\)](#) 06/19/2014
- [Illumina Experiment Manager v1.8](#) 06/19/2014
- [IEM TruSeq DNA, RNA, or ChIP Quick Reference Card](#) 04/30/2014

User Guide

- [Illumina Experiment Manager Guide \(15031335\)](#)

The Illumina Experiment Manager software helps you create and edit well-formed sample sheets for Illumina sequencers and analysis software.

Illumina recommends using the Experiment Manager before starting sample or library preparation. The Experiment Manager can detect and warn of sub-optimal index combinations. By creating the sample sheet prior to sample or library preparation, you can try a different index combination without risking your samples.

You can use the Illumina Experiment Manager to create sample sheets for any Illumina sequencer and for any Nextera or TruSeq libraries.

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illumina®

Log in to get personalized account information. Quick Order View Cart

Contact Us MyIllumina Tools

APPLICATIONS SYSTEMS INFORMATICS CLINICAL SERVICES SCIENCE SUPPORT COMPANY

Support » Sequencing » Sequencing Software » Experiment Manager » Downloads

Experiment Manager Support

- Overview
- Computing Requirements
- Downloads
- Documentation & Literature
- Training

Downloads

DESCRIPTION	FILE INFO	DATE
Illumina Experiment Manager v1.8 The Illumina Experiment Manager v1.8 is a desktop tool that guides you in building library plates and creating the sample sheet used for sequencing on the MiSeq, NextSeq, or HiSeq systems. Minimum System Requirements: Windows Vista or 7 Operating System; 32-bit system, 4 GB RAM or 64-bit system, 8 GB RAM; Microsoft.Net Framework 4.0 or above; Microsoft Office Excel (recommended).	ZIP (2 MB)	08/19/2014
Illumina Experiment Manager v1.8 Release Notes	PDF (< 1 MB)	08/19/2014

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- Support>Sequencing Supports>Downloads
>Illumina Experiment Manager v1.8

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EM Illumina Experiment Manager

Illumina Experiment Manager

Sample Sheet Wizard - Sample Selection

Samples to include in sample sheet * - required field Maximize

Sample ID*	Sample Name	Plate	Well	Index1 (I7)*	I7 Sequence	Index2 (I5)*	I5 Sequence	Genome Folder*
1	sample1							
2	sample2							

N701
N702
N703
N704
N705
N706
N707
N708
N709
N710
N711
N712

S501
S502
S503
S504
S505
S506
S507
S508
S517

Add

Sample Sheet Status: Invalid
Reason: Not all samples in this sample sheet have all the required fields

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EM Illumina Experiment Manager

Illumina Experiment Manager

Sample Sheet Wizard - Sample Selection

Samples to include in sample sheet * - required field Maximize

Sample ID*	Sample Name	Plate	Well	Index 1 (I7)*	I7 Sequence	Index 2 (I5)*	I5 Sequence	Genome Folder*
1	Sample1							Arabidopsis_thaliana#NCBI#build9.1#Sequence#WholeGenomeFasta
2	Sample2							Bos_taurus#Ensembl#UMD3.1#Sequence#WholeGenomeFasta

Escherichia_coli_K_12_DH10B#NCBI#2008-03-17#Sequence#WholeGenomeFasta
Homo_sapiens#UCSC#hg19#Sequence#WholeGenomeFasta
Mus_musculus#UCSC#mm9#Sequence#WholeGenomeFasta
PhiX#Illumina#RTA#Sequence#WholeGenomeFasta
Rattus_norvegicus#UCSC#rn4#Sequence#WholeGenomeFasta
Saccharomyces_cerevisiae#UCSC#sacCer2#Sequence#WholeGenomeFasta
Staphylococcus_aureus_NCTC_8325#NCBI#2006-02-16#Sequence#WholeGenomeFasta

Sample Sheet Status: **Invalid**
Reason: **Not all samples in this sample sheet have all the required fields**

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The screenshot shows the 'Illumina MiSeq' software window. At the top, there is a green checkmark icon and a question mark icon. Below them is a text box that says 'Please make any needed changes to the Folders.' The main area contains four rows of folder selection fields, each with a 'Browse...' button to its right:

- Recipe Folder: D:\Illumina\MiSeq Control Software\CustomRecipe
- Sample Sheet Folder: D:\Illumina\MiSeq Control Software\SampleSheets (highlighted with a red box)
- Manifest Folder: D:\Illumina\MiSeq Control Software\Manifests
- Output Folder: D:\Illumina\MiSeqOutput

At the bottom of the window, there are two large buttons: 'Back' on the left and 'Save and Continue' on the right. The Windows taskbar is visible at the very bottom, showing system icons and a temperature display.

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The screenshot shows a dialog box titled 'Illumina Experiment Manager' with a subtitle 'Sample Sheet Wizard - Review'. The main text inside the dialog asks: 'Would you like to view your sample sheet in Excel?'. At the bottom of the dialog, there are two buttons: 'Yes' and 'No'.

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[Header]												
IEMFileVersion									4			
Date									#####			
Workflow									Resequencing			
Application									Resequencing			
Assay									Nextera XT			
Description									Amplicon			
Chemistry												
[Reads]												
									151			
									151			
[Settings]												
FlagPCRDuplicates									1			
ReverseComplement									0			
VariantFilterQualityCutoff									30			
outputgenomevcf									FALSE			
Adapter									CTGTCTCTTATACATCT			
[Data]												
Sample_ID	Sample_Na	Sample_Pla	Sample_We	I7_Index_ID	index	I5_Index_ID	index2	GenomeFo	Sample_Prc	Description		
1	Sample1			N701	TAAGGCG.S501	TAGATCG		Homo_sapiens	UCSC	hg19	Sequence	WholeGenomeFasta
2	Sample2			N702	CGTACTA(S502	CTCTCTA		Homo_sapiens	UCSC	hg19	Sequence	WholeGenomeFasta

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Sample Sheet Status: **Warning**

Reason: There is not enough diversity in the index cycles. Every cycle in the index reads must have at least one base in the green channel (G,T) and one in the red channel (A,C) to ensure proper image processing

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Nextera Low Pooling Guidelines

http://support.illumina.com/content/dam/illumina-marketing/documents/products/technotes/technote_nextera_low_plex_pooling_guidelines.pdf

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
Additional Information

- ▶ **Illumina Experiment Manager Guide**

http://support.illumina.com/sequencing/sequencing_software/experiment_manager.html

- ▶ **MiSeq Sample Sheet Quick Reference Guide**

http://support.illumina.com/downloads/miseq_sample_sheet_quick_reference_guide_15028392.html

- ▶ **MiSeq System User Guide**

http://support.illumina.com/downloads/miseq_system_user_guide_15027617.html

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Illumina Experiment Manager

Sample Sheet Wizard - Instrument Selection

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
Eri Kibukawa
Bioinformatics Support Scientist
Illumina APAC

© 2013 Illumina, Inc. All rights reserved.
Biology, Biotechnology, Bioinformatics, Bioinformatics, Life Sciences, Computing, Data, Genomics, Energy, Genome Analysis, Genomics, Health, Microarrays, Sequencing, Systems, and Life Sciences. Illumina, Inc. and the Illumina logo are trademarks of Illumina, Inc. All other marks and names contained herein are the property of their respective owners.

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