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© 2015 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDX, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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
\* Sequence data performed in duplicate (50 ng) and Triplicate (25 ng)

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

Extraction Inhibitors	Inhibitory Effect	Likely Source	Methods to Minimize Inhibition
EDTA	Chelation of metal ions	TE buffer	Reduce EDTA concentration in TE buffer or simply use Tris-HCl (10 mM) or water as elution buffer
Alcohols	Enzyme denaturation	Ethanol, isopropanol, isoamyl alcohol	Dry pellet and resuspend or use silica-based purification
Excess Salts	Template blocking	KC <sub>l</sub> , NaCl, CsCl, NaAc	Wash with 70% ethanol or use silica-based purification
Chaotropic Salts	Enzyme denaturation	Guanidinium chloride; magnesium chloride, urea	Wash with 70% ethanol or use silica-based purification
Phenol:Chloroform	Enzyme denaturation	Organic carryover	Use PVP, PVP/ammonium acetate, incorporation of 1.2% citric acid at the DNA extraction step
Detergents/DDT	Enzyme denaturation	Sodium deoxycholate, SDS, Tween 20	Wash with 70% ethanol
Exogenous DNA/RNA	Template competition	Carryover	DNase I for DNA removal; RNase A for RNA, RNA:DNA hybrid removal
Carriers	Template competition/blocking	RNA, heparin, glycogen	Only use carriers that do not serve as template or block the template such as linear acrylamide, N- or P-carriers
Agarose	Template blocking	Gel extractions	Use a spin column with chaotropic salt buffer; dialysis
Excess Metal Ions	Reduce oligo specificity	Mg <sup>++</sup> from PCR buffer	Dialysis against PBS (pH 7.4), phenol: chloroform extraction followed by EtOH precipitation

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