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© 2014 Illumina, Inc. All rights reserved. Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeqFGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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| Nextera | Nextera XT |
|---|--|
| Large / complex genomes | Small genomes, amplicons, plasmids |
| Human genomes | PCR Amplicons (> 300 bp)* |
| non-human mammalian genomes (e.g. mouse, rat, bovine) | Plasmids |
| Plant genomes (e.g. arabidopsis, maize, rice) | Microbial Genomes (e.g. Prokaryotes, archea) |
| Invertebrates genomes (e.g. Drosophila) | Concatenated Amplicons |
| | double-stranded cDNA |

* Illumina recommends > 300 bp to ensure even coverage across the length of the DNA fragment. An expected drop off in sequencing coverage about 50 bp from each distal end of a fragment may be seen. This is because the tagmentation reaction cannot add an adapter right at the distal end of a fragment. For PCR amplicon sequencing this can be easily averted by simply designing your amplicons to be ~100 bases larger than the desired insert to be sequenced.

GTATCAATAGATTAATTGATCCACTGATTCAGCCTAACGGANDTATCAATDAGACTAAATAATACGGTATCCATTAAGAGCCTGACGGTCCAGCAGCGAAGGAATGATAACGATACACACTCTCTTTAAGCCTTAAGCGAAGTATCACCTTAAGATTACT
ACTTACACACTTCTGTTAAGGHTAAGATTAATTGATCCACTGATTCAGCCTAACGGANDTATCAATDAGACTAAATAATACGGTATCCATTAAGAGCCTGACGGTCCAGCAGCGAAGGAATGATAACGATACACACTCTCTTTAAGCCTTAAGCGAAGTATCACCTTAAGATTACT
CCTTGCACACCTGAGCAGCTTGTGTAACTTAAAGATTAATTGATCCACTGATTCAGCCTAACGGANDTATCAATDAGACTAAATAATACGGTATCCATTAAGAGCCTGACGGTCCAGCAGCGAAGGAATGATAACGATACACACTCTCTTTAAGCCTTAAGCGAAGTATCACCTTAAGATTACT
CTTCTAAGCCTTAAGATTAATTGATCCACTGATTCAGCCTAACGGANDTATCAATDAGACTAAATAATACGGTATCCATTAAGAGCCTGACGGTCCAGCAGCGAAGGAATGATAACGATACACACTCTCTTTAAGCCTTAAGCGAAGTATCACCTTAAGATTACT
CAAGCCTTAAGATTAATTGATCCACTGATTCAGCCTAACGGANDTATCAATDAGACTAAATAATACGGTATCCATTAAGAGCCTGACGGTCCAGCAGCGAAGGAATGATAACGATACACACTCTCTTTAAGCCTTAAGCGAAGTATCACCTTAAGATTACT
ATAAGCCTTAAGATTAATTGATCCACTGATTCAGCCTAACGGANDTATCAATDAGACTAAATAATACGGTATCCATTAAGAGCCTGACGGTCCAGCAGCGAAGGAATGATAACGATACACACTCTCTTTAAGCCTTAAGCGAAGTATCACCTTAAGATTACT
GTATCAATAGATTAATTGATCCACTGATTCAGCCTAACGGANDTATCAATDAGACTAAATAATACGGTATCCATTAAGAGCCTGACGGTCCAGCAGCGAAGGAATGATAACGATACACACTCTCTTTAAGCCTTAAGCGAAGTATCACCTTAAGATTACT

Nextera XT DNA sample prepに必要なキット

▶ Nextera XT DNA Sample Prep kit

- 24 samples FC-131-1024
- 96 samples FC-131-1096

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▶ Nextera XT Index Kit

- 24 Indices, 96 samples FC-131-1001
- 96 Indices, 384 samples FC-131-1002

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▶ TruSeq Dual Index Sequencing Primers

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- ▶ Nextera DNA Sample Prep kit
 - 24 samples FC-121-1030
 - 96 samples FC-121-1031

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- ▶ Nextera DNA Index Kit
 - 24 Indices, 96 samples FC-121-1011
 - 96 Indices, 384 samples FC-121-1012

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- ▶ TruSeq Dual Index Sequencing Primers
 - Single Read Run FC-121-1003
 - Paired End Run PE-121-1003

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| Carryover | Inhibitory Effect | Likely Source | Methods to minimize inhibition |
|---------------------|---------------------------------------|---|---|
| EDTA | Chelation of metal ions | TE buffer | Reduce the concentration of EDTA to 0.1 mM in the TE buffer or simply use Tris-HCl (10 mM) or water to bring DNA in solution. |
| Alcohols | Enzyme precipitation and denaturation | Ethanol, isopropanol, isoamyl alcohol | Dry pellet and resuspend |
| Excess Salts | Template blocking | KCl, NaCl, CsCl, NaAc | Wash with 70% ethanol or use silica-based purification |
| Chaotropic salts | Enzyme denaturation | Guanidinium chloride; magnesium chloride; urea | Wash with 70% ethanol or use silica-based purification |
| Phenol:Chloroform | Enzyme denaturation | Organic carryover | Use PVP, PVP/ammonium acetate. Incorporation of 1.2% citric acid at the DNA extraction step |
| Polysaccharides | Template blocking | Dextrans, glycans | High salt precipitation, CTAB buffer, chloroform extraction, pectinase, cellulase, hemicellulase and/or α -amylase digestion |
| Protein | Template blocking | BSA, immunoglobins, PEG | Use SDS, CTAB or guanidinium buffers, proteinase K |
| Fat | Template blocking | Glycerol | Lipase or hexane treatment and chloroform extraction |
| Detergents / DDT | Enzyme precipitation and denaturation | Sodium deoxycholate, sarkosyl, SDS, Nonidet P-40, Tween 20, Triton X-100, N-octyl glucoside | Wash with 70% ethanol |
| Proteases | Protein degradation | Proteinase K, sample handling | Phenol:chloroform extraction followed by EtOH precipitation |
| Nucleases | Template degradation | Sample handling, restriction enzymes, micrococcal nuclease, S1 | Use B-ME, EGTA or SDS during protein precipitation |
| Exogenous DNA / RNA | Template competition | Carryover | DNase I for DNA removal; RNase A for RNA, RNA:DNA hybrid removal |
| Carriers | Template competition and/or blocking | RNA, heparin, glycogen | Only use carriers that do not serve as template or block the template such as linear acrylamide, N- or P-carriers |
| Agarose | Template blocking | Gel extractions | Use a spic column with chaotropic salt buffer; dialysis |
| Excess metal ions | Reduce oligo specificity | Mg ⁺⁺ from PCR buffer | Dialysis against PBS (pH 7.4), phenol:chloroform extraction followed by EtOH precipitation |

<https://icom.illumina.com/MyIllumina/Bulletin/QqttpskaC0ywyI9co83qXg/considerations-for-dna-isolation-when-using-nexter>

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Technical Note: DNA Sequencing illumina®

Nextera® Library Validation and Cluster Density Optimization

Guidelines for generating high-quality data with Nextera sample preparation kits.

Introduction

Each kit in the Nextera family of sample preparation (DNA, XT DNA, Custom Enrichment, and Exome Enrichment) requires different considerations for final library validation and optimization of cluster density. Because each Nextera protocol is designed for a distinct application, it is important to tailor the library validation and clustering processes to each specific Nextera kit.

The Nextera "tagmentation" reaction simultaneously fragments and tags DNA with adapters. As with other enzymatic fragmentation protocols, the amount of input DNA can affect the outcome of the assay. Using too little DNA may result in libraries with shorter insert sizes, and conversely using too much DNA will result in libraries with larger insert sizes. This technical note describes best practices for optimizing the final libraries generated with each Nextera kit.

Starting Material

Nextera kits are ideal for high-throughput studies, providing fast sample preparation methods with low input requirements. Nextera DNA and Enrichment protocols are optimized for 50 ng of high-quality genomic DNA, and the Nextera XT DNA kit requires 1 ng of small genome or amplicon DNA. To ensure that nucleic acids aside from double-stranded DNA do not interfere with quantitation, Qubit² or PicoGreen² assays are recommended to quantify the input genomic DNA. To ensure accuracy, Illumina recommends quantifying a concentrated solution of each sample and subsequently diluting to the input concentration specified in each sample preparation user guide.

Organic contaminants (including ethanol) have been shown to interfere with the Nextera tagmentation reaction, and should be eliminated from the input DNA prior to quantitation using standard methods. EDTA in DNA elution or dilution buffers can also interfere with the tagmentation reaction, and should be avoided.

Figure 1: Distribution of Final Library Sizes

A

B

The final library traces for Nextera DNA libraries were generated using a BioAnalyzer. These example size distributions are also applicable to Nextera XT libraries that have not been normalized. To optimize clustering concentrations based on average library size of 500 bp, the sample shown in Panel A requires the conversion factor 1 ng/μL = 3 nM. Clustering concentrations for the sample shown in Panel B should be converted using the equation 1 ng/μL = 1.5 nM, based on an average library size of 1 Kb.

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
http://support.illumina.com/content/dam/illumina-marketing/documents/products/technotes/technote_nextera_library_validation.pdf

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
Index CV for
20-sample pools

| Pool A | Pool B |
|--------|--------|
| 15.8% | 18.2% |

Normalization
beadsで処理:
Bind, Wash,
Elute

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Consumables and Equipment

Check to make sure that you have all of the necessary user-supplied consumables and equipment before proceeding to library preparation. These consumables and equipment are Illumina recommended for the Nextera XT DNA Library Preparation protocols.

Table 3 User-Supplied Consumables


| Consumable | Supplier |
|-------------------------------|----------------------|
| 10 µl pipette tips | General lab supplier |
| 10 µl multichannel pipettes | General lab supplier |
| 10 µl single channel pipettes | General lab supplier |
| 1000 µl pipette tips | General lab supplier |
| 1000 µl multichannel pipettes | General lab supplier |

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
Nextera XT DNA Library Preparation Workflow

The following diagram illustrates the workflow using the Nextera XT DNA Library Preparation Kit. Safe stopping points are marked between steps.

Figure 1 Nextera XT DNA Library Preparation Workflow (For 8 samples)

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Technical Notes:

[Nextera Library Validation and Cluster Density Optimization](#)

[Nextera Low Plex Pooling Guidelines](#)

Support Bulletins:

[Nextera Bioanalyzer traces and sequenced insert size: a comparison](#)

[Considerations for DNA isolation when using Nextera kits](#)

[Decreasing DNA input amount leads to shorter inserts for Nextera DNA Libraries](#)

[Pool up to 384 Samples with Nextera XT](#)

[Nextera Tagmentation Troubleshooting](#)

[New Online Troubleshooting Tools for Nextera DNA and TruSeq DNA PCR-Free](#)

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