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Table 19 Covaris S220 Settings

Setting	350 bp Insert	550 bp Insert
Duty factor	5%	
Peak Incident Power	175 W	
Cycles per burst	200	
Duration	50 seconds	25 seconds
Mode	Frequency sweeping	
Temperature	5.5° to 6°C	

Table 20 Covaris M220 Settings

Setting	350 bp Insert	550 bp Insert
Duty factor	20%	
Peak Incident Power	50 W	
Cycles per burst	200	
Duration	65 seconds	45 seconds
Temperature	20°C	

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


What Inhibitors Are Found in Different Sample Types?			
Sample Source Inhibitors	Inhibitory Effect	Likely Source	Methods to Minimize Inhibition
Fat	Template blocking	Adipose Tissue; Glycerol	Upase or hexane treatment and chloroform extraction, silica-based purification
Bile Salts	Template blocking	Feces, Stool	Wash with 70% ethanol or use silica-based purification
Collagen	Template blocking	Skin, Connective Tissue	Use SDS, CTAB or guanidinium buffers, proteinase K, silica-based purification
Heme	Competition with MgCl ₂	Blood	Wash with 70% ethanol or use silica-based purification
Humic Acid	Chelation of metal ions	Soil, Plant Material	Wash with 70% ethanol or use silica-based purification
Melanin and Eumelanin	Enzyme Binding/ Template Blocking	Hair, Skin	Wash with 70% ethanol or use silica-based purification
Myoglobin	Chelation of metal ions	Muscle Tissue	Wash with 70% ethanol or use silica-based purification
Complex Polysaccharides	Template blocking	Feces, Plant Material	High salt precipitation, CTAB buffer, chloroform extraction, or use silica-based purification
Calcium Ions	Competition with MgCl ₂	Milk, Bone	Wash with 70% ethanol or use silica-based purification
Urea	Enzyme denaturation	Urine	Wash with 70% ethanol or use silica-based purification
Hemoglobin, Lactoferrin	Competition with MgCl ₂	Blood	Wash with 70% ethanol or use silica-based purification
Immunoglobulin G (IgG)	Template blocking	Blood	Use SDS, CTAB or guanidinium buffers, proteinase K, silica-based purification
Indigo Dye; Tannic Acids	Template blocking	Specific Plants	Wash with 70% ethanol or use silica-based purification
What Inhibitors Are Found in Different Sample Treatments/Extraction Methods?			
Extraction Inhibitors	Inhibitory Effect	Likely Source	Methods to Minimize Inhibition
EDTA	Chelation of metal ions	TE buffer	Reduce EDTA concentration in TE buffer or simply use Tris-HCl (10 mM) or water as elution buffer
Alcohols	Enzyme denaturation	Ethanol, isopropanol, isoamyl alcohol	Dry pellet and resuspend or use silica-based purification
Excess Salts	Template blocking	KC _l , NaCl, CsCl, NaAc	Wash with 70% ethanol or use silica-based purification
Chaotropic Salts	Enzyme denaturation	Guanidinium chloride; magnesium chloride, urea	Wash with 70% ethanol or use silica-based purification
Phenol:Chloroform	Enzyme denaturation	Organic carryover	Use PVP, PVP/ammonium acetate, incorporation of 1.2% citric acid at the DNA extraction step
Detergents/DDT	Enzyme denaturation	Sodium deoxycholate, SDS, Tween 20	Wash with 70% ethanol
Exogenous DNA/RNA	Template competition	Carryover	DNase I for DNA removal; RNase A for RNA, RNA:DNA hybrid removal
Carriers	Template competition/ blocking	RNA, heparin, glycogen	Only use carriers that do not serve as template or block the template such as linear acrylamide, N- or P-carriers
Agarose	Template blocking	Gel extractions	Use a spin column with chaotropic salt buffer; dialysis
Excess Metal Ions	Reduce oligo specificity	Mg ⁺⁺ from PCR buffer	Dialysis against PBS (pH 7.4), phenol: chloroform extraction followed by EtOH precipitation

<https://my.illumina.com/MyIllumina/Bulletin/pr9ujrEx7USo794O2k6G4Q/dnarna-isolation-considerations-when-using-truseq>

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




* Read lengths greater than 2 x 151 bp produce a significantly higher percentage of overlapping read-pairs.

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

Table 1: Stranded Parameters

RNA	TruSeq Sample Prep Method	Percent Aligned	Percent Stranded
Universal Human Reference	mRNA*	84.9%	99.6%
Human Reference Brain	mRNA	79.6%	99.0%
Universal Human Reference	Total RNA with Ribo-Zero**	79.0%	98.6%
Human Reference Brain	Total RNA with Ribo-Zero	73.0%	98.6%

Indexed, 2 x 75 cycle run on HiSeq® 2000; *Four mRNA samples per lane; ** Two Ribo-Zero samples per lane.

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TruSeq RNA V2
 TruSeq Stranded mRNA
 TruSeq Stranded total RNA
 TruSeq Small RNA
 TruSeq Targeted RNA
 TruSeq RNA Access

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




Purified mRNA Input

You can also use previously isolated mRNA as starting material. Use the entire fraction of mRNA purified from 0.1 μ g to 4 μ g of total RNA. If you start with isolated mRNA, follow the Illumina recommendations for isolated mRNA specified in the introduction of the Purify and Fragment mRNA procedures. Begin mRNA fragmentation with *Incubate RFP* on page 20 for LS processing or *Incubate RFP* on page 57 for HS processing.

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http://www.mbl.edu/jbpc/files/2014/05/Bioanalyzer_for_NGS_slideshow.pdf

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


Figure 15 Shortened Fragmentation Time Results

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TruSeq DNA PCR-Free Sample Prep LS Protocol

Experienced User Card

Adenylyate 3' Ends

A single 'A' nucleotide is added to the 3' ends of the blunt fragments to prevent them from ligating to one another during the adapter ligation reaction. A corresponding single 'T' nucleotide on the 3' end of the adapter provides a complementary overhang for ligating the adapter to the fragment. This strategy ensures a low rate of chimera (concatenated template) formation.

Consumables

Item	Quantity	Storage	Supplied By
(Optional) A-Tailing Control (CTA)	1 tube per 48 reactions	-15° to -25°C	Illumina
A-Tailing Mix (ATL)	LT kit - 1 tube per 24 reactions or HT kit - 1 tube per 48 reactions	-15° to -25°C	Illumina
Resuspension Buffer (RSB)	1 tube	2° to 8°C	Illumina
Microseal 'B' Adhesive Seal	1	15° to 30°C	User
RNase/DNase-free Reagent Reservoirs (if using multichannel pipettes)	3	15° to 30°C	User
RNase/DNase-free Strip Tubes and Caps (if using multichannel pipettes)	3	15° to 30°C	User

Add ATL

- 1 Centrifuge the thawed A-Tailing Control (if using A-Tailing Control) and A-Tailing Mix tubes to 600 xg for 5 seconds.
- 2 Do one of the following:
 - If using the in-line control reagent, add 2.5 µl of thawed A-Tailing Control to each well of the ALP plate.
 - If not using the in-line control reagent, add 2.5 µl of Resuspension Buffer to each well of the ALP plate.

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The screenshot shows the Illumina website's navigation menu. The 'SUPPORT' menu item is highlighted with a red box. A red callout box points to it with the text: "Support"プルダウンメニューから"Kit&Reagents"を選ぶ. Below the menu, the 'Kits & Reagents' sub-item is also highlighted with a red box. A red arrow points from this box to a red-bordered box containing two product support pages: 'TruSeq DNA PCR-Free HT Library Prep Kit Support' and 'TruSeq DNA PCR-Free LT Library Prep Kit Support'. A red callout box points to this area with the text: 購入したキットを選択. The background of the screenshot shows a 'HiSeq Series' banner with images of sequencing equipment and a 'Log In' button.

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The screenshot shows the support page for the TruSeq DNA PCR-Free HT Library Prep Kit. On the left, a navigation menu lists various resources, each highlighted with a colored box and connected to a callout box on the far left. The callout boxes contain the following text:

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The main content area of the page includes:


- Header: Support » Sequencing » Kits and Reagents » TruSeq DNA PCR-Free HT Library Prep Kit
- Product Name: TruSeq DNA PCR-Free HT Library Prep Kit
- Catalog IDs: FC-121-3003
- Latest Updates section with links to: Custom Protocol Selector 03/10/2015, Illumina Experiment Manager v1.9 03/09/2015, and Illumina Experiment Manager Guide 03/09/2015.
- Features section listing: Easy-to-use 96-well plate, pre-loaded with 96 unique index combinations; Robust dual index system; Optimized shearing for whole-genome resequencing with 350 bp and 550 bp insert size workflows; Bead-based size selection reagents included in kit; Designed for manual or automated preparation; Compatible with all Illumina sequencers; End-to-end solution from sample prep to data analysis.

http://support.illumina.com/sequencing/sequencing_kits/truseq-dna-pcr-free-ht-kit.html

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

Library Prep Workflow

The following figure illustrates the TruSeq DNA PCR-Free Library Prep LS workflow.

Figure 1 TruSeq DNA PCR-Free Library Prep LS Workflow

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[TruSeq Nano DNA 350 bp library](#)

[TruSeq Nano DNA 550 bp library](#)

[TruSeq Stranded mRNA Sample Prep data set](#)

Support Bulletins:

[DNA/RNA Isolation Considerations When Using TruSeq Library Prep Kits](#)

[Sequencing Kits for FFPE samples](#)

[Considerations for RNA-Seq Read Length and Coverage](#)

[How do I deplete bacterial ribosomal RNA for Illumina RNA sequencing?](#)

[Will a TruSeq Stranded Total RNA kit work with my organism?](#)

[Sequencing Solutions for Low Inputs of RNA](#)

[New Online Troubleshooting Tools for Nextera DNA and TruSeq DNA PCR-Free Sample Prep Kits](#)

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