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© 2014 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, GAIx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeqFGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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Application Note: Sequencing illumina

NextSeq™ 500 System Exome Sequencing Solution

A cost-effective, high-coverage exome sequencing solution that delivers the most accurate variant calling.

Highlights

- **Fast Sample-to-Data Exome Solution**
Simple workflow with the lowest hands-on time
- **Comprehensive Exome Coverage**
Interrogate more of the exome than ever before, even in challenging regions
- **Easy Data Analysis**
Walkaway sample-to-results solution with analysis performed locally or in the cloud
- **Most Accurate Variant Detection**
Accurate calls and low detection limit to identify common mutations and rare somatic events
- **End-to-End Burnina Support**
Burnina scientists and engineers are there every step of the way, providing installation, training, applications, and data analysis support


Figure 1: NextSeq 500 System Sample-to-Answer Exome Sequencing Workflow

Prep	Sequencing	Analysis	Share
1.5	29	5	1
hours	hours	hours	hours

The NextSeq 500 System's simple workflow delivers highly accurate sequencing data. Data analysis includes alignment and variant calling.

Base calls are generated on the NextSeq 500 System and data analysis (including alignment and variant calling) can be performed in the cloud or on the Burnina premises computing environment (Figure 2). With its intuitive interface accessed through a common web browser, BaseSpace provides researchers with access to a rich ecosystem of commercial and open-source DNA data analysis software tools designed primarily to analyze Burnina data.

Figure 2: NextSeq 500 Sequencing System

The NextSeq 500 System integrates the latest advances in 300x chemistry and the industry's simplest workflow.

Introduction

The NextSeq 500 System Exome Sequencing Solution enables researchers to investigate the protein encoding (exonic) regions of the genome, identifying variants for population genetics, genetic disease, and cancer studies efficiently and cost effectively. It leverages industry leading Burnina next generation sequencing (NGS) technology responsible for > 90% of global exome sequencing delivering the best data quality and highest accuracy to identify true coding variants. The NextSeq 500 System Exome Sequencing Solution includes integrated library preparation and exome enrichment, push-button sequencing, and simple data analysis. With minimal hands-on time, the NextSeq 500 System Exome Sequencing Solution is the most flexible, comprehensive tool for interrogating more of the exome quickly and efficiently.

A Simple, Efficient Exome Sequencing Workflow

The NextSeq 500 exome sequencing workflow simplifies exome sequencing, enabling researchers to maximize their productivity (Figure 1). It begins with library preparation and exome enrichment using a Nextera® Rapid Capture Exome Kit. The NextSeq 500 System's dual sequencing input mode enables researchers to pool their exome studies (20, 40, 60, or 120 Chx per run).

What is a "true coding variant"?

A true coding variant is an accurate base call that differs from the consensus sequence within a coding region. It is not a false positive (where a variant is called but does not truly exist) or a false negative (where a variant that truly exists is not called). A system with a high false positive call rate requires extensive downstream validation, increasing costs and experimental time. A system with a high false negative call rate is failing to detect potentially important findings, often in regions that are highly repetitive or that contain homopolymer stretches.

The NextSeq system leverages Burnina industry leading sequencing technology to deliver the highest number of true coding variants.

Technical Note: Sequencing illumina

Read Length and Nextera® Rapid Capture Exome Data

Nextera Rapid Capture Exome data sets with longer read lengths deliver higher mean coverage and more variant calls.

Highlights

- **Paired-End Data With Exome Sequencing**
PE sequencing provides better read depth and a higher number of variant calls
- **Longer Reads Produce Higher Coverage**
Longer read lengths increased mean coverage for all inbred data sets
- **Higher Coverage Increases True Variant Calls**
Longer read lengths generated more variant calls with high precision and sensitivity rates


Introduction

The protein coding portion of the human genome—the exome—represents approximately 1% of the human genome. At first glance this seems like a small percentage; however, studies have shown roughly 80% of all mutations associated with Mendelian diseases fall within exome boundaries.¹ The high concentration of disease-associated variants within the exome as well as the lower cost compared to whole genome sequencing (WGS) have led to a steady increase in the use of whole exome sequencing (WES) over the past several years.² The ability to focus on a targeted subset of the genome, and therefore decrease sequencing coverage, significantly increases the power of WES methods to identify, not only causative variants of common Mendelian genetic diseases, but to detect rare variants associated with complex disorders.³

The purpose of this technical note is to provide guidance regarding optimization of sequencing parameters for exome sequencing. We also review the data quality parameters affected by changes in read length and coverage depth. As sequencing depth is the cornerstone of successful variant discovery, any WES method that balances coverage levels will also increase experimental accuracy.⁴ To that end, Nextera Rapid Capture Exome sequencing, which utilizes a paired-end approach, imparts several key advantages (Figure 1). In addition to producing twice the number of reads for the same amount of hands-on time and effort, sequences aligned as read pairs enable more accurate read alignment, higher numbers of single-nucleotide variant (SNV) calls, and an ability to detect insertions and deletions (indels) that is not possible with single read data.⁵ Analysis of differential read pair spacing also allows removal of PCR duplicates.

Along with the paired-end approach, sequence read length is another means of tuning sensitivity and accuracy. To investigate these relationships, we sequenced and analyzed a well-characterized CEP110 to assess the impact of read length on mean coverage, coverage uniformity, and variant calling.

Figure 1: Nextera Rapid Capture Exome Workflow and Specifications

The workflow includes: Project, deconvoluted Nextera sequencing library; Hybridize biotinylated probes to exome targets; Wash away non-bound probes; Probes bind to oligonucleotide beads; Wash away unbound beads; Rapid capture using oligonucleotide beads; Sequencing library fragments from beads; Elute sequencing-ready fragments from beads.

Specification	Nextera Rapid Capture Exome
Target region size	27 Mb
Number of target genes	214,026
Capture DNA input	50 ng
Hybridize time	5 hours
Total time	1.5 days

The Nextera Rapid Capture Exome V2 provides a fast, simple method for isolating the human exome. The streamlined workflow combines library preparation and exome enrichment steps, and can be completed in 1.5 days with low Chx sample input.

<http://www.illumina.com/content/dam/illumina-marketing/documents/products/technotes/technote-nrc-exome-read-length.pdf>
<http://www.illumina.com/content/dam/illumina-marketing/documents/products/appnotes/appnote-nextseq-exome.pdf>
http://www.illumina.com/documents/products/datasheets/datasheet_nextera_rapid_capture_exome.pdf

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
BaseSpace
HiSeq2500: Nextera Rapid Capture Exome (12plex, CEPH Trio replicates) : BWA Enrichment (12plex, Manifest v1.2) 01/16/2014 9:15:22

Analysis Info

Inputs

Output Files

Analysis Reports

Aggregates Summary

PDF Reports

FC1_NA12878_01

FC1_NA12878_04

FC1_NA12878_07

FC1_NA12878_10

FC1_NA12891_02

FC1_NA12891_05

FC1_NA12891_08

FC1_NA12891_11

ALIGNMENT SUMMARY 1

Sample	Sample Name	Mean Coverage	Target Coverage at 1X	Target Coverage at 10X	Target Coverage at 20X	Target Coverage at 50X
1	FC1_NA12878_01	101.6 X	99.3%	95.5%	90.0%	68.6%
2	FC1_NA12878_04	115.1 X	99.5%	96.3%	91.9%	73.6%
3	FC1_NA12878_07	119.2 X	99.5%	96.5%	92.3%	75.0%
4	FC1_NA12878_10	77.1 X	99.2%	94.0%	86.0%	57.0%
5	FC1_NA12891_02	76.2 X	99.3%	93.2%	84.4%	55.5%
6	FC1_NA12891_05	121.7 X	99.6%	96.6%	92.1%	74.8%
7	FC1_NA12891_08	102.2 X	99.6%	95.9%	90.4%	68.8%
8	FC1_NA12891_11	83.9 X	99.4%	94.5%	87.1%	60.4%





Variants by Sequence Context

	SNVs	Insertions	Deletions
Number in Genes	32,157	1,046	1,228
Number in Exons	22,949	305	295
Number in Coding Regions	21,641	216	205
Number in UTR Regions	1,308	89	90
Number in Splice Site Regions	1,310	99	152


Variants by Consequence

	SNVs	Insertions	Deletions
Frameshift	-	104	82
Non-synonymous	10,464	112	122
Synonymous	11,076	-	-
Stop Gained	76	0	0
Stop Lost	24	0	0

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PDF Summary Report

SAMPLE INFORMATION

Total PF Reads	Percent Q30	Percent Duplicate Paired Reads	Fragment Length Median	Fragment Length Standard Deviation
65,262,238	82.3%	6.5%	227 bp	80 bp

ENRICHMENT SUMMARY

Total Length of Targeted Reference	Padding Size
44,296,481 bp	150 bp

Read Level Enrichment

Total Aligned Reads	Percent Aligned Reads	Target Aligned Reads	Read Enrichment	Padded Target Aligned Reads	Padded Read Enrichment
62,912,089	96.4%	39,472,997	62.7%	42,236,377	67.1%

Base Level Enrichment

Total Aligned Bases	Target Aligned Bases	Bases Enrichment	Padded Target Aligned Bases	Padded Base Enrichment
8,778,180,848	4,378,065,923	49.9%	5,837,362,539	66.5%

SMALL VARIANTS SUMMARY

	SNVs	Insertions	Deletions
Total Passing	31,080	989	1,185
Percent Found in dbSNP	99.0%	88.0%	78.3%

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The screenshot shows the BaseSpace Public Data interface. The top navigation bar includes 'BaseSpace', 'Dashboard', 'Prep', 'Runs', 'Projects', 'Apps', 'Public Data', and 'Help'. The user is logged in as 'Eri Kibukawa'. The main content area is titled 'Public Data' and features a search bar. Below the search bar, there are links for 'Nextera Rapid Capture Exome Data Sets' and 'Targeted Sequencing Exome'. A specific data set is highlighted: 'HiSeq 2500: Nextera Rapid Capture Exome (CEPH Trio)'. The description states: 'Data generated from 4 replicates each of the CEPH trio NA12878, NA12891 and NA12892 sequenced on HiSeq 2500. Samples prepared using Nextera Rapid Capture Exome reagent kit.' Below the description, there are tabs for 'Resequencing', 'Exome', and 'Targeted Sequencing'. A table lists the data set details:


Run	HiSeq2500_NRC_CEPH_Trio_2x150_FC1 (71.51 GB)	Import
Project	HiSeq 2500: Nextera Rapid Capture Exome (12plex, CEPH Trio replicates) ()	Import

The 'Import' button for the Project row is circled in red. To the right of the main content, there are two panels: 'Research Areas' with buttons for 'Cancer Research', 'Genetic Disease', 'Complex Disease', and 'Microbial Research'; and 'Categories' with buttons for 'Exome', 'Resequencing', 'Small RNA', 'Targeted Sequencing', 'De Novo Assembly', 'RNA-Seq', 'Gene Fusion', and 'ChIP-Seq'.

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
BWA Enrichment v2.1
illumina, inc

Analysis Name: BWA Enrichment v2.1 05/14/2015 8:02:4

Save Results To: Select Project(s):
HiSeq 2500: Nextera Rapid Capture Exome (12plex, CEPH Trio replicates)

Sample(s) [96 maximum]: Select Sample(s):

Reference Genome: Human (UCSC hg19)

Targeted Regions: Nextera Rapid Capture Exome v1.2

Custom Targeted Manifest: Select File(s):

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


**Isaac
Enrichment**

VariantStudio

http://support.illumina.com/downloads/basespace_core_apps_user_guides.html
http://support.illumina.com/sequencing/sequencing_software/variantstudio.html

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Appendix.

http://www.illumina.com/applications/sequencing/dna_sequencing/exome-sequencing.html

The image displays three overlapping screenshots of the Illumina website's Exome Sequencing section. The top-left screenshot shows the 'Exome Sequencing' overview page, highlighting the benefits of exome sequencing and providing a sign-up form for interested parties. The top-right screenshot shows the 'Exome Sequencing Data Analysis' page, detailing the BaseSpace environment and the Exome Data Analysis with BaseSpace Apps. The bottom screenshot shows a table titled 'Table: Optimal Samples per Run Across' which provides technical specifications for different sequencing systems.

	MI20	MI250 500	MI250 1500
	V3 REAGENT KIT	MD OUTPUT	HIGH OUTPUT
Rolling Retry	1	Up to 1	
	3	Up to 3	
	5		Up to 5
	9		Up to 9
	12		Up to 12

Number of exomes pooled before enrichment.
For Research Use Only.

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