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© 2015 Illumina, Inc. All rights reserved.
Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDX, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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登録: http://jp.illumina.com/events/ilmn_support_program/ilmn-ischool.html

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Supporting Information

Consumables and Equipment

Make sure that you have the required user-supplied consumables and equipment before starting the protocol.

The protocol has been optimized and validated using the items listed. Comparable performance is not guaranteed when using alternate consumables and equipment.

NOTE

- Use a dedicated set of consumables and equipment for pre-PCR and post-PCR procedures.
- The TruSeq Custom Amplicon Low Input library prep protocol requires different magnetic stands for pre-PCR and post-PCR procedures.

Consumables

Consumable	Supplier
10 N NaOH (prepare from tablets or use a standard solution)	General lab supplier
20 µl barrier pipette tips	General lab supplier
20 µl multichannel pipettes	General lab supplier
20 µl single channel pipettes	General lab supplier
200 µl barrier pipette tips	General lab supplier
200 µl multichannel pipettes	General lab supplier
200 µl single channel pipettes	General lab supplier
1000 µl barrier pipette tips	General lab supplier
1000 µl multichannel pipettes	General lab supplier
1000 µl single channel pipettes	General lab supplier
One of the following plate types:	

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Thermal Cyclers

Thermal Cyclers

Use the following recommended settings for selected thermal cycler models. Before performing library prep, validate any thermal cyclers not listed.

Thermal Cycler	Block Type	Ramp Rate	Lid Temp	Block Rate	Vessel Type
Bio-Rad DNA Engine Tetrad 2	Standard	0.1°C	Heated, Constant at 100°C	--	Bio-Rad Hard-Shell 96-Well Skirted PCR Plates, low-profile, skirted
Bio-Rad S1000	Standard	0.1°C	Heated, Constant at 100°C	--	Bio-Rad Hard-Shell 96-Well Skirted PCR Plates, low-profile, skirted
Bio-Rad C1000	Standard	0.1°C	Heated, Constant at 100°C	--	Bio-Rad Hard-Shell 96-Well Skirted PCR Plates, low-profile, skirted
Bio-Rad T100	Standard	0.1°C	Heated, Constant at 100°C	--	Eppendorf twin.tec PCR Plate 96, semi-skirted
Applied Biosystems GeneAmp PCR System 9700	Gold	4%	Heated, Constant at 100°C	9600 (for HYB and EXT_LIG programs) or MAX (for PCR program)	Eppendorf twin.tec PCR Plate 96, semi-skirted
Applied Biosystems Veriti 96-Well	Alloy	2%	Heated, Constant at 100°C	--	Eppendorf twin.tec PCR Plate 96, semi-skirted

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One of the following magnets, depending on the PCR plate type used:

- DynaMag-96 Side Skirted Magnet (use with 96-well full-skirted PCR plates)
- DynaMag-96 Side Magnet (use with Eppendorf 96-well twin.tec PCR plates)

Magnetic stand-96 (use with midi 96-well storage plates)

Life Technologies:

- catalog # 12027
- catalog # 12331D

Life Technologies, catalog # AM10027

Magnetic Stand

Magnetic Stand

Plate Shaker

One of the following:

- BioShake iQ high-speed thermal mixer
- BioShake XP high-speed lab shaker

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
TruSeq Custom Amplicon Low Input Library Prep Workflow

Pre-amplification

- 1 Quantify and Dilute DNA
- 2 Hybridize Oligo Pool
- 3 Remove Unbound Oligos
- 4 Extend-Ligate Bound Oligos
- 5 Amplify PCR

Post-amplification

- 6 Clean Up PCR
 - 7 Normalize Libraries
 - 8 Pool Libraries
- Cold Storage Option
- [Click each step to learn more.](#)

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http://support.illumina.com/downloads/prepare_libraries_for_sequencing_miseq_15039740.html

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TSCA Low Input Kit Training動画

● Training

https://support.illumina.com/sequencing/sequencing_kits/truseq-custom-amplicon-low-input/training.html

illumina® MyIllumina Quick Order View Cart Contact Us English

AREAS OF INTEREST ▾ TECHNIQUES ▾ SYSTEMS ▾ PRODUCTS & SERVICES ▾ INFORMATICS ▾ SCIENTIFIC CONTENT ▾ COMPANY ▾ SUPPORT ▾ Q SEARCH

Support > Sequencing > Kits and Reagents > TruSeq Custom Amplicon Low Input Library Prep Kit > Training

illumina® TruSeq Custom Amplicon Low Input Help

Menu Transcript

- 1. Welcome
- 2. Getting Started
- 3. TruSeq Custom Amplicon Low Input Overview
- 4. Place Order and Set Up
- 5. Download Manifest
- 6. Plan Index Setup and Prepare Libraries
 - 6.1. How Does TruSeq Custom Amplicon Low Input Library Prep Kit Work?
 - 6.2. Plan Index Setup
 - 6.3. TruSeq Custom Amplicon Low Input Library Prep Kit
 - 6.4. Summary of Best Practices
- 7. Summary

Remove Unbound Oligos

00:27 / 01:31

Search... [Q] [▶] [↺] [↻] [◀ BACK] [NEXT ▶]

Innovative... At Illumina, we don't just... not even in... our custom... highest level... advancements... For Research... © 2016 Illumina, Inc. All rights reserved.

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- **TruSeq Custom Amplicon Low Input Kit**
<http://jp.illumina.com/products/truseq-custom-amplicon-low-input.html>
- **DataSheet**
http://jp.support.illumina.com/content/dam/illumina-marketing/apac/japan/documents/pdf/datasheet_tru_seq_custom_amplicon_low_input.pdf
- **Training**
https://support.illumina.com/sequencing/sequencing_kits/truseq-custom-amplicon-low-input/training.html
- **User Guide**
<https://support.illumina.com/downloads/truseq-custom-amplicon-low-input-library-prep-reference-guide-100000002191.html>

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