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© 2016 Illumina, Inc. All rights reserved.

Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiniSeq, MiSeq, MiSeqDx, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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Local Run Manager

localhost/#/login

Apps For quick access, place your bookmarks here on the bookmarks bar. [Import bookmarks now...](#)

illumina

Local Run Manager

admin

.....

Login

[Forgot your user name or password?](#)

Instrument: MiniSeq ML-P2-14
Local Run Manager Version: 1.2.0

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The screenshot shows the Local Run Manager interface. At the top, there are statistics for run statuses: 0 Ready, 0 In Progress, 0 Stopped or Unsuccessful, 5 Complete, and 5 Total. A red callout box labeled "新規ランの開始" (Start New Run) points to the "Create Run" button. The dropdown menu for "Create Run" includes options: Amplicon DS, Amplicon, DNA Enrichment, GenerateFASTQ, Resequencing, Small RNA, and Targeted RNA. Below the statistics is a table of runs.

RUN NAME / ID	MODULE	STATUS	LAST MODIFIED
H00MWL_TargetedRNA_ELM4_Pool3 160111_ML-P2-14_0059_AH00MW...	TARGETED RNA	Analysis Completed	2016-01-11 21:12
H00G52_PoolT_TrusightTumor 160106_ML-P2-14_0056_AH00G5...	AMPLICON DS	Analysis Completed	2016-01-07 10:01
H00KNK_PoolT_TrusightTumor 160105_ML-P2-14_0055_AH00KN...	AMPLICON DS	Analysis Completed	2016-01-06 11:30
TruSeq_Amplicon_LowInput_Po... 151220_ML-P2-14_0051_ATruSe...	AMPLICON	Analysis Completed	2015-12-21 12:00
TruSight_Enrichment_PoolX 151216_ML-P2-14_0048_ATruSL...	DNA ENRICHMENT	Analysis Completed	2015-12-17 14:03

Showing 1-5 of 5 Active Runs

Instrument: MiniSeq ML-P2-14
Local Run Manager Version: 1.2.0

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Run Description

Run Settings

Library Prep Kit* TruSeq Targeted RNA Expression

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Module-Specific Settings

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No module-specific settings available for the selected Library Prep Kit.

Show advanced module settings...

Remember to upload the manifests. Run cannot be saved until manifests are uploaded.

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Local Run Manager

localhost/#/runupload/TargetedRNAWorkflow

Index Reads* 0 1 2

Read Lengths* 51 6 8 --

Custom Primers Read 1 Index 1 Index 2 --

Module-Specific Settings

No module-specific settings available for the selected Library Prep Kit.

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	SAMPLE ID*	SAMPLE NAME*	SAMPLE DESCRIPTION	INDEX 1 (I7)*	INDEX 2 (I5)*	MANIFEST*	GENE NORMALIZATION	
1	Library_1	Sample_A	A_biological_replicate	R701	A501	TruSeq CRT Manifest TC0028756-CRT.txt		✕
2	Library_2	Sample_A	A_biological_replicate	R702	A502	TruSeq CRT Manifest TC0028756-CRT.txt		✕
3	Library_3	Sample_B	B_biological_replicate	R703	A503	TruSeq CRT Manifest TC0028756-CRT.txt		✕
4	Library_4	Sample_B	B_biological_replicate	R704	A504	TruSeq CRT Manifest TC0028756-CRT.txt		✕

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Perform wash

ML-P2-10 | admin

Log In | **Select Run** | Load | Run Setup | Check | Sequence

Select a run from the list of runs from Local Run Manager. Select Next to show run details.

Search for Run Name

Targeted_RNA_160119

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Created by _____
Date _____

Read 1	Index 1	Index 2	Read 2
_____	_____	_____	_____

Run Description _____

Back | Exit | Next

ML-P2-10 | admin

Log In | **Select Run** | Load | Run Setup | Check | Sequence

Select a run from the list of runs from Local Run Manager. Select Next to show run details.

Search for Run Name

Targeted_RNA_160119

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Created by	Admin User		
Date	1/13/2016		
Read 1	Index 1	Index 2	Read 2
51	6	8	0

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ML-P2-10 admin

Log In Select Run **Load** Run Setup Check Sequence

i Load the flow cell and reagent cartridge, and empty the spent reagents bottle. Select Next.

Flow Cell
000H00MY7
✓ Flow cell clamp
✓ Flow cell door

Reagent

✗ Reagent cartridge
✓ Spent reagents
✓ Reagent door

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The screenshot shows the BaseSpace TruSeq Amplicon interface. The top navigation bar includes 'Dashboard', 'Prep', 'Runs', 'Projects', 'Apps', 'Public Data', and 'Help'. The user is logged in as 'Stephen Gross'. The breadcrumb trail indicates the current location: 'MiniSeq: TruSeq Custom Amplicon Low Input (replicates of Horizon samples and NA12878) > TruSeq Amplicon Horizon samples'. Below the navigation bar are buttons for 'Download Analysis', 'Rename analysis', and 'Move to Trash', along with a 'View Trash' link. The left sidebar contains 'Analysis Info', 'Inputs', 'Output Files', and 'Analysis Reports'. The main content area is titled 'Aggregate Summary' and 'Amplicon Summary'. A large red box highlights a table of 'Read Level Statistics'.

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Sample	Sample Name	Total Aligned Reads (R1/R2)	Percent Aligned Reads (R1/R2)	Overall Aligned Reads
1	HD701-rep2	1,752,951 / 1,763,108	93.59% / 94.13%	93.86%
2	HD701-rep4	1,736,843 / 1,746,096	94.17% / 94.67%	94.42%
3	HD200-rep4	2,177,036 / 2,197,990	91.50% / 92.38%	91.94%
4	HD200-rep3	2,035,933 / 2,054,105	93.53% / 94.37%	93.95%
5	HD200-rep1	2,221,159 / 2,238,305	92.41% / 93.12%	92.77%
6	HD701-rep1	2,208,827 / 2,223,237	94.42% / 95.03%	94.72%
7	HD200-rep2	1,950,528 / 1,963,673	93.33% / 93.96%	93.65%
8	HD701-rep3	1,912,418 / 1,926,364	93.35% / 94.03%	93.69%

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BaseSpace® Dashboard Prep Runs **Projects** Apps Public Data Help Stephen Gross illumina®

MiniSeq: TruSeq Custom Amplicon Low Input (replicates of Horizon samples and NA12878) : TruSeq Amplicon Horizon samples

Download Analysis Rename analysis Move to Trash View Trash

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Amplicon Summary


Number of Amplicon Regions	Total Length of Amplicon Regions
212	36,024 bp

Read Level Statistics

Read	Total Aligned Reads	Percent Aligned Reads
1	2,208,827	94.42%
2	2,223,237	95.03%

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BaseSpace Dashboard Prep Runs **Projects** Apps Public Data Help Stephen Gross illumina

MiniSeq: TruSeq Custom Amplicon Low Input (replicates of Horizon samples and NA12878) : TruSeq Amplicon Horizon samples

Download Analysis Rename analysis Move to Trash View Trash

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