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1. Rahman N. Integrating germline cancer predisposition gene testing into routine clinical practice - The Mainstreaming Cancer Genetics programme [AACR abstract 4705]. *Cancer Res.* 2014;74(19suppl):4705.
2. Münz M, Ruark E, Renwick A, et al. CSN and CAVA: variant annotation tools for rapid, robust next-generation sequencing analysis in the clinical setting. *Genome Med.* 2015;7(76):1-8.
3. Wimberger P, Benet-Pages A, Schubert S, et al. Identification of ERCC2 as novel susceptibility gene for hereditary breast and ovarian cancer [ASCO abstract 1512]. *J Clin Oncol.* 2015;33(suppl):1512.
4. Sykorova V, Dvorakova S, Vcelak J, et al. Search for new genetic biomarkers in poorly differentiated and anaplastic thyroid carcinomas using next generation sequencing. *Anticancer Res.* 2015;35(4):2029-2036.
5. Shendrue J, and Hanlee J. Next-generation DNA sequencing. *Nat Biotech.* 2008;26:1135-1145.
6. Schuster SC. Next-generation sequencing transforms today' s biology. *Nat Methods.* 2008;5(1):16-18.
7. Tuononen K, Maki-Nevala S, Sarhadi VK, et al. Comparison of targeted next-generation sequencing (NGS) and real-time PCR in the detection of EGFR, KRAS, and BRAF mutations on formalin-fixed, paraffin-embedded tumor material of non-small cell lung carcinoma- superiority of NGS. *Genes, Chromosomes and Cancer.* 2013;52(5):503-511.
8. Liu L, Li Y, Li S, et al. Comparison of next-generation sequencing systems. *J Biomed Biotechnol.* 2012;2012:251364.
9. Ross MG, Russ C, Costello M, et al. Characterizing and measuring bias in sequence data. *Genome Biol.* 2013;14:R51.
10. Sherry ST, Ward MH, Kholodov M, et al. dbSNP: the NCBI database of genetic variation. *Nucleic Acids Res.* 2001;29(1):308-311.

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