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

Figure 2: RNA Quality from FFPE Samples


RNA isolated from FFPE samples was examined using an Agilent Bioanalyzer. RNA Integrity Numbers (RINs)<sup>®</sup> were calculated from the Bioanalyzer traces.

Figure 3: RIN versus DV<sub>200</sub> and Library Yield

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
\*The " $DV_{200}$ " is the percentage of RNA fragments > 200 nucleotides. Although RIN values for these samples lie within a relatively narrow range (2.3–3.2), the size distribution of the RNA varies greatly among the samples.

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Table 2: Recommended RNA Input Based on DV<sub>200</sub>

Quality**††	DV200	Recommended Input Quantity
High	> 70 %	20 ng
Medium	50 – 70%	20 – 40 ng
Low	30 – 50%	40 – 100 ng
Too Degraded	< 30 %	Not Recommended

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
## Identification of a Novel *PARP14-TFE3* Gene Fusion from 10-Year-Old FFPE Tissue by RNA-Seq

Weihua Huang,<sup>1</sup> Michael Goldfischer,<sup>2</sup> Sabina Babyeva,<sup>1</sup> Yong Mao,<sup>3</sup> Konstantin Volyansky,<sup>3</sup> Nevenka Dimitrova,<sup>3</sup> John T. Fallon,<sup>1</sup> and Minghao Zhong<sup>1\*</sup>

<sup>1</sup>Department of Pathology, New York Medical College, Westchester Medical Center, Valhalla, NY

<sup>2</sup>Department of Pathology, Hackensack University Medical Center, Hackensack, NJ

<sup>3</sup>Clinical Informatics Solutions and Services (CISS), Philips Research North America, Briarcliff Manor, NY

Novel *PARP14-TFE3* gene fusion in RCC suggested by RNA Seq results.

### RNA-Seq

After xylene deparaffinization, total RNA was extracted using the AllPrep DNA/RNA FFPE Kit (Qiagen, Hilden, Germany). The quantity and integrity of the RNA was measured using the Agilent 2100 Bioanalyzer. One hundred nanograms of the RNA, with DV200 at 30%, was applied for sequencing library preparation using the TruSeq RNA Access Library Prep Kit (Illumina, San Diego) as per the manufacturer's protocol. Paired-end sequencing (75 × 2) was performed using the MiSeq Reagent V3 Kit (150 cycles) and the MiSeq sequencing system (Illumina). Bowtie2 was used for alignment and mapping of short sequence reads to the human genome reference hg19 or the fusion transcript *PARP14-TFE3*. Both STAR and TopHat algorithms were used for detection of any potential *TFE3* fusion. Integrative Genomics Viewer (IGV) was used for data visualization.

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The image shows a WebEx session interface. At the top, a toolbar contains several icons: a back arrow, a microphone icon labeled '音声プロ...', a group of people icon labeled '参加者', a speech bubble icon labeled 'チャット', and a question mark icon labeled 'Q&A'. The 'Q&A' icon is highlighted with a red box. Below the toolbar, a green notification bar reads 'Mio Tonouchiのアプリケーションを表示...'. The main content area displays a slide with the following text:

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illumina®

At the bottom right, a 'Q&A' window is open. It has a title bar with a question mark icon and the text 'Q&A'. Below the title bar, there is a tab labeled 'すべて (0)'. The main area of the window is empty. At the bottom, there is a '質問:' label followed by a dropdown menu showing 'すべてのパネリスト'. Below this, there is a text input field containing the question '質問です。Nexteraの場合' and a '送信' button. The entire Q&A window area is highlighted with a red box.

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©2015 Illumina, Inc. All right reserved.