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

Figure 2: RNA Quality from FFPE Samples

RNA isolated from FFPE samples was examined using an Agilent Bioanalyzer. RNA Integrity Numbers (RINs)⁵ were calculated from the Bioanalyzer traces.

Figure 3: RIN versus DV₂₀₀ and Library Yield

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
Figure 4 Example: Low Quality FFPE ($DV_{200} = 30\%$)


Figure 3 Example: Medium Quality FFPE ($DV_{200} = 55\%$)



Figure 5 Example: FFPE Quality Not Recommended for Use ($DV_{200} = 8\%$)


TruSeqRNA Access LibraryPrep Guide , Part # 15049525 Rev. B <Technical Note>

http://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/samplepreps_truseq/truseqrnaaccess/truseq-rna-access-library-prep-guide-15049525-b.pdf

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^{*}The "DV₂₀₀" is the percentage of RNA fragments > 200 nucleotides. Although RIN values for these samples lie within a relatively narrow range (2.3–3.2), the size distribution of the RNA varies greatly among the samples.

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
Table 1 FFPE RNA Input Recommendations

Quality	DV ₂₀₀	Input Requirement Per Reaction
High	> 70%	20 ng
Medium	50–70%	20–40 ng
Low	30–50%	40–100 ng
Too Degraded	< 30%	Not recommended

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	TruSight Tumor (26 genes)	TruSight Cancer Hotspot Panel	TruSeq Custom Amplicon Low Input	TruSight Tumor (15 genes)
用途	Somatic mutation detection in solid tumors	Somatic mutational hotspots in a broad spectrum of cancers	Somatic mutational detection in a broad spectrum of cancers	Somatic mutational hotspots in specific cancers
製品番号	FC-130-2001/ TG-130-2001	FC-135-2002, FC-135-2001 (& FC-121-9999 for FFPE QC)	FC-134-2001, FC-134-2002 (& FC-121-9999 for FFPE QC)	OP-101-1001, OP-101-1002
手法	Amplicon (double stranded)	Amplicon	Amplicon	Multiplex PCR
ターゲットサイズ	21kb (174 amplicons, 26 genes)	>35kb (212 amplicons, 48 genes)	(16 - 1536 amplicons)	44kb (250 amplicons, 15 genes)
DNA 量	30 - 300ng	10ng~	10ng~	20ng
FFPE 対応	Yes	Yes	Yes	Yes
リード長	2x121bp	2x150bp	2x150bp	2x150bp
カバレッジ	At least 1000x per amplicon	~1000x average coverage	~1000x average coverage	At least 500x per amplicon
キットサイズ	48 samples	16, 96 samples	16, 96 samples	24 samples
対応機種/ サンプル数	MiSeq: 6 samples NextSeq: 34 samples /48 samples HiSeq: 48 samples	MiSeq: up to 96 samples NextSeq: up to 96 samples HiSeq: up to 96 samples	MiSeq: up to 96 samples NextSeq: up to 96 samples HiSeq: up to 96 samples	MiSeq: 8 samples
アライメント/ バリエーションコール	MiSeq Reporter with Amplicon - DS workflow; AmpliconDS for HAS AmpliconDS	MiSeq Reporter with Somatic Variant Caller BaseSpace App TruSeq Amplicon	MiSeq Reporter with Somatic Variant Caller BaseSpace App TruSeq Amplicon	MiSeq Reporter Predefined Report
フィルタリング/ アノテーション	VariantStudio	VariantStudio	BaseSpace App / Onsite	VariantStudio
	固形腫瘍	腫瘍全般	カスタム	固形腫瘍

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
The screenshot shows the Illumina website's 'Oncology Breakthrough' page. At the top, there's a navigation bar with 'illumina' logo and various menu items. Below the navigation, the main heading is 'Oncology Breakthrough' with a sub-heading '次世代シーケンサーの活用が確。最新癌研究のブレイクスルー' (Next-generation sequencing is confirmed for use. Breakthrough in the latest cancer research). There's a search bar and social media links. A central image shows a woman in a lab coat. To the right, there's a '製品検索' (Product Search) section with a dropdown menu and a 'ログイン' (Login) form. Below the main content, there are several icons representing different research areas: '最新発表はここから' (Latest publications from here), 'Oncology Breakthrough', 'イルミナ サポートプログラム' (Illumina Support Program), and '生産性のパワーを拡大' (Expand productivity power). At the bottom, there's a section for 'イルミナ School' (Illumina School) with various topics like '癌ゲノム' (Cancer Genomics), '遺伝性疾患' (Genetic Diseases), '微生物ゲノム' (Microbial Genomics), '創薬におけるゲミクス' (Genomics in Drug Discovery), 'アグリゲノム' (Agricultural Genomics), 'リプロダクティブヘルス' (Reproductive Health), and '法医学ゲノム' (Forensic Genomics). A yellow arrow points from the 'Oncology Breakthrough' banner to the campaign details on the right.

The poster is for the '2015 Oncology Breakthrough RNA Access & TruSight Tumor FFPE 対応ライブラリー調製キット お試しキャンペーン' (2015 Oncology Breakthrough RNA Access & TruSight Tumor FFPE Compatible Library Preparation Kit Trial Campaign). It features a background image of a person in a lab coat. The main text reads: 'FFPE サンプルの mRNA 解析に最適な RNA Access ライブラリー調製キット、FFPE 処理によるアーチファクトに対応した TruSight Tumor を約 30%OFF でご提供' (Optimal RNA Access library preparation kit for mRNA analysis of FFPE samples, TruSight Tumor compatible with artifacts from FFPE processing, provided at approximately 30% OFF). Below this, there's a section for '対象製品' (Target Products) listing 'TruSeq RNA Access · TruSight Tumor'. The '内容' (Content) is '対象製品を約 30% OFF でご提供' (Provide target products at approximately 30% OFF). The '期間' (Period) is '2015年8月3日(月) から2015年12月11日(金) まで (弊社受注分まで)' (From August 3, 2015 (Monday) to December 11, 2015 (Friday) (up to our order volume)). The 'プロモーションコード' (Promotion Code) is 'OBFFPE (ご注文時に備考欄にお書きください)' (OBFFPE (please write in the remarks field when ordering)). The 'お問い合わせ先' (Contact) is '03-4578-2800'. At the bottom, there's a 'キャンペーンの最新情報はウェブから' (Get the latest campaign information from the web) and the website 'www.illumina.co.jp'. The Illumina logo is at the bottom right.

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Breast
Melanoma
Colon
Prostate
Lung
Ovary
Gastric

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


FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

TruSight® Tumor Panel (15 Genes)

AMHSminor

Sample ID: AMHSminor Report Definition: TST_15Genes-ReportDefinition-v1
Run Name: Report Generated On: 2015-06-07 20:59
Run ID: 141219_M00894_0164_000000000-ACK09

Variants identified as specified in the report definition

Detected SNVs, Insertions, and Deletions

Gene	Amino Acid Change	Variant Type	Nucleotide Change	Variant Frequency	Transcript
KRAS	p.Gly13Asp	missense variant	c.38G>A	0.541	ENST00000256078
PIK3CA	p.Glu545Lys	missense variant	c.1633G>A	0.495	ENST00000263967
EGFR	p.Leu861Gln	missense variant	c.2582T>A	0.532	ENST00000263967

No Calls

Gene	Chromosome	Coordinate	Failed Filter
NRAS	chr1	115252203	Low Coverage,Low Genotype Quality
NRAS	chr1	115252204	Low Coverage,Low Genotype Quality
TP53	chr17	7579651	Low Variant Frequency
TP53	chr17	7579654	Low Variant Frequency
TP53	chr17	7579660	Low Variant Frequency
TP53	chr17	7579666	Low Variant Frequency
TP53	chr17	7579671	Low Variant Frequency
TP53	chr17	7579679	Low Variant Frequency
TP53	chr17	7579684	Low Variant Frequency,Low Genotype Quality

illumina®

1 of 2

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

TruSight® Tumor Panel (15 Genes)

AMHSminor

Sample ID: AMHSminor Report Definition: TST_15Genes-ReportDefinition-v1
Run Name: Report Generated On: 2015-06-07 20:59
Run ID: 141219_M00894_0164_000000000-ACK09

Report Details

Library IDs: A04-AMHSminor-MxA-54056-1001,B04-AMHSminor-MxB-54056-1001
Library Index Pairs: R702/A501,R702/A502
Instrument Name: M00894
Sequencing Start Date: 2014-12-19 00:00
Sequenced By User: svc_cpuser
RTA Version: 1.18.54
MSR Version: 2.6.1.16
MCS Version: 2.4.1.3

Detected Variants

Please refer to the MiSeq Reporter TruSight Tumor Panel (15 Genes) Workflow reference guide for more information regarding the filtered variant results displayed in this report. Refer to the corresponding sample genome VCF file for additional unfiltered variant information.

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
Figure 1: Comprehensive Analysis and Interpretation of Variant Data

The VariantStudio tool is a powerful software application for analyzing and interpreting variant data. This tool aggregates information from a collection of databases to streamline annotation. It also provides flexible filtering options for analyzing variant data and tools to enable classification and reporting of actionable variants.

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


TruSeq Stranded Total RNA
TruSeq RNA Access

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
GENES, CHROMOSOMES & CANCER 54:500–505 (2015)

Identification of a Novel *PARP14-TFE3* Gene Fusion from 10-Year-Old FFPE Tissue by RNA-Seq

Weihua Huang,¹ Michael Goldfischer,² Sabina Babyeva,¹ Yong Mao,³ Konstantin Volyanskyy,³ Nevenka Dimitrova,³ John T. Fallon,¹ and Minghao Zhong^{1*}

¹Department of Pathology, New York Medical College, Westchester Medical Center, Valhalla, NY

²Department of Pathology, Hackensack University Medical Center, Hackensack, NJ

³Clinical Informatics Solutions and Services (CISS), Philips Research North America, Briarcliff Manor, NY

Xp11 (*TFE3*) translocation renal cell carcinoma (RCC) is officially recognized as a distinct subtype of RCC in the 2004 WHO classification. This neoplasm is characterized by several chromosomal translocations between the *TFE3*-involving Xp11.2 breakpoint and various fusion partners. To date, five partner genes have been identified, that is, *PRCC* in 1q21, *PSF* in 1q34, *ASPL* in 17q25, *CLTC* in 17q23, and *NONO* in Xq12; and three additional translocations have been reported with no partner gene being defined: t(X;3)(p11;q23), t(X;10)(p11;q23), and t(X;19)(p11;q13). Here, we report the identification of a novel *TFE3* fusion partner, *PARP14* in chromosome band3q21. We used RNA-seq on a 10-year-old FFPE (formalin-fixed, paraffin-embedded) tissue sample, which carried t(X;3)(p11;q23) as detected in the original cytogenetic study. The fusion transcript connected the 5'-end of the first two exons of *PARP14* to the 3'-end of five exons of *TFE3*, which was verified by reverse transcription PCR (RT-PCR) and Sanger sequencing. Similar to other *TFE3* fusions previously reported, the predicted *PARP14-TFE3* product retains the nuclear localization and DNA-binding domains of *TFE3*. This finding expands the list of *TFE3* translocation partner genes and re-emphasizes the essential oncogenic role of *TFE3* fusion proteins in this tumor. Our result also clearly demonstrated the feasibility of identifying chromosomal translocation by RNA-seq in clinical FFPE, which are easily accessible and associated with valuable clinical information. © 2015 Wiley Periodicals, Inc.

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