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TruSight HLA Assign 2.0 72-Sample Panel on MiSeq System Project File

TruSight HLA Assign 2.0 72-Sample Panel on MiniSeq System Project File

TruSight HLA v2 Sequencing Panel Reference Guide

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- Robinson J, Halliwell JA, Hayhurst JH, et al. The IPD and IMGT/HLA database: allele variant databases. *Nucleic Acids Res.* 2015;43:D423-431.
- Marsh SGE, Albert ED, Bodmer WF, et al. Nomenclature for factors of the HLA system, 2010. *Tissue Antigens.* 2010;75:291-455.
- Cough SCL and Simmonds MJ. The HLA region and Autoimmune Disease: Associations and Mechanisms of Action. *Curr Genomics.* 2007;8:453-465.
- Woodrow JC and Ilchysyn A. HLA antigens in psoriasis and psoriatic arthritis. *J Med Genetics.* 1985;22:492-495.
- Fogdell A, Hillert J, Sachs C, and Olerup O. The multiple sclerosis- and narcolepsy-associated HLA class II haplotype includes the DRB5*0101 allele. *Tissue Antigens.* 1995;46:333-336.
- Woodrow JC and Ilchysyn A. HLA antigens in psoriasis and psoriatic arthritis. *J Med Genetics.* 1985;22:492-495.
- Carrington M, Nelson GW, Martin MP, et al. HLA and HIV-1: Heterozygote Advantage and B*35-Cw*04 Disadvantage. *Sci.* 1999;283:1748-1752.
- Mallal S, Nolan D, Witt C, et al. Association between presence of HLA-B*5701, HLA-DR7, and HLA-DQ3 and hypersensitivity to HIV-1 reverse-transcriptase inhibitor abacavir. *Lancet.* 2002;359:727-732.
- Ferrell PB and McLeod HL. Carbamazepine, HLA-B*1502 and risk of Stevens-Johnson syndrome and toxic epidermal necrolysis: US FDA recommendations. *Pharmacogenomics.* 2008;9:1543-1546.
- Daly AK, Donaldson PT, Bhatnagar P, et al. HLA-B*5701 genotype is a major determinant of drug-induced liver injury due to flucloxacillin. *Nat Genetics.* 2009;41:816-819.
- Alfirevic A and Pirmohamed M. Predictive Genetic Testing for Drug-Induced Liver Injury: Considerations of Clinical Utility. *Clin Pharmacol Ther.* 2012;92:376-380.
- Romano A, De Santis A, Romito A, et al. Delayed Hypersensitivity to Aminopenicillins Is Related to Major Histocompatibility Complex Genes. *Ann Allergy Asthma Immunol.* 1998;80:433-437.
- Lonjou C, Borot N, Sekula P, et al. A European study of HLA-B in Stevens-Johnson syndrome and toxic epidermal necrolysis related to five high-risk drugs. *Pharmacogenet and Genomics.* 2008;18:99-107.
- Dekker JW, Nizankowska E, Schmitz-Schumann M, et al. Aspirin-induced asthma and HLA-DRB1 and HLA-DPB1 genotypes. *Clin Exp Allergy.* 1997;27:574-577.
- International Histocompatibility Working Group Website, Reference Panels. <http://ihwg.org/reference/index.html>. Accessed April 18, 2016.
- Mack SJ, Cano P, Hollenbach JA, et al. Common and well-documented HLA alleles: 2012 update to the CWD catalogue. *Tissue Antigens.* 2013;81:194-203.
- Wittig M, Anmarkrud JA, Kassens JC, et al. Development of a high-resolution NGS-based HLA-typing and analysis pipeline. *Nucleic Acids Res.* 2015;43:e70.

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© 2016 Illumina, Inc. All rights reserved.