# MiSeq Reporter (MSR)をはじめよう

癸生川絵里 (Eri Kibukawa)  
Bioinformatics Support Scientist

MiSeq Reporter

© 2012 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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Mozilla Firefox window showing the Illumina MiSeq Reporter homepage. The address bar displays [www.illumina.com/help/miseq\\_reporter/default.htm](http://www.illumina.com/help/miseq_reporter/default.htm). The page content includes a navigation menu on the left with links for Home, The Toolbar, Jobs, MiSeq Files, Technical Assistance, and Revision History. The main content area features a welcome message, a navigation bar with BACK and FORWARD buttons, and a table of contents on the left. The footer contains copyright information for Illumina, Inc. and contact details.

Summary Information page in the Illumina MiSeq Reporter. The page title is 'Summary Information'. It contains a table of contents on the left and a main content area with a 'Table: Job Information' and a 'Table: High Percentage Graph Information'.

Row	Description
Analysis Folder	The full path to the analysis output folder on-instrument
Completion Time	The time the experiment was completed
Copy Folder	The full path to the network folder to which analysis output is copied
Cycles	A string describing the number of cycles in each read and whether the read is an index read. For example, 151 [N] 8[Y] 151[N] indicates a first read of 151 cycles that is not an index read, a second read of 8 cycles that is not an index read, and a third read of 151 cycles that is not an index read.
Experiment Name	A user-supplied descriptive name of the experiment.
Investigator	The name of the investigator (optional).
Project	An optional user-supplied descriptive name of the project.
Run Folder	The output folder produced by RTA primary analysis software. This folder contains all primary and secondary analysis results for the run.
Start Time	The clock time that the secondary analysis was started.

y-axis	x-axis	Description
Percent	20/ 1 1	The ratio of intensities at Cycle 20 to the intensities at Cycle 1 for Read 1
	20/ 1 2	The ratio of intensities at Cycle 20 to the intensities at Cycle 1 for Read 2
Align 1		The <b>percentage</b> of clusters that aligned to the reference genome in Read 1. This field is not available for Assembly, Metagenomics, and SmallRNA workflows.
Align 2		The <b>percentage</b> of clusters that aligned to the reference in Read 2. This field is not available for Assembly, Metagenomics, and SmallRNA workflows.
PE Orientation		The <b>percentage</b> of paired-end alignments that have the expected orientation
PE Resynthesis		The ratio of first cycle intensities for Read 1 to first cycle intensities for Read 2
PF		The <b>percentage</b> of clusters passing filters

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




The screenshot displays the MiSeq Reporter interface. On the left, there are sections for 'Current Jobs' (showing an 'Idle' job) and 'Completed Jobs' (listing 'DenovoAssembly', 'Resequencing', and 'SmallRNA'). The main area shows a 'Status' tab selected, which is circled in red. This tab displays a table of job details. A 'STOP JOB' button is visible at the bottom right of the status window. In the top right corner of the interface, a wrench icon is also circled in red.

Status:	
<b>Experiment Name</b>	2x151
<b>Project</b>	MiSeq
<b>Investigator</b>	Shamili
<b>Run Folder</b>	111104_M00102_0074_AMS0003254-00300
<b>Start Time</b>	11/21/2011 4:23:17 PM
<b>Completion Time</b>	11/21/2011 6:04:21 PM
<b>Cycles</b>	151 [N], 151 [N]
<b>Analysis Folder</b>	
<b>Copy Folder</b>	

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The screenshot displays the MiSeq Reporter software interface. On the left, there are sections for 'Current Jobs' (showing an 'Idle' job) and 'Completed Jobs' (listing jobs like 'DenovoAssembly', 'Resequencing', and 'SmallRNA'). The main area shows a tabbed interface with 'Sample Sheet' selected and circled in red. The 'Sample Sheet' tab displays the following information:

[Header]	
Investigator Name	Shamili
Project Name	MISeq
Experiment Name	2x151
Date	9/21/2011
Workflow	Resequencing
Chemistry	Default

[Reads]	
151	
151	


[Manifests]	
A	ManifestNameHere

At the bottom right of the main area, there is a 'STOP JOB' button. The interface also includes a toolbar with icons for home, refresh, settings (circled in red), print, and help.

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
The screenshot displays the MiSeq Reporter software interface. On the left, there are sections for 'Current Jobs' (showing an 'Idle' job) and 'Completed Jobs' (listing 'DenovoAssembly', 'Resequencing', and 'SmallRNA'). The main area shows a job window with tabs for 'Status', 'Sample Sheet', 'Logs', and 'Errors'. The 'Logs' tab is selected and circled in red. The log text shows the start of a run for analysis D:\Illumina\MiSeqAnalysis\111104\_M00102\_0074\_AMS0003254-00300, including details about genome folders, lane 1, tile groups, and FASTQ file building. A 'STOP JOB' button is visible at the bottom right of the log window.

```
11/6/2011,04:01:12.007,MiSeq Reporter 1.1.1.1: Starting Logging for
111104_M00102_0074_AMS0003254-00300
11/6/2011,04:01:12.074,Starting Execution for Analysis D:\Illumina\MiSeqAnalysis
\111104_M00102_0074_AMS0003254-00300
11/6/2011,04:01:12.504,Genome folders: C:\Illumina\MiSeq Reporter\Genomes\PhiX
11/6/2011,04:01:12.927,Queueing FASTQ file(s) for lane 1, tile group 1 (tiles: 1,2,3,4,5,6,7)
11/6/2011,04:01:12.933,Queueing FASTQ file(s) for lane 1, tile group 2 (tiles:
8,9,10,11,12,13,14)
11/6/2011,04:01:13.083,Building FASTQ file for lane 1, tile group 2, read 1
11/6/2011,04:01:13.084,Building FASTQ file for lane 1, tile group 1, read 1
11/6/2011,04:06:24.563,Building FASTQ file for lane 1, tile group 1, read 2
11/6/2011,04:06:25.048,Building FASTQ file for lane 1, tile group 2, read 2
11/6/2011,04:11:26.892,Invoking alignment wrapper
```

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
The screenshot displays the MiSeq Reporter software interface. On the left, there are sections for 'Current Jobs' (showing an 'Idle' status) and 'Completed Jobs' (listing tasks like DenovoAssembly, Resequencing, and SmallRNA). The main area features a tabbed interface with 'Status', 'Sample Sheet', 'Logs', and 'Errors' tabs. The 'Errors' tab is highlighted with a red circle. Above the tabs, a toolbar contains icons for home, refresh, settings (a wrench icon circled in red), print, and help. The 'Errors' tab content shows a log entry for a fatal processing error: 'Fatal Processing Error for Analysis Runs\111104\_M00102\_0074\_AMS0003254-00300: Error - GenomeFolder 'PhiX' was not found. Please check the sample sheet.' Below the log text is a stack trace and a 'STOP JOB' button.

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
The screenshot displays the MiSeq Reporter software interface. On the left, the 'Current Jobs' section shows an 'Idle' job with a green progress bar. Below it, the 'Completed Jobs' section lists three jobs: 'DenovoAssembly' (2x151, 12/3/2011 5:13:21 PM), 'Resequencing' (2x151, 11/21/2011 6:04:21 PM), and another 'Resequencing' job (1x151 Run, 2/7/2012 3:43:52 PM). The main area shows a 'Status' window with tabs for 'Status', 'Sample Sheet', 'Logs', and 'Errors'. The 'Errors' tab is active, displaying a log entry: '11/21/2011,16:13:39.836,Fatal Processing Error for Analysis C:\Users\Runs\111104\_M00102\_0074\_AMS0003254-00300: Error - GenomeFolder 'PhiX' was not found. Please check the sample sheet. Stack Trace: at Illumina.SecondaryAnalysis.SampleSheet.GetGenomes(String GenomeRepository) in c:\Projects\Illumina.Bolt\Dev\1.1'. A 'STOP JOB' button is visible at the bottom right of the error log window. A red circle highlights a button in the top right corner of the interface, which is the 'Change Repository' button.

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
- ▶ D:\¥Illumina¥MiSeqAnalysis¥当該ランフォルダ名¥Data¥Intensities¥BaseCalls

Name	Size	Type
Alignment		File Folder
Alignment2		File Folder
L001		File Folder
Matrix		File Folder
Phasing		File Folder
config.xml	6 KB	XML Document
Sample1_S1_L001_R1_001.fastq.gz	38,666 KB	WinZip File
Sample1_S1_L001_R2_001.fastq.gz	40,738 KB	WinZip File
Sample2_S2_L001_R1_001.fastq.gz	31,584 KB	WinZip File
Sample2_S2_L001_R2_001.fastq.gz	33,492 KB	WinZip File
SampleSheet.csv	3 KB	Microsoft Office Exc...
TruSeq_Amplicon_Cancer_Panel_Manifest.txt	101 KB	Text Document
TruSeq_Custom_Amplicon_Control_Manifest.txt	144 KB	Text Document
TSCA-CONTROL_S17_L001_R1_001.fastq.gz	55,609 KB	WinZip File
TSCA-CONTROL_S17_L001_R2_001.fastq.gz	60,355 KB	WinZip File
Undetermined_S0_L001_R1_001.fastq.gz	20,533 KB	WinZip File
Undetermined_S0_L001_R2_001.fastq.gz	21,945 KB	WinZip File

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MiSeq Reporter

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The screenshot shows a Mozilla Firefox browser window with the address bar displaying [www.illumina.com/help/miseq\\_reporter/default.htm](http://www.illumina.com/help/miseq_reporter/default.htm). The page content includes a navigation menu on the left with items: Home, The Toolbar, Jobs, MiSeq Files, Technical Assistance, and Revision History. Below this menu are sections for Table of Contents, Search, and Favorites. The main content area features the heading "Welcome to Illumina MiSeq Reporter" in orange, followed by "Introduction" and "MiSeq Concepts" in blue. There are "BACK" and "FORWARD" navigation buttons. At the bottom, it states "© 2011 Illumina, Inc. All rights reserved" and provides the contact URL [www.illumina.com/contact](http://www.illumina.com/contact). The Illumina logo is visible in the bottom right corner of the page.

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- RTAComplete.txt
- SampleSheet.csv
- RunInfo.xml

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