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© 2012 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GALiX, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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
The screenshot displays the 'RUN CONFIGURATION' tab in the Illumina sequencing software. The 'Align Lanes' section is visible, showing a grid of 8 lanes, each with a checkbox. All checkboxes are checked, indicating that all lanes are selected for alignment. A red arrow points from the text on the left to the 'Align Lanes' section. To the right of the interface is a diagram of a flow cell with 8 vertical lanes.

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Level	Yield Total (G)	Projected Total Yield (G)	Yield Perfect (G)	Yield <=3 errors (G)	Aligned (%)	% Perfect [Num Cycles]	% <=3 errors [Num Cycles]	Error Rate (%)	Intensity Cycle 1	% Intensity Cycle 20	% >= Q30
Read 1	160.8	160.8	119.3	154.9	96.63	75.5 [100]	98.0 [100]	0.54	8271	76.2	93.3
Read 2	160.8	160.8	136.8	153.6	97.25	87.3 [100]	98.0 [100]	0.47	7254	76.4	91.4
Total	321.7	321.7	256.1	308.5	96.94	81.4	98.0	0.51	7762	76.3	92.3

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
RUN CONFIGURATION PRE-RUN SETUP INITIATE RUN

Scan Recipe Reagents Review

Please enter "Flowcell ID".

Flow Cell ID: Experiment: User Name: SBSUser

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Browse

Refresh

Analysis Imaging Summary Tile Status TruSeq Controls Indexing

Run Summary

Level	Yield Total (G)	Projected Total Yield (G)	Yield Perfect (G)	Yield <=3 errors (G)	Aligned (%)	% Perfect [Num Cycles]	% <=3 errors [Num Cycles]	Error Rate (%)	Intensity Cycle 1	% Intensity Cycle 20	% >= Q30
Read 1	68.5	68.5	6.6	8.0	97.23	82.3 [50]	99.6 [50]	0.47	5698	91.9	92.9
Read 2 (I)	8.2	8.2	0.0	0.0	0.00	0.0 [6]	0.0 [6]	0.00	5303	0.0	78.1
Total	76.7	76.7	6.6	8.0	97.23	82.3	99.6	0.47	5500	91.9	91.1

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▶ PhiX TechNote:

- http://www.illumina.com/documents/product/s/technotes/technote_phixcontrolv3.pdf

▶ Illumina homepage

- <http://www.illumina.com/>

▶ Sequencing Support

- <http://support.illumina.com/sequencing.ilmn>

Technical Note: Sequencing

illumina

Using a PhiX Control for HiSeq® Sequencing Runs

A low-concentration spike-in of Illumina PhiX Control v3 provides quality and calibration controls.

Highlights

- For most libraries, Illumina recommends using PhiX Control v3 (Catalog # FC-110-3001) in a low-concentration spike-in (1%) for improved sequencing quality control
- For samples with bias base composition, use a dedicated PhiX control lane to improve cross-talk and phasing calculations
- For samples with low diversity, use a dedicated control lane and a high-concentration spike-in (40% or higher) of PhiX to create a more diverse set of clusters

What are PhiX Control Libraries?

Control libraries generated from the PhiX virus serve as an effective control in sequencing runs. Characteristics of the PhiX genome provide several benefits:

- **Small**—PhiX is a small genome, which enables quick alignment and estimation of error rates.
- **Diverse**—The PhiX genome contains approximately 45% GC and 55% AT.
- **Well-Defined**—PhiX has a well-defined genome sequence.

Illumina cluster detection algorithms are optimized around a balanced representation of A, T, G, and C nucleotides. Illumina PhiX Control v3, with an insert size of 425–525 bp, is a balanced and diverse library that can help mitigate sequencing challenges in unbalanced and low diversity libraries.

Benefits of Using a PhiX Control

PhiX libraries provide a quality control for cluster generation, sequencing, and alignment, and a calibration control for cross-talk matrix generation, phasing, and prephasing.

Cluster Generation

Illumina PhiX Control v3 is shipped as a ready-to-use 10 nM library that can be used as a positive control in the clustering process. If a problem occurs in sample preparation, PhiX will still generate clusters. This helps you discern whether a lack of clusters is due to sample preparation failure or a failure in the cluster generation process.

Cross-Talk Matrix Generation

During an Illumina sequencing run, the cross-talk due to spectral overlap between the four fluorescently-labeled nucleotides is calculated during template generation (cluster identification) in cycles 1–4.

For proper cross-talk calculation, HiSeq Control Software (HCS) v1.3.8, or higher, requires approximately equal numbers and at least 50,000 clusters/mm² in each of the four bases. Therefore, it is imperative to have a balanced representation of bases at the beginning of each read, excluding the index read.

Phasing and Prephasing

During sequencing by synthesis, each DNA strand in a cluster extends by one base per cycle. A small proportion of strands may become out of phase with the current cycle, either falling a base behind (phasing) or jumping a base ahead (prephasing). The phasing and prephasing rates define the fraction of molecules that become phased or prephased per cycle.

Calculation of these rates requires a balanced and random base composition in cycles 2–12. Any library that does not comprise a balanced base composition (e.g., initial cycle indexing, restriction enzyme libraries) should use a control lane.

High GC samples typically show higher phasing rates. However, if the sample has good diversity (for example, whole-genome sequencing libraries with 40–60% GC), it does not require a control lane.

Alignment

Because PhiX has a small, well-defined genome sequence, it is an excellent alignment control. If a lane containing PhiX (control lane or low-concentration spike-in) is designated for your run, Real Time Analysis (RTA) software aligns complete sequences to the PhiX reference beginning after the 25th cycle is accumulated and calculates error rates, providing an indication of sequencing success during the run.

Using a PhiX Control for Unbalanced Samples

Most mammalian genomic or whole-transcriptome RNA samples have a balanced genomic composition (approximately equal proportions of A, T, G, and C). These samples do not require a dedicated control lane to generate accurate matrix and phasing estimations. Unbalanced samples contain genomes with high AT or GC content (less than 40% or greater than 60%). Some examples are Arabidopsis, Plasmodium, some bacteria, and bisulfite conversion studies.

Use a dedicated control lane with unbalanced samples to improve cross-talk and phasing calculations.

Using a PhiX Control for Low Diversity Libraries

Low diversity libraries are libraries where a significant number of the reads have the same sequence. This shifts the base composition because the reads are no longer random.

Low diversity can occur with some expression studies with > 25%

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