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© 2013 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners. The Genetic Energy streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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
1a. Abundance estimation (**genes and transcripts**), fusion detection, and variant calling workflow.


TopHat Alignment

1b. Creating and merging assemblies and differential expression calculation.

Cufflinks Assembly/
Diff. Exp

2. Fast Alignment and differential expression calculation of **genes**

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

http://support.illumina.com/sequencing/sequencing_software/basespace/documentation.ilmn

The screenshot shows the Illumina support website interface. At the top left is the Illumina logo. A navigation bar contains links for APPLICATIONS, SYSTEMS, CLINICAL, SERVICES, SCIENCE, SUPPORT (circled), and COMPANY. A search box is on the right. Below the navigation bar is a breadcrumb trail: Support > Sequencing > Sequencing Software > BaseSpace > Documentation & Literature (circled). On the left, a 'BaseSpace Support' sidebar lists Overview, Computing Requirements, and Documentation & Literature (circled). The main content area is titled 'Documentation & Literature' and includes links for Documentation, Product Literature, and Tech Notes. Under a 'DESCRIPTION' section, there is a list of user guides: BaseSpace User Guide (expanded), BaseSpace Core App User Guides (collapsed), and several specific app guides: Instructions for using the BWA Whole-Genome Sequencing App (15050952 A), Instructions for using the BWA Enrichment App (15050958 A), Instructions for using the Tumor Normal App (15050950 A), and Instructions for using the TopHat and Cufflinks RNA-seq apps (15050962 A).

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http://support.illumina.com/help/BaseSpace_App_enr_BWA_help/BWAEnrichmentHelp.htm
http://support.illumina.com/help/BaseSpace_App_enr_Isaac_help/IsaacEnrichmentHelp.htm

The screenshot shows the Illumina support website interface. The browser address bar displays the URL: support.illumina.com/help/BaseSpace_App_enr_BWA_help/BWAEnrichmentHelp.htm#Vault/Informatics/Sequencing. The page title is "Running BWA Enrichment".

Running BWA Enrichment

- 1 Navigate to the project or sample that you want to analyze.
- 2 Click the **Apps** button and select **BWA Enrichment** from the dropdown list.
- 3 Read the End User License Agreement and permissions, and click **Accept** if you agree.
- 4 Fill out the required fields in the BWA Enrichment input form:
 - a **Analysis Name:** Provide the analysis name. Default name is the app name with the date and time the analysis was started.
 - b **Save Results To:** Select the project that stores the app results.
 - c **Sample(s):** Browse to the sample you want to analyze, and select the check box. You can analyze multiple samples.
 - d **Reference Genome:** Select the reference genome. Currently, you can only use hg19.
 - e **Targeted Regions:** Select the targeted region of your enrichment.
 - f **Base Padding:** Select the padding you want. Padding defines the amount of sequence immediately upstream and downstream of the targeted regions that is also used in enrichment analysis.
 - g **Annotation:** Choose which gene and transcript annotation reference database to use.
- 5 If desired, fill out the advanced fields in the BWA Enrichment input form:
 - a **Depth Threshold:** The GATK variant caller filters variants if the coverage depth at that location is less than the specified threshold. Decreasing this value will increase variant calling sensitivity, but raise the risk of false positives. The variant caller reports the variants, but filters them in the VCF files by adding a LowDP flag. Default value for GATK: 5.
 - b **Trim Nextera Rapid Capture Adapters:** If selected, Nextera Rapid Capture adapters will be trimmed. Use this setting only if not already applied as a sample sheet setting.
 - c **Flag PCR Duplicates:** If selected, PCR duplicates are flagged in the BAM files and not used for variant calling. PCR duplicates are defined as two clusters from a paired-end run where both clusters have the exact same alignment positions for each read. Optical duplicates are already filtered out during RTA processing.
 - d **Generate Picard HS Metrics:** If selected, Picard HS metrics are generated. See [Picard Metrics](#) for more information.

Figure 2 BWA Enrichment Input Form

The input form shows the following fields:

- Analysis Name: BWA Enrichment 12/02/2013 1:49:23
- Save Results To: Select a Project
- Sample(s): Select Sample
- Reference Genome: Human (JCSC hg19)
- Targeted Regions: Nextera Rapid Capture Exom
- Base Padding: 150
- Annotation: RefSeq @ Ensembl

There is an "Advanced" section at the bottom of the form.

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- Accurate whole human genome sequencing using reversible terminator chemistry.
Nature 456: 53-59 [PMID: 18987734]

▶ RNA-Seq, RPKM

- Mapping and quantifying mammalian transcriptomes by RNA-Seq
Nature Methods, Volume 5, 621 – 628 [PMID: 18516045]

▶ BaseSpace RNA-Seq アプリ

- TopHat: discovering splice junctions with RNA-Seq.
Bioinformatics. 2009 May 1;25(9):1105-11 [PMID: 19289445]
- Cufflinks/Cuffdiff
Transcript assembly and quantification by RNA-Seq reveals unannotated transcripts and isoform switching during cell differentiation. Nat Biotechnol. 2010 May;28(5):511-5 [PMID: 20436464]
- Cuffdiff2 Differential analysis of gene regulation at transcript resolution with RNA-seq
Nature Biotechnology 31, 46–53 (2013) [PMID: 23222703]
- TopHat-Fusion: an algorithm for discovery of novel fusion transcripts.
Genome Biol. 2011 Aug 11;12(8):R72. [PMID: 21835007]
- **STAR: ultrafast universal RNA-seq aligner.**
Bioinformatics, 2013 Jan 1;29(1):15-21. [PMID: 23104886]
- **DESeq**
Differential expression analysis for sequence count data.
Genome Biol. 2010;11(10):R106. [PMID: 20979621]
DESeq2: www.bioconductor.org/packages/2.13/bioc/html/DESeq2.html

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<http://www.illumina.com/landing/basespace-core-apps-for-rna-sequencing/>
<http://res.illumina.com/documents/products/technotes/technote-basespace-rna-seq.pdf>
http://support.illumina.com/help/BaseSpace_App_RNAseq_help/RNAseq_Apps_Help.htm
<http://www.illumina.com/applications/sequencing/rna.ilmn>

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© 2014 Illumina, Inc. All rights reserved.

Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

illumina®