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© 2012 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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16S Metagenomic Sequencing Library Preparation

Preparing 16S Ribosomal RNA Gene Amplicons for the Illumina MiSeq System

Introduction	2
16S Library Preparation Workflow	5
Amplicon PCR	6
PCR Clean-Up	8
Index PCR	10
PCR Clean-Up 2	13
[Optional] Validate Library	15
Library Quantification, Normalization, and Pooling	16
Library Denaturing and MiSeq Sample Loading	17
MiSeq Reporter Metagenomics Workflow	20
Supporting Information	21

IMPORTANT NOTICE This document provides information for an application for Illumina technology that has been demonstrated internally and may be of interest to customers. This information is provided as-is and is not an Illumina product and is not accompanied by any rights or warranties. Customers using or adapting this information should obtain any licenses required and materials from authorized vendors. Illumina products mentioned herein are for research use only unless marked otherwise. While customer feedback is welcomed, this application is not supported by Illumina Technical Support and Field Application Scientists.

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
Noninvasive Identification and Monitoring of Cancer Mutations by Targeted Deep Sequencing of Plasma DNA

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Table S1: Target-specific primers. Each target-specific primer consists of a universal 5' end and a target-specific 3' end as follows:

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ForsheW T, et al., Noninvasive identification and monitoring of cancer mutations by targeted deep sequencing of plasma DNA. *Sci Transl Med.* 2012 May 30;4(136):136ra68.

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Table S2: Unique sequencing barcodes. Platform-specific adaptors and barcodes are attached through PCR following the single-plex amplification step. The primers consisted of the PE1 and PE2 sequences for Illumina cluster generation, a 10-bp barcode, and the CS1 and CS2 adaptors, used in pairs: PE1-CS1 with PE2-BC-CS2, and PE1-CS2 with PE2-BC-CS1.

P5配列に相当 : AATGATACGGCGACCACCGAGATCT. AACTGACGACATGGTTCTACA
 P7配列に相当 : CAAGCAGAAGACGGCATAACGAGAT. Barcode TACGGTAGCAGAGACTTGGTCT

Barcode name	Barcode sequence	FLD0060	GCACGTAGCT	FLD0121	CAGAGCTAGT
FLD0001	GTATCGTCGT	FLD0061	TCACGCTATG	FLD0122	CGCAGAGCAT
FLD0002	GTGTATGCGT	FLD0062	CGTACTACGT	FLD0123	TGTACAGCGA
FLD0003	TGCTCGTAGT	FLD0063	CAGCTGAGTA	FLD0124	ACGTCAGTAT
		FLD0064	GAGATCAGTC	FLD0125	TCACAGCATA

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Forshe T, et al., Noninvasive identification and monitoring of cancer mutations by targeted deep sequencing of plasma DNA. *Sci Transl Med.* 2012 May 30;4(136):136ra68.

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Open

The ISME Journal (2012), 1–4
© 2012 International Society for Microbial Ecology. All rights reserved. 1751-7362/12
www.nature.com/ismej

SHORT COMMUNICATION

Ultra-high-throughput microbial community analysis on the Illumina HiSeq and MiSeq platforms

J Gregory Caporaso¹, Christian I. Lauber², William A Walters³, Donna Berg-Lyons⁴, James Huntley⁴, Noah Fierer^{2,5}, Sarah M Owens⁶, Jason Betley⁷, Louise Fraser⁸, Markus Bauer⁹, Niall Gormley⁷, Jack A Gilbert¹⁰, Geoff Smith¹¹ and Rob Knight^{12,13}

¹Department of Computer Science, Northern Arizona University, Flagstaff, AZ, USA; ²Cooperative Institute for Research in Environmental Sciences, UCB 216, University of Colorado, Boulder, CO, USA; ³Department of Molecular, Cellular and Developmental Biology, UCB 347, University of Colorado, Boulder, CO, USA; ⁴Colorado Initiative in Molecular Biotechnology, UCB 347, University of Colorado, Boulder, CO, USA; ⁵Department of Ecology and Evolutionary Biology, UCB 334, University of Colorado, Boulder, Colorado, USA; ⁶Argonne National Laboratory, Argonne, IL, USA; ⁷Illumina Cambridge Ltd., Chesterford Research Park, Saffron Walden, Essex, UK; ⁸Department of Ecology and Evolution, University of Chicago, Chicago, IL, USA; ⁹Department of Chemistry and Biochemistry, UCB 215, University of Colorado, Boulder, CO, USA and ¹⁰Howard Hughes Medical Institute, University of Colorado at Boulder, UCB 215, Boulder, CO, USA

DNA sequencing continues to decrease in cost with the Illumina HiSeq2000 generating up to 600 Gb of paired-end 100 base reads in a ten-day run. Here we present a protocol for community amplicon sequencing on the HiSeq2000 and MiSeq Illumina platforms, and apply that protocol to sequence 24 microbial communities from host-associated and free-living environments. A critical question as more sequencing platforms become available is whether biological conclusions derived on one platform are consistent with what would be derived on a different platform. We show that the protocol developed for these instruments successfully recaptures known biological results, and additionally that biological conclusions are consistent across sequencing platforms (the HiSeq2000 versus the MiSeq) and across the sequenced regions of amplicons.

The ISME Journal advance online publication: 8 March 2012; doi:10.1038/ismej.2012.8

Subject Category: microbial ecology and functional diversity of natural habitats

Keywords: illumina; barcoded sequencing; QIIME

DNA sequencing cost continues to decline: a vast price per sequence decrease on Illumina HiSeq2000 and MiSeq platforms further supports democratization of sequencing (Trings and Gupenholzer, 2009). Interest in amplicon sequencing on Illumina is growing (Bartram *et al.*, 2011; Caporaso *et al.*, 2011; Zhou *et al.*, 2011), largely due to lower cost per sequence than other platforms, enabling high-throughput microbial ecology at the greatest coverage yet possible. Although some technical issues exist with community sequencing, such as PCR primer biases and differential DNA extraction efficiency from different organisms in complex communities, these techniques continue to vastly expand our understanding of the microbial world.

Here we present an amplicon sequencing protocol for the HiSeq2000 and MiSeq platforms, and apply this protocol to sequence host-associated and free-living microbial communities to verify that biological conclusions drawn from the data are consistent across platforms and sequence reads. The HiSeq and MiSeq platforms differ markedly in scale. The HiSeq2000 produces >50 Gb per day, and in the course of a 10.8 day run produces 1.6 billion 100-base paired-end reads. By contrast, the MiSeq is for single-day experiments, and generates 1.5 Gb per day from 5 million 150-base paired-end reads. Our results capture known differences between microbial communities on each platform; biological conclusions drawn are consistent across platforms and sequence reads. This protocol is therefore ready for widespread use in microbial community analysis, such as by the Earth Microbiome Project (Gilbert *et al.*, 2010), which has adopted it for amplicon sequencing. Details on the sequencing protocol are provided as Supplementary Methods.

Twenty-four samples were sequenced on three paired-end Illumina HiSeq2000 lanes, and in one paired-end MiSeq run. The samples represented soil (source: USA; n=8) and several host-associated environment types: human feces (source: USA;

Correspondence: R Knight, Howard Hughes Medical Institute, University of Colorado at Boulder, UCB 215, Boulder, CO 80309, USA.
E-mail: rsknight@colorado.edu
Received 12 September 2011; revised 13 January 2012; accepted 19 January 2012

Caporaso JG, et al., Ultra-high-throughput microbial community analysis on the Illumina HiSeq and MiSeq platforms. *ISME J.* 2012 Mar 8.

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Preparing Libraries
for Sequencing on the MiSeq®

FOR RESEARCH USE ONLY

Revision History	3
Introduction	4
Prepare HT1	6
Prepare a Fresh Dilution of NaOH	7
Denature and Dilute DNA	8
Prepare PhiX Control	11
Combine Sample Library and PhiX Control	13
Next Steps	14
Technical Assistance	

ILLUMINA PROPRIETARY
Part # 15039740 Rev. D
October 2013

illumina®

http://support.illumina.com/downloads/prepare_libraries_for_sequencing_miseq_15039740.ilmn

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