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© 2012 Illumina, Inc. All rights reserved.
Illumina, illuminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DASL, DesignStudio, Eco, GAIIX, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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© 2012 Illumina, Inc. All rights reserved.
Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPPro, DNAS, DesignStudio, Eco, GATK, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, Sentrix, SeqMonitor, Solexa, TruSeq, ViewCode, the jumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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
Nextera Mate Pair

TruSeq Custom
Amplicon

TruSeq Custom
Amplicon Cancer Panel

Nextera Rapid
Capture Exome

Nextera Rapid Capture
Custom Enrichment

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Below the grid, there is a section titled 'サンプル調製キット詳細' (Sample Preparation Kit Details) containing the following text:

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The screenshot shows the 'Sample Sheet Wizard - Workflow Parameters' window in the Illumina Experiment Manager. The window is divided into two main sections: 'TruSeq Amplicon Run Settings' and 'TruSeq Amplicon Workflow-Specific Settings'.

TruSeq Amplicon Run Settings:

- Reagent Cartridge Barcode*: [Empty text field]
- Sample Prep Kit: [TruSeq Amplicon dropdown menu]
- Index Reads: [Radio buttons for 0, 1, and 2; 2 is selected]
- Experiment Name: [Empty text field]
- Investigator Name: [Empty text field]
- Description: [Empty text field]
- Date: [2014/06/20 date picker]
- Read Type: [Radio buttons for Paired End and Single Read; Paired End is selected]
- Cycles Read 1: [151 spinner]
- Cycles Read 2: [151 spinner]
- * - required field

TruSeq Amplicon Workflow-Specific Settings:

- Use Somatic Variant Caller (Recommended for Cancer Panel)
- Variant Quality: [30 spinner]

At the bottom of the window, there are three buttons: 'Cancel', 'Back', and 'Next'.

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Illumina Experiment Manager

Sample Sheet Wizard - Sample Selection

Samples to include in sample sheet * - required field Maximize

Sample ID*	Sample Name	Plate	Well	Index1 (I7)*	I7 Sequence	Sample Project	Description
1	Test1			D701	ATTACTCG		
2	Test2			D703	CGCTCATT		
3	Test3			D705	ATTCAGAA		

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Sample Sheet Status: **Warning**

Reason: There is not enough diversity in the index cycles. Every cycle in the index reads must have at least one base in the green channel (G,T) and one in the red channel (A,C) to ensure proper image processing

Buttons: Add Blank Row, Remove Selected Rows, Cancel, Back, Finish

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EM Illumina Experiment Manager

Illumina Experiment Manager

Sample Sheet Wizard - Sample Selection

Samples to include in sample sheet * - required field Maximize

Sample ID*	Sample Name	Plate	Well	Index 1 (I7)*	I7 Sequence	Sample Project	Description
1	Test1			D701	ATTACTCG		
2	Test2			D702	TCGGGAGA		
3	Test3			D704	GAGATTCC		

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

Sequencing By Synthesis (SBS) 法 (サンガー法の改良法)

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



Q Score	Probability of error
Q30	1 in 1000
Q20	1 in 100

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© 2013 Illumina, Inc. All rights reserved.
Illumina, IlluminaDx, BaseSpace, BeadArray, BeadXpress, cBot, CSPRO, DASL, DesignStudio, Eco, GAlIx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, Infinium, iSelect, MiSeq, Nextera, NuPCR, SeqMonitor, Solexa, TruSeq, TruSight, VeraCode, the pumpkin orange color, and the Genetic Energy streaming bases design are trademarks or registered trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.

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


Analysis Type	Outputs	Outputs
<p>Sequencing (MCS)</p>		 <p>Images/TIFF files</p>
<p>Primary Analysis (RTA)</p>		 <p>Intensities Base Calling</p>
<p>Secondary Analysis (MSR or BaseSpace)</p>	 <p>MiSeq Reporter</p>	 <p>Alignments and Variant Detection</p>

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Reference Guide

http://support.illumina.com/sequencing/sequencing_software/misseq_reporter.ilmn

MiSeq Reporter Theory of Operation

http://support.illumina.com/sequencing/sequencing_software/misseq_reporter/documentation.ilmn

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The image displays a collection of report pages from MiSeq Reporter. On the left, there are two document thumbnails: one labeled 'PDF' and another labeled 'HTML', both showing the 'Resequencing Report' header for sample NA12878_S1. To the right, several overlapping report pages are shown, detailing various metrics:

- Sample Information:** Sample ID: NA12878_S1, Sample Name: NA12878_S1
- Read Statistics:** Table with columns: Read, Percent Aligned, Percent Q30, Mismatch Rate. Data: Read 1 (96.09%, 83.52%, 0.37%), Read 2 (95.42%, 81.01%, 0.46%).
- Coverage Histogram:** A line graph showing the distribution of read coverage.
- Small Variants Summary:** Table with columns: SNVs, Insertions, Deletions. Data: Total Passing (1,489,088 SNVs, 245,114 Insertions, 219,298 Deletions), Percent Found in dbSNP (99.27%, 83.09%, 87.87%), HetHom Ratio, TiTv Ratio.
- Structural Variants Summary:** Table with columns: SV Type, Total, Number in Genes. Data: CNV (74), Inversion (11), Tandem Duplication (3), Deletion (1), Inversion (1).
- Duplicate Information:** Percent Duplicate Pairs: 1.1%.
- Settings:** Table with columns: Setting Name, Value. Data: Reference Genome (Human (UCSC hg19)), Create SLCNV Calling (True), Annotation Source (Ensembl), Min GQ for Variants (30), Max Strand Bias for Variants (10), Flag PCR Duplicates (True).
- Software Versions:** Table with columns: Software, Version. Data: Isaac Whole Genome Sequencing (BaseSpace Workflow) (0.7.0.0), IAS (Analysis Software) (2.4.55.0), Isaac (Alignment) (01.13.10.21), SAMtools (0.1.10 (982.295)), Isaac Variant Caller (2.0.14), IAS (Annotation Service) (VEP 72.4).

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