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Cell

Resource

Expressed Pseudogenes in the Transcriptional Landscape of Human Cancers


Shanker Kalyana-Sundaram,^{1,2,6,7} Chandan Kumar-Sinha,^{1,2,7} Sunita Shankar,^{1,2} Dan R. Robinson,^{1,2} Yi-Mi Wu,^{1,2} Xuhong Cao,^{1,3} Irfan A. Asangani,^{1,2} Vishal Kothari,¹ John R. Prensner,^{1,2} Robert J. Lonigro,^{1,2} Matthew K. Iyer,¹ Terrence Barrette,^{1,2} Achiraman Shanmugam,⁶ Saravana M. Dhanasekaran,^{1,2} Nallasivam Palanisamy,^{1,2} and Arul M. Chinnaiyan^{1,2,3,4,5,*}

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<http://cancer.sanger.ac.uk/cancergenome/projects/classic>


The screenshot shows the COSMIC website interface. At the top left is the COSMIC logo with the tagline "Catalogue of somatic mutations in cancer". A navigation bar contains links for Home, About, Download, Publications, News, Contact, Help, and FAQ. The main heading is "COSMIC Curated Genes & Fusions". Below this, a paragraph explains that COSMIC contains curated somatic mutation data from scientific literature and the Cancer Genome Project's data. It notes that the first three tabs summarize important cancer genes, and the fourth tab links to the Gene Overview Page. The "Curated Fusions" tab is selected, showing a list of gene fusion entries. A dropdown menu is set to "All" entries. The list includes:

Genes
ACBD6/RRP15
ACSL3_ENST00000357430/ETV1
ACTB/GLI1
AGPAT5/MCPH1
AGTRAP/BRAF
AKAP9_ENST00000356239/BRAF
ARFIP1/FHDC1
ARID1A/MAST2_ENST00000361297
ASPSCR1/TFE3
ATG4C/FBXO38

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July 2008 - Vol 5 No 7
Nature Methods 5 (2008)
The beginning of the end for microarrays?

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
BaseSpace®

The screenshot displays the BaseSpace web interface. At the top, there's a navigation bar with 'Dashboard', 'Prep', 'Runs', 'Projects', 'Apps', 'Public Data', and 'Help'. Below this, a progress bar shows four steps: 1. Biological Samples, 2. Libraries, 3. Pools, and 4. Planned Runs. The main area is titled 'BaseSpace Applications' and features a grid of application icons including 'BWA Enrichment', 'BWA Whole Genome Sequencing v1', and 'Cufflinks Assembly & DE Beta'. A 'TopHat Alignment Beta' application window is open, showing fields for 'App Session Name', 'Save Results To', 'Samples', and 'Reference Genome' (set to 'Homo sapiens/hg19 (RefSeq)'). It also has 'Options' for 'Call Fusions' and 'Trim TruSeq Adapters'. To the right, a scatter plot shows 'log2(generator FPKM)' vs 'log2(control FPKM)'. Below the scatter plot, there are bar charts and a genomic track visualization with labels A through K.

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TopHat Alignment Beta [Details»](#)

A>T Illumina

App Session Name: TopHat Alignment Beta 12/18/2013 9:27:32 ⓘ

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The screenshot shows the BaseSpace web interface for a project named 'NextSeq_RNA-Seq_Demo'. The interface includes a navigation bar with 'Dashboard', 'Prep', 'Runs', 'Projects', 'Apps', 'Public Data', and 'Help'. The 'Projects' tab is active. Below the navigation bar, there are several action buttons: 'Launch app', 'Download Project', 'Import', 'Share project', 'Get link', 'Edit project', and 'Transfer Owner'. The main content area displays a table of samples. The table has columns for 'Sample Id', 'Date Created', 'Read Length', 'Origin', 'Genome', and 'Size'. There are 8 samples listed, with the first three being 'HumanBrain2x125', 'HumanLiver2x125', and 'HumanHeart2x125'. A red oval is overlaid on the table, containing the text 'BaseSpace プロジェクトとしてデータをシェアしています'.

Sample Id	Date Created	Read Length	Origin	Genome	Size
1 HumanBrain2x125	Jun 18 2014				
2 HumanLiver2x125					
3 HumanHeart2x125					
4 HumanHeartDiseased2x125					
5 HumanLiver					
6 HumanHeartDiseased	Jun 07 2014	2 x 151	NextSeq_RNA-...		17.86 GB
7 HumanHeart	Jun 07 2014	2 x 151	NextSeq_RNA-...		25.13 GB
8 HumanBrain	Jun 07 2014	2 x 151	NextSeq_RNA-...		16.41 GB

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Liu, S., Lin, L., Jiang, P., Wang, D. & Xing, Y. A comparison of RNA-Seq and high-density exon array for detecting differential gene expression between closely related species. *Nucleic Acids Res.* 39, 578–588 (2011).

Marioni, J.C., Mason, C.E., Mane, S.M., Stephens, M. & Gilad, Y. RNA-seq: an assessment of technical reproducibility and comparison with gene expression arrays. *Genome Res.* 18, 1509–1517 (2008).

Mortazavi, A., Williams, B.A., McCue, K., Schaeffer, L. & Wold, B. Mapping and quantifying mammalian transcriptomes by RNA-Seq. *Nat. Methods* 5, 621–628 (2008).

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Shi, L. et al. The MicroArray Quality Control (MAQC) project shows inter- and intraplatform reproducibility of gene expression measurements. *Nat. Biotechnol.* 24, 1151–1161 (2006).

SEQC/MAQC-III Consortium. A comprehensive assessment of RNA-seq accuracy, reproducibility and information content by the Sequencing Quality Control Consortium. *Nat. Biotechnol.* doi:10.1038/nbt.2957 (24 August 2014).

't Hoen, P.A. et al. Reproducibility of high-throughput mRNA and small RNA sequencing across laboratories. *Nat. Biotechnol.* 31, 1015–1022 (2013).

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