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BaseSpace®  
ONSITE

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**bcl2fastq2**

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
Categories	
Resequencing	Small RNA
Targeted Sequencing	De Novo Assembly
RNA-Seq	Gene Fusion Detection
ChIP-Seq	Methyl-Seq
Metagenomics	Tumor Normal
Variant Analysis	Differential Expression
Quality	Proteomics
Synthetic Long Reads	Visualization

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
## FASTQ Toolkit

- ▶ Sub-sample reads
- ▶ Trim Adapters
- ▶ Trim Bases
- ▶ Ploy A/T trimming
- ▶ Quality Trimming
- ▶ Read Filtering
- ▶ Reverse Complement

## FastQC

- ▶ Perform QC of raw sequencing data.
- ▶ Determine adapter contamination

## VCAT v2.3

- ▶ Compare Variant Call Sets to standards
- ▶ Intersect variant call sets.

## SRA Import v0.0.3

- ▶ Import up to 25GB of sequencing data from SRA

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### [Data]

Sample_ID	Sample_N	Sample_P	Sample_V	Sample_P	Description
test	test		test		

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## MaskAdapter, MaskAdapterRead2

[Header]	
IEMFileVersion	4
Investigator Name	Mr.X
Experiment Name	Example
Date	#####
Workflow	GenerateFastq
Application	GenerateFastq
Assay	TruSeq LT
Description	Example
Chemistry	Default
[Reads]	
	151
	151
[Settings]	
Adapter	AGATCGG
AdapterRead2	AGATCGG
[Data]	
Sample_ID	Sample_1
test	test

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▶ Modify Reads

BaseSpace Labs:

☐ I acknowledge that I am using it AS IS. Illumina has no obligation to provide support, training, or other resources. (iv) Illumina has no liability for any damages, including consequential damages, arising from the use of the software.

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Adapter trim stringency (0.01-0.99):  
0.9

Select an adapter to trim or specify the adapter sequence(s) below:  
None selected

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Adapter sequence(s) to trim from the 3'-end:  
CTGTCTCTTATACACATCTCCGAG

Trim Ns from the 3'-end before identifying adapters:

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BWA  
(backtrace)

Sample	Sample Name	Total Aligned Reads	Percent Aligned Reads
1	NA12892	354,882	77.4%
2	NA12892-trim	450,007	98.2%

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BWA  
(mem)

Sample	Sample Name	Total Aligned Reads	Percent Aligned Reads
1	NA12892	457,611	99.8%
2	NA12892-trim	457,542	99.8%

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Assembly metrics	Before adapter trimming	After adapter trimming
N50	21	29,791
Maximum contig	553	174,326
Assembly length	18,497,207	4,876,437
Number of contigs	1,387,508	1,115

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

- ▶ TruSeq Small RNA

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- ▶ Human HG19
- ▶ Mus musculus
- ▶ Rattus norvegicus

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- ▶ Alignment
- ▶ Classification of miRNAs, isomiRs, and piRNAs
- ▶ Novel miRNA discovery
- ▶ miRNA Precursor discovery
- ▶ Differential Expression of miRNAs, precursor groups, miRNA families, and piRNAs

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“Adapter trimming” > “Adapter sequence(s) to trim from the 3’ end”:

“TGGAATTCTCGGGTGCCAAGG” (This is the TruSeq smallRNA adapter)

## ▼ Adapter Trimming

Use these parameters if you want to trim adapter sequences. The adapter sequence can be specified separately for the 5'- and 3'-end and is a required input for adapter trimming. In addition to the adapter sequence, you can specify a stringency value that determines how similar a sequence has to be in order to be identified as an adapter sequence.

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Select an adapter to trim or specify the adapter sequence(s) below:

None selected

Adapter sequence(s) to trim from the 5'-end:

Adapter sequence(s) to trim from the 3'-end:

TGGAATTCTCGGGTGCCAAGG

Trim from the 5'-end before identifying adapters:

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4

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“Read Filter” > “Minimum Read length: 15” (変更可能)

▶ Quality Trimming

▼ Read Filtering

Use these parameters if you want to filter reads. Paired-end reads are only filtered (removed from the sample) if both reads are filtered out. Otherwise, the filtered mate is replaced by a sequence of Ns (number of Ns will be the minimum read length) to keep the order of pairs in the FASTQ files (necessary for many secondary analysis tools).

Minimum read length:

15

Maximum read length:

Minimum mean quality score:

Maximum mean quality score:

*Note, that leaving as default will result in conversion of sequences <32bp to “N” strings*

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► Modify Reads

BaseSpace Labs:

- ☒ I acknowledge and agree that (i) this is a BaseSpace Labs App, (ii) I am using it AS-IS without any warranty of any kind, (iii) Illumina has no obligation to provide any technical support for this App, and (iv) Illumina has no liability for my use of this App, including without limitation, any loss of data, incorrect results, or any costs, liabilities, or damages that result from use of this App.

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
Read Length Distributions

Note: Reads that only contain Ns (no-calls) in the output sample are represented with length zero in the charts above.

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# BaseSpace Small RNAアプリ

**Small RNA v1.0**  
ILLUMINA  
[Launch](#)

Free

Precursor: **hsa-let-7a-1**

[Save as PNG](#)

**hsa-let-7a-1 in miRBase**  
5' Counts: 203771  
3' Counts: 129  
5'/3' ratio: 1579.62

```
tggaTGAGGTAGTAGTTGTATAGTTtaggggtcacaccaccactggagataaCTATACAATCTACTGCTTTCcta
((((((((((((((((((((((((((((((((((((((((((((((((((((((((((((((((((((((((
...TCTGAGGTAGTAGTTGTATAGTT..... 6
....TGAGGTAGTAGTTGTATAGTT..... 196698
....TGAGGTAGTAGTTGTATAGT..... 7061
....TGAGGTAGTAGTTGTGTTT..... 6
.....CTATACAATCTACTGCTTTC... 33
```

[Overview](#)

### Description

The BaseSpace Small RNA v1.0 app analyzes small RNA samples. The Small RNA app supports TruSeq Small RNA Sample Preparation Kits and assumes that reads have been properly adapter-trimmed. The app aligns reads against four reference databases (abundant, mature miRNA, other RNA, and genomic) and outputs hits to mature miRNAs, isomiRs, and piRNAs. Optionally, the app performs novel precursor discovery and pairwise differential expression analysis. Pairwise differential expression analysis identifies differentially expressed miRNAs, precursor groups, miRNA families, and piRNAs for each pair of sample groups.


This application is for research use only.

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The full list is available at Hits.txt.

## IsomiRs (Known Precursor)

A read is counted as isomiR if it is a subsequence on the same strand of a pre-miRNA sequence. The precursor sequences were obtained from miRBase. An isomiR hit is labeled as "start position}\_{sequence read}\_{precursor ID}", where the start position is 0-based. Mismatches are allowed. IsomiR hits are filtered to remove artifacts of PCR amplification.



Number of IsomiRs (Known Precursor) with Reads	Total Number of Reads
464	22,967

## ces

Reads are listed for each type. If there are less than 10 sequences with a certain length, they are not listed.


Records are counted. A read must align to the start of a reference sequence and have the same length of the reference sequence. No mismatches are allowed.

Reads	Total Number of Reads
281	22,967

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


## GenerateFastq in MSR/ BaseSpace /bcl2fastq2)

[Header]			
IEMFileVersion	4		
Investigator Name	Mr.X		
Experiment Name	Example		
Date	4/05/2015		
Workflow	Resequencing		
Application	Resequencing		
Assay	TruSeq LT		
Description	Example		
Chemistry	Default		
[Reads]			
	151		
	151		
[Settings]			
Adapter	AGATCGGAAGAGCACACGTCTGAA		
Adapter_Read2	AGATCGGAAGAGCGGTCTGTAGGGA		
QualityScoreTrim	20		
[Data]			
Sample_ID	Sample_Name	Sample_Plate	
test	test		test

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Save Processed Sample(s) to:

Select Project(s):

Add this string to the output sample name(s):

▶ Sub-sampling

▶ Adapter Trimming

▶ Base Trimming

▶ Poly-A/T Trimming

▶ Quality Trimming

▶ Read Filtering

▶ Modify Reads

Trim the 3'-end of reads with quality level:

BaseSpace Labs:

☐ I acknowledge and agree that (i) this is a BaseSpace Labs App, (ii) I am using it AS-IS without any warranty of any kind, (iii) Illumina has no obligation to provide any technical support for this App, and (iv) Illumina has no liability for my use of this App, including without limitation, any loss of data, incorrect results, or any costs, liabilities, or damages that result from use of this App.

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
Save Processed Sample(s) to:

Select Project(s)

Add this string to the output sample name(s):

▶ Sub-sampling

▶ Adapter Trimming

▶ Base Trimming

▶ Poly-A/T Trimming

▶ Quality Trimming

▶ Read Filtering

▶ Modify Reads

BaseSpace Labs:

☐ I acknowledge am using it AS-is. Illumina has no obligation to provide support without limitation of liabilities, or damages.

▼ Sub-sampling

Sub-sampling is required when only a subset of the sample can be processed by an application (e.g. de novo assembly with memory constraints) or it is not necessary to process a full sample (e.g. for validating an approach at varying levels of genomic coverage). Samples can be sub-sampled to a user defined number of reads or number of bases, or to a specified fraction of the input sample (e.g. 50% of input sample) - again, based on reads or bases.

Maximum number of FASTQ entries to keep:

Maximum percentage of FASTQ entries to keep:

Maximum number of bases to keep:

Maximum percentage of bases to keep:

Choose reads at random when sub-sampling instead of taking the first n reads:

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The screenshot shows the BaseSpace web interface for importing FASTQ files. The browser address bar displays <https://basespace.illumina.com/projects/16424410/import>. The navigation bar includes links for Dashboard, Prep, Runs, Projects (active), Apps, Public Data, and Help. The user is logged in as Gabriel Kolle.

The main content area is titled "Fastq upload test : Import". Under "Import type:", the "Sample" option is selected. The sample name is "Sample1" and the reference is "Homo Sapiens - UCSC (hg19)".

Two files are being uploaded:

- `Sample1_S1_L001_R1_001.fastq.gz` (231.70 KB) - 100% complete (green bar with checkmark)
- `Sample1_S1_L001_R2_001.fastq.gz` (253.25 KB) - 50% complete (blue bar)

A purple box highlights the "Complete Import" button. Below the button, a summary of the upload statistics is shown:

Read 1:	151
Read 2:	151
Paired End:	Yes
Reads Passing Filter	100.00%
Total Reads	5,000
Total File Size	484.94 KB

The "Illumina FASTQ file standard" link is also visible.

On the right, there is a "Drag and drop or select files" area with a file type filter set to ".fastq.gz". A "contact us" button is located on the far right.

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FASTQ Toolkit v2.0.0  
BaseSpace Labs

Analysis Name: FASTQ Toolkit\_09/04/2015 12:36:21

Input Sample(s) to Process: **Select Sample(s):**

Save Processed Sample(s) to: **Select Project(s):**

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Default Settings

Minimum read length: 32

The image shows the FASTQ Toolkit v2.0.0 interface. A dashed orange box highlights the 'Select Sample(s):' button under 'Input Sample(s) to Process:' and the 'Select Project(s):' button under 'Save Processed Sample(s) to:'. The 'Analysis Name' field contains 'FASTQ Toolkit\_09/04/2015 12:36:21'. The 'Minimum read length' is set to 32.

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[https://my.illumina.com/MyIllumina/Bulletin/qFYNf9hn\\_kW5SyEZwGOUrA/adapter-sequences-for-use-with-casava-or-bcl2fastq](https://my.illumina.com/MyIllumina/Bulletin/qFYNf9hn_kW5SyEZwGOUrA/adapter-sequences-for-use-with-casava-or-bcl2fastq)

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[http://support.illumina.com/content/dam/illumina-support/documents/documentation/software\\_documentation/miseqreporter/miseq-reporter-generatefastq-workflow-guide-15042322-b.pdf](http://support.illumina.com/content/dam/illumina-support/documents/documentation/software_documentation/miseqreporter/miseq-reporter-generatefastq-workflow-guide-15042322-b.pdf)

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