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© 2014 Illumina, Inc. All rights reserved. Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPro, CytoChip, DesignStudio, Epicentre, GAllx, Genetic Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, ForenSeq, MiSeq, MiSeqDx, MiSeqFGx, NeoPrep, Nextera, NextBio, NextSeq, Powered by Illumina, SeqMonitor, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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© 2014 Illumina, Inc. All rights reserved.  
Illumina, iScan, BaseSpace, BaseAnalyzer, Bluefish, BlueFuse, BlueGenome, eShot, CDP+, CytoChip, DesignStudio, EpiCentre, GATK, Genetix Energy, Genome Analyzer, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSeq, NextSeq, NextSeq 500, MiSeq, MiSeq Dx, NextSeq, Nextera, NextSeq, powered by Illumina, SingleShot, SureSelect, TruSight, TruSight, TruSight, Understand Your Genome, UVO, Veracode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the U.S. and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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


## Library Prep Workflow

The following diagram illustrates the workflow using a Nextera Rapid Capture Enrichment kit. Safe stopping points are marked between steps.

Figure 1 Nextera Rapid Capture Enrichment Workflow

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## Supporting Information

### Consumables and Equipment

Check to make sure that you have all of the necessary user-supplied consumables and equipment before proceeding to the library preparation and enrichment procedures.

**NOTE**  
The Nextera Rapid Capture Enrichment protocol has been optimized and validated using the items listed. Comparable performance is not guaranteed when using alternate consumable and equipment.

Table 6 User-Supplied Consumables

| Consumable                                                                                                                  | Supplier                             |
|-----------------------------------------------------------------------------------------------------------------------------|--------------------------------------|
| 1.7 ml microcentrifuge tubes                                                                                                | General lab supplier                 |
| 20 µl barrier pipette tips                                                                                                  | General lab supplier                 |
| 20 µl multichannel pipettes                                                                                                 | General lab supplier                 |
| 20 µl single channel pipettes                                                                                               | General lab supplier                 |
| 200 µl barrier pipette tips                                                                                                 | General lab supplier                 |
| 200 µl multichannel pipettes                                                                                                | General lab supplier                 |
| 200 µl single channel pipettes                                                                                              | General lab supplier                 |
| 1000 µl barrier pipette tips                                                                                                | General lab supplier                 |
| 1000 µl multichannel pipettes                                                                                               | General lab supplier                 |
| 1000 µl single channel pipettes                                                                                             | General lab supplier                 |
| Adhesive seal roller                                                                                                        | General lab supplier                 |
| 96-well flat clear bottom black microplates<br>Note: Used when quantifying samples with a SpectraMax M5 spectrofluorometer. | Corning, part # 3904                 |
| 96-well storage plates, round well,<br>0.8 ml ("MIDI" plate)                                                                | Fisher Scientific,<br>part # AB-0859 |

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Copyright © 2015 Thermo Fisher Scientific Inc. Used under permission.

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