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© 2015 Illumina, Inc. All rights reserved.  
Illumina, 24sure, BaseSpace, BeadArray, BlueFish, BlueFuse, BlueGnome, cBot, CSPPro, CytoChip, DesignStudio, Epicentre, ForenSeq, Genetic Energy, GenomeStudio, GoldenGate, HiScan, HiSeq, HiSeq X, Infinium, iScan, iSelect, MiSeq, MiSeqDX, MiSeq FGx, NeoPrep, NextBio, Nextera, NextSeq, Powered by Illumina, SureMDA, TruGenome, TruSeq, TruSight, Understand Your Genome, UYG, VeraCode, verifi, VeriSeq, the pumpkin orange color, and the streaming bases design are trademarks of Illumina, Inc. and/or its affiliate(s) in the US and/or other countries. All other names, logos, and other trademarks are the property of their respective owners.

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
<http://jp.illumina.com/products/truseq-custom-amplicon-low-input.html>

# TruSeq Custom Amplicon Design Workflow

Assay Type      Configure Design      Add Targets      Submit Design      Review Design

1 ..... 2 ..... 3 ..... 4 ..... 5

### Configure Design

Assay Version ?

Comparison Chart

Variant source ?

Amplicon size ?

**TruSeq Custom Amplicon Low Input:**  
Used for low input and/or FFPE applications. Not to be used with applications requiring 425bp.

Design dual pools ?

**Population:**


AFR  AMR  ASN  EUR

1000 Genomes: common variants detected in at least 1% of the population(s) studied.

Next

# DesignStudio Definition of Amplicon

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

The screenshot displays the 'Configure Design' step (step 2) of the workflow. At the top, a progress bar shows five steps: 1. Assay Type, 2. Configure Design (highlighted), 3. Add Targets, 4. Submit Design, and 5. Review Design. The main content area is titled 'Configure Design' and includes several configuration options:

- Assay Version:** Two options are shown: '1.5' and 'Low Input'. The 'Low Input' option is highlighted in orange. To its right, text reads: 'TruSeq Custom Amplicon Low Input: Used for low input and/or FFPE applications. Not to be used with applications requiring 425bp.'
- Comparison Chart:** A checkbox is checked.
- Variant source:** Two options are shown: 'dbSNP 138' and '1000 Genomes'. The '1000 Genomes' option is highlighted in orange. To its right, text reads: 'Population:  AFR  AMR  ASN  EUR.   
1000 Genomes: common variants detected in at least 1% of the population(s) studied.'
- Amplicon size:** Three options are shown: '150 bp', '175 bp', and '250 bp'. The '175 bp' option is highlighted in orange.

A 'Next' button is located at the bottom right of the interface.

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