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1. Caporaso JG, Lauber CL, Walters WA, Berg-Lyons D, Lozupone CA, et al. (2011) Global patterns of 16S rRNA diversity at a depth of millions of sequences per sample. Proc Natl Acad Sci USA 108:4516–4522.
2. Caporaso JG, Lauber CL, Walters WA, Berg-Lyons D, Huntley J, et al. (2011) Manuscript in preparation.
3. Caporaso JG, Kuczynski J, Stombaugh J, Bittinger K, Bushman FD, et al. (2010) QIIME allows analysis of high-throughput community sequencing data. Nat Methods 7(5):335–336.
4. Desantis TZ, Hugenholtz P, Larsen N, Rojas M, Brodie EL, et al. (2006) Greengenes, a chimera-checked 16S rRNA gene database and workbench compatible with ARB. Appl Environ Microbiol 72(7): 5069–5072.
5. Cole JR, Wang Q, Cardenas E, Fish J, Chai B, et al. (2009) The Ribosomal Database Project: improved alignments and new tools for rRNA analysis. Nucl Acids Res 37:D141–D145.
6. Liu Z, Lozupone C, Hamady M, Bushman FD, Knight R. (2007) Short pyrosequencing reads suffice for accurate microbial community analysis. Nucl Acids Res 35:18.
7. Liu Z, DeSantis TZ, Andersen GL, Knight R. (2008) Accurate taxonomy assignments from 16S rRNA sequences produced by highly parallel pyrosequencers. Nucl Acids Res 36:18.
8. http://qiime.org/tutorials/processing_illumina_data.html

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© 2014 Illumina, Inc. All rights reserved.