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1. Sie D, Snijders PJ, Meijer GA, et al. Performance of amplicon-based next generation DNA sequencing for diagnostic gene mutation profiling in oncopathology. Cell Oncol. 2014;37(5):353-361.
2. Luthra R, Patel KP, Reddy NG, et al. Next generation sequencing based multigene mutational screen for acute myeloid leukemia using Miseq: applicability for diagnostics and disease monitoring. Haematologica. 2014;99(3):465-473.
3. Thomas M, Sukhai M, Zhang T, et al. Clinical testing and implementation of the TruSight Myeloid Next Generation Sequencing (NGS) panel for identification of clinically relevant variants in hematological malignancies. Cancer Res. 2015;75(15):4260.
4. Thiede C, Bullinger L, Hernandez-Rivas JM, et al. Results of the "Evaluation of NGS in AML-Diagnostics (ELAN)" Study-an Inter-Laboratory Comparison Performed in 10 European Laboratories. Blood. 2014;124(21): 2374.
5. Shendrue J, and Hanlee J. Next-generation DNA sequencing. Nat Biotech. 2008;26:1135-1145.
6. Schuster SC. Next-generation sequencing transforms today's biology. Nat Methods. 2008;5(1):16-18.
7. Tuononen K, Maki-Nevala S, Sarhadi VK, et al. Comparison of targeted next-generation sequencing (NGS) and real-time PCR in the detection of EGFR, KRAS, and BRAF mutations on formalin-fixed, paraffin-embedded tumor material of non-small cell lung carcinoma-superiority of NGS. Genes, Chromosomes and Cancer. 2013;52(5):503-511.
8. Ross MG, Russ C, Costello M, et al. Characterizing and measuring bias in sequence data. Genome Biol. 2013;14:R51.
9. Liu L, Li Y, Li S, et al. Comparison of next-generation sequencing systems. J Biomed Biotechnol. 2012;2012:251364.
10. Cross D, Burmester JK. The promise of molecular profiling for cancer identification and treatment. Clin Med Res. 2004;2(3):147-150.
11. Lindeman NI, Cagle PT, Beasley MB, et al. Molecular testing guideline for selection of lung cancer patients for EGFR and ALK tyrosine kinase inhibitors: guideline from the College of American Pathologists, International Association for the Study of Lung Cancer, and Association for Molecular Pathology. J Thorac Oncol. 2013;8:823-859.
12. NCCN Clinical Practice Guidelines in Oncology
13. WHO/IARC Cancer Pathology and Genetics. www.iarc.fr/en/publications/pdfs-online/pat-gen/
14. Staudt LM. Molecular diagnosis of the hematological cancers. N Eng J Med. 2003;348:1777-1785.
15. Gardner JA, de Abreu F, Peterson J, et al. Identification of somatic mutations in acute myeloid leukemia patients using the TruSight Myeloid Sequencing Panel. A J Clin Path. 2015;144(2):A238.
16. Mertens F, Johansson B, Fioretos T, et al. The emerging complexity of gene fusions in cancer. Nat Rev Cancer. 2015;15:371-381.
17. Parker BC, Zhang W. Fusion genes in solid tumors: an emerging target for cancer diagnosis and treatment. Chin J of Cancer. 2013;32(11):594-603.

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18. Mardis ER. The translation of cancer genomics: time for a revolution in clinical cancer care. Genome Med. 2014;6(3):22.
19. Bentley DR, Balasubramanian S, Swerdlow HP, et al. Accurate whole human genome sequencing using reversible terminator chemistry. Nat. 2008;456:53-59.
20. Illumina. TruSight Tumor 15 Data Sheet. www.illumina.com/content/dam/illumina-marketing/documents/products/datasheets/trusight-tumor-15-data-sheet-1170-2015-003.pdf

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© 2016 Illumina, Inc. All rights reserved.