Three scientists, two men and one woman, are standing in a row, wearing white lab coats. The man on the left is smiling, the woman in the middle is looking directly at the camera, and the woman on the right is holding a clipboard and smiling. The background is white with a wavy, colorful band of DNA sequence motifs (A, T, C, G) in orange, green, and blue. The text 'Illumina Sequencing Systems – Power, Simplicity, Versatility.' is written in a large, black, sans-serif font on the right side of the image.

Illumina Sequencing Systems – Power, Simplicity, Versatility.

Jeremy Preston, Ph.D.
Director, Product Marketing
Systems & Consumables
San Diego

Illumina Sequencing

It keeps getting easier – do more with less!

MiSeq

liSeq 2

HiSeq 2500

Circa Sep 2009



**Single Instrument Workflow for
MiSeq and HiSeq**


iNTEGRATION

Streamlined end-to-end solution



The Pace of Change in Sequencing

Making sequencing cheaper and faster




Optimizing Cluster Size and Spacing

Random cluster generation to patterned flow cell

- ▶ New patterned substrate technology
 - Transforms random cluster generation into densely packed ordered arrays
- ▶ Universal technology that can benefit multiple Illumina platforms
- ▶ Increased performance
 - Potential increase in reads
 - Faster cluster detection timing

Random Cluster Generation




HiSeq 2500 Sequencing System

Speed and versatility



HIGH OUTPUT

**50 – 600Gb
2 – 11 days
2 x 100bp max**




**Larger projects,
fewer runs**



RAPID RUN

**10 – 180Gb
7 – 40 hours
2 x 150bp max**



**Smaller projects,
quick results**



HiSeq 2500 WGS – from DNA to Disease in ~50hrs

Rapid sequencing and annotation of neonatal intensive care patients

Neonate at 5 mo -
developmental
regression &
seizures

Samples prepared,
sequenced and
analyzed
(HiSeq 2500)

Identified novel
variant in gene
linked to copper
metabolism

Confirmed
diagnosis of
Menkes disease



RESEARCH /

DIAGNOSTICS

Rapid W
Diagnosis


Carol Jean Sau
Darrell Lee Din
Melanie LeAnn
Zoya Kingsbury
Elliott Margulie
Joshua Erin Pet

The New York Times



Sample to Variant Call in Less than 24 Hours!

Blinded re-run of Menkes disease sample using Illumina pipeline




HiSeq 2500 – Extending Performance in 2013

Commercial path to 2x250bp reads in rapid run mode

- ▶ 2x250bp in rapid mode
- ▶ Up to 300Gb in ~60 hours
 - ~150Gb per flow cell
- ▶ Increased application breadth
 - Transcript isoforms/gene fusions
 - Metagenomics & complex de novo assembly
 - Resolve difficult regions of the genome
- ▶ Improved accuracy for all read lengths

New: reagents & software

Sample	Yield per Flow Cell	Depth (x)	% >Q30	Reads per Flow Cell
HCC2218 Normal	129Gb	38	78.6	258M
HCC2218 Tumor	156Gb	45	76.3	311M
Current R&D performance	147Gb	42	87.5	294M





The Longer Read Advantage

Synthetic gut Metagenomics

- ▶ *In vitro* model of large intestine inoculated with human gut microbiota
- ▶ 2x250 bp on HiSeq 2500
- ▶ Examine bacterial community & assess metabolic pathways (MG-RAST)
- ▶ Reads mapping to known bacterial genes **highly dependent on read length**
- ▶ **Better resolution** of closely related sequences with longer reads (species/strain)

"We were able to see a very large number of those functional genes that we were not observing before."





- ▶ Reads mapping to genes
- ▶ 100nt 3,963,158
- ▶ 150nt 7,170,312
- ▶ 200nt 23,763,325
- ▶ 251nt 25,160,962

Sample and analysis courtesy of Bas Dutilh,
Radboud University Medical Centre Nijmegen, Holland

HiSeq 2500 – R&D Demonstrated Scalability

Overlapping 2 x 400 bp Reads (In Research)



Run metrics

Read length	2 x 400
Output (single flow cell)	>180Gb
Quality	70% Q30 (All bases)


Read Length # of reads

>= 500	64.8M
>= 600	42.5M
>= 700	5.3M
>= 750	739K
>= 775	128K

- ▶ Human NA12882 (CEPH Child) library – 650 bp median insert size

MiSeq – Continuous Performance Improvements

Path towards 15Gb per run; enabling broader range of applications



"This information is intended to outline general product direction and it should not be relied on in making a purchasing decision. This material is for information purposes only and may not be incorporated into any contract. This information is not a commitment, promise, or legal obligation to deliver this functionality. The development, release, and timing of any features or functionality described for our products remains at our sole discretion."


MiSeq – R&D Demonstrated Scalability

Early Development data

Output	15.2 Gb
Clusters	25.1M
Read length	2 x 300
Quality	90% \geq Q30
MMR	0.17% (R1), 0.37% (R2)
Run time	72 hours

**New: 150 & 600
cycle PE kits at
launch**

“Research” *longer read runs in excess of 20Gb*





Output	22 Gb
Clusters	26.5M
Read length	450 (R1), 375 (R2)
Quality	71% \geq Q30
MMR	1.7% (R1), 2.4% (R2)

MiSeq – Greater than 85% of Desktop Data

Analysis of data submissions to the NCBI Sequence Read Archive (SRA)

Number of Benchtop Sequencer Experiments




**n = 80+ publications*

NIH's primary archive of NGS data and part of the international partnership of archives at the NCBI, the European Bioinformatics Institute and the DNA Database of Japan.

MiSeq – Transforming Clinical Microbiology

Rapid bacterial genome sequencing



MiSeq – Adopted by Leading Public Health Organizations



Identification of food-borne pathogen outbreaks



Cost effective, public health genomic epidemiology

Novel Avian-Origin Influenza A (H7N9) Virus Genome Sequencing at China CDC Using MiSeq

The NEW ENGLAND JOURNAL of MEDICINE

ORIGINAL ARTICLE

Preliminary Report: Epidemiology of the Avian Influenza A (H7N9) Outbreak in China

Qun Li, M.D., Lei Zhou, M.D., Minghao Zhou, Ph.D., Zhiping Chen, M.D.,
Furong Li, M.D., Huanyu Wu, M.D., Nijuan Xiang, M.D., Enfu Chen, M.P.H.,
Fenyang Tang, M.D., Dayan Wang, M.D., Ling Meng, M.D., Zhiheng Hong, M.D.,
Wenxiao Tu, M.D., Yang Cao, M.D., Leilei Li, Ph.D., Fan Ding, M.D., Bo Liu, M.D.,
Mei Wang, M.D., Rongheng Xie, M.D., Rongbao Gao, M.D., Xiaodan Li, M.D.,
Tian Bai, M.D., Shumei Zou, M.D., Jun He, M.D., Jiayu Hu, M.D., Yangting Xu, M.D.,
Chengliang Chai, M.D., Shiwen Wang, M.D., Yongjun Gao, M.D., Lianmei Jin, M.D.,
Yanping Zhang, M.D., Huiming Luo, M.D., Hongjie Yu, M.D., M.P.H.,
Lidong Gao, M.D., Xinghuo Pang, M.D., Guohua Liu, M.D., Yuelong Shu, Ph.D.,
Weizhong Yang, M.D., Timothy M. Uyeki, M.D., M.P.H., M.P.P., Yu Wang, M.D.,
Fan Wu, M.D., and Zijian Feng, M.D., M.P.H.

MiSeq – Phased HLA sequencing

Hosomichi *et al. BMC Genomics* 2013, **14**:355
<http://www.biomedcentral.com/1471-2164/14/355>




METHODOLOGY ARTICLE

Open Access

Phase-defined complete sequencing of the HLA genes by next-generation sequencing

Kazuyoshi Hosomichi¹, Timothy A Jinam¹, Shigeki Mitsunaga², Hirofumi Nakaoka¹ and Ituro Inoue^{1*}

Division of Human Genetics, National Institute of Genetics, Shizuoka, Japan



Illumina's Suite of Sample Prep Solutions



DNA Sample Prep

- Nextera & Nextera XT
- Nextera Mate Pair
- TruSeq DNA PCR-free
- TruSeq Nano DNA

Targeted Resequencing



- TruSeq Custom Amplicon
- Nextera Rapid Capture Custom & Exome
- TruSight Panels

RNA & Regulation Sample Prep

- TruSeq Small RNA
- TruSeq RNA
- TruSeq Stranded RNA (FFPE)
- TruSeq ChIP
- TruSeq Targeted RNA Expression

New: TruSeq PCR Free and Nano

Improved coverage of challenging genomic regions




- ▶ **TruSeq PCR Free – 1ug input**
 - Fast, no PCR duplicates
- ▶ **TruSeq Nano – 100 ng input**
 - Low input, precious samples

- ✓ High library diversity
- ✓ Excellent coverage uniformity
- ✓ All-inclusive kit w/gel-free size-selection
- ✓ 24 & 96 sample options

Nextera Rapid Capture Exome

Fastest exome sequencing - sample to data in less than 2.5 days!

- ▶ **Fast assay**
 - Only 1.5 days to prep 96 exomes
 - 12-plex pre-enrichment pooling
- ▶ **Lowest DNA input**
 - Only 50ng
- ▶ **Efficient**
 - 37Mb exonic content requires only 4Gb total sequence*
- ▶ *Also available:*
 - *Nextera Rapid Capture Expanded Exome*
 - *Nextera Rapid Capture Custom Enrichment*



TruSeq Targeted RNA

Rapid & economical RNA profiling

- ▶ Accurate targeting for Human, Mouse & Rat transcriptomes
 - Alternative isoforms
 - Individual exons and splice junctions
 - cSNP detection for allele specific expression
 - Non-coding RNA transcripts
- ▶ Add custom content to Fixed Panels

- ✓ **Rapid – sample to answer in 1.5 days; <4hrs HOT**
- ✓ **Low RNA input – 50 ng or less**
- ✓ **Low price per sample; 48-384 spls/run**

- ✓ **Validation of over 10,000 assay designs**
- ✓ **Custom and pre-validated fixed panels available**



Pre-Validated Fixed Panels

Immune Response

Lung Cancer

Breast Cancer

Stem Cell

P53 Pathway

Cytochrome P450

NFκB Pathway

Cardiotox

Apoptosis

Neuro Panel

Prostate Cancer

Wnt Pathway

Cell Cycle

Hedgehog

TruSeq Targeted RNA

Rapid RNA profiling

Single MiSeq run



=

Equivalent to 15,000 qPCR
reactions
(40 x 384 well plates)




Molecuro Technology Enables Synthetic Long Reads

Up to 10Kb from Illumina short reads

- ✓ Synthetic long reads 8 – 10kb
- ✓ Enables fully phased genomes
- ✓ Accurate de novo assembly of large, complex genomes


Available:

- ✓ Illumina services 2H13
- ✓ Kit format early 2014



Platform Synergy

Leveraging synergy for fast development



Thank You!

Jeremy Preston
jpreston@illumina.com