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
Workflow	Library Prep	Sequencer Fit	Analysis	Data Sharing
Shotgun Metagenomics (RNA)	 <p>Requires RNA extraction followed by TruSeq stranded total RNA</p>	 <p>NextSeq, NovaSeq</p>	 <ul style="list-style-type: none"> • DRAGEN Metagenomics 	
Enrichment	 <p>Requires cDNA followed by Nextera DNA Flex for Enrichment + Respiratory Virus Oligo Panel</p>	 <p>iSeq, MiniSeq, MiSeq</p>	 <ul style="list-style-type: none"> • DRAGEN Enrichment • DRAGEN RNA Pathogen Detection 	 <p>GISAID Submission</p>
Amplicon	 <p>AmpliSeq for Illumina SARS-CoV-2 Research Panel</p>	 <p>iSeq, MiniSeq, MiSeq</p>	 <ul style="list-style-type: none"> • DNA Amplicon 	

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

Table 2: Basic metrics using LRM Resequencing Module

Sample	Aligned reads	Mean coverage	Small variants
■ CoVOC43CtrlPat	5.5%	1786.8	47
■ CoV229EctrlPat	6.1%	2266.5	24
■ CoVOC43CtrlCult	78.8%	25,971.5	47
■ CoV229EctrlCult	71.5%	33,833.3	24

Alignment-Based Coverage Maps, Viral Genome Assembly

5% CoV
+95%UHR

100%
CoV

Comprehensive workflow for detecting coronavirus using Illumina benchtop systems

<https://jp.illumina.com/content/dam/illumina-marketing/documents/products/appnotes/ngs-coronavirus-app-note-1270-2020-001.pdf>

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カタログ番号	製品名	①ライブラリ調製試薬 (サンプル数)	③濃縮試薬 (Enrichment反応量)
20025523	旧) Nextera DNA Flex Pre-Enrichment Library Prep and Enrichment Reagents – 16 samples	16 sample	16 Enrichment
	新) Illumina DNA Prep with Enrichment (S) Tagmentation 16 Samples		
20025524	旧) Nextera DNA Flex Pre-Enrichment Library Prep and Enrichment Reagents – 96 samples	96 sample	8 Enrichment
	新) Illumina DNA Prep with Enrichment (S) Tagmentation 96 Samples		
20025519	旧) Nextera DNA Flex Pre-Enrichment Library Prep Reagents (16 samples)	16 sample	(なし)
	新) Illumina DNA Prep (S) Tagmentation 16 Samples		
20025520	旧) Nextera DNA Flex Pre-Enrichment Library Prep Reagents (96 samples)	96 sample	(なし)
	新) Illumina DNA Prep (S) Tagmentation 96 Samples		

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https://support.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/samplepreps_nextera/nextera-flex-enrichment/nextera-flex-for-enrichment-rna-probes-demonstrated-protocol-1000000070581-01.pdf

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




Enriched and Indexed Library Ready for Sequencing

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Illumina Workflows

Workflow	Library Prep	Sequencer Fit	Analysis	Data Sharing
<p>Shotgun Metagenomics (RNA)</p>	 <p>Requires RNA extraction followed by TruSeq stranded total RNA</p>	 <p>NextSeq, NovaSeq</p>	<p>NEW</p>  <ul style="list-style-type: none"> • DRAGEN Metagenomics 	
<p>Enrichment</p>	 <p>Requires cDNA followed by Nextera DNA Flex for Enrichment + Respiratory Virus Oligo Panel</p>	 <p>iSeq, MiniSeq, MiSeq</p>	<p>NEW</p>  <ul style="list-style-type: none"> • DRAGEN Enrichment • DRAGEN RNA Pathogen Detection 	<p>NEW</p>  <p>GISAIID Submission</p> 
<p>Amplicon</p>	 <p>AmpliSeq for Illumina SARS-CoV-2 Research Panel</p>	 <p>iSeq, MiniSeq, MiSeq</p>	 <ul style="list-style-type: none"> • DNA Amplicon 	

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
The image shows a screenshot of the DRAGEN configuration form. Several sections are highlighted with red dashed boxes. The "Reference Database" section includes a dropdown menu set to "MiniKraken2 (March 2020)". The "Custom Reference Database" section has a "Select Dataset File(s)" button. The "Host Reference" section has a dropdown menu set to "Human (hg38 Alt-Aware, with HLAs)". The "Dehosting" section has a checked checkbox for "Enable Dehosting". The "Organism Detection Report" section has three radio buttons, with "From text input" selected. The "Organism Identifiers" section has a text input field containing "Severe acute respiratory syndrome coronavirus 2", "Influenza A", and "Influenza B".

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Organism Detection report, Organism Identifier

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結果サマリー CoV2_Enriched

AGGREGATE SUMMARY

Aggregate Organism Detection Report [Export \(CSV\)](#)

Color Palette: 

[Copy](#) [CSV](#) [Excel](#)

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Search:

Sample ID	Human coronavirus OC43	Influenza A	Influenza B
CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	High Detected (1740576)	ND (0)	ND (0)
CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep2	High Detected (1308243)	ND (0)	ND (0)
CoVOC43_Thermo_cDNA_Synthesis_Enriched_Rep1	High Detected (5235615)	ND (0)	ND (0)
CoVOC43_Thermo_cDNA_Synthesis_Enriched_Rep2	High Detected (5576985)	ND (0)	ND (0)

Showing 1 to 4 of 4 entries

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Sample identifierの指定

<https://www.ncbi.nlm.nih.gov/Taxonomy/Browser/wwwtax.cgi>

The "Token set" option returns longer names that include the search terms, e.g., hybrid taxa. See what happens if you query "Bos taurus" using the "Complete match" option versus the "Set of tokens" option. The "Phonetic search" option can be used when you are not sure about the exact spelling of an organism name. It tries to find the phonetically closest match.

This is the top level of the taxonomy database maintained by NCBI/GenBank. You can explore any of the taxa listed below by clicking on the link.

- [Archaea](#)
- [Bacteria](#)
- [Eukaryota](#)
- **[Viruses](#)**
- [Other](#)
- [Unclassified](#)

These are direct links to some of the organisms commonly used in molecular research projects:

Arabidopsis thaliana	Escherichia coli
Bos taurus	Hepatitis C virus
Caenorhabditis elegans	Homo sapiens
Chlamydomonas reinhardtii	Mus musculus
Danio rerio (zebrafish)	Mycoplasma pneumoniae
Dictyostelium discoideum	Oryza sativa
Drosophila melanogaster	Plasmodium falciparum

Disclaimer: The NCBI taxonomy database is not an authoritative source for nomenclature or classification - please consult the literature.

NCBI Taxonomy Browserで参照できる名称

Sample identifierの指定

Kraken 2 Taxonomic Sequence Classification System

Home Manual Downloads CCB » Software » Kraken2

2020/03 UPDATE

All links are currently being updated.
ftp://ftp.ccb.jhu.edu/pub/data/kraken2

名前	サイズ	更新日
16S_Greengenes13.5_20200326.tgz	72.1 MB	2020/03/26 20:19:00
16S_RDP11.5_20200326.tgz	168 MB	2020/03/26 20:21:00
16S_Silva132_20200326.tgz	117 MB	2020/03/26 18:28:00
16S_Silva138_20200326.tgz	112 MB	2020/03/27 13:09:00
README.md	1.9 kB	2020/03/27 15:29:00
README.txt		
minikraken_8GB_20200326.tgz		
minikraken_8GB_20200312.tgz		
old/		

```
>kraken:taxid|273360|NC_043646.1 Potosi virus strain IL94-1899 segment L, complete sequence
>kraken:taxid|273360|NC_043645.1 Potosi virus strain IL94-1899 segment M, complete sequence
>kraken:taxid|273360|NC_043647.1 Potosi virus strain IL94-1899 segment S, complete sequence
>kraken:taxid|204928|NC_021926.1 Calibrachoa mottle virus isolate California, complete genome
>kraken:taxid|21698|NC_006213.1 Human coronavirus OC43 strain ATCC VR-759, complete genome
>kraken:taxid|31631|NC_006213.1 Human coronavirus OC43 strain ATCC VR-759, complete genome
>kraken:taxid|75734|NC_006213.1 Human coronavirus OC43 strain ATCC VR-759, complete genome
>kraken:taxid|2169971|NC_006213.1 Human coronavirus OC43 strain ATCC VR-759, complete genome
>kraken:taxid|11020|NC_006213.1 Human coronavirus OC43 strain ATCC VR-759, complete genome
>kraken:taxid|1076255|NC_006213.1 Human coronavirus OC43 strain ATCC VR-759, complete genome
>kraken:taxid|10869|NC_006213.1 Human coronavirus OC43 strain ATCC VR-759, complete genome
```

Human coronavirus OC43 strain ATCC VR-759, complete genome
NCBI Reference Sequence: NC_006213.1
[FASTA](#) [Graphics](#)

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SAMPLE INFORMATION

Fragments Total	Fragments After De-hosting	Fragments Classified	Fragments Unclassified
6948512		3200309	373947

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DRAGEN METAGENOMICS RESULTS FOR SAMPLE [COVOC43_ILLUMINA_CDNA_SYNTHESIS_ENRICHED_REP1](#)

- Kraken2 Classification Report
- Organism Detection Report

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SAMPLE INFORMATION

Fragments Total	Fragments After De-hosting	Fragments
6948512		3200309

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







MiniKraken2: **15min**


Fragments Classified	Fragments Unclassified
2699548	295099


Extended Kraken2: **52min**

Fragments Classified	Fragments Unclassified
2994126	521

De-hostingの有無 CoV2+UHR

De-hosting ON **12min**
1,136 / 3,580 reads after De-hosting

De-hosting OFF **29min**
1,423 / 4,613,308 reads


DRAGEN RNA Pathogen Detection Pipeline

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Configuration

Reference ⓘ

Human (hg38) and Illumina Respiratory Virus Panel

Reference

- Human (hg38) and Illumina Respiratory Virus Panel
- Human (hg38) and Respiratory Viruses from the Seattle Flu Study

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Accession list for coverage plots and FASTA generation ⓘ

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Percent of pathogen genome with > 5x coverage

Color Palette:

Copy CSV Excel

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サンプル名	Human coronavirus OC43	Human Respiratory syncytial virus 9320 (type B)	Human metapneumovirus (CAN97-83)	Human parainfluenza virus 1	Human coronavirus NL63	Human coronavirus 229E	Human parainfluenza virus 2	Influenza A virus (A/Hong Kong/1073/99(H9N2))
CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Detected (99.97%)	ND (4.91%)	ND (3.25%)	ND (1.54%)	ND (0.56%)	ND (1.21%)	ND (1.20%)	ND (1.16%)
CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep2	Detected (99.96%)	ND (3.98%)	ND (2.43%)	ND (1.38%)	ND (0.62%)	ND (1.22%)	ND (1.37%)	ND (0.85%)
CoVOC43_Thermo_cDNA_Synthesis_Enriched_Rep1	Detected (100.00%)	Detected (10.15%)	Detected (7.73%)	ND (3.34%)	ND (1.76%)	ND (2.40%)	ND (3.21%)	ND (2.89%)
CoVOC43_Thermo_cDNA_Synthesis_Enriched_Rep2	Detected (100.00%)	Detected (10.46%)	Detected (8.43%)	ND (4.03%)	ND (2.01%)	ND (2.87%)	ND (3.74%)	ND (2.64%)

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
Sample Information

Type	Read 1	Read 2	Total Reads	Mapped Reads
Paired	151	151	6948512	3473135

Percentage of viral genome with coverage >1x, >5x

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Severe acute respiratory syndrome coronavirus 2 isolate Wuhan-Hu-1, complete genome

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Map/Align Metrics

マッピング/アライメントメトリクス (一部分のみ表示)		
	Total input reads	6948512
	Number of duplicate marked reads	5496499
	Number of duplicate marked and mate reads removed	NA
	Number of unique reads (excl. duplicate marked reads)	1452013
	Reads with mate sequenced	6948512
	Reads without mate sequenced	0
	QC-failed reads	0
	Mapped reads	6947188
	Mapped reads R1	3474053
	Mapped reads R2	3473135
	Number of unique & mapped reads (excl. duplicate marked reads)	1450689
	Unmapped reads	1324

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Total Input reads	6948512
Number of duplicate marked reads	5496499
Number of duplicate marked and mate reads removed	NA
Number of unique reads (excl. duplicate marked reads)	1452013
Reads with mate sequenced	6948512

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
VARIANT CALLER SUMMARY		Number of samples	1
VARIANT CALLER SUMMARY		Reads Processed	2403780
VARIANT CALLER SUMMARY		Child Sample	NA
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Total	483
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Biallelic	483
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Multiallelic	0
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	SNPs	476
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Insertions (Hom)	0
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Insertions (Het)	2
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Deletions (Hom)	0
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Deletions (Het)	5
VARIANT CALLER PREFILTER	CoVOC43_Illumina_cDNA_Synthesis_Enriched_Rep1	Indels (Het)	0

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
Target Regions, Custom Target BED File

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SAMPLE INFORMATION ⁱ

Total PF Reads	Unique PF Reads	Percent Q30 Bases	Mean Target Coverage Depth
1,025,784	580,990	95.16%	2390.1

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ENRICHMENT SUMMARY ⁱ

Target Manifest	Total Length of Targeted Reference	Number of Targeted Regions	Padding Size
RespiratoryVirus.20200409_used.bed	596,713 bp	118	150 bp

Note: All enrichment values are calculated without padding unless otherwise stated. If any targeted region overlaps another region, the region positions will be adjusted to remove overlaps (check manifest file in the output directory for differences). The total length of the targeted reference excludes no-call positions (Ns) in the reference sequence and might result in a smaller value than expected from the targeted regions coordinates alone.

Read Level Enrichment ⁱ

Total Aligned Reads	Percent Aligned Reads	Targeted Aligned Reads	Read Enrichment	Padded Target Aligned Reads	Padded Read Enrichment
580,388	99.90%	639,293	110.15%	639,293	110.15%

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Base Level Enrichment ⁱ

Total Aligned Bases	Percent Aligned Bases	Targeted Aligned Bases	Base Enrichment	Padded Target Aligned Bases	Padded Base Enrichment
71,668,731	99.87%	71,472,294	99.73%	71,472,294	99.73%

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SAMPLE INFORMATION ⁱ

Total PF Reads	Unique PF Reads	Percent Q30 Bases	Mean Target Coverage Depth
1,025,784	580,990	95.16%	2390.1

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ENRICHMENT SUMMARY ⁱ

Target Manifest	Total Length of Targeted Reference	Number of Targeted Regions	Padding Size
RespiratoryVirus.20200409_used.bed	596,713 bp	118	150 bp

Note: All enrichment values are calculated without padding unless otherwise stated. If any targeted region overlaps another region, the region positions will be adjusted to remove overlaps (check manifest file in the output directory for differences). The total length of the targeted reference excludes no-call positions (Ns) in the reference sequence and might result in a smaller value than expected from the targeted regions coordinates alone.

Read Level Enrichment ⁱ

Total Aligned Reads	Percent Aligned Reads	Targeted
580,388	99.90%	639,293

Base Level Enrichment ⁱ

Total Aligned Bases	Percent Aligned Bases	Targeted
71,668,731	99.87%	71,472,293

SMALL VARIANTS SUMMARY ⁱ

	SNVs	Insertions	Deletions
Total Passing	38	24	5
Het/Hom Ratio	0	0	0
Ts/Tv Ratio	1.11	-	-


VARIANTS TABLE ⁱ

📄 Load Variants Table

Please note that for samples with a large number of variants, loading the variants table may result in slow performance. You can alternatively view the variants table content outside of the browser by downloading the "varianttable.txt" file available under "Output Files".

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
Mean Target Coverage Depth

2390.1 X

Depth of Coverage in Targeted Regions

Target Coverage at 1X

99.91%

Target Coverage at 10X

99.90%

Note: Reads flagged as duplicates are excluded from the coverage statistics. The target coverage is based on the number of unique reads. No-calls (Ns) do not add to the target coverage.

FRAGMENT LENGTH SUMMARY ⁱ

Fragment Length Median

132 bp

Minimum

35 bp

Maximum

477 bp


Standard Deviation

69 bp

Note: The minimum and maximum are calculated from values within approximately three standard deviations (excluding the lower and upper 0.15% of the data) to account for potential outliers.

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
AmpliSeq for Illumina SARS-CoV-2 Research Panel

This sequencing panel offers a quick, accurate, and cost-effective method to research the SARS-CoV-2 virus. [Read More...](#)

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Panel

0 AmpliSeq™ SARS-CoV-2 Research Panel for Illumina® ⓘ
20020496 [Sign in to see pricing and favorite products.](#)

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- AmpliSeq for Illumina SARS-CoV-2 Panel Manifest File
product_file | ZIP< 1 MB
- AmpliSeq for Illumina SARS-CoV-2 Panel BED File
product_file | BED< 1 MB
- AmpliSeq for Illumina SARS-CoV-2 Custom Reference Genome File
product_file | ZIP< 1 MB

Upload Files

Select the type of files you want to upload

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The screenshot shows the 'DNA Amplicon v2.1.1' configuration form. It includes sections for 'Analysis Name', 'Save Results To', 'Biosample(s)', 'Custom Genome (Optional)', 'Targeted Amplicons', and 'Custom Manifest File'. Red dashed boxes and arrows highlight specific fields: 'Analysis Name' (containing 'DNA Amplicon 05/18/2020 7:09:30'), 'Custom Genome (Optional)' (containing 'MN908947v3.fasta'), 'Targeted Amplicons' (containing 'Custom Panels (select manifest file below)'), and 'Custom Manifest File' (containing 'SARS-CoV-2.dna_manifest.20200408.txt').

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SAMPLE INFORMATION ⁱ			
Total PF Reads	Percent Q30 Bases	Percent On-target Aligned Reads	Autosome Call Rate
1,025,784	95.16%	100%	0%

AMPLICON SUMMARY ⁱ	
Number of Amplicon Regions	Total Length of Target Regions
242	30,253 bp

Read Level Statistics ⁱ		
Read	Total Aligned Reads	Percent Aligned Reads
1	5,289	1.03%
2	5,275	1.03%

Base Level Statistics ⁱ				
Read	Percent Q30	Total Aligned Bases	Percent Aligned Bases	Mismatch Rate
1	96.08	690,747	0.89%	0.48%
2	94.23	689,647	0.89%	0.51%

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SAMPLE INFORMATION ⁱ

Total PF Reads	Percent Q30 Bases	Percent
1,025,784	95.16%	100%

AMPLICON SUMMARY ⁱ

Number of Amplicon Regions
242

Read Level Statistics ⁱ

Read	Total Aligned Reads
1	5,289
2	5,275

Base Level Statistics ⁱ

Read	Percent Q30	Total Aligned Bases
1	96.08	690,747
2	94.23	689,647


SMALL VARIANTS SUMMARY ⁱ

	SNVs	Insertions	Deletions
Total Passing	15	7	3
Het/Hom Ratio	2.75	6	-
Ts/Tv Ratio	0.667	-	-

COVERAGE SUMMARY ⁱ


Amplicon Mean Coverage	Uniformity of Coverage (Pct > 0.2*mean)
43.7	45.87

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
Easy Sharing & Collaboration


Push-button submission to public infectious disease databases

Share accepted

3 minutes ago

You accepted John Duddy's invitation to the Project [HiSeq 2000: NeoPre...](#)

Share pending

4 minutes ago

julie baher invited you to the Project **Study 123**

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Analysis Name
SRA Import 05/07/2020 10:34:41

SRA Accession ⓘ
SRX7857604

Save Results To
Select Project
Covid19_analysis

BaseSpace Labs Disclaimer
BaseSpace Labs
 I acknowledge and agree that (i) this is a BaseSpace Labs App, (ii) I am using it AS-IS without any warranty of any kind, (iii) Illumina has no obligation to provide any technical support for this App, and (iv) Illumina has no liability for my use of this App, including without limitation, any loss of data, incorrect results, or any costs, liabilities, or damages that result from use of this App.

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SRX7857604: RNA-Seq of Severe acute respiratory syndrome coronavirus 2
1 ILLUMINA (Illumina MiSeq) run: 61,329 spots, 11.1M bases, 2.5Mb downloads

Design: random hexamer first and second strand synthesis followed by Nextera XT
Submitted by: UNIVERSITY OF WASHINGTON
Study: Washington SARS-CoV-2 isolate sequences
[PRJNA610428](#) • [SRP251618](#) • [All experiments](#) • [All runs](#)
[show Abstract](#)

Sample:
[SAMN14308028](#) • [SRS6266706](#) • [All experiments](#) • [All runs](#)
Organism: [Severe acute respiratory syndrome coronavirus 2](#)

Library:
Name: WAB_UW5
Instrument: Illumina MiSeq
Strategy: RNA-Seq
Source: VIRAL RNA
Selection: unspecified
Layout: SINGLE

Runs: 1 run, 61,329 spots, 11.1M bases, [2.5Mb](#)

Run	# of Spots	# of Bases	Size	Published
SRR11247076	61,329	11.1M	2.5Mb	2020-03-05

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

The image shows the main configuration area of the GISAID Submission v1.0.0 application. It contains several sections, each highlighted with a red dashed box and connected to explanatory text by a red dashed arrow. The sections are: "Metadata file (CSV)", "Choose Source of Consensus FASTAs", "Consensus sequences for upload (FASTA)", "Variant files (VCF)", and "Reference genome that was used to call variants (hCoV-19/Wuhan/WIV04/2019|EPI_ISL_402124 recommended)".

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© 2019 Illumina, Inc. All rights reserved.