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
The screenshot shows the Illumina website interface. At the top, there is a navigation bar with links for SIGN IN, 予算申請書類リクエスト (Request for budget application documents), 論文要旨集リクエスト (Request for paper abstract collection), お問い合わせ (Contact us), and ログイン (Login). Below the navigation bar, there is a main content area with several sections:


- Support for Your System**: Includes a section for MiSeq with an image of the MiSeq sequencer.
- MiSeq Support Links**: Includes links for MiSeq System Guide, Questions and Answers, Bulletins, and a link to view all MiSeq support.
- Popular Kits for MiSeq**: Lists various sequencing kits like Nextera DNA Flex Kit, Nextera XT DNA Library Prep Kit, AmpliSeq for Illumina Custom DNA Panel, AmpliSeq for Illumina BRCA Panel, and AmpliSeq for Illumina Focus Panel, with a link to view all supported kits.
- Support for Popular Kits**: Includes links for Nextera DNA Library Prep Kit Support, MiSeq Reagent Kit Support, and links to view all sequencing and array kits.
- その他のリソース** (Other resources): Includes links for desktop sharing, support notifications, webinars, safety data sheets, Japanese protocol guides, and technical tips.
- Support for your organization**: Includes links to view all systems in your organization and reorder kits with Quick Order.

Below the main content area, there is a blue banner with the text "お客様のワークフローをサポート" (Support your workflow) and "お役立ちツールおよびヒント" (Helpful tools and hints). Underneath this banner, there is a workflow diagram with steps: Plan > Prep > Run > **Analyze** (highlighted with a red circle and a mouse cursor). To the right of the workflow, there are two promotional cards: "充実したカスタムプロトコールセレクトター" (Comprehensive custom protocol selector) and "ライブラリー調製キットおよびアレイキットセレクトター" (Library preparation kit and array kit selector). A blue button on the right says "すべてのサポートオプションを" (All support options).

At the bottom of the page, there are four small images: a book, a person's face, a person's face, and a person's face.

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
Support » Sequencing » Sequencing Software » Illumina Experiment Manager (IEM)

**Illumina Experiment Manager Support**

- Overview >
- Computing Requirements
- Downloads
- Documentation & Literature
- Training

## Illumina Experiment Manager

### Latest Updates

- Illumina Experiment Manager v1.15 01/28/2018**
- Illumina Experiment Manager User Guide 01/08/2018
- Illumina Experiment Manager v1.14 10/02/2017

### User Guide

#### Illumina Experiment Manager User Guide

The Illumina Experiment Manager software helps you create and edit well-formed sample sheets for Illumina sequencers and analysis software.


Illumina recommends using the Experiment Manager before starting sample or library preparation. The Experiment Manager can detect and warn of sub-optimal index combinations. By creating the sample sheet prior to sample or library preparation, you can try a different index combination without risking your samples.

You can use the Illumina Experiment Manager to create sample sheets for any Illumina sequencer and for any Nextera or TruSeq libraries.

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The screenshot shows the Illumina website's download page for Illumina Experiment Manager v1.15. The page includes a navigation bar with links like 'SIGN IN', 'お問い合わせ', and 'ロケーション選択'. Below the navigation bar, there's a breadcrumb trail: 'Support » Downloads » Illumina Experiment Manager v1.15'. The main heading is 'Illumina Experiment Manager v1.15', followed by a description: 'The Illumina Experiment Manager v1.15 is a desktop tool that guides you in building library plates and creating the sample sheet used with the MiniSeq, NextSeq, HiSeq, and NovaSeq systems.' Below this, there's a table of files for download. The first row, 'Illumina Experiment Manager v1.15.1', is circled in red. The table also lists software guides and release notes for v1.15.1 and v1.15.0.

FILE NAME	FILE INFO	DATE POSTED
<a href="#">Illumina Experiment Manager v1.15.1</a>	ZIP (2 MB)	01/29/2018
<a href="#">Illumina Experiment Manager v1.15 Software Guide (15031335 v07)</a>	PDF (< 1 MB)	01/09/2018
<a href="#">Illumina Experiment Manager Software v1.15.1 Release Notes</a>	PDF (< 1 MB)	01/28/2018
<a href="#">Illumina Experiment Manager Software v1.15.0 Release Notes</a>	PDF (< 1 MB)	01/28/2018

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




\*TruSight Panel製品, TruSeq Amplicon, TruSeq Exome, TruSeq Targeted RNA Expression, TruSeq Bovineなど

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	A	B	C	D	E	F	G	H	I	J	K	L	M	
1	[Header]													
2	Customer Name	Illumina												
3	Product Type	15025138												
4	Date Manufactured													
5	Lot													
6	DesignStudio ID	TSCA_Control												
7	Target Plexity	341												
8														
9	[Probes]													
10	Target Region Name	Target Region ID	Target ID	Species	Build ID	Chromosome	Start Position	End Position	Submitter	ULSO Sequence	ULSO Genotype	DLSO Sequence	DLSO Genotype	Probe
11	snp_rs10138962	snp_rs10138962	snp_rs101	Homo sap	hg19	chr14	1.05E+08	1.05E+08	+	AAGCCAA	282	TGCGCAC	1200	+
12	snp_rs1149580	snp_rs1149580	snp_rs114	Homo sap	hg19	chr11	76549431	76549489	+	GGATCAA	32	ACTGGGG	68	+
13	snp_rs10988298	snp_rs10988298	snp_rs109	Homo sap	hg19	chr9	1.32E+08	1.32E+08	+	GAGGCTG	1200	ATTGTGC	75	+
14	snp_rs10086550	snp_rs10086550	snp_rs100	Homo sap	hg19	chr8	1.43E+08	1.43E+08	+	GTGTGAC	51	CATAGAG	73	+
15	snp_rs10409978	snp_rs10409978	snp_rs104	Homo sap	hg19	chr19	57796594	57796646	+	TGTCTTAC	50	TTATGGG	115	+

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Illumina Experiment Manager

## Illumina Experiment Manager Settings

Sample Plate Repository

C:\Users\y\_aga\AppData\Roaming\Illumina Inc\Illumina Experiment Manager\Plates

Sample Sheet Repository

D:\Illumina\MiSeq Control Software\SampleSheets

Manifest Repository

D:\Illumina\MiSeq Control Software\Manifests

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


EM Illumina Experiment Manager

Illumina Experiment Manager

Sample Sheet Wizard - MiSeq Application Selection

Select Category

Targeted Resequencing

Small Genome Sequencing

RNA Sequencing

Other

Select Application

Resequencing

Plasmids

Assembly

Cancel

Back

Next

# Workflow Parameters の設定

Illumina Experiment Manager

## Sample Sheet Wizard - Workflow Parameters

**Resequencing Run Settings**

Reagent Cartridge Barcode\*

Library Prep Workflow

Index Adapters

Index Reads ☐ 0 (None) ☐ 1 (Single) ☒ 2 (Dual)

Experiment Name

Investigator Name

Description

Date

Read Type ☒ Paired End ☐ Single Read

Cycles Read 1

Cycles Read 2

\* - required field

**Resequencing Workflow-Specific Settings**

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# Workflow-Specific Settings の比較

## Resequencing

### Resequencing Workflow-Specific Settings

- ☐ Custom Primer for Read 1
- ☐ Custom Primer for Index
- ☐ Custom Primer for Read 2
- ☐ Use Somatic Variant Caller
- ☒ Flag PCR Duplicates
- ☐ Reverse Complement
- ☐ Indel Realignment GATK
- Variant Quality Filter
- ☐ Export to gVCF
- ☐ BWA-backtrack
- ☒ Use Adapter Trimming

## Enrichment

### Enrichment Workflow-Specific Settings

- ☒ Indel Realignment GATK
- ☐ Flag PCR Duplicates
- ☐ Run Picard HsMetrics
- Variant Quality Filter
- ☐ Export to gVCF
- ☐ BWA-backtrack
- ☒ Use Adapter Trimming
- ☒ Use Adapter Trimming Read 2

# Workflow-specific Setting: Resequencing

Illumina Experiment Manager Software Guide より

## Resequencing Workflow-Specific Settings

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# Workflow-specific Setting

Illumina Experiment Manager Software Guide より

## Enrichment Workflow-Specific Settings

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# Sample Selection ①

Illumina Experiment Manager

## Sample Sheet Wizard - Sample Selection

Sample Plate

Select Plate... New Plate...

Table View Plate View

Sample ID	Sample Name	Index Well	Index 1 (17)
-----------	-------------	------------	--------------

Samples to include in sample sheet

\* - required field ☐ Maximize

Sample ID*	Sample Name	Plate	Well	Index Well*	Index 1 (17)*	17 Seq
------------	-------------	-------	------	-------------	---------------	--------

# Sample Selection ①

Illumina Experiment Manager

## Sample Sheet Wizard - Sample Selection

Samples to include in sample sheet

\* - required field ☒ Maximize

Sample ID*	Sample Name	Plate	Well	Index1 (I7)*	I7 Sequence	Index2 (I5)*	I5 Sequence	Genome Folder*

< >

**Add Blank Row** Remove Selected Rows ?

Sample Sheet Status: Invalid  
Reason: Not all Samples in this sample sheet have all the required fields

Cancel Back Finish

# Sample Selection ①

Illumina Experiment Manager

## Sample Sheet Wizard - Sample Selection

Sample ID 1000 Index情報 \* - required field ☒ Maximize

Sample ID*	Sample Name	Plate	Well	Index Well*	Index1 (I7)*	I7 Sequence	Index2 (I5)*	I5 Sequence	Genome Folder*
PhiX			A01		D701	ATTACTCG	D501	TATAGCCT	
PhiV			A02		D702	TCCGGAGA	D501	TATAGCCT	
PhiV			A03		D703	CGCTCATT	D501	TATAGCCT	
PhiT			A04		D704	GAGATTCC	D501	TATAGCCT	
PhiS			A05		D705	ATTCAGAA	D501	TATAGCCT	
PhiR			A06		D706	GAATTCGT	D501	TATAGCCT	
PhiQ			A07		D707	CTGAAGCT	D501	TATAGCCT	
PhiP			A08		D708	TAATGCGC	D501	TATAGCCT	
PhiO			A09		D701	ATTACTCG	D502	ATAGAGGC	
PhiN			A10		D702	TCCGGAGA	D502	ATAGAGGC	
PhiM			A11		D703	CGCTCATT	D502	ATAGAGGC	
PhiL			A12		D704	GAGATTCC	D502	ATAGAGGC	
PhiK			B01		D705	ATTCAGAA	D502	ATAGAGGC	
PhiJ			B02		D706	GAATTCGT	D502	ATAGAGGC	
PhiI			B03		D707	CTGAAGCT	D502	ATAGAGGC	
PhiH			B04		D708	TAATGCGC	D502	ATAGAGGC	
PhiG			B05		D701	ATTACTCG	D503	CCTATCCT	
PhiF			B06		D702	TCCGGAGA	D503	CCTATCCT	
			B07						
			B08						
			B09						
			B10						
			B11						
			B12						

Sample Sheet Status: Invalid

Reason: Not all Samples in this sample sheet have all the required fields

# Sample Selection ②

Sample Sheet Wizard - Sample Selection

Samples to include in sample sheet

\* - required field ☒ Maximize

Sample ID*	Sample Name	Plate	Well	Index Well*	Index1 (I7)*	I7 Sequence	Index2 (I5)*	I5 Sequence	Genome Folder*
PhiX				A01	D701	Arabidopsis_thaliana#NCBI#build9.1#Sequence#WholeGenomeFasta			
PhiV				A02	D702	Bos_taurus#Ensembl#UMD3.1#Sequence#WholeGenomeFasta			
PhiT				A03	D703	Escherichia_coli_K_12_DH10B#NCBI#2008-03-17#Sequence#WholeGenomeFasta			
PhiS				A05	D705	Homo_sapiens#UCSC#hg19#Sequence#WholeGenomeFasta			
PhiR				A06	D706	Mus_musculus#UCSC#mm9#Sequence#WholeGenomeFasta			
PhiQ				A07	D707	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta			
PhiP				A08	D708	Rattus_norvegicus#UCSC#rn4#Sequence#WholeGenomeFasta			
PhiO				B01	D701	Saccharomyces_cerevisiae#UCSC#sacCer2#Sequence#WholeGenomeFasta			
PhiN				B02	D702	Staphylococcus_aureus_NCTC_8325#NCBI#2006-02-13#Sequence#WholeGenomeFasta			
PhiM				B03	D703	CGTCATT D502 ATAGAGGC			
PhiL				B04	D704	GAGATTCC D502 ATAGAGGC			
PhiK				B05	D705	ATTCAGAA D502 ATAGAGGC			
PhiJ				B06	D706	GAATTCGT D502 ATAGAGGC			
PhiI				B07	D707	CTGAAGCT D502 ATAGAGGC			
PhiH				B08	D708	TAATGCGC D502 ATAGAGGC			
PhiG				C01	D701	ATTACTCG D503 CCTATCCT			
PhiF				C02	D702	TCCGAGGA D503 CCTATCCT			

Sample Sheet Status: **Invalid**

Reason: **Not all Samples in this sample sheet have all the required fields**

## Sample Selection ③

Sample Sheet Wizard - Sample Selection

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
[Header]											
ITEMFileVersion	5										
Date	2018/7/										
Workflow	Resequencing										
Application	Resequencing										
Instrument Type	MiSeq										
Assay	TruSeq DNA PCR-Free										
Index Adapters	TruSeq DNA CD Indexes (96 Indexes)										
Description											
Chemistry	Amplicon										
[Reads]											
151											
151											
[Settings]											
FlagPCRDuplicates	1										
ReverseComplement	0										
VariantFilterQualityCut	30										
outputgenomevcf	FALSE										
RunBwaAln	0										
Adapter	AGATCGGAAGAGCACACGTCTGAACTCCAGTCA										
AdapterRead2	AGATCGGAAGAGCGTCGTGTAGGGAAAGAGTGT										
[Data]											
Sample_ID	Sample_Name	Sample_Plate	Sample_Well	Index_Plate	I7_Index_ID	index	I5_Index_ID	index2	GenomeFolder	Sample_Project	Description
PhiX				A01	D701	ATTACTCG	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiW				A02	D702	TCCGGAGA	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiV				A03	D703	CGCTCATT	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiT				A04	D704	GAGATTCC	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiS				A05	D705	ATTCAGAA	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiR				A06	D706	GAATTCGT	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiQ				A07	D707	CTGAAGCT	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiP				A08	D708	TAATGCGC	D501	TATAGCCT	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		
PhiO				B01	D701	ATTACTCG	D502	ATAGAGGC	PhiX#Illumina#RTA#Sequence#WholeGenomeFasta		

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
- ▶ **Illumina Experiment Manager Software Guide**

[http://support.illumina.com/sequencing/sequencing\\_software/experiment\\_manager.html](http://support.illumina.com/sequencing/sequencing_software/experiment_manager.html)

- ▶ **MiSeq Sample Sheet Quick Reference Guide**

[http://support.illumina.com/downloads/miseq\\_sample\\_sheet\\_quick\\_reference\\_guide\\_15028392.html](http://support.illumina.com/downloads/miseq_sample_sheet_quick_reference_guide_15028392.html)

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## IEM Quick Reference Card も参考になります

The screenshot shows the Illumina support website for the TruSight One Sequencing Panel. The sidebar on the left contains links to Downloads, Documentation & Literature (highlighted), Training, Questions & Answers, Best Practices, Bulletins, Webinars, and Product Ordering Information. The main content area features a welcome message and links to the TruSight One Sequencing Panel Protocol Guide, Checklist, and Consumables and Equipment. Below this is a section titled 'Experimental Design Documentation' which contains a table with columns for DESCRIPTION, FILE INFO, and DATE POSTED. The table lists the 'IEM TruSight Quick Reference Card' as a PDF document posted on 02/03/2016. The page also includes a 'Supporting Documentation' section at the bottom.

DESCRIPTION	FILE INFO	DATE POSTED
<a href="#">IEM TruSight Quick Reference Card</a> Instructions for using IEM with TruSight kits.	PDF (< 1 MB)	02/03/2016

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