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**nature biotechnology** ARTICLES

### Full-length transcriptome assembly from RNA-Seq data without a reference genome

Manfred G Grabherr<sup>1,2</sup>, Brian J Haas<sup>1,3,4</sup>, Moran Yassour<sup>1,3,4</sup>, Joshua Z Levin<sup>1</sup>, Dawn A Thompson<sup>1</sup>, Ido Amit<sup>1</sup>, Xian Adiconis<sup>1</sup>, Lin Fan<sup>1</sup>, Raktima Raychowdhury<sup>1</sup>, Qiangdong Zeng<sup>1</sup>, Zehua Chen<sup>1</sup>, Evan Muehler<sup>1</sup>, Nir Hacohen<sup>1</sup>, Andreas Gnirke<sup>1</sup>, Nicholas Rhind<sup>1</sup>, Federica Di Palma<sup>1</sup>, Bruce W Birren<sup>1</sup>, Chad Nusbaum<sup>1</sup>, Kerstin Lindblad-Toh<sup>1,5</sup>, Nir Friedman<sup>1,6</sup> & Aviv Regev<sup>1,2,7</sup>

Massively parallel sequencing of cDNA has enabled deep and efficient probing of transcriptomes. Current approaches for transcript reconstruction from such data often rely on aligning reads to a reference genome, and are thus unsuitable for samples with a partial or missing reference genome. Here we present the Trinity method for *de novo* assembly of full-length transcripts and evaluate it on samples from fission yeast, mouse and whitefly, whose reference genome is not yet available. By efficiently constructing and analyzing sets of de Bruijn graphs, Trinity fully reconstructs a large fraction of transcripts, including alternatively spliced isoforms and transcripts from recently duplicated genes. Compared with other *de novo* transcriptome assemblers, Trinity recovers more full-length transcripts across a broad range of expression levels, with a sensitivity similar to methods that rely on genome alignments. Our approach provides a unified solution for transcriptome reconstruction in any sample, especially in the absence of a reference genome.

Recent advances in massively parallel cDNA sequencing (RNA-Seq) provide a cost-effective way to obtain large amounts of transcriptomic data from many organisms and tissue types<sup>1,2</sup>. In principle, such data can allow us to identify all expressed transcripts<sup>3,4</sup>, as complete and contiguous mRNA sequences from the transcription start site to the transcription end, for multiple alternatively spliced isoforms. However, reconstruction of all full-length transcripts from short reads with considerable sequencing error rates poses substantial computational challenges<sup>5,6</sup>: (i) some transcripts have low coverage, whereas others are highly expressed; (ii) read coverage may be uneven across the transcript's length, owing to sequencing biases; (iii) reads with sequencing errors derived from a highly expressed transcript may be more abundant than correct reads from a transcript that is not highly expressed; (iv) transcripts encoded by adjacent loci can overlap and thus can be erroneously fused to form a chimeric transcript; (v) data structures need to accommodate multiple transcripts per locus, owing to alternative splicing; and (vi) sequences that are repeated in different genes introduce ambiguity. A successful method should address each challenge, be applicable to both complex mammalian genomes and gene-dense microbial genomes, and be able to reconstruct transcripts of variable sizes, expression levels and protein-coding capacity.

There are two alternative computational strategies for transcriptome reconstruction<sup>7</sup>. Mapping-first approaches<sup>8</sup>, such as Scripture<sup>9</sup> and Cufflinks<sup>10</sup>, first align all the reads to a reference (unannotated) genome and then merge sequences with overlapping alignment, spanning splice junctions with reads and paired ends. Assembly-first (*de novo*) methods, such as ABySS<sup>11</sup>, SOAPdenovo<sup>12</sup> or Oases (E. Birney, European Bioinformatics Institute, personal communication), use the reads to assemble transcripts directly, which can be mapped subsequently to a reference genome, if available. Mapping-first approaches promise, in principle, maximum sensitivity but depend on correct read-to-reference alignment, a task that is complicated by splicing, sequencing errors and the lack or incompleteness of many reference genomes. Conversely, assembly-first approaches do not require any read-reference alignments, important when the genomic sequence is not available, is spliced, highly fragmented or substantially altered, as in cancer cells.

Successful mapping-first methods were developed in the past year<sup>8</sup>, but substantially less progress was made to date in developing effective assembly-first approaches. As the number of reads grows, it is increasingly difficult to determine which reads should be joined into contiguous sequence contigs. An elegant computational solution is provided by the *de Bruijn* graph<sup>13</sup>, the basis for several whole-genome assembly programs<sup>14,15</sup>. In this graph, a node is defined by a sequence of a fixed length *k* nucleotides (a *k*-mer with *k* considerably shorter than the read length), and nodes are connected by edges, if they perfectly overlap by *k* - 1 nucleotides, and the sequence data support this connection. This compact representation allows for enumerating all possible solutions by which linear sequences can be reconstructed given overlaps of *k* - 1.

<sup>1</sup>Broad Institute of Massachusetts Institute of Technology and Harvard, Cambridge, Massachusetts, USA. <sup>2</sup>School of Computer Science, Brown University, Providence, Rhode Island, USA. <sup>3</sup>Department of Biology, Massachusetts Institute of Technology, Cambridge, Massachusetts, USA. <sup>4</sup>Department of Biochemistry and Molecular Pharmacology, University of Massachusetts Medical School, Worcester, Massachusetts, USA. <sup>5</sup>Science to Life Laboratory, Department of Medical Biochemistry and Microbiology, Uppsala University, Uppsala, Sweden. <sup>6</sup>Maximilian Ruberman Institute of Life Sciences, Hebrew University, Jerusalem, Israel. <sup>7</sup>Howard Hughes Medical Institute, Department of Biology, Massachusetts Institute of Technology, Cambridge, Massachusetts, USA. <sup>8</sup>These authors contributed equally to this work. Correspondence should be addressed to M.G. (grabherr@mit.edu) or A.R. (regev@mit.edu).

Received 3 December 2010; accepted 28 April 2011; published online 15 May 2011; doi:10.1038/nbt.1883

**NATURE BIOTECHNOLOGY** ADVANCE ONLINE PUBLICATION

Grabherr MG, Haas BJ, Yassour M et al. (2011) Full-length transcriptome assembly from RNA-Seq data without a reference genome. *Nat Biotechnol.* 2011 May 15

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
July 2008 - Vol 5 No 7

*Nature Methods* 5 (2008)

The beginning of the end for microarrays?

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### Review

*Nature Reviews Genetics* 12, 87-98 (February 2011) | doi:10.1038/nrg2934

ARTICLE SERIES: [Applications of next-generation sequencing](#)


**RNA sequencing: advances, challenges and opportunities**

Fatih Ozsolak<sup>1</sup> & Patrice M. Milos<sup>1</sup> [About the authors](#)


PMID: 21191423

## RNA-seq for detection of alternative splicing events

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

## Identification of the transforming *EML4-ALK* fusion gene in non-small-cell lung cancer

Manabu Soda<sup>1,2</sup>, Young Lim Choi<sup>1</sup>, Munehiro Enomoto<sup>1,2</sup>, Shuji Takada<sup>1</sup>, Yoshihiro Yamashita<sup>1</sup>, Shunpei Ishikawa<sup>5</sup>, Shin-ichiro Fujiwara<sup>1</sup>, Hideki Watanabe<sup>1</sup>, Kentaro Kurashina<sup>1</sup>, Hisashi Hatanaka<sup>1</sup>, Masashi Bando<sup>2</sup>, Shoji Ohno<sup>2</sup>, Yuichi Ishikawa<sup>6</sup>, Hiroyuki Aburatani<sup>3,7</sup>, Toshiro Niki<sup>3</sup>, Yasunori Sohara<sup>4</sup>, Yukihiko Sugiyama<sup>2</sup> & Hiroyuki Mano<sup>1,7</sup>

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Vol. 25 no. 9 2009, pages 1105–1111  
doi:10.1093/bioinformatics/btp120

Sequence analysis

### TopHat: discovering splice junctions with RNA-Seq

Cole Trapnell<sup>1,\*</sup>, Lior Pachter<sup>2</sup> and Steven L. Salzberg<sup>1</sup>

<sup>1</sup>Center for Bioinformatics and Computational Biology, University of Maryland, College Park, MD 20742 and

<sup>2</sup>Department of Mathematics, University of California, Berkeley, CA 94720, USA

Received on October 23, 2008; revised on February 24, 2009; accepted on February 26, 2009

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
Transcript assembly and quantification by RNA-Seq reveals unannotated transcripts and isoform switching during cell differentiation


Cole Trapnell<sup>1-3</sup>, Brian A Williams<sup>4</sup>, Geo Pertea<sup>2</sup>, Ali Mortazavi<sup>4</sup>, Gordon Kwan<sup>4</sup>, Marijke J van Baren<sup>5</sup>, Steven L Salzberg<sup>1,2</sup>, Barbara J Wold<sup>4</sup> & Lior Pachter<sup>3,6,7</sup>

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
### Resource

#### Systematic identification of long noncoding RNAs expressed during zebrafish embryogenesis

Andrea Pauli,<sup>1,7,8</sup> Eivind Valen,<sup>2,7</sup> Michael F. Lin,<sup>3,4</sup> Manuel Garber,<sup>4</sup> Nadine L. Vastenhouw,<sup>1</sup> Joshua Z. Levin,<sup>4</sup> Lin Fan,<sup>4</sup> Albin Sandelin,<sup>2</sup> John L. Rinn,<sup>4,5</sup> Aviv Regev,<sup>3,4,6,8</sup> and Alexander F. Schier<sup>1,4,8</sup>

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PMID:19188922

nature  
REVIEWS GENETICS

Long non-coding RNAs:  
insights into functions

Tim R. Mercer, Marcel E. Dinger and John S. Mattick

PROGRESS

multigene transcripts that encompass several genes or even the whole chromosome, further confound efforts for systematic classification<sup>23</sup>.

**Widespread functionality of long ncRNAs**  
Given their unexpected abundance, long ncRNAs were initially thought to be spuri-


#### The Xist lncRNA Exploits Three-Dimensional Genome Architecture to Spread Across the X Chromosome

Jesse M. Engreitz,<sup>1,2</sup> Amy Pandya-Jones,<sup>3</sup> Patrick McDonel,<sup>1</sup> Alexander Shishkin,<sup>1</sup> Klara Sirokman,<sup>1</sup> Christine Surka,<sup>1</sup> Sabah Kadri,<sup>1</sup> Jeffrey Xing,<sup>1</sup> Alon Goren,<sup>1</sup> Eric S. Lander,<sup>1,4,5\*</sup> Kathrin Plath,<sup>3\*</sup> Mitchell Guttman<sup>1\*†</sup>

PMID: 23828888

#### RNA Antisense Purification (RAP):


A Method to Map lncRNA Interactions with Chromatin

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SCIENTIFIC REPORTS


OPEN

## Technical Variations in Low-Input RNA-seq Methodologies

SUBJECT AREAS: TRANSCRIPTOMICS

Vipul Bhargava<sup>1</sup>, Steven R. Head<sup>2</sup>, Phill

: PMID: 24419370



Smart-Seq、DP-Seq、CEL-Seqの比較

illumina®

# Single Cell Transcriptome解析

ddSEQ™ Single-Cell Isolator  
Bio-Rad

SureCell™ WTA 3' Library Prep Kit  
illumina



Fig. 2. Unbiased Cluster Analysis of PBMCs in Seurat—Nine cell clusters were identified with 3354 cells down-sampled to approximately 70,000 reads per cell, with a resolution setting of 0.80 and 100 genes as a cut off.

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
The screenshot shows the Illumina website's navigation bar with the logo and links for Customer Login, MyIllumina, Quick Order, Contact Us, and Select Language. Below the navigation bar are menu items for AREAS OF INTEREST, TECHNIQUES, SYSTEMS, PRODUCTS & SERVICES, INFORMATICS, SCIENTIFIC CONTENT, COMPANY, and SUPPORT, along with a search icon. The main content area features a decorative graphic of four overlapping diamond shapes containing icons for DNA, a molecular structure, a DNA double helix, and a microarray chip. Below this graphic is the title "Library Prep and Array Kit Selector" in orange text, followed by a subtitle: "This tool will help you determine the best kit for your needs based on your project type, starting material, and the method or application." A form section titled "Please select your project type :" contains two radio button options: "Research Use Only" and "Molecular Diagnostics", each with an information icon.

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Low Input

Ribosomal Profiling

	TruSeq RNA v2	TruSeq Stranded mRNA	TruSeq RNA Access
Intended Use	Simple, cost-effective solution for analysis of the coding transcriptome.	Clearest, most complete view of the coding transcriptome with precise strand information.	Clearest, most complete view of the coding transcriptome with precise strand information and reduced input requirement.
Available for NeoPrep System	No	Yes	No
Input Amount	0.1 – 1ug High Quality Total RNA 10-400ng previously isolated mRNA	NeoPrep: 25-100ng LT/HT: 0.1 – 1ug Total RNA or 10-100ng previously isolated mRNA	10ng High Quality Total RNA 20ng Degraded Total RNA
Hands-on Time <sup>a</sup>	~4.5 hours	NeoPrep: ~30 minutes LT/HT: ~4.5 hours	~11 hours
Assay Time <sup>b</sup>	~10.5 hours	NeoPrep: ~10.5 hours LT/HT: ~10.5 hours	~2 days
FFPE Compatible	No	No	Yes
Capture Method	Oligo dT beads capture poly-A tail	Oligo dT beads capture poly-A tail	Capture probes targeting coding RNA sequence
Capture Content	Coding Transcriptome	Coding Transcriptome	Coding Transcriptome
Stranded	No	Yes	Yes

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TruSeq Stranded mRNA Sample Prep Kit

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TruSeq RNA Access/ TruSight Pan-Cancer Library Prep Kit

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TruSeq Stranded Total RNA Ribo-Zero H/M/R

TruSeq Stranded Total RNA Ribo-Zero Gold

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TruSeq small RNA Sample Prep Kit

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TruSeq Stranded Total RNA RiboZero Plant kit

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Save Results To: Select Project(s):

Samples: Select Sample(s):

Select All	Stranded?	
	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	
BrainmRNA-20131201	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	x
BrainTotalRNA-20131117	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	x
UHRRmRNA-20131205	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	x
BrainTotalRNA-20131118	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	x
UHRRTotalRNA-20131122	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	x
BrainmRNA-20131202	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	x
UHRRmRNA-20131206	<input type="radio"/> No <input checked="" type="radio"/> First Strand <input type="radio"/> Second Strand	x

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QC Reads (Read Pairs): 100000

Novel Transcript Assembly:

Adjust Transcript Assembly for Samples Without PolyA Selection:

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