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
Platform	PhiX Aligned (%)†
iSeq 100	minimum 5%
MiniSeq	10-50%*
MiSeq (MCS 2.2 or higher)	minimum 5%
NextSeq	10-50%*
HiSeq 2500 (HCS 2.2.38 or higher)	minimum 10%
HiSeq 3000/4000 (HCS 3.3.76 or lower)	10-50%*
HiSeq 3000/4000 (HCS 3.4.0 or higher)	5-20%*
NovaSeq	minimum 10%

How much PhiX spike-in is recommended when sequencing low diversity libraries on Illumina platforms?  
[https://jp.support.illumina.com/content/illumina-support/apac/ja\\_JP/bulletins/2017/02/how-much-phix-spike-in-is-recommended-when-sequencing-low-divers.html](https://jp.support.illumina.com/content/illumina-support/apac/ja_JP/bulletins/2017/02/how-much-phix-spike-in-is-recommended-when-sequencing-low-divers.html)

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

Table 1: Optimal Raw Densities for Illumina sequencing systems

	MiniSeq™	MiSeq		NextSeq®	HiSeq 2500 Rapid Run (RR)	HiSeq 2500 High Output (HO)	
Versions	High and Mid Output	v2	v3	v2 High and Mid Output	v1 and v2	v3	v4
Raw Density (K/mm <sup>2</sup> )	170–220	1000–1200	1200–1400	170–220	850–1000	750–850	950–1050

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

Cluster Optimization Overview Guide

[https://support.illumina.com/content/dam/illumina-support/documents/documentation/system\\_documentation/cluster-optimization-overview-guide-1000000071511-00.pdf](https://support.illumina.com/content/dam/illumina-support/documents/documentation/system_documentation/cluster-optimization-overview-guide-1000000071511-00.pdf)

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
- **Technical Notes:**

[Low-Diversity Sequencing on the Illumina MiSeq Platform](#)

[Low-Diversity Sequencing on the Illumina HiSeq Platform](#)

[Cluster Optimization Overview Guide](#)

[Using a PhiX Control for HiSeq Sequencing Runs](#)

- **Application Notes:**

[Whole-Genome Bisulfite Sequencing on the HiSeq 3000/HiSeq 4000 Systems](#)

[16S Metagenomics Studies with the MiSeq System](#)

[16S metagenomics sequencing with the iSeq™ 100 System](#)

- **Support Bulletins:**

[Best Practices for Low-Diversity Sequencing on the NextSeq and MiniSeq Systems](#)

[How much PhiX spike-in is recommended when sequencing low diversity libraries on Illumina platforms?](#)

[Cluster density guidelines for Illumina sequencing platforms](#)

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- **Support Bulletins:**

[Library quantification and quality control quick reference guide](#)

[Cluster Density Specifications for Illumina Sequencing Platforms](#)

[Considerations when migrating non-Illumina libraries between sequencing platforms](#)

[Cluster density considerations when migrating Illumina libraries between sequencing platforms](#)

[How to achieve more consistent cluster density on Illumina sequencing platforms](#)

[Nextera Library Validation and Cluster Density Optimization](#) – Lists some guidelines for loading concentrations based on average size of library

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© Illumina Sequencing Analysis Viewer 2.4.7 - 180305\_M04219\_0169\_000000000-BNDF6