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

A shotgun metagenomics NGS workflow for assessing microbial populations in complex samples

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https://support.illumina.com/sequencing/sequencing_instruments/nextseq-1000-2000/downloads.html

DRAGEN BCL Convert Workflow Installer

Installation and release note files for DRAGEN BCL Convert (Generate FastQ) Workflow.

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Reference Builder (Instruments)
Illumina, Inc.
[Bookmark this app](#) [? Help](#)

FOR RESEARCH USE ONLY

This version of the Reference Builder (Instruments) App includes the following features:

Pricing
Compute cost Free

Version
2.0.0

[Launch Application](#)

Additional Information
[Release Notes](#)
Version 2.0.0

New Features

- Upgraded to DRAGEN 4.2.7.
- Generates v9 references.

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[Product Documentation]: [Denature and Dilute Libraries for the NextSeq 1000/2000 Sequencing Systems \(Onboard\)](#)

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
Library Type	XLEAP-SBS Final Loading Concentration (pM)			
	XLEAP-SBS P1	XLEAP-SBS P2	XLEAP-SBS P3	XLEAP-SBS P4
Illumina DNA Prep with Enrichment	1000	1000	1000	1000
Illumina Stranded Total RNA with Ribo-Zero Plus	750	750	750	750
TruSeq DNA Nano 550 (600-cycle kit only)	1500	1500	Not applicable	Not applicable
100% PhiX (300-cycle kits and below)	650	650	450	450
100% PhiX (600-cycle kit only)	450	450	Not applicable	Not applicable

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Standard

RTA3

2-Channel Chemistry				
		G		
	A		T	C
Image 1				
Image 2				

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RTA4

2-Channel Chemistry				
		G		
	C		T	A
Image 1				
Image 2				

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- [Illumina Knowledge](#)

The screenshot shows the Illumina Knowledge website. At the top left is the Illumina logo. To the right of the logo are the links "Illumina.com" and "Illumina Support". Below the logo is a navigation menu with the following items: "Illumina Knowledge", "What's New", "Announcements", "INSTRUMENTATION", "General", "iSeq 100", "MiniSeq", "MiSeq", "NextSeq 500/550", "NextSeq 1000/2000" (highlighted with a yellow box), "cBot", "HiSeq 1000/2500", "HiSeq X", and "NovaSeq X/X Plus". The "NextSeq 1000/2000" dropdown menu is open, showing "FAQ", "Reference Material", and "Troubleshooting". The main content area is titled "Illumina Knowledge" and contains a large orange banner. Below the banner is a paragraph: "Illumina Knowledge is a repository of 1310 FAQs, troubleshooting articles and reference material for Illumina products and workflows, covering all Illumina Microarrays and every stage of Sequencing. Use the navigation bar on the left of the page to find areas of interest and content types, use the search function at the top right to query keywords." Below this is a section titled "What's New" with four bullet points: "Instructions to back up DRAGEN license files when upgrading the operating system on a DRAGEN server (March 25, 2024)", "Proper care of IAPS with ILASS robot fluidics (March 25, 2024)", "Troubleshooting Illumina Service Connector not connected for Window computer (March 25, 2024)", and "What is the required iScan Control Software version and scan setting for Infinium EX BeadChips? (March 25, 2024)".

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- [Common Questions Related to Sequencing Analysis Viewer \(SAV\) software and the NextSeq 1000/2000](#)
- [How to create new user accounts on the NextSeq 1000/2000](#)
- [How to obtain a NextSeq 1000/2000 v2 sample sheet?](#)
- [How To Power Cycle NextSeq 1000/2000 Video](#)
- [Understanding loading and percent loading concentration on the NextSeq 1000/2000](#)
- [Troubleshooting Failed to execute connectivity check input parameters on the NextSeq 1000/2000](#)
- [Troubleshooting 'mount' Errors when Networking NextSeq 1000/2000 Instruments](#)

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Illumina Knowledgeより

[Library Preparation]

- [What is nucleotide diversity and why is it important?](#)

[Run Planning]

- [Why sequencing 26 or more cycles in Read 1 is recommended](#)
- [How many cycles of SBS chemistry are in my kit?](#)
- [Maximum read length for Illumina sequencing platforms](#)
- [How to thaw and store sequencing reagents for optimal performance](#)
- [Run time estimates for each sequencing step on the Illumina sequencing platforms](#)
- [Approximate sizes of sequencing run output folders](#)

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- [NextSeq 1000/2000: Introduction](#)
- [NextSeq 1000/2000: Introduction to Analysis with DRAGEN](#)
- [Illumina sequencing innovations powered by XLEAP-SBS chemistry](#)
- [Advancing epigenetic breakthroughs with the Illumina NextSeq 1000/2000](#)

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- [Index Adapters Pooling Guide](#)
- [Index Sequencing Overview Guide](#)

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© 2024 Illumina, Inc. All rights reserved.