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The screenshot displays the Illumina Connected Analytics interface. The top navigation bar is blue and contains the Illumina logo, the text 'Connected Analytics', and a user profile icon with the text '<User名>'. The left sidebar is white and contains a menu with the following items: 'Projects' (highlighted with an orange border), 'Reference Data', 'Bundles', 'Event log', 'System Settings' (expanded), 'Notifications', 'Metadata Models', 'Docker Repository', 'Tool repository', 'Credentials', and 'Storage'. Below the sidebar is an orange button labeled 'Main Menu'. The main content area is titled 'Projects' and features a '+ Create' button. It includes filters for 'Show only favorites' (checked) and 'Show hidden' (unchecked), a 'Filter by workgroup' dropdown, and a search bar labeled 'Search specific projects'. Two project cards are displayed:

Project Name	Location	Size	Analyses	My role
Training_CLS	United States - North Virginia (US)	321.2 MB	2	Administrator
TrainingProject	United States - North Virginia (US)	22.9 KB		Administrator

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The screenshot displays the 'Team' management interface in the Illumina Connected Analytics application. The left sidebar shows the navigation menu with 'Team' highlighted. The main content area is divided into 'Project owner' and 'Members' sections. An orange callout box highlights the 'Add' dropdown menu, which lists three options: 'Workgroup', 'User of your tenant', and 'User by email'. The 'Add' button is also highlighted with an orange arrow.

Team

Project owner

Name

Email

Members

+ Add

- Workgroup
- User of your tenant
- User by email

Remove Resend invitation Export

Search specific members

status Tenant Email Download Upload Project access Flow access Base access Bench access

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Name User2

Email User2@illumina.com

Tenant global-support-training , United States

Members

+ Add Remove Resend invitation Export

☐	Name	Status	Tenant	Email	Download	Upload	Project access	Flow access	Base access	Bench access
☐	AMR Training	✓	n/a	n/a	✓	-	Viewer	Viewer	No access	No access
☐	User	?	n/a	User@illumina.com	-	✓	Contributor	Contributor	Viewer	Contributor
☐	User2	✓	global-suppor...	User2@illumina.com	✓	✓	Administrator	Contributor	Contributor	Contributor

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Data management

2-4

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The screenshot displays the Illumina Connected Analytics web interface. The left sidebar shows a navigation menu with 'Data' highlighted under the 'Training_CLS' project. The main content area is titled 'Data' and shows a list of files. At the top right, it indicates 'Total Data : 321.16 MB' and 'Total Analyses : 2'. Below the title, there are tabs for 'Data' and 'Data History', and a row of action buttons: 'Schedule for Download', 'View', 'Load to', 'Manage', and 'Export'. A large grey box with a dashed border contains the text 'Choose a file or drag it here' and a note about upload speed. Below this, a filter bar shows 'Showing Data in top level' and various filter options (Type, Name, Status, Format) along with 'Add filter' and 'Clear filters' buttons. A search bar on the right contains the text 'Training_CLS'. The main table lists files with columns for checkboxes, Type, Path, Status, Format, Tenant, Size, and Owning Project. All files are listed as 'Available' and are owned by 'Training_CLS'.

	Type	Path	Status	Format	Tenant	Size	Owning Project
☐	File	/Sample-B_R1.fastq.gz	Available	FASTQ	global-support-training , United States	8.61 MB	Training_CLS
☐	File	/sampleX.final.count.tsv	Available	TSV	global-support-training , United States	22.9 KB	Training_CLS
☐	File	/Sample-A_R1.fastq.gz	Available	FASTQ	global-support-training , United States	8.61 MB	Training_CLS
☐	File	/forward_reads.fastq	Available	FASTQ	global-support-training , United States	14 B	Training_CLS
☐	File	/reverse_reads.fastq	Available	FASTQ	global-support-training , United States	14 B	Training_CLS
☐	File	/sampleY.expression.txt	Available	TXT	global-support-training , United States	22.9 KB	Training_CLS
☐	File	/Sample-C_R1.fastq.gz	Available	FASTQ	global-support-training , United States	8.61 MB	Training_CLS

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Start analysis

General

User reference *
Gemline_test
Enter a reference to be able to identify this analysis

Pipeline
DRAGEN_Germline_Whole_Genome_4-4-6

User tags
Type to add a tag

Notify me when analysis is completed
Type or select an email address
Leave blank if you don't want to be notified

Output folder
/
Select a folder in which the output folder of the analysis should be located. When no folder is selected, the output folder will be located in the root of the project.

Logs folder
/
Select a folder in which the logs of the analysis should be located. When no folder is selected, the logs folder will be located in the output folder.

Pricing

Subscription *
ICA_Ent-JP_Pipeline_Entitlement - ICA_lap-apn1_IAPEntSubscription - HPK2-A1ML-6K1I-ZHHJ

Input files

FASTQs
No files selected

Input FASTQ files. Multiple samples may be provided. If a supplementary FASTQ list is not provided, the pipeline will attempt to determine which FASTQ files belong together via the naming scheme. The RGSM will be taken as the name up to "_SX" and the suffix after "_RX." "_LXXX" in the filename denotes the lane number and "_RX" to denotes the read number. If either is omitted, lane 1 and read 1 will be used in

Settings

Map Align Options

Enable Map/Align
☒ True ☐ False
Corresponds to DRAGEN argument "--enable-map-align".

Map/Align Output
BAM
Corresponds to DRAGEN argument "--output-format".

Enable Duplicate Marking
☒ True ☐ False
Corresponds to DRAGEN argument "--enable-duplicate-marking".

Downsample Fraction

Keep a fraction of input reads to perform DRAGEN analysis. If set, runs DRAGEN with --enable-fractional-down-sampler true and --down-sampler-normal-subsample [fraction]. This is applied directly to all the input reads prior to any filtering or trimming.

Small Variant Caller Options

Enable Variant Caller
☒ True ☐ False
Corresponds to DRAGEN argument "--enable-variant-caller".

Emit Ref Confidence

Corresponds to DRAGEN argument "--vc-emit-ref-confidence". With GVCF, reference calls are grouped into blocks. Typically, BP_RESOLUTION is only used for debugging, as it takes a longer time to run, and will create very large GVCF files (around 2 hours for a 30x sample, generating 20 GB GVCF).

Enable Personalization
☒ True ☐ False
Enable computing a personalized 2-haplotype reference to impute variants and boost accuracy of variant calling (Default=false). Personalization is only supported when used with a pangenome reference with map align enabled. Runtime will increase when enabling this option. Corresponds to DRAGEN argument "--enable-personalization".

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User reference *

Somatic_test

Enter a reference to be able to identify this analysis

Pipeline

DRAGEN_Somatic_4-4-6

User tags

Type to add a tag

Notify me when analysis is completed

Type or select an email address

Leave blank if you don't want to be notified

Output folder

/

Select a folder in which the output folder of the analysis should be located. When no folder is selected, the output folder will be located in the root of the project.

Logs folder

/

Select a folder in which the logs of the analysis should be located. When no folder is selected, the logs folder will be located in the output folder.

Pricing

Subscription *

ICA_Ent-JP_Pipeline_Entitlement - ICA_iap-apn1_JAPEntSubscription - HPK2-A1ML-6K1I-ZHHJ

Resources

Storage size *

Medium

2.4TB

Input

Input Type

FASTQ GZ

Choose an input type

Tumor FASTQs

No files selected

Input FASTQ files. Multiple samples may be provided. The RGSM will be taken as the name up to ".SX" and the suffix after ".RX" ".LXXX" in the filename denotes the lane number and ".RX" to denotes the read number. If either is omitted, lane 1 and read 1 will be used in the FASTQ list. The tool will automatically write a FASTQ list from all files provided and process each sample.

Normal FASTQs

No files selected

Optional input normal FASTQ files. Only one tumor-normal pair may be provided. If a supplementary FASTQ list is not provided, the pipeline will attempt to determine which FASTQ files belong together via the naming scheme using the method described above. Detection of multiple normal samples will result in an error.

Tumor FASTQ List

No file selected

A supplementary FASTQ list file can be provided to specify the groups of tumor FASTQ files. For the file path, use ONLY the filename of the file. FASTQs still need to be selected via the FASTQs field.

Normal FASTQ List

No file selected

A supplementary FASTQ list file can be provided to specify the groups of normal FASTQ files. Only one sample in the normal FASTQ list as the pipeline can only analyze one tumor-normal pair. For the file path, use ONLY the filename of the file. FASTQs still need to be selected via the FASTQs field.

Reference

Homo sapiens [1000 Genomes] hg38 v5

Select the reference genome.

Target BED File

No file selected

Within the target BED file, the contig names must match those of the chosen reference, and the BED file must be sorted. This file is applicable only to the variant calling stage.

Sample Sex

Auto-Detect

Library Specific Settings

Library preparation kit used to process the samples. The Illumina FFPE DNA Prep with Exome 2.5 Enrichment library is only applicable in Tumor-Normal mode.

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☐ rnaTest001 Preparing Inputs Apr, 9 2025 15:28:26 DRAGEN RNA 4-2-4-v2

Total data : 129.6 KB Total analyses : 1

☐ Show only analyses I have started Search specific analyses

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← Projects
Training03
Data
Samples
Activity
▼ Flow
Reference Data
Pipelines
Analyses
▼ Base
▼ Bench
Workspaces
▼ Variant Viewer
▼ Cohorts
▼ Project Settings
Details
Team
Connectivity
Notifications

Data

Total data : 156.49 MB Total analyses : 1

Choose a file or drag it here

While uploads are ongoing, please keep this browser tab open to ensure that the uploads succeed. The upload speed is dependent on your internet connection. If you want to upload a large number of files, or continuously synchronize a folder on your computer, we recommend using the connector instead.

Download View Load to Manage Export

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The screenshot displays the Illumina Connected Analytics web interface. The left sidebar shows a navigation menu with 'Data' selected. The main panel shows a list of analysis output files for the project 'Training03'. An orange box with the text '解析出力フォルダを開いて、解析出力ファイルを確認' (Open the output folder and check the output files) is overlaid on the file list.

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The screenshot shows the 'Analysis run information' page in the Illumina software interface. The page is divided into several tabs: 'Details', 'Output files', 'Steps', and 'Execution report'. The 'Details' tab is currently selected. The page is further divided into two main sections: 'General' and 'Settings'.

General Section:

- User reference:** testRnaRun002
- ID:** 3064cb8c-eb95-4275-9cb8-63c4e169de85 (Annotated with '解析ID (Analysis ID)')
- Project:** iw-temp-project, US , North Virginia , United States - global
- Pipeline:** DRAGEN RNA 4-2-4-v2
- Status:** Failed (Annotated with 'Show failed step(s)')
- Progress:** 100%
- Error:** EXECUTOR_ERROR {"error_cause": "Task execution failed with \".ica/user/bpe1-3064cb8c-eb95-4275-9cb8- ... (Annotated with 'Show more')
- Duration:** 21m 17s
- Logs:** ica_logs (Annotated with 'Data ページ内のログファイルへのリンク')
- Request date:** Oct, 31 2025 05:48:15
- Requestor:** Ching Aeng Lim (at 'global-support-training')

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
GWAS analysis

Builder for Python and R - 2.0.0

Demo of GWAS analysis on subset of 1000 genomes data using HAIL

Medium

5 GB

Stopping

ICA Benchの使い方

The screenshot displays the ICA Bench web interface. On the left is a sidebar with a menu including 'Connected Analytics', 'Data', 'Samples', 'Activity', 'Flow', 'Reference data', 'Pipelines', 'Analyses', 'Base', 'Bench', 'Cohorts', 'Project set', 'Details', 'Team', 'Connectivity', and 'Notifications'. The 'Workspaces' item under the 'Bench' category is highlighted with an orange box. An orange arrow points from this box to the 'Workspaces' header in the main content area, which also has an orange box around its '+ Create' button. Below this, the 'Create workspace' form is shown. It includes a 'General' tab with a warning message: 'Once the workspace is started, it will be restarted every 30 days. An automatic shutdown exceeding these 30 days will not prevent from restarting the workspace at the planned interval.' The form contains several input fields: 'Name *' (filled with 'test-for-webinar'), 'Automatic restart reminder' (filled with '5h'), 'Automatic shutdown reminder' (empty), 'Storage size *' (filled with '50'), 'Description' (empty), 'Resource model *' (filled with 'standard-small'), and 'Storage size *' (filled with '50'). The 'Save and run' and 'Save' buttons at the top right of the form are also highlighted with an orange box.

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[Details](#) [Docker builds](#) [Activity](#)

General

Once the workspace is started, it will be restarted every 30 days. An automatic shutdown exceeds

Status	Starting
Automatic restart reminder	5h The time (days/hours) prior to an automatic restart at which an email reminder for this event is sent out. For example: 1d 2h
Automatic shutdown	10h The time (days/hours) between start of the workspace and automatic shutdown. For example: 5d 12h
Automatic shutdown reminder	The time (days/hours) prior to the automatic shutdown at which an email reminder for this event is sent out. For example: 1d 2h

Notebook

Python 3 (ipykernel) R

Console

Python 3 (ipykernel) R

Other

Terminal Text File Markdown File Python File R File Show Contextual Help

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

Unlock deeper biological insights

Proteomics

- Illumina Protein Prep
- SomaLogic Proteomics

Epigenomics

- Illumina 5-Base DNA Prep
- Illumina methylation array
- scATAC-Seq

Genomics

- VCF

Bulk Transcriptomics

- Illumina Stranded Total and mRNA Prep
- Illumina miRNA Prep
- Data in Salmon or count matrix formats

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




Illumina Single Cell 3' RNA Prep (Custom)

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
ABOUT

 **Interactive Demos**

Explore multiomics data with an interactive click-through demo of Connected Multiomics. Click  to expand the selection for easier viewing.

The screenshot displays the Illumina Connected Multiomics web application. The top navigation bar includes 'Studies', 'Protein - Cancer/CVD Demo', 'Analyses', and 'Protein - Cancer/CVD'. The main content area shows a 'Spliceosome' pathway diagram with various protein complexes and RNA molecules. An orange overlay box with white text reads: 'Thanks for touring Illumina Connected Multiomics! Ready to try it yourself? Request a trial'.

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