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[Bookmark this app](#) [Help](#)

Pricing

Compute cost ..... 4.80 iCredits  
per node hour

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RNA-Seq Alignment  
Illumina, Inc.

[Bookmark this app](#) [Help](#)

The RNA-Seq Alignment workflow performs the following main functions

- Read mapping using the STAR aligner
- Quantification of reference genes and transcripts using salmon

Pricing  
Compute cost ..... 3.00 iCredits per node hour

Version  
2.0.2

[LAUNCH APPLICATION](#)

## What's New

- **This app will be obsoleted by June 30, 2021. Illumina recommends usage of the DRAGEN RNA Pipeline app, whose outputs can be used with DRAGEN Differential Expression for differential expression analysis.**

Please visit these pages for more information:

- <https://basespace.illumina.com/apps/10500490/DRAGEN-RNA-Pipeline>
- <https://basespace.illumina.com/apps/10502492/DRAGEN-Differential-Expression>

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DRAGEN v3.7: Single Cell RNA, PrecisionFDA Accuracy Gains, and More

<https://blog.software.illumina.com/2020/11/10/dragen-v3-7-single-cell-rna-precisionfda-accuracy-gains-and-more/>

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Upload Files

Projects: 2021-FEB Test Project for webinar

### FASTQ Dataset

Default project  
2021-FEB Test Project for webinar  
(Not editable)

Save upload to\*  
**SELECT BIOSAMPLE**

Library name\*  
(No special characters, 60 max)

Library prep kit\*  
SELECT

FILE NAME	FILE TYPE	SIZE	UPLOAD STATUS
ForTest01_S1_L001_R1_001.fastq.gz	FASTQ	1.39 GB	
ForTest01_S1_L001_R2_001.fastq.gz	FASTQ	1.39 GB	

Drop more files or select files to add to this session

2 Files Total (16 max) | 3 GB Total (250 GB max) | [Additional upload guidance](#) | [Illumina FASTQ Guidance](#)

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Upload Files

Projects: 2021-FEB Test Project for webinar

---

FASTQ Dataset

Default project  
2021-FEB Test Project for  
(Not editable)

Save upload to \*  
**SELECT BIOSAMPLE**

Library name \*  
(No special characters, 60 max)

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Upload Files  
Projects: 2021-FEB Test Project for webinar

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## File Upload Error

- Error: unable to extract metadata from reference file. Please ensure the FASTQ file adheres to the [Illumina FASTQ Guidance](#).

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× **Normal\_01**\_S1\_L001\_R1\_001.fastq.gz

FASTQ File Upload Requirements

[https://support.illumina.com/help/BaseSpace\\_Sequence\\_Hub/Source/Informatics/BS/UploadFastqReq\\_swBS.htm](https://support.illumina.com/help/BaseSpace_Sequence_Hub/Source/Informatics/BS/UploadFastqReq_swBS.htm)

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BasespaceDRAGEN RNA Pipelineで予め用意のあるHuman Reference

- Human (UCSC hg19 Alt-Aware)
- Human (UCSC hg19 No Alts, No decoys)
- Human (Ensembl GRCh37)
- Human (hg38 with HLAs)
- Human (hg38 without HLAs)
- Human (hg38 No Alts, with decoys)
- Human (Ensembl hs37d5: GRCh37+decoy)

# DRAGEN Reference Builder Appの使い方

Reference

HOME RUNS PROJECTS ANALYSES BIOSAMPLES **APPS** DEMO DATA

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Bowtie  
David Sachs  
Map reads to a reference using Bowtie.

BWA Aligner  
BaseSpace Labs  
Align samples with the BWA-MEM aligner to a reference genome, including custom references created from imported FASTA ...More

DRAGEN Reference Builder  
Edico Genome Inc.  
This app accepts FASTA files, and builds the proprietary reference used by the DRAGEN apps.

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DRAGEN Reference Builder  
Edico Genome Inc.  
Bookmark this app Help

FOR RESEARCH USE ONLY

Overview  
This app accepts a FASTA file and builds the proprietary DRAGEN Reference (also referred to as a Hash Table). This version of the app creates a V8 Hash Table. This reference can be used as input by any DRAGEN 3.5.x, 3.6.x, and 3.7.x apps which accept a Custom Reference File. The DRAGEN app typically provide built-in support for hg19, hg38 (HLA and non-HLA), GRCh37, and hg3745. The Reference Builder app is meant for users who need to generate a custom, non-human or non-standard reference. [READ MORE](#)

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# DRAGEN Reference Builder App 設定 1/2

Configuration

Analysis Name ⓘ  
DRAGEN Reference Builder 01/20/2021 3:01:22

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Liftover Validation

Seed Length ⓘ

Check this, to build the RNA hash table (in addition to the DNA hash table)

**RNA** ⓘ

Include RNA Data in Reference

**CNV** ⓘ

Include CNV Data in Reference

**Methylation** ⓘ

Include Methylation Data in Reference

Advanced Settings (WARNING: for experienced users only)

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
Identification of gene fusion transcripts by transcriptome sequencing in BRCA1-mutated breast cancers and cell lines

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
The screenshot shows a search interface with the text 'DRAGEN RNA' in the search bar. Below the search bar, two application cards are visible. The first card is for 'DRAGEN RNA Pathogen Detection' by Illumina Inc. The second card is for 'DRAGEN RNA Pipeline' by Edico Genome Inc. An arrow points from the first instruction to the search bar, and another arrow points from the second instruction to the 'DRAGEN RNA Pipeline' card.

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The screenshot shows the application page for 'DRAGEN RNA Pipeline' by Edico Genome Inc. The page includes a 'FOR RESEARCH USE ONLY' warning, an 'Overview' section, and a 'Pricing' section. The 'Pricing' section shows a compute cost of 5.00 iCredits per node hour. The 'Version' section shows a dropdown menu with '3.7.5' selected. A large blue button labeled 'LAUNCH APPLICATION' is prominently displayed at the bottom right. An arrow points from the third instruction to the 'LAUNCH APPLICATION' button.

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RNA-seq library type:

- A - Auto-Detect, the default. For this value, DRAGEN examines the first reads/pairs in the dataset to automatically determine the correct library type.
- IU - Paired-end unstranded library.
- ISR - Paired-end stranded library in which read2 matches the transcript strand (eg, TruSeq RNA).
- ISF - Paired-end stranded library in which read1 matches the transcript strand.
- U - Single-end unstranded library.
- SR - Single-end stranded library in which reads are in reverse orientation to the transcript strand (eg, TruSeq RNA).
- SF - Single-end stranded library in which reads match the transcript strand.

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The screenshot shows the DRAGEN RNA Pipeline configuration interface. It includes sections for Reference, Custom Reference File, Gene Fusion, RNA Quantification, Annotation-Assisted Alignment, Custom Gene Annotation File, and Map/Align Output. Red boxes and arrows highlight specific configuration options, with corresponding numbered instructions in Japanese text boxes.

**Reference**  
Human (hg38 No Alts, with decoys)

**Custom Reference File**  
SELECT DATASET FILE(S)

**Gene Fusion**  
 Enable Gene Fusion Detection

**RNA Quantification**  
 Enable RNA Quantification

**Annotation-Assisted Alignment**  
 Enable Annotation-Assisted Alignment

**Custom Gene Annotation File**  
SELECT DATASET FILE(S)

**Map/Align Output**  
BAM

**Reference**  
Human (hg38 No Alts, with decoys)  
Human (UCSC hg19 Alt-Aware)  
Human (UCSC hg19 No Alts, No decoys)  
Human (Ensembl GRCh37)  
Human (hg38 Alt-Aware, with HLAs)  
Human (hg38 Alt-Aware, No HLAs)  
**Human (hg38 No Alts, with decoys)**  
Human (Ensembl hs37d5: GRCh37+decoy)  
Custom

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Use Only. Not for use in diagnostic procedures.

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Gene Fusion ⓘ  
 Enable Gene Fusion Detection

RNA Quantification ⓘ  
 Enable RNA Quantification

Annotation-Assisted Alignment ⓘ  
 Enable Annotation-Assisted Alignment

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Advanced Settings (General)

Advanced Settings (Gene Fusion, Quantification)

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
LAUNCH APPLICATION

Detailed description: This block shows a screenshot of the 'Advanced Settings' menu in BaseSpace. Three menu items are listed: 'Advanced Settings (General)', 'Advanced Settings (Gene Fusion, Quantification)', and 'Additional Arguments'. A red box highlights the first two items, with an arrow pointing to them from a text box labeled '7. Advanced Option' which states '基本的に変更する必要はなし' (Basically, no need to change). Below the menu items is a blue 'LAUNCH APPLICATION' button, which is circled in red. An arrow points from a text box labeled '8. Launch Applicationをクリックし、アプリをスタートさせる' (Click Launch Application and start the app) to this button.

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General Info


Name	DRAGEN RNA Pipeline Test for Webinar
Application	DRAGEN RNA Pipeline   Version: 3.7.5
Date Started	2021-01-05 23:22
Date Completed	N/A
Duration	2 minutes 40 seconds
Compute Charge	0.00 iCredits
Session Type	Multi-Node
Node Count	4 Nodes Running
Status	Running
Delivery	None

Logs [View Files](#)

Please view the Multi-Node details page to see logs for this analysis.

Detailed description: This block shows a screenshot of the 'General Info' section of the BaseSpace interface. It contains a table with the following data: Name: DRAGEN RNA Pipeline Test for Webinar; Application: DRAGEN RNA Pipeline | Version: 3.7.5; Date Started: 2021-01-05 23:22; Date Completed: N/A; Duration: 2 minutes 40 seconds; Compute Charge: 0.00 iCredits; Session Type: Multi-Node; Node Count: 4 Nodes Running; Status: Running; Delivery: None. The 'Node Count' and 'Status' rows are highlighted with a red box. Below the table, there is a 'Logs' section with a 'View Files' link and a message: 'Please view the Multi-Node details page to see logs for this analysis.'

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General Info

Name	DRAGEN RNA Pipeline Test for Webinar
Application	DRAGEN RNA Pipeline   Version: 3.7.5
Date Started	2021-01-05 23:22
Date Completed	2021-01-05 23:39
Duration	16 minutes 45 seconds
Compute Charge	5.00 iCredits
Session Type	Multi-Node
Node Count	4 Nodes Complete
Size	31.78 GB
Status	Complete
Delivery	None

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The screenshot displays the DRAGEN RNA Pipeline web interface. At the top, there are navigation tabs: SUMMARY, REPORTS, INPUTS, and FILES. Below these are five circular icons representing different actions: a globe, a document, a timer, a share icon, and a trash can. The main content area is divided into two columns. The left column shows a 'Summary' section with a list of sample entries. The first entry, 'BrainmRNA20131201\_5...', is circled in red. Below it are other entries: 'BrainmRNA20131202\_5...', 'UHRmRNA20131205\_59...', and 'UHRmRNA20131206\_59...'. Each entry has links for 'Report' and 'Metrics'. The right column shows a 'Sample Report' for the selected sample: 'Sample: BrainmRNA20131201\_5900916 (ID=244607401)'. It includes version information: 'DRAGEN RNA Pipeline 3.7.5' and 'DRAGEN Host Software Version 05.021.595.3.7.5 and Bio-IT Processor Version 0x04261818'. Below this is a 'Multi-QC Report' section with a link to 'Open in new window' and a 'MultiQC v1.9.dev0' logo. A large red box with Japanese text is overlaid on the right side, pointing to the 'Metrics' link of the first sample entry.

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Sample: BrainmRNA20131201\_5900916 (ID=244607401)

DRAGEN RNA Pipeline 3.7.5

DRAGEN Host Software Version 05.021.595.3.7.5 and Bio-IT Processor Version 0x04261818

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**MultiQC**

A modular tool to aggregate results from bioinformatics analyses across many samples into a single report.

Report generated on 2021-01-05, 14:37 based on data in: `/data/output/appresults/220174962/BrainmRNA20131201_5900916`

📘 Welcome! Not sure where to start? [Watch a tutorial video \(6:06\)](#)

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## General Statistics

Copy table    Configure Columns    Plot    Showing 1/1 rows and 7/24 columns.

Sample Name	M Input reads	Unmap	Dup	Prop pair
BrainmRNA20131201	128.3	2.0%	35.1%	97.3%

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DRAGEN Differential Expression  
BaseSpace Labs  
Bookmark this app Help

FOR RESEARCH USE ONLY

Overview  
The Illumina® DRAGEN RNA Differential Expression application performs Next Generation Sequencing (NGS) secondary analysis of RNA transcripts. It runs the DESeq2 algorithm on Salmon quantification files to output genes and transcripts that are differentially expressed between two sample groups.

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Group Label

UHR

Control Salmon quantification files

**SELECT DATASET FILE(S)**

BrainmRNA20131201.quant.sf

BrainmRNA20131202.quant.sf

Comparison Salmon quantification files

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The screenshot shows the configuration interface for the DRAGEN Differential Expression application. It includes a dropdown menu for 'Reference Annotation' set to 'Human (hg38)', a section for 'Custom Gene Annotation File' with a 'SELECT DATASET FILE(S)' link, and a 'BaseSpace Labs Disclaimer' section with a checked checkbox. At the bottom, there is a blue 'LAUNCH APPLICATION' button. Three red-bordered boxes with arrows point to these elements, containing the following instructions:

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


# DRAGEN Differential Expressionの結果の見方

## Differential Expression Metrics

Annotation Gene Count	60609	<a href="#">Download</a>
Expressed Gene Count	22997	<a href="#">Download</a>
Differentially Expressed Gene Count	17674	<a href="#">Download</a>
Merged Gene Counts		<a href="#">Download</a>
Differentially Expressed Genes		<a href="#">Download</a>
Differentially Expressed Transcripts		<a href="#">Download</a>

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# DRAGEN Differential Expressionの結果の見方

## Differential Expression Metrics

Annotation Gene Count	60609
Expressed Gene Count	22997
Differentially Expressed Gene Count	17674
Merged Gene Counts	<a href="#">Download</a>
Differentially Expressed Genes	<a href="#">Download</a>
Differentially Expressed Transcripts	<a href="#">Download</a>

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
# Appendix 資料

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BaseSpace Sequence Hub Help

[https://support.illumina.com/help/BaseSpace\\_Sequence\\_Hub/Source/HomePages/Home\\_Page\\_BaseSpace\\_Sequence\\_Hub.htm](https://support.illumina.com/help/BaseSpace_Sequence_Hub/Source/HomePages/Home_Page_BaseSpace_Sequence_Hub.htm)

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