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Local Run Manager **Instrument Name** RUN DASHBOARD TOOLS Username | illumina

43 Ready 0 In Progress 1 Stopped or Unsuccessful 1 Complete 45 Total + Create Run

RUN NAME / ID	MODULE	STATUS	LAST MODIFIED	
Training2 --	DNA AMPLICON	Ready for Sequencing	2018-05-29 08:46	Actions
training --	GENERATEFASTQ	Ready for Sequencing	2018-05-29 07:39	Actions

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
---

# Local Run Manager (LRM) v3の新機能:

- Library Prep KitsとIndex Kitsの区別

The screenshot shows a 'Run Settings' section with two dropdown menus. The first is labeled 'Library Prep Kit \*' and is set to 'AmpliSeq Library PLUS for Illumina'. The second is labeled 'Index Kit \*' and is set to 'AmpliSeq UD Indexes'.

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## UPDATE SECURITY QUESTIONS

Your security questions will be used to verify your identity in case you forget your password.

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The screenshot shows the 'Local Run Manager' interface. At the top, there are navigation tabs for 'RUN DASHBOARD' and 'TOOLS'. The main area features a dashboard with five status indicators: 'Ready' (0), 'In Progress' (0), 'Stopped or Unsuccessful' (2), 'Complete' (1), and 'Total' (3). A 'Create Run' button is highlighted in a green box, with a dropdown menu showing 'GenerateFASTQ' and 'Resequencing'. Below the dashboard is a table with columns for 'RUN NAME / ID', 'MODULE', 'STATUS', and 'LAST MODIFIED'. The table contains three rows of data, each with an 'Actions' button.

RUN NAME / ID	MODULE	STATUS	LAST MODIFIED
GenerateFastQRun 201029_M99999_0010_FC1234567-ABBBB	GENERATEFASTQ	Sequencing Stopped	2020-10-29 09:22
Test 201009_M99999_0007_FC1234567-ABBBB	GENERATEFASTQ	Sequencing Stopped	2020-10-09 01:50
Demo1 200721_M99999_0005_FC1234567-ABBBB	GENERATEFASTQ	Analysis Completed	2020-10-09 01:24

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[Show Advanced Setting](#)

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[Show Index Sequence](#)

	SAMPLE ID*	SAMPLE DESCRIPTION	INDEX 1 (I7)*	INDEX 2 (I5)*	SAMPLE PROJECT
1					

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The screenshot shows the 'Create Run' interface with the 'Index Kit' dropdown menu open, displaying a list of compatible index kits. The 'Index Reads' table is also visible, showing the number of reads for each index type.

READ 1	INDEX 1	INDEX 2	READ 2
151	0	0	151

Index Reads table details:

- Read 1:  Read 1
- Index 1:  Index
- Index 2:  Index
- Read 2:  Read 2

Index Kit options:

- Ampliseq Library PLUS for Illumina
- Nextera DNA Flex
- Nextera DNA Flex PCR-Free
- Nextera Flex for Enrichment
- Nextera Mate Pair
- Nextera Rapid Capture Enrichment
- Nextera XT
- ScriptSeq Complete
- ScriptSeq v2
- SureCell Single Cell RNA
- SureCell WTA 3'
- TruSeq Amplicon Cancer Panel

Index Reads table:

Index Reads
151

Read Type:  Single Read  Paired End

Library Prep Kit: TruSeq DNA PCR-Free

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
The screenshot shows a software interface for creating a run. At the top, it says 'Create Run GENERATEFASTQ'. Below this, there is a blue button with a white border labeled 'Import Sample Sheet'. Underneath the button is a text input field labeled 'Run Name\*' with the placeholder text 'Run Name'. At the bottom of the interface, there is a large blue button with a white border and a white floppy disk icon, labeled 'Save Run'. A large pink arrow points from the 'Import Sample Sheet' button down to the 'Save Run' button.

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Test  
GENERATEFASTQ  
Ready for Sequencing  
2018-05-10  
04:23  
Actions

1  
**Actions**  
を選択  
Actions

2  
**Edit**  
を選択  
Edit  
Hide  
Import

3  
**Continue**  
を選択

EDIT RUN CONFIRMATION  
Are you sure you would like to edit the run 'GenerateFastQRun'?  
Cancel Continue

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
## Index & Library Prep Kits

LIBRARY PREP KIT	INDEX KIT		
		<a href="#">Download Template</a>	<a href="#">Add Library Prep Kit</a>
NAME	DESCRIPTION	SUPPORTED MODULES	VERSI...
Illumina DNA Prep with Enrichment	Illumina DNA Prep with Enrichment Library Prep Kit	GenerateFASTQ, Library QC	1
Illumina DNA Prep	Illumina DNA Prep Library Prep Kit	GenerateFASTQ, Library QC, Resequencing	1
Illumina DNA PCR-Free 2x151	Illumina DNA PCR-Free 2x151 Library Prep Kit	GenerateFASTQ, Library QC	1
Illumina DNA PCR-Free 2x149	Illumina DNA PCR-Free 2x149 Library Prep Kit	GenerateFASTQ, Library QC	1


```
<remove anything inside the chevrons prior to upload>

<please reference the user guide for more template information>

[LibraryPrepKit]
Name <Enter the name of the Prep Kit (Must be unique)>
DisplayName <Enter the display name of the prep kit>
Alias <Comma-separated list of prep kit names to overwritten by this prep kit>
Version <Enter Prep Kit version>
Description <Enter Description of the Prep Kit Here>
ReadType <Options are: All, Single, and Paired>
DefaultReadLength1 <Enter Default Read Length for Read 1>
DefaultReadLength2 <Enter Default Read Length for Read 2>

[Resources]
Name Type Format Value IndexKits
<Custom parameters can be set here. This section is optional. If resource applies to all IndexKits, leave IndexKits>
Read1UMILength UMISetting string 7
Read2UMILength UMISetting string 7
Read1StartFromCycle UMISetting string 9
Read2StartFromCycle UMISetting string 9
showReverseComplement kitSettings bool true

[SupportedModules]
<Select Modules for the Prep Kit to be associated with. Options available (multiple allowed): >
AmpliconDSWorkflow
AssemblyWorkflow
DNAAmpliconWorkflow
DNAEnrichmentWorkflow
GenerateFastQWorkflow
LibraryQCWorkflow
```

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

Import Manifests Show Index Sequence

	NAME	INDEX 1 (I7)	INDEX 2 (I5)	MANIFEST*	SAMPLE PROJECT	
1	Sample1	N701	E501	Example_manifest.txt		x

+ 1 Rows

Cancel Export Sample Sheet Save Run

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
Local Run Manager >

Select a run prepared with Local Run Manager software.

Sample Sheet >

Create a run from a sample sheet.

Manual >

Enter run information manually. No on-instrument secondary analysis is performed.

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The screenshot shows the BaseSpace interface for a run. At the top, there is a blue header with a menu icon and the text 'M99999'. Below the header is a navigation bar with 'BaseSpace' and several menu items: 'BaseSpace', 'Run Settings', 'Load Flow Cell', 'Load Reagents', 'Review', 'Pre-Run Check', 'Sequence', and 'Post-Run Wash'. A green checkmark icon is visible in the top right corner. The main content area features a cloud icon with a padlock, followed by the heading 'Run Monitoring, Analysis, and Collaboration'. Below this, there is a paragraph: 'Monitor runs, perform analysis, and collaborate anywhere with colleagues by using **multilayered security and encryption** to send run data to BaseSpace Sequence Hub. Illumina Proactive Support is automatically included when using this feature.' There are two radio button options: 'Use BaseSpace™ Sequence Hub for this run' (which is selected) and 'Run Monitoring from anywhere'. At the bottom, there are three buttons: '< Back', '× Exit', and 'Next >'.

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M99999

Run Setup

BaseSpace Run Settings Load Flow Cell Load Reagents Review Pre-Run Check Sequence Post-Run Wash

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
PhiX Verification Run - 20pm GENERATEFASTQ

**RUN OVERVIEW** | SEQUENCING INFORMATION | SAMPLES & RESULTS

Run Name	PhiX Verification Run - 20pm
Run ID	
Created By	System User
Description	Training Cert
Output Run Folder	<input type="text"/> <span>Copy</span> <span>Edit</span>

Analysis Completed  
Successfully completed in 2.38 minutes

Total Clusters	29,108,048
% Clusters PF	90.2%
% ≥ Q30 (Read 1)	96.6%
% ≥ Q30 (Read 2)	94.9%
Last Scored Cycle	150

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PhiX Verification Run - 20pm GENERATEFASTQ

**RUN OVERVIEW** | **SEQUENCING INFORMATION** | SAMPLES & RESULTS

Instrument Name: MISEQ  
Sequenced By: System User  
Sequencing Start: 2020-10-21 15:45

Chemistry Type: Default  
Chemistry Version: 3  
Module Version: 3.0.0  
(Created With)

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




RUN NAME / ID	MODULE	STATUS	LAST MODIFIED	
Enrichment Example	DNA ENRICHMENT	Ready for Sequencing	2020-12-16 13:49	Actions

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IMPORT RUN

Provide run folder path by copying or typing it into field below.

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[View All](#)

**Product Information**

- Technical Bulletins
- Specifications
- Proactive Instrument Monitoring
- Services and Warranties
- Regulatory and Quality Information

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- [https://jp.support.illumina.com/sequencing/sequencing\\_software/local-run-manager.html](https://jp.support.illumina.com/sequencing/sequencing_software/local-run-manager.html)

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**Local Run Manager**  
Support Resources

Support Center:  
Local Run Manager Software Support

This page supports on- and off-instrument versions of Local Run Manager.

Documentation	Software	Product Information
<a href="#">Local Run Manager v3 Software Guide</a>	<a href="#">Local Run Manager Generate FASTQ Analysis Module v3.0</a>	<a href="#">Technical Bulletins</a>
<a href="#">Generate FASTQ Analysis Module v3.0 Guide</a>	<a href="#">Local Run Manager Library QC Analysis Module v2.1</a>	<a href="#">Specifications</a>
<a href="#">Library QC Analysis Module v2.1 Guide</a>	<a href="#">Local Run Manager Resequencing Analysis Module v3.0</a>	
<a href="#">Resequencing Analysis Module v3.0 Guide</a>		

[View All](#)

**Local Run Manager Module Selector**

Select the version of the Local Run Manager module for your instrument and control software.

Instrument
<a href="#">iSeq 100</a>
<a href="#">MiniSeq</a>
<a href="#">MiSeq</a>
<a href="#">MiSeq Dx</a>
<a href="#">NextSeq 500/550</a>
<a href="#">NextSeq Dx</a>
<a href="#">Off-instrument</a>

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Start Run with MCS

Create run Prepare libraries Prepare reagent cartridge and load libraries Start run with MCS

MCS v4 and later Local Run Manager v3 and later

MCS v3.1 Local Run Manager v2

MCS v2.6 and earlier MiSeq Reporter

Select your version of MCS to view a demo of the steps.

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
- MiSeq how to start a run

- [https://support.illumina.com/content/dam/illumina-support/courses/MiSeq\\_How\\_to\\_Start\\_a\\_Run/story.html](https://support.illumina.com/content/dam/illumina-support/courses/MiSeq_How_to_Start_a_Run/story.html)

- Local Run Manager v3

- <https://support.illumina.com/content/dam/illumina-support/courses/local-run-manager-v3/story.html>

MENU TRANSCRIPT

- 1. Welcome
  - 1.1. Welcome
  - 1.2. Course Objectives
- 2. Overview
- 3. Create Run
- 4. Prepare Libraries
- 5. Prepare Reagent Cartridge and Load Libraries
- 6. Start Run with MCS
- 7. Summary

MiSeq: How to Start a Run


**MiSeq: How to Start a Run**

🕒 20 min 🔊 Contains audio

[Start Course](#)

illumina®


Document # 1000000149878 v00  
For Research Use Only. Not for use in diagnostic procedures.

MENU TRANSCRIPT

- 1. Welcome
  - 1.1. Local Run Manager v3
  - 1.2. Course Objectives
- 2. Overview
- 3. Creating and Monitoring Runs
- 4. Viewing Analysis Results
- 5. Managing Input Files
- 6. Summary

Local Run Manager v3

**Local Run Manager v3**

🕒 20 min 🔊 Contains audio

[Start Course](#)

illumina®

Document # 1000000149880 v00  
For Research Use Only. Not for use in diagnostic procedures.

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