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BaseSpace Help center : <https://help.basespace.illumina.com/>

The screenshot shows the BaseSpace Sequence Hub interface. At the top, there is a navigation bar with the 'illumina' logo on the left and 'SEQUENCE HUB' on the right. Below the logo is a search icon. The navigation bar contains several tabs: 'HOME', 'RUNS', 'PROJECTS', 'ANALYSES', 'BIOSAMPLES', 'APPS', and 'DEMO DATA'. On the right side of the navigation bar, there are three icons: a notification bell, a help icon (a question mark inside a circle), and a trash can icon. A red box highlights the help icon, with a red callout box pointing to it that says 'ここからもBSSHのヘルプにアクセスできます' (You can also access BSSH help from here). A grey callout box points to the 'HOME' tab, saying '各タブを選択してBSSH上の各情報、解析App、デモデータにアクセスできます' (You can access various information, analysis apps, and demo data on BSSH by selecting each tab). Another grey callout box points to the notification bell icon, saying '通知を閲覧できます' (You can view notifications). A grey callout box points to the user profile icon, saying 'アカウント設定' (Account settings). The main content area of the interface shows a 'Hello' message and a status message: 'Everything is looking good so far today. Check back later for any incoming items.'

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

The image shows a screenshot of the Illumina BSSH registration page. The page has a light blue background with the Illumina logo at the top. Below the logo, there is a 'Register' section with a link for 'Already have an account? Sign in'. The form includes fields for 'Email address', 'First name', 'Last name', 'Password' (with a 'Guidelines' link), and 'Confirm password'. There is also a 'Location' dropdown menu currently set to 'United States'. At the bottom of the form, there is a checkbox for 'I'm not a robot' and a reCAPTCHA logo. Two buttons, 'Create Account' and 'Cancel', are located at the bottom of the form.

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16S Metagenomics
Illumina, Inc.
Bookmark this app Help

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Home > [Runs](#) >

illumina SEQUENCE HUB

HOME **RUNS** PROJECTS ANALYSES BIOSAMPLES APPS DEMO DATA

Runs

ACTIVE PLANNED [NEW RUN](#)

There are no urgent actions. Well done. Return in a few hours to check incoming runs.

STATUS	RUN NAME	AVG%Q30	%PF	INSTRUMENT	CREATED
<input type="checkbox"/> Complete	[Redacted]	93.70%	84.82%	[Redacted]	2022/6/
<input type="checkbox"/> Complete	[Redacted]	92.64%	75.95%	[Redacted]	2022/6/
<input type="checkbox"/> Complete	[Redacted]	91.86%	91.49%	[Redacted]	2022/3/

1 - 7 of 7 S

Project

- Biosample
 - Library
 - FASTQ Dataset
 - Lane
 - Run
- Analysis
- Analysis

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
Instrument iSeq100-simulat...	92.64 AVG %Q30	75.95 %PF	Created 2022-06-10 11...	Instrument Type iSeq100
Run Status Complete	Lane QC Status QcPassed	Flow Cell Status QcPassed	File Count/Size 1,006 files (86...	File Status: Active
Latest Analysis FASTQ Generati...	Cycles 151 8 8 151	Yield 1.83 Gbp	Owner	User
			Flow Cell ID	Run ID --

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
Home > Projects > [Project name] > Analyses

illumina SEQUENCE HUB

HOME RUNS **PROJECTS** ANALYSES BIOSAMPLES APPS DEMO DATA

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Showing 4 of 4


Project hierarchy diagram:

- Project
 - Biosample
 - Library
 - FASTQ Dataset
 - Lane
 - Run
 - Analysis
 - Analysis

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Home > Biosamples > [Biosample name]

illumina SEQUENCE HUB

HOME RUNS PROJECTS ANALYSES BIOSAMPLES APPS DEMO DATA

Biosample: [Name]


SUMMARY LIBRARIES RUNS REQUEUES ANALYSES FASTQS OTHER DATASETS

STATUS	RUN NAME	FLOWCELL ID	ANALYSIS	UPDATED
Complete	[Name]	[ID]	Generate FAS...	2022-06-10

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

illumina SEQUENCE HUB

HOME RUNS PROJECTS ANALYSES BIOSAMPLES **APPS** DEMO DATA

Search Apps

ALL CATEGORIES

All Sequence Hub Applications

- 16S Metagenomics (Illumina, Inc.)
- 16S Metagenomics (BaseSpace Labs)
- 16S Metagenomics (Illumina, Inc.)
- Amplicon DS (Illumina, Inc.)
- BCL Convert (Illumina, Inc.)
- BCL Convert (Early Ac... (Illumina, Inc.)
- BCL Convert for ICA (B... (Illumina, Inc.)
- BWA Aligner (BaseSpace Labs)
- BWA Enrichment (Illumina, Inc.)
- Cancer Variant Caller (Samsung SDS)
- ChIPSeq (BaseSpace Labs)
- CosmosID Metagenomi... (CosmosID)

Bookmarked Apps
DRAGEN RNA Pipeline

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The screenshot displays the Illumina Sequence Hub interface. At the top, the 'illumina' logo and 'SEQUENCE HUB' are visible. A navigation bar includes 'HOME', 'RUNS', 'PROJECTS', 'ANALYSES', 'BIOSAMPLES', 'APPS', and 'DEMO DATA'. A search bar contains the text 'Metagenomics'. A red box highlights this search bar with the text 'Apps名の一部を入力して検索' (Enter part of the app name to search). Below the search bar, a dropdown menu shows 'ALL CATEGORIES'. The search results are displayed in a grid of six app cards. A red box highlights the first card, '16S Metagenomics' by Illumina, Inc., with the text '選択' (Select) and an arrow pointing to it. The app cards include icons, titles, descriptions, and category lists.

App Name	Developer	Description	Categories
16S Metagenomics	Illumina, Inc.	Analyzes DNA from amplicon sequencing of prokaryotic 16S small subunit rRNA genes. Provides visualizations of taxonomic composition.	Targeted Sequencing, Metagenomics, HIPAA
16S Metagenomics Database Creator	BaseSpace Labs	The 16S Metagenomics Database Creator app creates custom databases for the 16S Metagenomics app (v1.1.0 and above).	Targeted Sequencing, Metagenomics, HIPAA
16S Metagenomics Labs	Illumina, Inc.	Analyzes DNA from amplicon sequencing of prokaryotic 16S small subunit rRNA genes. Provides visualizations of taxonomic composition.	Targeted Sequencing, Metagenomics, HIPAA
CosmosID Metagenomics	CosmosID	A cloud-based software enables fast, accurate, and comprehensive analyses of unassembled next-generation sequencing data.	Metagenomics
DRAGEN Metagenomics Pipeline	Illumina Inc.	The DRAGEN Metagenomics pipeline performs taxonomic classification of reads and provides single sample and aggregate analysis.	Metagenomics

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The screenshot shows the Illumina Sequence Hub interface. At the top, there is a navigation bar with the Illumina logo, the text 'SEQUENCE HUB', and a user profile icon. Below this is a secondary navigation bar with tabs for 'HOME', 'RUNS', 'PROJECTS', 'ANALYSES', 'BIOSAMPLES', 'APPS', and 'DEMO DATA'. A search bar is located on the left side of the main content area. The main content area displays a grid of application cards under the heading 'All Sequence Hub Applications'. A dropdown menu is open on the right side of the page, showing a list of categories with their respective counts. The 'Metagenomics (15)' category is highlighted with a red box. A red arrow points from a text box on the right to the dropdown menu.

Search Apps

ALL CATEGORIES

- All Categories
- ChIP-Seq (4)
- De Novo Assembly (8)
- Differential Expression (18)
- Gene Fusion Detection (11)
- HiPAA (36)
- Metagenomics (15)
- Methyl-Seq (6)
- Proteomics (10)
- Quality (12)
- Resequencing (37)
- RNA-Seq (31)
- Small RNA (9)
- Synthetic Long Reads (1)
- Targeted Sequencing (40)
- Tumor Normal (9)
- Variant Analysis (32)

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The screenshot shows the Illumina Sequence Hub interface. At the top, there is a navigation bar with the Illumina logo and the text 'SEQUENCE HUB'. Below this is a search bar with the text 'Search Apps' and a search icon. To the right of the search bar, there is a filter dropdown menu showing 'METAGENOMICS (15)'. Below the search bar, there is a section titled '15 Results Found'. The first result is highlighted with a red box and a red arrow pointing to it from the left. This result is the '16S Metagenomics' app by Illumina, Inc. The app description reads: 'Analyzes DNA from amplicon sequencing of prokaryotic 16S sm all subunit rRNA genes. Provides visuals of taxonomic ...More'. The categories listed are 'Targeted Sequencing, Metagenomics, HIPAA'. Other visible results include '16S Metagenomics Database Creator' by BaseSpace Labs, '16S Metagenomics Labs' by Illumina, Inc., 'CosmosID Metagenomics' by CosmosID, 'DRAGEN Metagenomics Pipeline' by Illumina Inc., and 'Explify RPIP Analysis' by IDbyDNA.

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What's New
Users may find a lower % of species classified using some newer taxonomy databases such as the RDP-based database included with the v1.1.0 app vs-using the older GreenGenes based taxonomy database, but the RDP-based classifications will have less false positives. Species level classification from 16S amplicon data is generally less reliable vs genus-level classification.

Screenshots

Pricing
Compute cost 3.00 iCredits per node hour

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The 16S Metagenomics app performs taxonomic classification of 16S rRNA targeted amplicon reads using a taxonomic database. The app provides interactive visualizations and raw classification output for per-sample and aggregate analyses.

Classification is performed using the Illumina 16S Metagenomics workflow, which is also available in the MiSeq Reporter software.

The algorithm is a high-performance implementation of the Ribosomal Database Project (RDP) Classifier described in Wang, Q. et al (<http://dx.doi.org/10.1128%2FAEM.00062-07>).

1.1.0

LAUNCH APPLICATION

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Current limitations:

- Read lengths of at least 100bp
- Each dataset less than 50 Gigabases

Demo data showing analysis results from the app is here:

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The Illumina-curated versions of the Greengenes database are available here: [16S Metagenomics Database](#).
[Creator Demo](#).

The RefSeq RDP 16S v3 database is based on FASTA from: <https://benjjneb.github.io/dada2/training.html>. Citation: Ali Alishum. (2019). DADA2 formatted 16S rRNA gene sequences for both bacteria & archaea (Version Version 2)

The UNITE Fungal ITS Database v7.2 is based on FASTA from: UNITE Community (2017): UNITE general FASTA release. Version 01.12.2017. UNITE Community. <https://doi.org/10.15156/BIO/587475> Includes singletons set as RefS (in dynamic files).

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The screenshot shows a software application page with several sections and callouts:

- Additional Information**
 - Release Notes
 - Version 1.1.0
 - New Features** (highlighted with a red box and callout):
 - Added RefSeq RDP 16S v3 May 2018 reference taxonomy database using data from <https://benjjneb.github.io/dada2/training.html>.
 - Added UNITE Fungal ITS Database v7.2 reference taxonomy database.
 - Added support for custom reference taxonomy database.
 - Defect Repairs: N/A
 - Know Issues: N/A
 - [READ LESS](#)
- Developer Information**
 - Published: 2019-10-19 08:01
 - Developer: [Visit Website](#) (highlighted with a green box and callout: "Developer情報のリンク (Help)")
 - Terms
 - EULA
 - Privacy
- Categories** (highlighted with a blue box and callout: "Appに関連付けられたカテゴリー"):
 - Targeted Sequencing
 - Metagenomics
 - HIPAA
- Version and Action** (highlighted with a white box):
 - 1.1.0 (with a green dot)
 - [LAUNCH APPLICATION](#) (highlighted with a blue box and callout: "Appに関連付けられたカテゴリー")

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Biosample(s) ⓘ

SELECT BIOSAMPLE(S):

Unknown	x
Unknown	x
Unknown	x
Unknown	x
Unknown	x

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illumina SEQUENCE HUB

HOME RUNS PROJECTS **ANALYSES** BIOSAMPLES APPS DEMO DATA

Analysis: 16S Metagenomics 08/30/2022 11:05:49

Project

SUMMARY REPORTS INPUTS FILES

General Info

Name	16S Metagenomics 08/30/2022 11:05:49
Application	16S Metagenomics Version: 1.1.0
Date Started	2022-08-30 11:09
Date Completed	N/A
Duration	N/A
Compute Charge	0.00 iCredits
Session Type	N/A
Status	Initializing
Delivery	None

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Home > Analyses > [Analysis name] > Reports

HOME RUNS PROJECTS ANALYSES

Analysis: 16S Metagenomics 07/26/2

Project MiSeq: 16S Metagenomics 2x300 (ATCC Microbiome Standard MSA-1003)

SUMMARY **REPORTS** INPUTS FILES

Aggregate Sum... Report

ATCCMiSeq2x3... Report

ATCCMiSeq2x3... Report

ATCCMiSeq2x3... Report

ATCCMiSeq2x3... Report

ATCCMiSeq2x3... Report

ATCCMiSeq2x3... Report



AGGREGATE RESULTS

- PDF Aggregate Summary Report
- Kingdom-Level Aggregate Results
- Phylum-Level Aggregate Results
- Class-Level Aggregate Results
- Order-Level Aggregate Results
- Family-Level Aggregate Results
- Genus-Level Aggregate Results
- Species-Level Aggregate Results

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Home > Analyses > [Analysis name] > Reports

Analysis: 16S Metagenomics 08/30/2022 11:05:49
Project 16S metagenomics demo

SUMMARY **REPORTS** INPUTS FILES

Aggregate Su... Report
ATCCMiSeq2x... Report
ATCCMiSeq2x... Report
ATCCMiSeq2x... Report
ATCCMiSeq2x... Report
ATCCMiSeq2x... Report

RESULTS FOR SAMPLE ATCCMISEQ2X300R1



PDF Summary Report
CSV Classification Summary

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
Home > Analyses > [Analysis name]

The screenshot displays the BaseSpace web interface for an analysis titled "Analysis: 16S Metagenomics 07/26/2018 12:22:19". The project is identified as "Project MiSeq: 16S Metagenomics 2x300 (ATCC Microbiome Standard MSA-1003)". The navigation tabs include SUMMARY, REPORTS, INPUTS, and FILES, with the FILES tab selected and highlighted by a red box. A red callout box with white text points to the FILES tab, stating "データの確認・個別ダウンロードをすることができます" (You can check the data and download it individually). Below the navigation tabs, there are icons for various actions: a document, a folder, a clock, a share icon, and a trash can. The main content area is divided into two sections: "AGGREGATE RESULTS" and "SAMPLE INFORMATION". The "AGGREGATE RESULTS" section lists several report types, each with a green 'X' icon indicating it is available for download: PDF Aggregate Summary Report, Kingdom-Level Aggregate Results, Phylum-Level Aggregate Results, Class-Level Aggregate Results, Order-Level Aggregate Results, Family-Level Aggregate Results, Genus-Level Aggregate Results, and Species-Level Aggregate Results. The "SAMPLE INFORMATION" section is partially visible at the bottom.

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Home > Analyses > [Analysis name] > Files

The screenshot shows the 'FILES' tab of an analysis titled 'Analysis: 16S Metagenomics 08/30/2022 11:05:49'. The interface includes a navigation bar with 'SUMMARY', 'REPORTS', 'INPUTS', and 'FILES'. Below the navigation bar are several icons. A table lists the files with columns for NAME, CREATED, TYPE, and PATH. The table contains the following entries:

NAME	CREATED	TYPE	PATH
ATCCMiSeq2x300R9_S10.b...	2022-08-30...	json	ATCCMiSeq2x300R9_S10.barct...
ATCCMiSeq2x300R9_S10.re...	2022-08-30...	html	ATCCMiSeq2x300R9_S10.repor...
ATCCMiSeq2x300R9_S10.re...	2022-08-30...	pdf	ATCCMiSeq2x300R9_S10.repor...
ATCCMiSeq2x300R9_S10.st...	2022-08-30...	xml	ATCCMiSeq2x300R9_S10.stats...

A red callout box with the text 'ダウンロードしたいファイルを選択' (Select the file you want to download) has a red arrow pointing to the PDF file row in the table. A second screenshot, overlaid on the right, shows a modal window for the selected PDF file 'ATCCMiSeq2x300R9_S10.report.pdf'. The modal displays a preview of the report, which includes the Illumina logo and the title '16S Metagenomics Report'. At the bottom of the modal, there are two buttons: 'CLOSE' and 'DOWNLOAD (332.81 KB)'. A red callout box with the text 'ダウンロード' (Download) has a red arrow pointing to the 'DOWNLOAD' button.

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Home > Analyses > [Analysis name] > 任意のタブ

The screenshot shows the 'FILES' tab of an analysis titled 'Analysis: 16S Metagenomics 08/30/2022 11:05:49'. A red box highlights the 'FILE' icon in the top navigation bar. A dropdown menu is open, with 'DOWNLOAD' highlighted by another red box. A red arrow points to the 'ANALYSIS' link in the file list. A red banner at the bottom of the screenshot contains the text: 'FILE > DOWNLOAD > ANALYSIS と進みます。' (Proceed with FILE > DOWNLOAD > ANALYSIS).

NAME	CREATED	TYPE	PATH	SIZE
ANALYSIS			TCCMiSeq2×300R9_S10.barchart.json	10.33 KB
DATASET			TCCMiSeq2×300R9_S10.report.html	595.08 KB
PROJECT			TCCMiSeq2×300R9_S10.report.pdf	332.81 KB
RUN			TCCMiSeq2×300R9_S10.stats.xml	625 B
ATCCMiSeq2×300R9_S10.s...	2022-08-30...	csv	ATCCMiSeq2×300R9_S10.summary.csv	127.12 KB

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Install the BaseSpace Sequence Hub Downloader to download files. It's a one-time installation, is required, and provides fast and secure downloads via SSL.

Select the file types to be downloaded:

- All file types including VCF, BAM, & FASTQ
- VCF
- BAM

DOWNLOAD CLOSE

NAME	SIZE
ATCCMISeq2*300R9_S10.b...	10.33 KB
ATCCMISeq2*300R9_S10.re...	595.08 KB
ATCCMISeq2*300R9_S10.rs...	332.81 KB
ATCCMISeq2*300R9_S10.st...	625 B
ATCCMISeq2*300R9_S10.s...	127.12 KB
ATCCMISeq2*300R9_S10.s...	5.80 KB
ATCCMISeq2*300R9_S10.tx...	8.31 MB

Confirm Download

Ready to download:

Directory: c:\BaseSpace

START DOWNLOAD CANCEL

BaseSpace Sequence Hub Downloader

Avg speed: 3.10 MB/s

Time remaining: 48 seconds

Data downloaded: 24.73 MB/176.38 MB

Status: 14.02% done, Total: 182 Completed: 24

16S metagenomics demo-36...

Files: 182, Progress: 14.02%, Completed: 24.73 MB/176.38 MB, TransferRate: 0.00 MB/s

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Showing 1 of 1

<input type="checkbox"/>	NAME	LAST MODIFIED	APPLICATION	SIZE	COMMENTS	DELIVERY	STATUS
<input type="checkbox"/>	FASTQ Generatio...	2022-04-12	FASTQ Generation	695 MB		None	Complete

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The screenshot displays the 'Project PhiX_220411_Training' interface. At the top, there are navigation tabs: SUMMARY, ANALYSES, BIOSAMPLES, FASTQS, and OTHER DATASETS. Below these are five circular icons: a globe, a document, a timer, a share icon, and a trash can. The share icon is highlighted, and a 'SHARE' dropdown menu is open, containing three options: 'SEND INVITATION', 'GET LINK', and 'TRANSFER'. Three red callout boxes with white text point to these options: 'データの共有 (Invitation)' points to 'SEND INVITATION', '共有リンクの有効化、無効化' points to 'GET LINK', and '所有権 (Ownership) の移譲' points to 'TRANSFER'. Below the menu, a table is partially visible with columns for NAME, SIZE, COMMENTS, DELIVERY, and STATUS. One row is shown with the name 'FASTQ Generatio...' and a status of 'Complete'.

NAME	SIZE	COMMENTS	DELIVERY	STATUS
FASTQ Generatio...				Complete

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The image illustrates the process of sharing a project in BaseSpace Informatics Suite. It consists of several overlapping screenshots:

- Project Overview:** Shows the 'Project PhiX_220411_Training' interface with tabs for 'SUMMARY', 'ANALYSES', and 'BIOSAMPLES'. A share icon is highlighted.
- Share Dialog (Left):** A 'Share this project' dialog box with the project name 'PhiX_220411_Training' and an 'Email address' field highlighted by a red box. A red callout box points to this field with the text: **共有先のBSSHアカウントのメールアドレスを入力** (Enter the email address of the BSSH account to be shared with).
- Share Dialog (Right):** A larger view of the 'Share this project' dialog, showing the 'Email address' field, an 'Optional message' field, and an 'ADD COLLABORATOR' button. A red arrow points to this button.
- Email Notification:** A screenshot of an email notification from BaseSpace Informatics Suite. The subject is 'PhiX_220411_Training' and the body says: 'would like to share the project PhiX_220411_Training with you. Go to PhiX_220411_Training to accept the share.' The footer includes the Illumina logo and address: 'This email was sent by: Illumina, Inc. 5200 Illumina Way San Diego, CA 92122 USA'.

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Share By Link

PhiX_220411_Training

Click 'Activate' to generate a new URL for sharing.

CLOSE ACTIVATE

Share By Link

PhiX_220411_Training

Copy the below URL and send it to anyone you want to share with.

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The screenshot displays the 'Project PhiX_220411_Training' interface. At the top, there are navigation tabs: SUMMARY, ANALYSES, BIOSAMPLES, FASTQS, and OTHER DATASETS. Below these are five circular icons: a globe, a document, an hourglass, a share icon, and a trash can. The share icon is highlighted, and a dropdown menu is open with the following options: SHARE, SEND INVITATION, GET LINK, and TRANSFER. The 'TRANSFER' option is highlighted with a red box. A red arrow points from this box to a red text box containing the Japanese text '所有権 (Ownership) の移譲'. Below the menu, a table is partially visible with columns for NAME, SIZE, COMMENTS, DELIVERY, and STATUS. One row is shown with the name 'FASTQ Generatio...', a status of 'Complete', and a 'Complete' link.

NAME	SIZE	COMMENTS	DELIVERY	STATUS
FASTQ Generatio...				Complete

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Project Ph

SUMMARY ANALY

Showing 1 of 1

NAME

FASTQ Generati

Transfer Ownership
PhiX_220411_Training

New Owner

Email address

Optional Message =

Add a personal message (optional)

(max 140 characters)

ⓘ Apps in your account that currently use this data may not function properly after the data is transferred out of your account. Please review the [BaseSpace User Terms of Use](#) for more information.

CANCEL CONTINUE



Transfer Ownership
PhiX_220411_Training

Please confirm you want to transfer ownership to

Once the transfer has been completed it cannot be undone. The recipient of the transfer will be considered the legal owner of this project and its associated data. Please make sure you are authorized by your organization to make this transfer before proceeding.

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- [BaseSpace Sequence Hub Online Help](#)
- [The Many Benefits of BaseSpace Sequence Hub](#)
- [Basespace Sequencing Hub Apps Quick Guide](#)

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