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illumina

Microbial Whole-Genome Sequencing with the iSeq™ 100 Sequencing System

Fast and efficient sequencing that provides uniform coverage and genome assembly for microbial species.

Highlights

- Streamlined Workflow
Access a comprehensive workflow from DNA to data
- Optimized Library Prep
Obtain robust, consistent results over a wide range of DNA input, even at low DNA input amounts (1 ng)
- Comprehensive Coverage
Produce sequencing data with uniform coverage for viruses, bacteria, and other microbes

Introduction

Next-generation sequencing (NGS) has been established as an important tool in microbiology research for analysis of small genomes (<5 Mb), including bacteria, viruses, and other microbes. Microbial NGS, including whole-genome sequencing (WGS) and targeted resequencing, enables mapping and de novo assembly of novel organisms, comparing genomes of known organisms, and comparing genomes across samples.

The development of Nextera™ chemistry shortened and simplified library preparation by consolidating DNA fragmentation and adapter tagging steps into a single reaction (called fragmentation) and eliminating the need for library quantification before pooling and sequencing.¹ The Nextera DNA Flex Library Preparation Kit represents the next step in the evolution of Illumina library prep. In addition to speed and efficiency gains in the workflow, Nextera DNA Flex offers exceptional flexibility for sample input type and amount and robust, consistent preparation of sequencing-ready libraries.

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The latest innovation in NGS is poised to advance microbiology genomics research. The compact iSeq 100 System (Figure 1) combines complementary metal-oxide-semiconductor (CMOS) technology with proven Illumina sequencing-by-synthesis (SBS) chemistry to deliver high-accuracy data with fast time to results. The iSeq 100 System is part of a streamlined NGS workflow for targeted and whole-genome microbial sequencing (Figure 2).

	Prepare DNA	Prepare Library	Sequence	Analyze Data
Microbial Sequencing with the iSeq 100 System				
	~ 1-2 hours DNA Extraction Target Enrichment	~ 2.5 hours Nextera DNA Flex Library Prep Kit	~ 17.5 hours iSeq 100 Sequencing System	~ 1 hour IllumiSpace Core App Integrative Genomics Viewer (DNA Aligner)

Figure 2: Microbial Sequencing Workflow—Microbial sequencing on the iSeq 100 System is part of a streamlined, comprehensive NGS workflow that includes Nextera DNA Flex library preparation, sequencing, and data analysis.

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Prepare Library | Sequence | Analyze Data

illumina

Mitochondrial DNA Sequencing on the iSeq™ 100 Sequencing System

Uncover a comprehensive picture of the mitochondrial genome.

Highlights

- **Cost-Efficient Mitochondrial Sequencing**
Rapid, affordable analysis of the entire mitochondrial genome
- **Simple, DNA-to-Data Solution**
Supported workflow from library prep to data analysis
- **Exceptional Data Accuracy**
Highly accurate data for the detection of new variants and transcripts compared to qPCR or Sanger sequencing

Introduction

The mitochondrial genome comprises all of the DNA located within mitochondria, organelles important for energy conversion and other cellular functions. Mitochondrial DNA (mtDNA) is organized as a circular, double-stranded DNA molecule that is ~16,500 base pairs in length and codes for 37 genes. In humans (and most animals), mitochondria and mtDNA are maternally inherited during sexual reproduction.^{1,2}

Mitochondrial diseases are a heterogeneous group of disorders caused by mitochondrial dysfunction, resulting from inherited or somatic mutations in mtDNA or nuclear genes that code for mitochondrial components. Genetic defects in mtDNA include point mutations, deletions, and copy number variations (CNVs).² Studying mutations associated with mitochondrial diseases remains challenging due to phenotypic variability and genetic heterogeneity among individuals. Adding to this complexity is the heteroplasmic nature of mtDNA mutations, which is similar to genetic mosaicism, but on a subcellular, organelle level.

Figure 1: The iSeq 100 System—The iSeq 100 System harnesses the power of NGS in the most affordable, compact benchtop sequencing system in the Illumina portfolio.

Conventional analysis of the mitochondrial genome often starts with PCR-based screens for common point mutations and large deletions. If negative, Sanger sequencing of the entire mitochondrial genome is used to attempt to identify less common variants. However, conventional Sanger sequencing is not sensitive enough to detect low-level heteroplasmy and is not reliable in quantifying the level of heteroplasmy.³ Mitochondrial genome sequencing with next-generation sequencing (NGS) technology addresses these challenges, enabling comprehensive detection and analysis of mtDNA mutations.

Figure 2: Mitochondrial Sequencing Workflow—Mitochondrial sequencing on the iSeq 100 System is part of a streamlined, comprehensive NGS workflow that includes Nextera DNA Flex library preparation, sequencing, and data analysis.

Prepare DNA	Prepare Library	Sequence	Analyze Data
 ~ 1-2 hours DNA Extraction Target Amplification	 ~ 2.5 hours Nextera DNA Flex Library Prep Kit	 ~ 17.5 hours iSeq 100 Sequencing System	 ~ 1 hour BaseSpace One App mtDNA Variant Analyzer mtDNA Variant Processor

For Research Use Only. Not for use in diagnostic procedures.

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