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Description	List Price	1-plex Price/Sample	8-plex Price/Sample	12-plex Price/Sample
TruSeq® DNA Library Prep for Enrichment (24 Samples)	¥162,000	¥6,750	¥6,750	¥6,750
IDT for Illumina – TruSeq DNA UD Indexes (24 Indexes, 96 Samples)	¥121,000	¥1,260	¥1,260	¥1,260
IDT Lockdown Reagents (16 rxn)	¥8,000	¥500	¥63	¥42
IDT Exome Panel (16 rxn)	¥480,000	¥30,000	¥3,750	¥2,500
IDT Universal Blockers (16 rxn)	¥92,500	¥5,781	¥723	¥482
Total	¥863,500	¥44,292	¥12,546	¥11,034

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

Today, IDT is the leading brand in DNA and RNA synthesis

- 9 locations and >1,000 employees
- >100,000 active customers
- >50,000 oligonucleotides & >650 gBlocks® Fragments and Genes synthesized per day (over 55 million bases per month)
- >3,800 orders per day
- >400,000 website visits per month
- >94% of orders from website

IDT strategic advantages

Vertical integration

Custom synthesis platform
(hardware and software)

Synthesis inputs are
produced and controlled
internally

Quality

The molecular weight of
every oligo is verified by
mass spectrometry

Process traceability based
on custom lab information
management system (LIMS)

Flexibility and scalability

Transition from R&D to
commercialization with GMP
and B2B support

World-class custom
formulation and mixing

Support

Regionalized customer and
technical support

Raw materials

Creation of synthesis inputs

Production

QC

Kitting, packaging,
and formulation

Distribution


Superior oligos are made from high quality raw materials

Most reagents are acquired in bulk to reduce lot-to-lot variability

All raw materials undergo QC analysis for impurities known to negatively impact oligo synthesis

All synthesis inputs used at IDT are created or formulated in-house

Every synthesis material is QCed before use in production

IDT manufactures DNA and RNA phosphoramidites, as well as chemical modifications for oligo synthesis

IDT produces Nucleoside Controlled Pore Glass (CPG) that is used as a solid support for synthesis

Raw materials

Creation of synthesis inputs

Production

QC


Kitting, packaging,
and formulation


Distribution


IDT oligos have the highest coupling efficiency in the industry

$$\text{Percent full-length product} = (\text{Coupling efficiency})^{(\text{Total bases}-1)}$$




Electrospray ionization mass spectrometry (ESI-MS) is the superior method for confirming the identity of oligonucleotides

IDT last used MALDI-TOF Mass Spec in 2012

MALDI-TOF cannot identify a 120 nt oligo

MW = 36,279.3 Da

Voyager MALDI-TOF

Allowed error: 72.6 Da

Measured error: 843.1 Da

Error: 2.324 %

Bruker MALDI-TOF

Allowed error: 72.6 Da

Measured error: N/A

Error: N/A

Substitution	MW Difference
dA<-->dC	24
dA<-->dG	16
dA<-->dT	9
dC<-->dG	40
dC<-->dT	15
dG<-->dT	25

ESI

Allowed error: 7.3 Da

Measured error: 0.5 Da

Error: 0.0014 %

Raw materials

Creation of synthesis inputs

Production

QC

Kitting, packaging,
and formulation

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The screenshot shows the product page for the xGen Exome Research Panel. The page title is "xGen® Exome Research Panel" with the tagline "Deep, uniform coverage of the human exome with best-in-class capture panel". Below this, there is a detailed description of the panel and a list of key features:

- Obtain consistent results with individually synthesized target capture probes manufactured to GMP standards
- Achieve deep and uniform coverage even across difficult GC-rich regions
- Integrate with common sequencing platforms

The "Resources" dropdown menu is highlighted with a mouse cursor. The page also includes an "Ordering" section with a table of products and prices, and an "ORDER NOW" button.

Product	Pricing
16 rxn xGen® Exome Research Panel v1.0	¥480,000 JPY
96 rxn xGen® Exome Research Panel v1.0	¥2,400,000 JPY

The screenshot shows the "Resources" section of the product page. It lists various downloadable files and guides:

- Panel information**
 - xGen Exome Research Panel probes (bed) - Download (2,396.4 KB)
 - xGen Exome Research Panel targets (bed) - Download (1,095.4 KB)
 - xGen Exome Research Panel gene list (txt) - Download (1.44 KB)
- Files for BaseSpace analysis**
 - xGen Exome Research Panel probe manifest (txt) - Download (3,488.5 KB)
 - xGen Exome Research Panel targeted regions manifest (txt) - Download (6,900 KB)
- User guides and protocols**
 - Hybridization capture of Nextera DNA libraries using xGen Lockdown Probes and Panels (490 KB)
 - Hybridization capture with xGen Lockdown Probes and SeqCap reagents (96.2 KB)
 - Hybridization capture with xGen Lockdown Probes and Reagents (950 KB)
 - Hybridization capture of TruSeq libraries using xGen Lockdown Probes and Reagents (530 KB)
- Supplementing xGen Lockdown Panels**
- Application notes and case studies**
 - xGen Exome Research Panel—exceptional uniformity and complete coverage (755 KB)
- Posters**
 - Targeted vs. whole genomes—high quality oligonucleotide capture baits shift the target enrichment landscape (2,134 KB)
 - Development of a novel exome capture kit—ACMG 2015—Iowa Institute of Human Genetics (1,166 KB)
- Certificates of analysis (COAs)**
 - Find COAs by batch or lot number