

Sequencing

Libraries were run in parallel on the MiSeq and MiniSeq Systems. Paired-end 150-base pair sequencing (2 × 150 bp) was used on both systems. On the MiSeq System, 24 samples were pooled into each run using the MiSeq Reagent Kit v2 Micro (300 cycles). On the MiniSeq System, all 72 samples were pooled onto a single MiniSeq High Output Reagent Kit (300 cycles). The MiSeq System runs averaged 19 hours per run and the MiniSeq System runs averaged 24 hours per run (Table 3).

Table 3: MiSeq and MiniSeq Systems Specifications for the TruSight HLA v2 Sequencing Panel

MiSeq System (4-channel chemistry)				
Reagent Kit	Max. Reads	2 × 150 Output	2 × 150 Run Time	Max. Samples
v2 Nano	1M	300 Mb	17 hours	6
v2 Micro	4M	1.2 Gb	19 hours	24
V2 Standard	15M	4.5 Gb	24 hours	96
V3 Standard	25M	7.5 Gb	39 hours	144
MiniSeq System (2-channel chemistry)				
Reagent Kit	Max. Reads	2 × 150 Output	2 × 150 Run Time	Max. Samples
Mid Output	8M	2.4 Gb	17 hours	48
High Output	25M	7.5 Gb	24 hours	144

Data Analysis

The compressed fastq.gz files (2 per sample; Read 1 and Read 2) output from the MiSeq and MiniSeq Systems were imported directly into TruSight HLA Assign 2.0 (v2.0.0.920) Software. IMGT/HLA version 3.23 and CWD version 2.0.0 software were used for analysis and typing of the results. Assign 2.0 Software performed the alignment of individual reads to a consensus reference sequence (2 consensus reference sequences for HLA-DQB1 and 4 consensus reference sequences for HLA-DRB1).

The software then phased all heterozygous positions using base positions under individual reads and in paired reads. Often, multiple pairs of reads need to be layered to align phase directly. This method allows for phasing of heterozygous positions as far apart as the longest fragment sequenced (usually 1000–1300 bases). These phased alignments generate a consensus sequence for each locus comprising a mean depth of coverage per locus of more than 100x. Following sample multiplexing recommendations, mean depth of coverage is usually over 200x per locus. These consensus sequences were compared to the IMGT/HLA database housed within the Assign 2.0 Software to generate an HLA type for each locus.

Results: MiSeq System

Evaluation of Concordance with Reference Typing

Sequencing of the 72 IHWG samples on the MiSeq System resulted in consensus sequences for 1294 alleles, 959 of which had references for comparison. 907 (94.58%) alleles were concordant with the reference typing. The 52 typing results that did not match the reference sequence fell into 3 categories: 42 (4.38% of total referenced alleles) had incorrect or incomplete reference typings that prompted updates to the reference typings (Table 4), 5 (0.52% of total referenced alleles) had novel exonic variants, and 5 (0.52% of total referenced alleles) were not assigned a type and required manual editing.

All alleles with reference typings for which the results were concordant, alleles with reference typings that appeared to be incorrect, and novel alleles were used to calculate accuracy for the TruSight HLA v2 Sequencing Panel (Figure 2). 5 alleles were considered inaccurate and required a total of 9 base edits for concordance. These 5 alleles were not typed rather than being assigned an incorrect typing.




Figure 2: Accuracy of the TruSight HLA v2 Sequencing Panel on the MiSeq System—The accuracy is calculated using 959 alleles with references. A TruSight HLA v2 result was considered accurate when a result was fully concordant (907 alleles, 94.58%), when the result matched an update to the reference (42 alleles, 4.38%), and when the allele contained a novel variant (5 alleles, 0.52%). This resulted in a total of 954 alleles called accurately (99.48%).

Resolving HLA Typing Ambiguities

There were 41 ambiguities found in the 1292 sequenced alleles (Figure 5). There are 3 types of ambiguities present in these data (Table 5).

- Amplicon Ambiguities:** These ambiguities were resolved outside the amplicon and include DPB1*13:01:01/107:01 (13:01:01G group), DRB1*12:01:01/12:10 (12:01:01G group), and DRB1*08:01:01/03. All 3 of these ambiguities were resolved in exon 1 of their respective loci and exons 1 of DPB1 and DRB1 are not covered by TruSight HLA v2. Amplicon ambiguities accounted for 8 of the 41 ambiguities (19.51%). These amplicon ambiguities are highlighted in blue.
- Conditional Ambiguities:** These ambiguities appear based on the combination of alleles and usually occur because the homology between the alleles leads to large gaps in heterozygous base positions that are used to phase. With gaps of 1Kb or more, phase will be lost resulting in a phase ambiguity. This most commonly occurs in DPB1 in which intron 2 is over 4500 bases in length with certain combinations having het gaps of more than 4Kb (eg DPB1*04:01:01 paired with DPB1*04:02:01). These conditional ambiguities will appear every time these alleles are combined with one another, but not every time the allele is present in a different combination. Conditional ambiguities occurred in 28 of the 41 (68.29%) ambiguities and account for all the DQB1 ambiguities and all the DPB1 ambiguities except for the 13:01:01/107:01. Conditional ambiguities are highlighted in green.
- Phase Ambiguities:** These ambiguities rarely occur and are specific to a sequencing run. In other words, the same sample and allele combination will usually sequence unambiguously, but on occasion may not have the appropriate read structure to resolve phase. There are 3 loci (5 of 41 alleles, 12.20%) in which this occurred in these runs and these are highlighted in pink.




Figure 5: Ambiguity Rates of the TruSight HLA v2 Sequencing Panel on the MiSeq System—The ambiguity rate was calculated using the 1292 alleles that were sequenced and these data were calculated independently of the references. There were 41 total alleles with ambiguities.

Table 5: TruSight HLA v2 Sequencing Panel Ambiguities on the MiSeq System

Sample ID	Locus	Ambiguity
IHW01018	DRB1	08:01:01/03
IHW01040	DPB1	03:01:01/124:01
IHW01040	DRB1	04:01:01/350:01
IHW01093	DPB1	08:01:01/03
IHW01137	DPB1	04:01:01/105:01
IHW01141	DPB1	04:02:01/126:01
IHW01141	DRB1	04:01:01/105:01
IHW01174	DPB1	04:02:01/126:01
IHW09044	DPB1	03:01:01/124:01
IHW09045	DPB1	04:01:01/350:01
IHW09045	DRB1	02:01:02/105:01
IHW09056	DPB1	04:02:01/416:01
IHW09056	DRB1	02:01:02/105:01
IHW09013	DPB1	04:02:01/416:01
IHW09114	DRB1	02:01:02/106:01
IHW09267	DQB1	19:01/414:01
IHW09273	DRB1	13:01:01/107:01
IHW09366	DQB1	06:03:01/06:39
IHW09378	DQB1	06:04:01/06:41
IHW09388	DPB1	12:01:01/12:10
IHW09398	DRB1	06:02:01/06:39
IHW09417	DPB1	06:04:01/06:41
IHW09501	DPB1	06:02:01/06:84
IHW09501	DPB1	06:03:01/06:41
IHW09502	A	06:03:01/06:41
IHW09502	DPB1	06:09:01/06:88
IHW09502	DPB1	06:02:01/06:39
IHW09502	DPB1	06:04:01/06:84
IHW09502	DPB1	06:03:01/06:41
IHW09502	DPB1	06:09:01/06:88
IHW09502	DPB1	04:01:01/105:01
IHW09502	DPB1	04:02:01/126:01
IHW09502	DPB1	04:01:01/105:01
IHW09502	DPB1	04:02:01/126:01
IHW09502	DPB1	04:04:01/133:01
IHW09502	DPB1	13:01:01/107:01/350:01
IHW09502	A	02:05:01/02:22:01
IHW09502	DPB1	02:01:01/02:14
IHW09502	DPB1	03:01:01/124:01
IHW09502	DPB1	04:01:01/350:01

